

Surface disinfectants



Surface disinfectant Meliseptol®

B. Braun.

Ready-to-use alcohol-based fast-acting disinfectant without additional active substances. Also suitable for hard-to-reach surfaces and objects, for wiping or spraying.

Spectrum of activity: Bactericidal (incl. MRSA, TbB and mycobacteria), fungicidal, limited virucidal efficacy (incl. HBV, HIV, HCV) and effective against rotavirus, adenovirus and noroviruses.

Listing: VAH, IHO and DGHM.

Short reaction time of 1 min.

Use disinfectants with care. Always read the labelling and product information before use.

Technical specifications:

Art. No.	EAA6.1	EAA6.2	EAA6.3
Delivery form	Ready-to-use liquid		
Ingredients	Alcohol		
Bactericide	Yes		
Partially virucidal	Yes		
Fungicide	Yes		
CE marking	Yes		
Listing	VAH, IHO, DGHM		

Warning H226

Medical device class IIa MDD



Surface disinfectant Meliseptol® rapid

B. Braun.

Aldehyde- and alkylamine-free alcohol-based fast-acting disinfectant. Also suitable for hard-to-reach surfaces and objects, for wiping or spraying.

Spectrum of activity: Bactericidal (incl. MRSA and TbB), fungicidal, limited virucidal efficacy (incl. HBV, HIV, HCV) and virucidal (incl. noroviruses).

Listing: VAH, IHO and DGHM. Ready to use. Very fast acting, reaction time of 1 minute; 10 minutes for virucidal efficacy.

Use disinfectants with care. Always read the labelling and product information before use.

Technical specifications:

Art. No.	LK31.1	LK31.2	LK31.3
Delivery form	Ready-to-use liquid		
Ingredients	Alcohol		
Bactericide	Yes		
Partially virucidal	Yes		
Fungicide	Yes		
CE marking	Yes		
Listing	VAH, IHO, DGHM		

Danger H226-H318-H336

Medical device class IIa MDD