



**WUHU RUIJIN MEDICAL INSTRUMENT
& DEVICE CO.,LTD**

Address/Tel/Fax: No.33, EastWanchun Road, Wuhu Economic
& Technology Development Zone, Wuhu Anhui China;
Tel: 86-553-2672511 Fax: 86-553-2672510
Website: www.whruijin.com



***Medical electric drill RJ-PD (MZ-100)
Micro medical electric spine drill***

***For application in spinal surgery and plastic
surgery,neuro.***

» Feature

- 01 *Refined design*
- 02 *Imported motor, autoclavable*
- 03 *Low noise, high speed*
- 04 *Long service life.*



Standard configuration »


Sterilization Case×1
Power line×1
Foot pedal×1
Motor handpiece×1
Straight handpiece×1
Angled handpiece×1
Cutting burr: sizeΦ1-Φ5
Diamond burr: sizeΦ1-Φ5
Power console×1



» Product performance and application

Straight handle: different length, applied in neursugery, spine surgery and microsurgery.

Angled handle: different length, applied in neursugery, spine surgery and microsurgery.



» Main Technical Data

Motor handpiece	Output power	≥50W
	Temperature rise	≤25°C
	Input voltage	AC 220V/50HZ
Non-loading speed	0-40000 RPM	
Non-loading noise	≤50dB (A)	
Axial run out	≤0.03mm	
Maximum torque	0.5 N*m	

Notes:


1. Power console cannot be placed in high temperature and high pressure environment.

2. After the instrument is cleaned, it can be used only after it has been sterilized and sterilized according to the hospital's specifications. The host, handle can be gauze soaked in water or wipe clean, then put them away.

3. This product is a dedicated medical device. The user of this product must be a medical staff with appropriate skills or appropriate training.



Aluminum case packing



Sterilization case for power tool