



Authorized Partner Thermal Analysis and Chromatography

Complete Laboratory Solutions

Planning | Designing | Laboratory Equipment Delivery | Installation | Personnel Training | Consumables Delivery | Maintenance | After Sales Service

CLIENT: UNIVERSITATEA DE STAT

| NR CRT | COD PRODUS | DESCRIERE | UM | QTY |
|-----------|---------------|---|-----|-----|
| | EUROMEX | | | |
| 1 | NZ.1703-S | Trinocular stereo zoom microscope NexiusZoom EVO, 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with rack and pinion stand. Incident and transmitted 3 W LED illuminations | buc | 1 |
| | | Pair wide field eyepieces HWF10x/23mm without eyecups | | |
| | | Trinocular head with 45° inclined tubes | | |
| | | Diopter adjustment ± 5° on two tubes. Interpupillary distance adjustable between 54 and | | |
| | | 75 mm | | |
| | | 0.65x to 5.5x zoom objective, magnification from 6.5x to 55x with click-stops | | |
| | | Field of view from 35 mm to 4.2 mm. Working distance 110 mm | | |
| | | Ergonomically rack and pinion stand. Maximum object height 58 mm | | |
| | | Supplied with opaque stage plate, white/black stage plate and 2 object clamps Incident and transmitted 3W LED illumination with internal power supply | | |
| | | Both illumination intensities can be adjusted separately | | |
| | | Paint with Antibacterial Protection Layer | | |
| | | Supplied with dust cover in Styrofoam case | | |
| 2 | VC.3024-HDS | HD-Mini camera, color High definition 2 MP high speed camera with HD screen, 1280 x | buc | 1 |
| | | 1080p, stand-alone usage with standard 4GB SD card, HDMI output. With built- in mousedriven capture software | | |
| | | Sensor: Sony 1/2.8" HD 1.2 MP CMOS sensor, HDTV 1080p (50/60 Hz) | | |
| | | Pixels: 1920 x 1080 (viewing and capture mode) | | |
| | | Pixel size: 2.9 x 2.9 µm | | |
| | | 1080p video recording mode | | |
| | | HDMI 60 fps @ 1980 x 1080 pixels HDMI | | |
| | | Manual and automaric exposure time 0.01-1000 ms | | |
| | | Scan: Progressive scan | | |
| | | Noise 3D Noise Reduction function. Gain, Saturation Sharpness, Gamma | | |
| | | White balance: Automatic and manual, ROI | | |
| | | Stand alone usage with mouse driven built-in embedded software | | |
| | | Video playback and image browser. Measurement function | | |





Authorized Partner
Thermal Analysis and Chromatography

Complete Laboratory Solutions

Planning | Designing | Laboratory Equipment Delivery | Installation | Personnel Training | Consumables Delivery | Maintenance | After Sales Service

| 3 | NZ.9850 | C-mount adapter with 0.5x lens for 1/2 inch cameras buc | 1 |
|---|---------|---|---|
| | | Dimensions: 50 x 50 x 61 mm | |
| | | Operation: -10 to 40° C, 30-85 % humidity | |
| | | Power supply: external main adapter 100-240 Vac to12Vdc/1A | |
| | | adapter rings for stereo microscopes and computer mouse | |
| | | Supplied with 0.37x objective, HD screen, HDMI cable, 23.2 to 30 mm and 23.2 to 30.5 mm | |
| | | Mounting: C-mount interface | |
| | | Storage: 4 Gb SD memory card | |
| | | Interface: Digital HDMI 1.4 | |
| | | Sensitivity 1.3 V/lux/s , dynamic range 85.3 dB and S/N ratio 44 dB | |
| | | Image format 2MP JPEG, up to 10x digital zoom function | |
| | | Video format 1920 x 1080p, H264 encoded MP4 file (50-60 fps) | |