

Connecting Tubings



INTRODUCTION

The product is intended to connect pipelines and devices which provides pathway of the liquid-injection or/and invasive blood pressure monitoring in interventional operation.

PRODUCT FEATURES

- High pressure resistance: 500psi/1200psi
- Transparent material makes the flow visible
- Customized sizes are available

ORDERING INFORMATION

500 PSI female to male luer lock

| External OD (mm) | Length (cm) | | | | | | |
|------------------|--------------|--------------|--------------|--------------|--------------|---------------|---------------|
| | 10 | 20 | 30 | 60 | 90 | 120 | 150 |
| 3.00 | BM-CT3.00-10 | BM-CT3.00-20 | BM-CT3.00-30 | BM-CT3.00-60 | BM-CT3.00-90 | BM-CT3.00-120 | BM-CT3.00-150 |

BRAIDED TUBING

1200 PSI female luer lock to male rotator

| External OD (mm) | Length (cm) | | | | | | |
|------------------|---------------|---------------|---------------|----------------|----------------|----------------|----------------|
| | 20 | 60 | 90 | 120 | 150 | 180 | 200 |
| 3.65 | BM-CT3.65-20B | BM-CT3.65-60B | BM-CT3.65-90B | BM-CT3.65-120B | BM-CT3.65-150B | BM-CT3.65-180B | BM-CT3.65-200B |

*CE not available for length of 180/200cm

NON-BRAIDED TUBING

1200 PSI female to male luer lock

| External OD (mm) | Length (cm) | | | | | | |
|------------------|---------------|---------------|---------------|----------------|----------------|----------------|----------------|
| | 20 | 60 | 90 | 120 | 150 | 180 | 200 |
| 3.20 | BM-CT3.20-20A | BM-CT3.20-60A | BM-CT3.20-90A | BM-CT3.20-120A | BM-CT3.20-150A | BM-CT3.20-180A | BM-CT3.20-200A |

*CE not available for length of 180/200cm

Angiographic Syringes and Accessories



PRODUCT FEATURES

- Excellent sealing of plunger and unique push button facilitates handling
- Clear scale and accurate control
- Broad variety of angiographic syringes and accessories

SET CONTENTS

- Single Syringe Set for CT/MRI
 - 1* Syringe
 - 1* Coiled Tubing
 - 1* Quick Fill Tube (J Tube)
- Single Syringe Set for DSA
 - 1* Syringe
 - 1* Quick Fill Tube (J Tube)
- Dual Syringe Set
 - 2* Syringe
 - 1* Y tube with Valve
 - 2* Quick Fill Tube (J Tube)
 - 2* Spike

MEDRAD(A)

| Injector | Compatible Models | Syringe Set | Volume (ml) | Catalogue No. |
|----------|--|-------------|-------------|----------------------|
| CT | MCT & PLUS; OP-100; VCT610; ECT710, Vistron, EnVison | Single | 200 | BM-AS-CT-200-A01 |
| | | Dual | 200/200 | BM-AS-CT-200/200-A02 |
| | STELLANT | Single | 200 | BM-AS-CT-200-A03 |
| | | Dual | 200/200 | BM-AS-CT-200/200-A04 |
| | SALIENT | Single | 190 | BM-AS-CT-190-A05 |
| | | Dual | 190/190 | BM-AS-CT-190/190-A06 |
| | Multi-Patient syringe system Stellant (12 hours) | Dual | 200/200 | BM-AS-CT-200/200-A08 |

| | | | | |
|-----|-----------------------------|--------|---------|-----------------------|
| DSA | Mark V & Mark V Provis | Single | 150 | BM-AS-DSA-150-A01 |
| | | Single | 200 | BM-AS-DSA-200-A03 |
| | | Single | 60 | BM-AS-DSA-60-A05 |
| | Mark III & Mark IV & Mark V | Single | 130 | BM-AS-DSA-130-A07 |
| | Mark 7 | Single | 150 | BM-AS-DSA-150-A09 |
| MRI | Spectris & Solaris | Dual | 65/65 | BM-AS-MRI-65/65-A02 |
| | | Dual | 115/115 | BM-AS-MRI-115/115-A04 |
| | | Dual | 65/115 | BM-AS-MRI-65/115-A06 |

LIBEL-FLARSHEIM(B)

| Injector | Compatible Models | Syringe Set | Volume (ml) | Catalogue No. |
|----------|------------------------------------|-------------|-------------|----------------------|
| CT | CT 9000 & CT 9000 ADV, Optivantage | Single | 200 | BM-AS-CT-200-B01 |
| | | Dual | 200/200 | BM-AS-CT-200/200-B02 |
| DSA | Angiomat 6000 | Single | 150 | BM-AS-DSA-150-B01 |
| | Angiomat Illumena | Single | 150 | BM-AS-DSA-150-B03 |
| MRI | Optistar LE, Elite | Single | 60 | BM-AS-MRI-60-B01 |
| | | Dual | 60/60 | BM-AS-MRI-60/60-B02 |

EZEM(C)

| Injector | Compatible Models | Syringe Set | Volume (ml) | Catalogue No. |
|----------|--------------------------------|-------------|-------------|----------------------|
| CT | Empower, EmpowerCT, EmpowerCTA | Single | 200 | BM-AS-CT-200-C01 |
| | | Dual | 200/200 | BM-AS-CT-200/200-C02 |

NEMOTO KYORINDO(D)

| Injector | Compatible Models | Syringe Set | Volume (ml) | Catalogue No. |
|----------|---------------------------------|-------------|-------------|----------------------|
| CT | A-25 | Dual | 100/200 | BM-AS-CT-100/200-D02 |
| | A-25, A-60 | Single | 100 | BM-AS-CT-100-D01 |
| | | Single | 200 | BM-AS-CT-200-D03 |
| | | Dual | 100/100 | BM-AS-CT-100/100-D04 |
| | | Dual | 200/200 | BM-AS-CT-200/200-D06 |
| | | Dual | 60/200 | BM-AS-CT-60/200-D08 |
| | DSA | Nemoto 120S | Single | 120 |
| MRI | Nemoto Dual Shot, Sonic Shot 50 | Single | 60 | BM-AS-MRI-60-D01 |
| | | Dual | 60/60 | BM-AS-MRI-60/60-D02 |

MEDTRON(E)

| Injector | Compatible Models | Syringe Set | Volume (ml) | Catalogue No. |
|----------|-------------------|-------------|-------------|-----------------------|
| CT | ACCUTRON | Single | 200 | BM-AS-CT-200-E01 |
| | | Dual | 200/200 | BM-AS-CT-200/200-E02 |
| MRI | ACCUTRON | Single | 65 | BM-AS-MRI-65-E01 |
| | | Dual | 65/65 | BM-AS-MRI-65/65-E02 |
| | | Dual | 65/200 | BM-AS-MRI-65/200-E-04 |



| Description | Length(cm) | Catalogue No. |
|---|------------|-----------------|
| · Coiled tube | 50 | BM-CT3.2-50SC0 |
| · Coiled tube | 100 | BM-CT3.2-100SC0 |
| · Coiled tube | 150 | BM-CT3.2-150SC0 |
| · Straight tube | 20 | BM-CT3.2-20SS0 |
| · Straight tube | 150 | BM-CT3.2-150SS0 |
| *Straight tube with one check valve | 12 | BM-CT3.2-12SS1 |
| *Straight tube with two male check valves | 150 | BM-CT3.2-150SS2 |
| *Straight tube with two male check valves | 250 | BM-CT3.2-250SS2 |
| *Straight tube with one needless adapter and spike | 25 | BM-CT3.2-25SN0 |
| *Y tube with two check valves | 15 | BM-CT3.2-15SY2 |
| *Y tube with two check valves | 150 | BM-CT3.2-150SY2 |
| *Y tube with two check valves | 250 | BM-CT3.2-250SY2 |
| *Y tube with one check valve | 150 | BM-CT3.2-150SY1 |
| *Y tube with one check valve | 250 | BM-CT3.2-250SY2 |
| *Double infusion line with two check valves | 150 | BM-CT3.2-150SD2 |
| *T tube with one coiled tube 150cm and 100cm spike tube | 150 | BM-CT3.2-150SQ0 |

· With CE mark, *Without CE Mark