



Electrophysiology Catheters
PRODUCT CATALOG


Index

Fixed Diagnostic Catheter Curve Reference Chart.....	5
Fixed Curve Diagnostic Catheters	6–7
Fixed Curve Diagnostic Catheters w/Guidewire Tip	8
Steerable Catheter Curve Reference Chart.....	9–11
Steerable Catheters	12–13
PV Catheters.....	13
Diagnostic Catheter Extension Cables.....	14
Ablation Catheters Curve Reference Chart.....	15–16
Ablation Catheters.....	17
Bi-Directional and Irrigated Ablation Catheters.....	18
Ablation Catheter Extension Cables	19

Map-iT™ Electrophysiology Catheters

CURVE REFERENCE CHART

Diagnostic Fixed Catheters



MAP-IT™ FIXED CURVE DIAGNOSTIC CATHETERS


Map-iT™ Fixed Quadripolar Catheter



Reference Number	French Size	Number of Electrodes	Inter-electrode spacing (±0.5 mm)	Curve Type	Length (±0.5 cm)	Ring Electrode width (±0.05mm)	Tip electrode width (±0.2mm)	Cable Type
901470	4F	4	10-10-10.	MPA	120	1	1.6	901950
901474	4F	4	10-10-10.	HIS	120	1	1.6	901950
901475	4F	4	10-10-10.	CS	120	1	1.6	901950
901477	4F	4	5-5-5.	MPA	120	1	1.6	901950
901478	4F	4	5-5-5.	MPB	120	1	1.6	901950
901479	4F	4	5-5-5.	RVA	120	1	1.6	901950
901480	4F	4	5-5-5.	RVB	120	1	1.6	901950
901482	4F	4	2-5-2.	MPA	120	1	1.6	901950
901484	4F	4	2-5-2.	HIS	120	1	1.6	901950
901485	4F	4	2-2-2.	MPA	120	1	1.6	901950
901486	4F	4	2-2-2.	MPB	120	1	1.6	901950
901423	5F	4	10-10-10.	CS	120	2	2	901950
901430	5F	4	10-10-10.	MPB	120	2	2	901950
901431	5F	4	10-10-10.	MPA	120	2	2	901950
901435	5F	4	5-5-5.	MPB	120	2	2	901950
901436	5F	4	5-5-5.	MPA	120	2	2	901950
901440	5F	4	2-5-2.	MPB	120	2	2	901950
901441	5F	4	2-5-2.	MPA	120	2	2	901950
901442	5F	4	2-5-2.	RVA	120	2	2	901950
901443	5F	4	2-5-2.	RVB	120	2	2	901950
901444	5F	4	2-5-2.	HIS	120	2	2	901950
901450	6F	4	5-5-5.	MPB	120	2	2	901950
901451	6F	4	5-5-5.	MPA	120	2	2	901950
901452	6F	4	5-5-5.	RVA	120	2	2	901950
901453	6F	4	5-5-5.	RVB	120	2	2	901950
901454	6F	4	5-5-5.	HIS	120	2	2	901950
901455	6F	4	2-5-2.	MPB	120	2	2	901950
901456	6F	4	2-5-2.	MPA	120	2	2	901950
901457	6F	4	2-5-2.	RVA	120	2	2	901950
901458	6F	4	2-5-2.	RVB	120	2	2	901950
901459	6F	4	2-5-2.	HIS	120	2	2	901950
901460	6F	4	2-2-2.	MPB	120	2	2	901950
901461	6F	4	2-2-2.	MPA	120	2	2	901950


MAP-IT™ FIXED CURVE DIAGNOSTIC CATHETERS

Map-iT™ Fixed Hexapolar Catheter



Reference Number	French Size	Number of Electrodes	Inter-electrode spacing (±0.5 mm)	Curve Type	Length (±0.5 cm)	Ring Electrode width (±0.05mm)	Tip electrode width (±0.2mm)	Cable Type
901130	6F	6	5-5-5-5-5.	MPA	120	2	2	901957
901131	6F	6	5-5-5-5-5.	MPB	120	2	2	901957
901132	6F	6	5-5-5-5-5.	HIS	120	2	2	901957
901133	6F	6	2-5-2-5-2.	MPA	120	2	2	901957
901134	6F	6	2-5-2-5-2.	MPB	120	2	2	901957
901135	6F	6	2-5-2-5-2.	HIS	120	2	2	901957
901139	6F	6	5-5-5-175-175.	MPA	120	2	2	901957
901140	6F	6	5-5-5-175-175.	MPB	120	2	2	901957

Map-iT™ Fixed Decapolar Catheter



901518	5F	10	2-5-2.	MPA	120	1	2	901957
901519	5F	10	2-5-2.	MPB	120	1	2	901957
901520	5F	10	2-5-2.	HIS	120	1	2	901957
901521	5F	10	2-2-2.	MPA	120	1	2	901957
901522	5F	10	2-2-2.	MPB	120	1	2	901957
901523	5F	10	2-2-2.	HIS	120	1	2	901957
901530	5F	10	2-8-2.	CS	120	1	2	901957
901579	5F	10	4-4-4-23-4-4-4-4-4	MPA	120	1	2	901957
901580	5F	10	4-4-4-29-4-4-4-4-4	MPA	120	1	2	901957
901532	6F	10	2-8-2.	CS	120	1	2	901957
901534	6F	10	2-5-2.	MPA	120	1	2	901957
901535	6F	10	2-5-2.	MPB	120	1	2	901957
901536	6F	10	2-5-2.	HIS	120	1	2	901957
901537	6F	10	2-2-2.	MPA	120	1	2	901957
901538	6F	10	2-2-2.	MPB	120	1	2	901957
901539	6F	10	2-2-2.	HIS	120	1	2	901957
901540	6F	10	5-5-5.	CS	120	1	2	901957
901823	6F	10	4-4-4-23-4-4-4-4-4	MPA	120	1	2	901957
901824	6F	10	4-4-4-29-4-4-4-4-4	MPA	120	1	2	901957

MAP-IT™ FIXED CURVE DIAGNOSTIC CATHETERS


Map-iT™ 4/3.3F Fixed Diagnostic Catheters with Guidewire Tip

Reference Number	French Size	Number of Electrodes	Inter-electrode spacing (±0.5 mm)	Curve Type	Length (±0.5 cm)	Ring Electrode width (±0.05mm)	Tip electrode width (±0.2mm)	Cable Type
903215	4/3.3F	10	2-8-2.	S	120	1	Guidewire	901957
903217	4/3.3F	10	5-5-5.	S	120	1	Guidewire	901957
903219	4/3.3F	10	2-5-2.	S	120	1	Guidewire	901957
903221	4/3.3F	10	2-10-2.	S	120	1	Guidewire	901957
903226	4/3.3F	10	2-2-2.	S	120	1	Guidewire	901957
903256	4/3.3F	20	2-8-2.	S	120	1	Guidewire	901957 x2
903257	4/3.3F	20	2-10-2.	S	120	1	Guidewire	901957 x2
903269	4/3.3F	20	2-2-2.	S	120	1	Guidewire	901957 x2
903271	4/3.3F	20	5-5-5.	S	120	1	Guidewire	901957 x2
903273	4/3.3F	20	2-5-2.	S	120	1	Guidewire	901957 x2

Map-iT™ 4/3.3F - Decapolar Fixed Curve Catheters with Guidewire Tip







Map-iT™ 4/3.3F - Duo-decapolar Fixed Curve Catheters with Guidewire Tip







Map-iT™ Steerable Catheters

CURVE REFERENCE CHART







Map-iT™ Steerable Catheters

CURVE REFERENCE CHART



Map-iT™ Steerable Catheters

CURVE REFERENCE CHART



MAP-IT™ STEERABLE DIAGNOSTIC CATHETERS

Reference Number	French Size	Number of Electrodes	Inter-electrode spacing (±0.5 mm)	Curve Type	Length (±0.5 cm)	Ring Electrode width (±0.05mm)	Tip electrode width (±0.2mm)	Cable Type
901601	5F	4	5-5-5.	MPD-D	110	1	2	901955
901604	5F	4	5-5-5.	MPD-S	110	1	2	901955
901700	5F	4	5-5-5.	MPD-DL	110	1	2	901955
901703	5F	4	5-5-5.	MPD-SL	110	1	2	901955
903703	5F	4	5-5-5.	MPD-CS-S	110	1	2	901955
903707	5F	4	5-5-5.	MPD-CS-M	110	1	2	901955
903711	5F	4	5-5-5.	MPD-CS-L	110	1	2	901955
901611	6F	4	5-5-5.	MPD-D	110	1	2	901955
901614	6F	4	5-5-5.	MPD-S	110	1	2	901955
901709	6F	4	5-5-5.	MPD-DL	110	1	2	901955
901713	6F	4	5-5-5.	MPD-SL	110	1	2	901955
903803	6F	4	5-5-5.	MPD-CS-S	110	1	2	901955
903807	6F	4	5-5-5.	MPD-CS-M	110	1	2	901955
903810	6F	4	2-2-2.	MPD-CS-L	110	1	2	901955
901651	5F	10	2-5-2.	MPD-D	110	1	2	901956
901654	5F	10	2-5-2.	MPD-G	110	1	2	901956
901656	5F	10	2-5-2.	MPD-S	110	1	2	901956
901705	5F	10	2-5-2.	MPD-DL	110	1	2	901956
901707	5F	10	2-5-2.	MPD-SL	110	1	2	901956
903730	5F	10	2-5-2.	MPD-CS-S	110	1	2	901956
903734	5F	10	2-5-2.	MPD-CS-M	110	1	2	901956
903738	5F	10	2-5-2.	MPD-CS-L	110	1	2	901956
901661	6F	10	2-5-2.	MPD-D	110	1	2	901956
901664	6F	10	2-5-2.	MPD-G	110	1	2	901956
901666	6F	10	2-5-2.	MPD-S	110	1	2	901956
901715	6F	10	2-5-2.	MPD-DL	110	1	2	901956
901717	6F	10	2-5-2.	MPD-SL	110	1	2	901956
901718	6F	10	2-8-2.	MPD-SL	110	1	2	901956
903830	6F	10	2-5-2.	MPD-CS-S	110	1	2	901956
903835	6F	10	2-5-2.	MPD-CS-M	110	1	2	901956
903839	6F	10	2-5-2.	MPD-CS-L	110	1	2	901956
903901	7F	10	2-5-2.	MPD-D	110	1	2	901956
903904	7F	10	2-5-2.	MPD-DL	110	1	2	901956
903907	7F	10	2-5-2.	MPD-G	110	1	2	901956
903909	7F	10	2-5-2.	MPD-S	110	1	2	901956
903911	7F	10	2-5-2.	MPD-SL	110	1	2	901956

MAP-IT™ STEERABLE DIAGNOSTIC CATHETERS

Reference Number	French Size	Number of Electrodes	Inter-electrode spacing (±0.5 mm)	Curve Type	Length (±0.5 cm)	Ring Electrode width (±0.05mm)	Tip electrode width (±0.2mm)	Cable Type
903913	7F	10	2-5-2.	MPD-CS-S	110	1	2	901956
903917	7F	10	2-5-2.	MPD-CS-M	110	1	2	901956
903922	7F	10	2-5-2.	MPD-CS-L	110	1	2	901956
901670	6F	20	2-10-2.	MPD-H	100	1	2	901957 x2
901671	6F	20	2-5-2.	MPD-H	100	1	2	901957 x2
901673	6F	20	2-8-2.	MPD-H	100	1	2	901957 x2
901674	6F	20	2-8-2-...60-...2-8-2.	MPD-H	100	1	2	901957 x2
901675	6F	20	2-8-2.	MPD-L	100	1	2	901957 x2
901680	7F	20	2-10-2.	MPD-H	100	1	2	901957 x2
901681	7F	20	2-5-2.	MPD-H	100	1	2	901957 x2
901683	7F	20	2-8-2.	MPD-H	100	1	2	901957 x2
901684	7F	20	2-8-2-...60-...2-8-2.	MPD-H	100	1	2	901957 x2
901685	7F	20	2-8-2.	MPD-L	100	1	2	901957 x2

MAP-IT™ STEERABLE PV CATHETERS (FIXED LOOP)

Steerable 7F Decapolar PV Catheters

Reference Number	French Size	Number of Electrodes	Inter-electrode spacing (±0.5 mm)	Curve Type	Length (±0.5 cm)	Ring Electrode width (±0.05mm)	Tip electrode width (±0.2mm)	Cable Type
901641	7F	10	3-3-3.	PV12	110	1	-	901956
901642	7F	10	4-4-4.	PV15	110	1	-	901956
901643	7F	10	5-5-5.	PV18	110	1	-	901956
901644	7F	10	6-6-6.	PV21	110	1	-	901956
901645	7F	10	7-7-7.	PV24	110	1	-	901956
903950	7F	20	1.25-1.25-1.25.	PV12	110	0.6	-	901957 x2
903951	7F	20	1.5-1.5-1.5.	PV15	110	0.6	-	901957 x2
903952	7F	20	2-2-2.	PV18	110	0.6	-	901957 x2
903953	7F	20	2.5-2.5-2.5.	PV21	110	0.6	-	901957 x2
903954	7F	20	3-3-3.	PV24	110	0.6	-	901957 x2


MAP-IT™ DIAGNOSTIC CATHETER EXTENSION CABLES



Reference Number	Cable Name	Connects with	Usable Length
901950	Shrouded Quad EXT CBL, black, 150cm	Fixed curve Quadripolar catheter	150cm
901955	Shrouded Quad EXT CBL, Steerable 4 pole, black, 150cm,	Steerable curve Quadripolar catheter	150cm
901956	Shrouded Deca EXT CBL, Steerable 10 pole, black, 150cm	Steerable curve Decapolar catheter	150cm
901957	Shrouded Deca EXT CBL, Fixed/20 pole, black, 210cm	Fixed Hexapolar, Octapolar, Decapolar and Duo-Decapolar and Steerable Duo-Decapolar	210cm




Map-iT™ Ablation Catheters

CURVE REFERENCE CHART



Map-iT™ Ablation Catheters

CURVE REFERENCE CHART



MAP-IT™ STANDARD ABLATION CATHETERS*

Reference Number	French Size	Number of Electrodes	Inter-electrode spacing (±0.5 mm)	Curve Type	Length (±0.5 cm)	Ring Electrode width (±0.05mm)	Tip electrode width (±0.2mm)	Temperature Type
902120	7F	4	2-5-2.	B	85	1	4	TC
902121	7F	4	2-5-2.	C	85	1	4	TC
902122	7F	4	2-5-2.	D	85	1	4	TC
902123	7F	4	2-5-2.	D	85	1	8	TC
902124	7F	4	2-5-2.	E	85	1	4	TC
902125	7F	4	2-5-2.	E	85	1	8	TC
902126	7F	4	2-5-2.	F	85	1	4	TC
902127	7F	4	2-5-2.	F	85	1	8	TC
902128	7F	4	2-5-2.	A	110	1	4	TC
902129	7F	4	2-5-2.	B	110	1	4	TC
902130	7F	4	2-5-2.	B	110	1	8	TC
902131	7F	4	2-5-2.	C	110	1	4	TC
902132	7F	4	2-5-2.	D	110	1	4	TC
902133	7F	4	2-5-2.	D	110	1	8	TC
902134	7F	4	2-5-2.	E	110	1	4	TC
902135	7F	4	2-5-2.	E	110	1	8	TC
902136	7F	4	2-5-2.	F	110	1	4	TC
902137	7F	4	2-5-2.	F	110	1	8	TC
902138	7F	4	2-5-2.	DL	85	1	4	TC
902139	7F	4	2-5-2.	DL	110	1	4	TC

* Please note Ablation catheters are also available in Thermistor (TH).

* Ablation catheters not currently available for sale in the United States, Japan or Korea.

MAP-IT™ BI-DIRECTIONAL STANDARD ABLATION CATHETERS*

Reference Number	French Size	Number of Electrodes	Inter-electrode spacing (±0.5 mm)	Curve Type	Length (±0.5 cm)	Ring Electrode width (±0.05mm)	Tip electrode width (±0.2mm)	Temperature Type
902200	7F	4	2-5-2.	AA	110	1	4	TC
902201	7F	4	2-5-2.	BB	110	1	4	TC
902202	7F	4	2-5-2.	CC	110	1	4	TC
902203	7F	4	2-5-2.	DD	110	1	4	TC
902204	7F	4	2-5-2.	DD	110	1	8	TC
902205	7F	4	2-5-2.	EE	110	1	4	TC
902206	7F	4	2-5-2.	FF	110	1	4	TC
902207	7F	4	2-5-2.	FF	110	1	8	TC

* Curves are symmetric to the Standard Ablation Catheter Curves

* Ablation catheters not currently available for sale in the United States, Japan or Korea.

MAP-IT™ IRRIGATED ABLATION CATHETERS*

Reference Number	French Size	Number of Electrodes	Inter-electrode spacing (±0.5 mm)	Curve Type	Length (±0.5 cm)	Ring Electrode width (±0.05mm)	Tip electrode width (±0.2mm)	Temperature Type
902300	7F	4	2-5-2.	A	110cm	1	4	TC
902301	7F	4	2-5-2.	B	110cm	1	4	TC
902302	7F	4	2-5-2.	D	110cm	1	4	TC
902303	7F	4	2-5-2.	F	110cm	1	4	TC
902316	7F	4	2-5-2.	E	110cm	1	4	TC
902338	7F	4	2-5-2.	DL	110cm	1	4	TC

* Ablation catheters not currently available for sale in the United States, Japan or Korea.

MAP-IT™ ABLATION CATHETER EXTENSION CABLES		
Reference Number	Product Description	Length (±0.5 cm)
902402	ABL Catheter EXT CBL for EPT-1000	300cm
902404	ABL Catheter EXT CBL for Stockert	300cm
902406	ABL Catheter EXT CBL for IBI	300cm
902408	ABL Catheter EXT CBL for Generic	300cm
902412	ABL Catheter EXT CBL for Medtronic	300cm

* Ablation catheters not currently available for sale in the United States, Japan or Korea.