

DOSE AREA PRODUCT (DAP) METER VACUDAP

VacuDAP - OEM

The dose area product measuring chamber contains the complete measurement electronics and has a serial RS485 interface. The measuring chamber can be controlled and operated through an ASCII protocol. The measuring results are shown on the control panel of the X-ray generator or a separate display unit. An RS485 - RS232 converter is available to connect the chamber to any host computer.



Downloads

Downloads - Registered User

CAN - Converter for VacuDAP - OEM



BROCHURE



FLYER VACUDAP-OEM