

SPINE CATALOGUE



VERTEBROPLASTIC

VERTEBROPLASTY NEEDLE – PLASTIC HANDLE

VERTEBROPLASTY NEEDLE WITH PLASTIC HANDLE

VERTEBROPLASTIC is used in Vertebroplasty procedures to inject cement into the vertebral body. Vertebroplasty is a very well established technique for treating back pain caused by severe osteoporosis or tumors with consequent loss of height or fracture of the vertebral body. Coaxial biopsy can also be performed.

FEATURES:

1. ERGONOMIC GRIP WITH LOCKING SYSTEM
2. THE ENHANCED DESIGN OF THE PLASTIC HANDLE ALLOWS EASY HANDLING OF THE NEEDLE
3. TIP DIRECTION INDICATOR
4. AISI 304 STAINLESS STEEL CANNULA AND STYLET

AGO PER VERTEBROPLASTICA CON IMPUGNATURA IN PLASTICA

VERTEBROPLASTIC viene utilizzato nelle procedure di vertebroplastica per iniettare cemento nel corpo vertebrale. La vertebroplastica è una tecnica consolidata per il trattamento delle fratture vertebrali causate da grave osteoporosi o da patologie tumorali con conseguente perdita di altezza del corpo vertebrale. È possibile effettuare una biopsia coassiale.

CARATTERISTICHE:

1. IMPUGNATURA ERGONOMICA CON SISTEMA DI BLOCCAGGIO
2. IL DESIGN MIGLIORATO DELL'IMPUGNATURA IN PLASTICA AGEVOLA L'UTILIZZO DELL'AGO
3. INDICATORE DI DIREZIONE PUNTA
4. CANNULA E STILETTO IN ACCIAIO AISI 304



Stylet with bevel tip

Cannula with sharp bevel plus tip to reduce trauma

GAUGE / CALIBRO	DIAMETER / DIAMETRO	100 mm	120 mm	150 mm
11G	3,00	KRVT1110200CBP-A	KRVT1112200CBP-A	KRVT1115200CBP-A
13G	2,50	KRVT1310200CBP-A	KRVT1312200CBP-A	KRVT1315200CBP-A
15G	1,80	KRVT1510200CBP-A	KRVT1512200CBP-A	KRVT1515200CBP-A



BIOPSYBELL S.R.L.

Via Aldo Manuzio 24
41037 Mirandola (MO) - Italy

Tel. +39 0535 27850

Fax. +39 0535 33526

international1@biopsybell.it

www.biopsybell.it

FDA (Establishment Registration Number) 9617616

©Copyright 2019



Rev. 03 del 03-09-19