

Super Sheath™


Introducer Sheath

Smooth Transitions – Total Performance

Focused design for ease of entry
and reduced trauma*

Excellent Device Passage

- Silicone-coated valve and smooth inner surface designed to device passage
- Tricuspid valve design promotes hemostasis
- Radiopaque markers facilitate accurate device placement



Efficient Operation

- Rotating suture wing allows for quick sheath repositioning
- Translucent hub enhances visibility
- Innovative dilator/hub twist-lock feature for security and safety during insertion

Smooth Transitions

- Smooth transitions facilitate ease of entry with low insertion force*
- Improved kink resistance and lubricity promote pushability while maintaining device integrity**
- Long, tapered tip aims to reduce potential for vessel trauma



* Bench testing performed by Togo Medikit Co., Ltd. Data on file.

** Bench testing performed by Boston Scientific Corporation.

Data on file. Bench test results not necessarily indicative of clinical performance.

Super Sheath™ Introducer sheath are products of Togo Medikit Co, Ltd and are distributed by Boston Scientific

Ordering Information

Introducer Sheath and Guidewire

Order number	UPN Code	Size (F) (mm)		Sheath Length	Guide Wire Size (in) (mm)		Units Per Box
15-704B	M00115704B1	4	1.33	7cm	0.035	0.89	10
15-705B	M00115705B1	5	1.67	7cm	0.035	0.89	10
15-706B	M00115706B1	6	2.0	7cm	0.035	0.89	10
15-707B	M00115707B1	7	2.33	7cm	0.035	0.89	10
15-708B	M00115708B1	8	2.67	7cm	0.035	0.89	10
15-710B	M00115710B1	4	1.33	11cm	0.035	0.89	10
15-711B	M00115711B1	5	1.67	11cm	0.035	0.89	10
15-712B	M00115712B1	6	2.0	11cm	0.035	0.89	10
15-713B	M00115713B1	7	2.33	11cm	0.035	0.89	10
15-714B	M00115714B1	8	2.67	11cm	0.035	0.89	10
15-715B	M00115715B1	9	3.0	11cm	0.035	0.89	10
15-716B	M00115716B1	10	3.33	11cm	0.035	0.89	10
15-717B	M00115717B1	11	3.67	11cm	0.035	0.89	10
15-718B	M00115718B1	12	4.0	11cm	0.035	0.89	10
15-741B	M00115741B1	14	4.37	11cm	0.035	0.89	10



Introducer Sheath (Guidewire Not Included)

Order number	UPN Code	Size (F) (mm)		Sheath Length	Max Guide Wire (in) (mm)		Units Per Box
15-720B	M00115720B1	4	1.33	11cm	0.038	0.97	10
15-721B	M00115721B1	5	1.67	11cm	0.038	0.97	10
15-722B	M00115722B1	6	2.0	11cm	0.038	0.97	10
15-723B	M00115723B1	7	2.33	11cm	0.038	0.97	10
15-724B	M00115724B1	8	2.67	11cm	0.038	0.97	10
15-725B	M00115725B1	9	3.0	11cm	0.038	0.97	10
15-726B	M00115726B1	10	3.33	11cm	0.038	0.97	10
15-727B	M00115727B1	11	3.67	11cm	0.038	0.97	10
15-728B	M00115728B1	12	4.0	11cm	0.038	0.97	10
15-740B	M00115740B1	14	4.37	11cm	0.038	0.97	10
15-730B	wM00115730B1	5	1.67	25cm	0.038	0.97	10
15-731B	M00115731B1	6	2.0	25cm	0.038	0.97	10
15-732B	M00115732B1	7	2.33	25cm	0.038	0.97	10
15-733B	M00115733B1	8	2.67	25cm	0.038	0.97	10
15-734B	M00115734B1	9	3.0	25cm	0.038	0.97	10
15-735B	M00115735B1	10	3.33	25cm	0.038	0.97	10
15-736B	M00115736B1	11	3.67	25cm	0.038	0.97	10
15-737B	M00115737B1	12	4.0	25cm	0.038	0.97	10
15-739B	M00115739B1	14	4.37	25cm	0.038	0.97	10

RO Sheath (Guidewire Not Included)

Order number	UPN Code	Size (F) (mm)		Sheath Length	Max Guide Wire (in) (mm)		Units Per Box
15-956B	M00115956B1	5	1.67	7cm	0.038	0.97	10
15-957B	M00115957B1	6	2.0	7cm	0.038	0.97	10
15-958B	M00115958B1	7	2.33	7cm	0.038	0.97	10
15-959B	M00115959B1	8	2.67	7cm	0.038	0.97	10
15-955B	M00115955B1	9	3.0	7cm	0.038	0.97	10
15-960B	M00115960B1	4	1.33	11cm	0.038	0.97	10
15-961B	M00115961B1	5	1.67	11cm	0.038	0.97	10
15-962B	M00115962B1	6	2.0	11cm	0.038	0.97	10
15-963B	M00115963B1	7	2.33	11cm	0.038	0.97	10
15-964B	M00115964B1	8	2.67	11cm	0.038	0.97	10
15-969B	M00115969B1	9	3.0	11cm	0.038	0.97	10
15-970B	M00115970B1	4	1.33	25cm	0.038	0.97	10
15-971B	M00115971B1	5	1.67	25cm	0.038	0.97	10
15-965B	M00115965B1	6	2.0	25cm	0.038	0.97	10
15-966B	M00115966B1	7	2.33	25cm	0.038	0.97	10
15-967B	M00115967B1	8	2.67	25cm	0.038	0.97	10
15-968B	M00115968B1	9	3.0	25cm	0.038	0.97	10

**Boston
Scientific**
Advancing science for life™

www.bostonscientific-international.com

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations.

PI-198004-AA NOV 2013 Printed in Germany by medicalvision.

© 2013 Boston Scientific Corporation
or its affiliates. All rights reserved.
DINPER2129EC