



CT Scanner Accessories





1. Arm holders

► Arm up holder

Used to support the patient's arms when the arms are elevated.



CAAH-014A

Dimensions	Materials
W: 530 mm L: 378 mm H: 288 mm	Poly dichloropentadiene

The metal-free arm up holder is made of a material that minimizes artifacts. It is also light weight and washable. The metal-free arm up holder simplifies patient positioning.


► Handle

This handle is grasped by the patient to hold up his/her arms.




CTA-AHC

Dimensions	Materials
W: 350 mm L: 75 mm H: 165 mm	Iron Chrome-plated



► Arm down holder (one side)

These holders are mounted to the couch top to support the patient arm (one side only).




CAAH010A

For 470-mm-wide couch tops.

Dimensions	Materials
W: 335 mm L: 400 mm H: 48 mm	Polyvinyl chloride

► Arm down holder (both sides)

These holders are mounted to the couch top to support the patient's arms (both sides).




CAAH012A

For 470-mm-wide couch tops.

Dimensions	Materials
W: 700 mm L: 400 mm H: 53 mm	ABS resin

2. Baby bed / Child holder

► **Baby bed**




CTA-TMB

Placed on the couch top to restrain infant movement.

Dimensions	Materials
W: 230 mm L: 510 mm H: 80 mm	Acrylic resin Nylon Protective felt

► **Child holder**

This holder immobilizes infants for a wide range of body and head examinations.



CACH-001A

Dimensions	Materials
W: 370 mm L: 1000 mm H: 200 mm	Base: PET resin, Polyethylene Pillow: Polyethylene foam Mat, cushions: Urethane leather, Flexible urethane foam Bands, Body cover: Urethane leather, Polyester



3. Head rests

Used to immobilize the patient’s head during scanning.

► Head rests



CAHR024A

The tilt angle can be set manually.

Dimensions	Materials
W: 200 mm L: 413 mm H: 148 mm	Carbon fiber Iron



CAHR-026A

For 3D subtraction.

Dimensions	Materials
W: 358 mm L: 629 mm H: 374 mm Forward/backward stroke: 160 mm Vertical stroke: 60 mm	Carbon fiber Polyvinyl chloride Iron



CAHR-028A

Dimensions	Materials
W: 190 mm L: 374 mm H: 208 mm	Carbon fiber reinforced plastic



CAHR-029A

Dimensions	Materials
W: 221 mm L: 478 mm H: 207 mm	Carbon fiber reinforced plastic

3. Head rests

Used to immobilize the patient's head during scanning.

► Coronal head rests



CTA-HRJ (prone)	
Used when the patient is prone for coronal image acquisition.	
Dimensions	Materials
W: 230 mm	Carbon fiber Iron
L: 375 mm	
H: 120 mm	



CTA-HRO (supine)	
Used when the patient is supine for coronal image acquisition.	
Dimensions	Materials
W: 190 mm	Acrylic resin Nylon Iron
L: 360 mm	
H: 170 mm	

► Head positioner for Xe-study



CTA-HRP	
Used to immobilize the patient head in Xe-study.	
Dimensions	Materials
W: 370 mm	Acrylic resin Nylon Iron
L: 405 mm	
H: 278 mm	

► Head rest for infant




CTA-HRN	
Used to immobilize the heads of infants.	
Dimensions	Materials
W: 150 mm	Acrylic resin Nylon Iron
L: 310 mm	
H: 137 mm	

4. Flat board

Placed on the couch top to make the couch-top upper surface flat.

► Flat board



CABM021A	
For 470-mm-wide couch tops.	
Dimensions	Materials
W: 440 mm	Acrylic foam Polyester Cotton
L: 1905 mm	
H: 40 mm	



5. Foot supports

Mounted to the end of the couch top to provide additional table length for feet support.

► Foot rests



CAEX-008A

400 mm extension
(according to foot mat sizes)

Dimensions	Materials
W: 470 mm	Carbon fiber
L: 450 mm	
H: 80 mm	

(To be exchanged with the head rest.)
(With foot-mat fastener)

CAEX-012A

310 mm extension
(according to foot mat sizes)

Dimensions	Materials
W: 465 mm	Carbon fiber reinforced plastic
L: 409 mm	
H: 74 mm	

6. Head mats

► Head mat (cushion)

Placed on the head rest.



CAHM-003A

Dimensions	Materials
W: 450 mm	Polyethylene foam
L: 470 mm	
H: 30 mm	



CAHM-004A

Dimensions	Materials
W: 390 mm	Polyethylene
L: 230 mm	
H: 20 mm	



CAHM-005A

Dimensions	Materials
W: 226 mm	Polyethylene
L: 281 mm	
H: 168 mm	

6. Head mats

► Wedge mats



CAWM-001A/1

10 degrees

Dimensions

W: 186 mm
L: 180 mm
H: 78 mm

Materials

Polyethylene
foam

CAWM-002A/1

10 degrees

Dimensions

W: 189 mm
L: 210 mm
H: 83 mm

Materials

Polyethylene
foam

CAWM-001A/2

20 degrees

Dimensions

W: 186 mm
L: 180 mm
H: 108 mm

Materials

Polyethylene
foam

CAWM-002A/2

20 degrees

Dimensions

W: 190 mm
L: 210 mm
H: 121 mm

Materials

Polyethylene
foam

CAWM-001A/3

30 degrees

Dimensions

W: 186 mm
L: 180 mm
H: 128 mm

Materials

Polyethylene
foam

CAWM-002A/3

30 degrees

Dimensions

W: 186 mm
L: 210 mm
H: 149 mm

Materials

Polyethylene
foam

► Table top pillow



CATP-001A

Dimensions

W: 188 mm
L: 274 mm
H: 115 mm

Materials

Urethane foam
Urethane resin

► Side mats

Placed between the side of the patient's head and the head rest to immobilize the head.



CTA-HMC

Dimensions

W: 140 mm
L: 80 mm
H: 30 mm

Materials

Urethane
Polyester
Cotton

CASM-001A

Dimensions

W: 140 mm
L: 80 mm
H: 30 mm

Materials

Polyurethane
Polyester
Cotton



7. Body mats

Placed on the couch top.

► Body mats



CABM-031A/1

For 2000-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 2360 mm H: 20 mm	Polyurethane	Zenith Blue

CABM-031A/4

For 300 kg patient couch with a 1500-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 1860 mm H: 20 mm	Polyurethane	Zenith Blue



CABM-031A/5

For 2000-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 2360 mm H: 20 mm	Polyurethane	Terra Gray

CABM-031A/6

For 1800-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 2090 mm H: 20 mm	Polyurethane	Terra Gray

CABM-031A/7

For 1500-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 1790 mm H: 20 mm	Polyurethane	Terra Gray

CABM-031A/8

For 300 kg patient couch with a 1500-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 1860 mm H: 20 mm	Polyurethane	Terra Gray

CABM-032A/2

For 1700-mm-long couch top

Dimensions	Materials	Color
W: 470 mm L: 2070 mm H: 20 mm	Polyurethane	Terra Gray

► Body mat covers



CAMC-001A/1

Dimensions	Materials
W: 1020 mm L: 2900 mm	Polyvinyl chloride

CAMC-001A/2

Dimensions	Materials
W: 1020 mm L: 2600 mm	Polyvinyl chloride

7. Body mats

Placed on the couch top.

► Triangular mats

Placed under the patient's knees during scanning.



CAFM-002A/1

Dimensions	Materials	Color
W: 350 mm L: 570 mm H: 220 mm	Polyvinyl chloride Urethane	Win Blue



CAFM-002A/2

Dimensions	Materials	Color
W: 350 mm L: 570 mm H: 220 mm	Polyurethane	Terra Gray

► Foot mats



CAFM-001B For CAEX-008A

Dimensions	Materials	Color
W: 470 mm L: 400 mm H: 20 mm	Polyvinyl chloride Urethane	Zenith Blue



CAFM-003A/2 For CAEX-008A

Dimensions	Materials	Color
W: 470 mm L: 400 mm H: 20 mm	Polyurethane	Terra Gray


CAFM-005A For CAEX-012A

Dimensions	Materials	Color
W: 470 mm L: 310 mm H: 20 mm	Polyurethane	Terra Gray



8. Head straps

► Head straps (fastener)



CAHB-003A/1

Used to immobilize the head.

Dimensions	Materials	Color
W: 50 mm L: 750 mm	Polyurethane	Blue



CAHB-003A/2

Used to immobilize the head.

Dimensions	Materials	Color
W: 50 mm L: 750 mm	Polyurethane	Terra Gray

► Chin straps (fastener)



CAHB-002A/1

Used to immobilize the chin.

Dimensions	Materials	Color
W: 50 mm L: 750 mm	Polyurethane	Blue



CAHB-002A/2


Used to immobilize the chin.

Dimensions	Materials	Color
W: 50 mm L: 750 mm	Polyurethane	Terra Gray

9. Body straps

► Body straps (fastener)

Sliding type bands to be used to restrain patient movement.



CABB-008A/1, /2*

For 2000-mm-long couch top

CABB-008B/1, /2*

For 1800-mm-long couch top

CABB-008C/1, /2*

For 300 kg patient couch with a 1500-mm-long couch top

CABB-008D/1, /2*

For 1500-mm-long couch top

CABB-010A/1

For Aquilion Precision with a 2000-mm-long couch top

CABB-010A/2

For 1700-mm-long couch top



• Body band rails

For 2000-mm-long couch top

Dimensions	Materials
W: 143 mm L: 2355 mm	Polyethylene

For 1800-mm-long couch top

Dimensions	Materials
W: 143 mm L: 2085 mm	Polyethylene

For 300 kg patient couch with a 1500-mm-long couch top

Dimensions	Materials
W: 143 mm L: 1855 mm	Polyethylene

For 1500-mm-long couch top

Dimensions	Materials
W: 143 mm L: 1785 mm	Polyethylene

For Aquilion Precision with a 2000-mm-long couch top

Dimensions	Materials
W: 143 mm L: 2305 mm	Polyethylene

For 1700-mm-long couch top

Dimensions	Materials
W: 143 mm L: 2015 mm	Polyethylene

• Body band 1

Dimensions	Materials
W: 300 mm L: 513 mm	Polyethylene

• Body band 2

Dimensions	Materials
W: 200 mm L: 600 mm	Polyethylene

• Body band 3

Dimensions	Materials
W: 200 mm L: 513 mm	Polyethylene

*: The color of the body straps for /1 differs from that for /2.
The color of the body straps for /2 is the same as that of body mats CABM-031A/5-8.



CANON MEDICAL SYSTEMS CORPORATION

<https://global.medical.canon>

©Canon Medical Systems Corporation 2003-2022. All rights reserved.
Design and specifications are subject to change without notice.
MCACT0093EAI 2022-02 CMSC/Produced in Japan.

Canon Medical Systems Corporation meets internationally recognized standards for Quality Management System ISO 9001, ISO 13485. Canon Medical Systems Corporation meets the Environmental Management System standard ISO 14001.

Aquilion ONE, Aquilion ONE PRISM, Aquilion ONE GENESIS, Aquilion, Aquilion Lightning, Aquilion Exceed LB and Made for Life are trademarks of Canon Medical Systems Corporation.

Made For life