



## FSTERMED Twain

- For the use of the patients carried by ambulance and first aid.
- Brings the two wrist together and enables the hands at the front of the body .
- Skin touch cotton surface absorbs the sweat.



- In case of need with the velcro lock at the bottom surface can be stable by fixing to the bad or litter.
- In case of long term use does not cause allergy and irritation.

CE



## FSTERMED Classic

- For the use of inpatient both bed or litter.
- Skin touch cotton surface absorbs the sweat.

- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.

 CE<sub>1783</sub>

### Fstermed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Deluxe	34	1	25	20 x 30 x 20	1215 gr.	ALH 2006
Deluxe / Çift	34	2	12	20 x 30 x 20	1175 gr.	ALH 2007
Twain	32	1	25	20 x 30 x 20	1275 gr.	TC 2009
QRB	27	1	25	20 x 30 x 20	1790 gr.	QRB 2013
QRB / Çift	27	2	12	20 x 30 x 20	1726 gr.	QRB 2014
Klasik El	25	1	25	20 x 30 x 20	850 gr.	ALH 2001
Klasik El / Çift	25	2	12	20 x 30 x 20	815 gr.	ALH 2008
Klasik Ayak	40	1	25	20 x 30 x 20	910 gr.	ALH 2002
Klasik Ayak / Çift	40	2	12	20 x 30 x 20	885 gr.	ALH 2010
Klasik Bacak	65	1	25	20 x 30 x 20	1100 gr.	ALH 2003
Metal El	25	1	25	20 x 30 x 20	1215 gr.	ALH 2011
Metal El / Çift	25	2	12	20 x 30 x 20	1130 gr.	ALH 2012

# BINDERMED

## Baby Limb Holder



- Easy to use and pratique.
- Both use for hand and leg wrist.
- Skin touch surface absorbs the sweat.

- Covers the wrist safely and stable with its connection ribbons.
- Does not irritate the skin and cause allergy.
- Gamma sterilized.

Bindermed Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Bindermed	3 x 11	1	25	11 x 20 x 25	263 gr.	BLH 3001

CE 1783

# NASOMED

## Baby Nasal Cannula Fixation Band



- Pratique and easy to use with its velcro bands.
- Can be used with every dimension and model of cannula.
- Covers the head smoothly, skin touch surface absorbs the sweat.
- Keeps the cannula safely at the babies nose getting through the hoops of the nasal cannula.
- Does not irritate the skin and cause allergy.
- Gamma sterilized.

Nasomed Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Nasomed	43	1	25	11 x 20 x 25	245 gr.	NTH 5001


CE 1783

# STAENDERMED


## Endotracheal Cannula Holders



CLP



CLP TP



### STAENDERMED - CLP

Clamp supported holder bands


Holder band fixes the cannula by gripping it with its locked clasp (CLP). There are 2 different models.

**CLP - Single piece, velcro connection is on the back of the head model.**

It grips and fixes the endotracheal tube with its clasp, it wraps around the neck and velcro s connects at the back of the head. There are 2 different sizes.

**CLP TP - Two pieces, velcro connection is on the side of the head model.**

It grips and fixes the endotracheal tube with its clasp, it wraps around the neck and velcro s connects at the side of the head. There is only one standard single size.




V1



V2



V3 TP



### STAENDERMED - V

Velcro Supported Holders

Holds the canula with the agency of velcro (V) lock, 3 different models

**V1 Single piece mid-supported**

The middle portion of the tube wraps, 2 sizes.

**V2 Single piece side-supported**

The part close to the tube-to-end wraps 3 sizes.

**V3 Two pieces model**

One size.

- Cotton surfaces touching the skin absorbs the sweat.
- Does not cause allergy and irritation in case of long term use.
- Cannula can be placed and removed without moving the patient's neck

- by the usage of single piece side-supported and two pieces models.
- Three different types of Endotracheal Cannula Holders: Velcro supported, silicon supported and clemp supported.
- Gamma sterilized.




S1



S2



S3 TP



## STAENDERMED - SC

Silicone Supported Holder Bands

Holder band fixes the cannula by gripping it with its silicone surface (S). There are 3 different models.

### **S1** Single piece, supported at middle part model.

The middle part wraps and grips the cannula. There are 2 sizes.

### **S2** Single piece, supported by side part model.

The part close to one end wraps and grips the cannula. There are 3 sizes.

### **S3 TP** Two pieces model

There is only one standard single size.

### Staendermed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
V1	55	1	25	11 x 20 x 25	330 gr.	ETH 101
	65	1	25	11 x 20 x 25	345 gr.	ETH 102
V2	45	1	25	11 x 20 x 25	330 gr.	ETH 201
	55	1	25	11 x 20 x 25	345 gr.	ETH 202
	65	1	25	11 x 20 x 25	330 gr.	ETH 203
V3 TP	55	1	25	11 x 20 x 25	330 gr.	ETH TP
S1	55	1	25	11 x 20 x 25	330 gr.	ETH-SC 101
	65	1	25	11 x 20 x 25	345 gr.	ETH-SC 102
S2	55	1	25	11 x 20 x 25	330 gr.	ETH-SC 201
	65	1	25	11 x 20 x 25	345 gr.	ETH-SC 202
	55	1	25	11 x 20 x 25	330 gr.	ETH-SC 203
S3 TP	50 + 28	1	25	11 x 20 x 25	330 gr.	ETH-SC TP
CLP	55	1	25	11 x 20 x 25	330 gr.	ETH CLP 101
	65	1	25	11 x 20 x 25	345 gr.	ETH CLP 102
CLP TP	50 + 25	1	25	11 x 20 x 25	330 gr.	ETH CLP TP