

Rubicon™
Support Catheters

- Enhances lesion access and entry
- Facilitates crossing of tight occlusions & challenging vasculature
- Enables delivery with and exchange of 0.014" (0.36mm), 0.018" (0.46mm) & 0.035" (0.89mm) guidewires

Rubicon 14


Order Number	Guidewire Compatibility (in/mm)	Sheath Compatibility	Length (cm)	Tip Shape	Tip Entry Profile	Proximal Shaft Profile	Marker Band Spacing
H74939212014131	0.014" (0.36 mm)	4F*	135	Straight	1.4F*	3.4F*	15mm
H74939212014151	0.014" (0.36 mm)	4F*	150	Straight	1.4F*	3.4F*	15mm

Rubicon 18

Order Number	Guidewire Compatibility (in/mm)	Sheath Compatibility	Length (cm)	Tip Shape	Tip Entry Profile	Proximal Shaft Profile	Marker Band Spacing
H74939239018091	0.018" (0.46 mm)	4F*	90	Straight	1.7F*	4F*	15mm
H74939239018131	0.018" (0.46 mm)	4F*	135	Straight	1.7F*	4F*	15mm
H74939239018151	0.018" (0.46 mm)	4F*	150	Straight	1.7F*	4F*	15mm

Rubicon 35

Order Number	Guidewire Compatibility (in/mm)	Sheath Compatibility	Length (cm)	Tip Shape	Tip Entry Profile	Proximal Shaft Profile	Marker Band Spacing
H74939240035061	0.035" (0.89 mm)	5F*	65	Straight	2.84F*	5.1F*	50mm
H74939240035091	0.035" (0.89 mm)	5F*	90	Straight	2.84F*	5.1F*	50mm
H74939240035131	0.035" (0.89 mm)	5F*	135	Straight	2.84F*	5.1F*	50mm
H74939240035151	0.035" (0.89 mm)	5F*	150	Straight	2.84F*	5.1F*	50mm



Rubicon 14 support catheter



* 1.4F = 0.47 mm (0.019")
1.7F = 0.57 mm (0.022")
2.84F = 0.95 mm (0.037")
3.4F = 1.13 mm (0.044")
4F = 1.33 mm (0.052")
5F = 1.67 mm (0.066")
5.1F = 1.7 mm (0.067")