

# NASAL EPISTAXIS BALLOON

Safe and Fast Nasal Bleeding Control



OTOLARYNGOLOGY  
SERIES




# NASAL EPISTAXIS BALLOON

## Fortune® Silicone Nasal Epistaxis Balloon Features:

- 100% medical grade silicone for superior biocompatibility.
- Inflatable flexible balloon provides direct uniform pressure.
- Soft and beveled tip facilitates insertion and minimizes trauma.
- Larger balloon can be inflated up to 30 cc to control anterior nasal bleeding.
- Smaller balloon can be inflated up to 10 cc to control posterior nasal bleeding.
- Airway channel design allows nasal breathing.
- Unique pilot balloon design allows cuff monitoring to prevent over-pressure.
- Ideal for nasal packing for septoplasty, rhinoplasty and intra-nasal surgical procedure.

- A** Nasal Epistaxis Balloon, Intra-Pack
- B** Nasal Epistaxis Balloon, Post-Pack
- C** Nasal Epistaxis Balloon, Double Balloon




# NASAL EPISTAXIS BALLOON

## SPECIFICATIONS

### Nasal Epistaxis Balloon, Intra-Pack:


| REF. No.  | Size  | Application                     |
|-----------|-------|---------------------------------|
| 3101-0830 | 30 cc | Intra-anterior<br>nasal chamber |



3101-0830

### Nasal Epistaxis Balloon, Post-Pack:


| REF. No.  | Size  | Application                |
|-----------|-------|----------------------------|
| 3102-2010 | 10 cc | Posterior<br>nasal chamber |



3102-2010

### Nasal Epistaxis Balloon, Double Balloon:

| REF. No.  | Size     | Application                             |
|-----------|----------|-----------------------------------------|
| 3103-0930 | 30/10 cc | anterior and posterior<br>nasal chamber |



3103-0930



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,  
Danshuei Dist, New Taipei City. 251 Taiwan  
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266  
www.fortunemed.com Mail:info@fortunemed.com

*Distribuito autorizat in Moldova*



[www.ericon.md](http://www.ericon.md)

**Tel:0-22-520108**

**mob.079410042**

**FONDAT 1992**

15,Zelinski str.  
Chişinău,MOLDOVA