

Disposable Hemostasis Clips and Delivery Systems


Intended use

Clamping the mucosa tissue of alimentary cavity under the guide of endoscope.

A) Hemostasis for

- 1) Mucosal / sub-mucosal defects <3cm
- 2) Bleeding ulcers
- 3) Arteries <2mm
- 4) Surgery sites
- 5) Diverticula in colon

B) Closure of GI luminal perforations



Features

The clip is pre-load in disposable delivery system, which is convenient to use.

- Accurate sync-rotation offers exact clip deployment.
- Reinforced design makes the clip performance more steady.
- The biggest opening is 14mm after maximize the clip.

Specifications

Model	Working length	Jaw angle	Working Channel
SX-HC-1650-090L	1650mm	90°	≥2.8mm
SX-HC-1650-135L	1650mm	135°	≥2.8mm
SX-HC-2300-090L	2300mm	90°	≥2.8mm
SX-HC-2300-135L	2300mm	135°	≥2.8mm