



## Inspire nCPAP™ Circuits

Non-invasive respiratory support



# Inspire nCPAP<sup>™</sup> Circuits

Gentle non-invasive respiratory care

## Profile

Nasal CPAP for treating infants has developed a great deal in the past 20 years with many devices now being available. Inspire nCPAP<sup>™</sup> from Inspiration Healthcare builds on the experience of the past to allow care-givers a cost effective solution for supporting the respiratory compromised infant.

nCPAP works well when the Added Work Of Breathing is reduced ensuring the infants do not tire so quickly and allowing them to use energy to grow. Using established fluidic principles, the Inspire nCPAP<sup>™</sup> range is more effective than traditional methods of delivering nCPAP.



The Inspire nCPAP<sup>™</sup> range of breathing circuits are compatible with the majority of flow drivers and humidifier systems found in the neonatal intensive care area. Curly heated wire technology helps reduce condensation (rainout) delivering optimum humidification to the infant.



IHC4064 / 20N

# Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit (Flow Drivers)

Compatible with the MR850 series humidifiers complete with nasal prongs






## Ordering Information

Part Number	Description	Case quantity
IHC4064/20N	Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit (Flow Drivers), compatible with the MR850 series humidifiers complete with nasal prongs	20

# Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit (Flow Drivers)

IHC4112 / 20  
IHC4113 / 10

Compatible with the MR850 series humidifiers



Autofeed Humidification Chamber



Nasal sizer



## Ordering Information

Part Number	Description	Case quantity
IHC4112/20	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Flow Drivers), compatible with the MR850 series humidifiers</b>	20
IHC4113/10	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Flow Drivers), compatible with the MR850 series humidifiers complete with autofeed chamber</b>	10

# Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit (Fabian Ventilator)

IHC4066 / 20N

Compatible with the MR850 series humidifiers complete with nasal prongs





## Ordering Information

Part Number	Description	Case quantity
IHC4066/20N	Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit, (Fabian Ventilator), compatible with the MR850 series humidifiers complete with nasal prongs	20

# Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit (Fabian Ventilator)

IHC4114 / 20  
IHC4115 / 10

Compatible with the MR850 series humidifiers



Autofeed Humidification Chamber



Nasal sizer



## Ordering Information

Part Number	Description	Case quantity
IHC4114/20	Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit, (Fabian Ventilator), compatible with the MR850 series humidifiers	20
IHC4115/10	Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit (Fabian Ventilator), compatible with the MR850 series humidifiers complete with autofeed chamber	10

# Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit (Flow Drivers and Fabian)

IHC4067 / 20N

Compatible with the MR850 series humidifiers complete with nasal prongs




## Ordering Information

Part Number	Description	Case quantity
IHC4067/20N	Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit, (Flow Drivers and Fabian), compatible with the MR850 series humidifiers complete with nasal prongs	20

# Inspire nCPAP™ Generator & Cloverleaf Curly Wire Breathing Circuit (Flow Drivers and Fabian)

IHC4116/20  
IHC4117/10

Compatible with the MR850 series humidifiers



## Ordering Information

Part Number	Description	Case quantity
IHC4116/20	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit, (Flow Drivers and Fabian), compatible with the MR850 series humidifiers</b>	20
IHC4117/10	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Flow Drivers and Fabian), compatible with the MR850 series humidifiers complete with autofeed chamber</b>	10

## Ordering Information

Part Number	Description	Case quantity	Prongs inc (S.M.L)	Autofeed Chamber inc	Connects to Fabian	Connects to Standard drivers
IHC4064/20N	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Flow Drivers), compatible with the MR850 series humidifiers complete with nasal prongs</b>	20	•			•
IHC4112/20	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Flow Drivers), compatible with the MR850 series humidifiers</b>	20				•
IHC4113/10	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit, compatible with the MR850 series humidifiers complete with autofeed chamber</b>	10		•		•
IHC4066/20N	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Fabian Ventilator), compatible with the MR850 series humidifiers complete with nasal prongs</b>	20	•		•	
IHC4114/20	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Fabian Ventilator), compatible with the MR850 series humidifiers</b>	20			•	
IHC4115/10	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Fabian Ventilator), compatible with the MR850 series humidifiers complete with autofeed chamber</b>	10		•	•	
IHC4067/20N	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Flow Drivers and Fabian), compatible with the MR850 series humidifiers complete with nasal prongs</b>	20	•		•	•
IHC4116/20	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Flow Drivers and Fabian), compatible with the MR850 series humidifiers</b>	20			•	•
IHC4117/10	<b>Inspire nCPAP™ Generator &amp; Cloverleaf Curly Wire Breathing Circuit (Flow Drivers and Fabian), compatible with the MR850 series humidifiers complete with autofeed chamber</b>	10		•	•	•
IHC602/20	<b>Inspire nCPAP™ Silencer</b>	20				

The details given in this leaflet are correct at time of going to press. The company reserves the right to improve the products shown.

500-116(2)



Tel: +44 (0) 1455 840 555  
Email: [info@inspiration-healthcare.com](mailto:info@inspiration-healthcare.com)  
Web: [inspiration-healthcare.com](http://inspiration-healthcare.com)