
isleeve Expandable Introducer Set



Included with the ACURATE *neo*[™] Aortic Valve System



iSLEEVE Introducer Set Specifications

Recommended Guidewire	0.035" (0.89 mm)
Sheath Working Length	32 cm
Dilator Length	48.7 cm
Valve Compatibility	For use with all ACURATE neo Aortic Valves
Use	Single Use
Sterilization	Ethylene Oxide Gas
Set Contents	 Introducer Sheath (1) Dilator (1) Product Mandrel (1) Flushing Tube (1)
Carton Size	171.45 mm × 28.58 mm × 711.21 mm
Shelf Life	2 Years
Indications for Use	The 14F iSLEEVE Introducer Set is intended to facilitate femoral access to the vascular system for introduction and removal of the ACURATE ^{M} TF Valve System and ancillary devices. The iSLEEVE Introducer Set is suitable for use in vessels \geq 5.5 mm diameter.