

Paketleme Ürünleri Packaging Products

4A Sterilizasyon Ruloları, steril olacak malzemelerin içerisine konulduğu paketleme poşetleri olup 6 ay boyunca ürünün steril olarak saklanması sağlayan, poşet içerisine hiçbir mikroorganizmanın girmesine izin vermeyen ürünlerdir.

4A Sterilizasyon Rulolarının kenarında bulunan işlem indikatörleri su bazlı olup ISO 11140-1 standardını karşılar.

4A Sterilizasyon Ruloları buhar, etilen oksit, formaldehit sterilizasyon yöntemleri için kullanılmaktadır. Sterilizasyon koşulları, sterilize edilecek malzeme ile ilgili olarak son kullanıcı tarafından belirlenmelidir.

4A Sterilizasyon Ruloları ISO 11607 standartlarını karşılayan ürünler olup, en yüksek kalite ve son teknoloji ile üretilmektedir.

OEM üretim yapılabilmektedir.

4A Sterilization Reels provide sterile storage in packages and keep package sterile for 6 months and do not permit microorganisms to enter inside the package after sterilization cycle.

The process indicators that are printed on the side of the Sterilization Reels are water based and conforms with ISO 11140.

4A Sterilization Reels are used for steam, ethylene oxide, formaldehyde sterilization methods. The sterilization conditions should be determined by the end user regarding to material to be sterilized.

4A Sterilization Reels meet ISO 11607 standards and they are manufactured using state-of-art technology with the highest quality.

Products can be produced with your trademark

Sterilizasyon Düz Ruloları Sterilization Flat Reels

4A Sterilizasyon Düz Ruloları medikal kâğıttan 60gr/m² veya 70gr /m² ve 52 mikron, PP+PET yeşil renkli filmden üretilmektedir. 5 cm ile 60 cm arası herhangi bir ölçü için özel üretim yapılabilmektedir. Talep üzerine 100 m uzunlukta da üretilebilmektedir.

4A Sterilization Flat Reels can be produced from medical paper as 60gr/m² veya 70gr /m² and 52 micron, PP+PET green film. Any dimension between 5cm - 60cm can be made on request. On request, 100 m length can be produced as well.



4AFR

4a medical®
producing health for the world

Sterilizasyon Körüklü Ruloları *Sterilization Gusset Reels*



4AGR

4A Sterilizasyon Körüklü Ruloları medikal kâğıttan 60gr/m² veya 70gr /m² ve 52 mikron, PP+PET yeşil renkli filmde üretilmektedir. 7,5 cm ile 50 cm arası herhangi bir ölçü için özel üretim yapılabilmektedir. Körüklü Rulolara hacimsel malzemeler konulabilmektedir ve körüklü rulolar düz rulolara kıyasla daha fazla iç hacme sahiptir.

4A Sterilization Gusset Reels can be produced from medical paper as 60gr/m² or 70gr /m² and 52 micron, PP+PET green film. Any dimension between 7,5cm - 50cm can be made on request. Gusset Reels has increased filling volume and can hold an increased amount of volume in comparison with Flat Reels.

Sterilizasyon Zarfları *Sterilization Pouches*

4A Sterilizasyon Zarfları tıbbi kâğıttan 60gr/m² veya 70gr / m² ve 52 mikron, PP+PET yeşil renkli filmde üretilmektedir. 4A Sterilizasyon Zarfları kullanıcının isteğine bağlı olarak, istenilen ölçülerde en boy olarak üretilip kullanıma hazır bir şekilde sunulmaktadır.

4A Sterilizasyon Zarfları ISO 11607 standartlarına uygun olup işlem indikatörleri ISO 11140-1 standartlarını karşılamaktadır.

OEM üretim yapılabilmektedir

4A Sterilization Flat Pouches can be produced from medical paper as 60gr/m² or 70gr /m² and 52 micron, PP+PET green film. 4A Sterilization Pouches are customized with regards to their width and length and they are ready-to-use.

4A products comply with ISO 11607 standards while process indicators meet ISO 11140 standards.

Products can be produced with your trademark.

4a medical[®]
producing health for the world

Sterilization Flat Reels

Product Features

- Medical Paper provides maximum safety.
- All imprint are located outside the packaging area to prevent ink pigment migration to the product.
- Transparent , multilayer co-polymer film allows easy identification.
- Reinforced film to avoid tear during opening.
- Superior barrier and quality with 60gsm or 70gsm medical kraft paper.
- Indicators which are water based, non toxic and give accurate result for steam, EO Gas and Formaldehyde (FO).
- Strong seal-strength, visible sealing lines.
- Clean fibre-free opening.
- All materials comply with international standards.
- Opening direction marked.
- Conformity with EN 868-5.




Note: Sterilization Rolls can be produced with your trademark. Any dimension between 5cm – 60cm can be made on request. On request, 100 m length can be produce as well.

Product Types and Sizes

| Product Code | Product Name | Product Sizes |
|--------------|-----------------------------|---------------|
| 4AFR050200 | 4A Sterilization Flat Reels | 50mmx200m |
| 4AFR075200 | 4A Sterilization Flat Reels | 75mmx200m |
| 4AFR100200 | 4A Sterilization Flat Reels | 100mmx200m |
| 4AFR150200 | 4A Sterilization Flat Reels | 150mmx200m |
| 4AFR200200 | 4A Sterilization Flat Reels | 200mmx200m |
| 4AFR250200 | 4A Sterilization Flat Reels | 250mmx200m |
| 4AFR300200 | 4A Sterilization Flat Reels | 300mmx200m |
| 4AFR350200 | 4A Sterilization Flat Reels | 350mmx200m |
| 4AFR400200 | 4A Sterilization Flat Reels | 400mmx200m |
| 4AFR450200 | 4A Sterilization Flat Reels | 450mmx200m |
| 4AFR500200 | 4A Sterilization Flat Reels | 500mmx200m |

FLAT REEL



Indication for Use

The reels are cut with 5-7cm part from top and bottom. After cutting process, the product is placed into the reel. Top and bottom part will be opened. It is closed with suitable heat and sealing machine. If necessary, the indicator can be put into the reel. After that it will be ready for sterilization process. There are 2 or 3 indicators on the side of reels. Indicators are for steam, ethylene oxide and formaldehyde.

Color change

After sterilization color change as;
(Steam) Blue to Dark Green
(EO) Pink to Orange
(FO) Purple to Green

Interpretation

Reels are 5, 7.5, 10, 15, 20, 25, 30, 35 and 40 cm (+/- 1 cm) width and 200m. length (+/-1 m) . Gusset reels should be preferred for larger volumetric products. The sterilization process is made at 134°C for 3,5min or 7 min 121°C for 15min.

Storage

Keep products at room temperature 25°C

Expiry Date

The expiry date is printed on the product. Do not use after expiration date listed on the packs.


Lot Number

Unique identification code, LOT is printed on each reel.

Interfering Substances or Conditions

There are no interfering substances or conditions that could affect the intended use of the indicator or adversely affect the reel performance.

FLAT REEL



Release of Toxic Substances

The indicator releases no toxic substances in sufficient quantities to cause either a health hazard or a hazard to the intended properties of the product being sterilized during or after the sterilization process.

Additional Information

Please pay attention to the deformation and tear on the packs.

Do not overload the packs.

Please pay attention to the seal strength of both ends.

Do not use damaged products.

Flat Reel Packaging

| Dimension | Weight of Box (kg) | Pcs in a Box | Box Dimension(cm) | | |
|---------------|--------------------|--------------|-------------------|----|----|
| 5 cm x 200m | 12 | 10 | 42 | 23 | 27 |
| 7,5 cm x 200m | 14.5 | 8 | 42 | 23 | 32 |
| 10 cm x 200m | 19 | 8 | 42 | 23 | 42 |
| 12cm x 200m | 12.5 | 4 | 42 | 23 | 27 |
| 15 cm x 200m | 15 | 4 | 42 | 23 | 32 |
| 20 cm x 200m | 19.5 | 4 | 42 | 23 | 42 |
| 25 cm x 200m | 12 | 2 | 42 | 23 | 27 |
| 30 cm x 200m | 16 | 2 | 42 | 23 | 32 |
| 35 cm x 200m | 18 | 2 | 42 | 23 | 37 |
| 40 cm x 200m | 22 | 2 | 42 | 23 | 42 |

DECLARATION OF CONFORMITY

MANUFACTURER: DÖRT-A TIP MALZEMELERİ SANAYİ İTHALAT İHRACAT TİCARET LİMİTED ŞİRKETİ

Address: Balıkhisar Mahallesi Köyiçi Serpmeleri No795A Akyurt, Ankara, Turkey

Products: Bowie Dick Test Pack, Indicator Strips (H2O2 Indicator Strip, Formaldehyde Indicator Strip, Ethylene Oxide Indicator Strip, Dry heat Indicator Strip, Type 4 Indicator Strip 4ASI, Type 5 Indicator Strip, Type 6 Indicator Strip), Type 5 Integrator, Rapid Steam Biological Indicator, Longtime Steam Biological Indicator, Ethylene Oxide Biological Indicator, H2O2 (Plasma) Biological Indicator, Helix Group Tests, PCD Group Tests, Ethylene Oxide load control test, Autoclave Tapes (Steam, Ethylene Oxide, Plasma, Formaldehyde), Sterile Container System, Container Label, Container Seal, Container Filter, Documentation Labels with Indicator, Reel Barcode Labels with Indicator; Label Gun, Washer Disinfectors, Ultrasonic Devices and Washing Control Tests of Surgical Instruments (Pro Test, Hemo Tests, Washer Test, Cannula control Test, Sonicontrol Test.

Above described products complied with below norms.

| DOCUMENT NO | TITLE | EDITION / DATE OF ISSUE |
|-------------------|---|-------------------------|
| TS EN ISO 11140-1 | Chemical Indicators | 18.02.2015 |
| TS EN ISO 11140-4 | Bowie Dick Test Pack | 31.01.2008 |
| TS EN 11138-1 | Biological Indicators | 24.04.2017 |
| TS EN 15883-1 | Washing Machine Disinfectant Residue Test | 30.10.2014 |

Additional information: The development, production and the distribution is supported with a Quality Management System according to the requirements of the ISO 9001:2015 and ISO 13485:2016. The Quality Management System is certificated through the notified body ROYALCERT (Certificate No: 108 / DOR09B and 108 / DOR13A)

If no changes made Declaration of Conformity will remain valid.

On 24th February, 2022 at Ankara, Turkey

Canan Öktem

General Manager

4a medical[®]
producing health for the world

**Dört-A Tıp Malzemeleri
Sanayi İth. İhr. Tic. Ltd. Şti.**

ADRES : Balıkhisar Mah. Köyiçi Serpmeleri No: 795 / A
Akyurt - ANKARA / TÜRKİYE
TEL : 0312.363 50 52-53
FAX : 0312.363 52 10
WEB : http://www.4amedical.com
E-MAIL : ankara@4amedical.com

DECLARATION OF CONFORMITY

We, DÖRT-A TIP MALZEMELERİ SANAYİ İTHALAT İHRACAT TİCARET LİMİTED ŞİRKETİ, declare that the following products are in compliance with the provisions of the Medical Device Regulation 2017/745 for medical devices. All supporting documents are stored within the manufacturer. We declare that we will fulfill the obligations imposed by the quality system and that we will maintain the quality system adequately and effectively.

MANUFACTURER: DÖRT-A TIP MALZEMELERİ SANAYİ İTHALAT İHRACAT TİCARET LİMİTED ŞİRKETİ

ADDRESS: Balıkhisar Mah. Köyiçi Serpmeleri No:795/A AKYURT/ ANKARA Turkey

Products: Sterilization Reels and Pouches (Flat, Gusset, Tyvek), Flat and Self Sealing Sterilization Pouches, SMS Paper, Crepe Paper Sheet

Class: Class-I Nonsteril Products-ANNEX VIII

BASIC UDI-DI : 868082970SG3

BASIC UDI-DI : 869744275SJD

BASIC UDI-DI : 86811764SW2

The products depicted above comply with the following directives.

| DOCUMENT NO | TITLE | PRINT DATE OF REGULATION |
|--------------------|---|---------------------------------|
| 2017/745/MDR | Annex IV of the 2017/745/ Medical Device Regulation | 25.05.2017 |
| TS EN 868-5 | Packaging Materials | 28.01.2019 |
| TS EN 868-9 | Packaging Materials | 28.01.2019 |
| TS EN 11607-1 | Packaging Materials | 09.11.2020 |
| TS EN 11607-2 | Packaging Materials | 09.11.2020 |

| Primary Product Number | Product Name | BASIC UDI-DI |
|-------------------------------|--|---------------------|
| 8680829700019 | SMS Paper 500 x 500 mm | 868082970SG3 |
| 8680829700026 | SMS Paper 600 x 600 mm | 868082970SG3 |
| 8680829700040 | SMS Paper 900 x 900 mm | 868082970SG3 |
| 8680829700057 | SMS Paper 1000 x 1000 mm | 868082970SG3 |
| 8680829700064 | SMS Paper 1200 x 1200 mm | 868082970SG3 |
| 8680829700095 | Sterilization Gusset Reel 200 mm x 60 mm x 100 m, 70g 2 Indicators | 868082970SG3 |
| 8680829700101 | Tyvek Reel 300 mm x 70 m | 868082970SG3 |
| 8680829700118 | Sterilization Flat Reel 250 mm x 200 m, 60g 2 Indicators | 868082970SG3 |
| 8680829700125 | Sterilization Pouch 100 mm x 300 mm | 868082970SG3 |
| 8680829700132 | Sterilization Pouch 100 mm x 150 mm | 868082970SG3 |
| 8680829700149 | Sterilization Pouch 150 mm x 300 mm Gusset | 868082970SG3 |
| 8680829700156 | Sterilization Pouch 100 mm x 300 mm Gusset | 868082970SG3 |
| 8680829700163 | Sterilization Pouches 70 mm x 255 mm | 868082970SG3 |
| 8680829700170 | Sterilization Pouches 90 mm x 255 mm | 868082970SG3 |
| 8680829700187 | Sterilization Reel Without Print 75 mm x 200 m | 868082970SG3 |
| 8680829700194 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 70g 2 Indicators | 868082970SG3 |
| 8680829700200 | Sterilization Flat Reel 400 mm x 200 m, 70g 2 Indicators | 868082970SG3 |
| 8680829700217 | Sterilization Pouch 90 mm x 230 mm | 868082970SG3 |
| 8680829700224 | Tyvek Reel 200 mm x 70 m | 868082970SG3 |
| 8680829700231 | Sterilization Flat Reel 420 mm x 200 m, 70g 2 Indicators | 868082970SG3 |

| | | |
|---------------|--|--------------|
| 8680829700262 | Sterilization Flat Reel 450 mm x 200 m, 60g 3 Indicators | 868082970SG3 |
| 8680829700279 | Sterilization Gusset Reel 500 mm x 85 mm x 100 m, 60g 2 Indicators | 868082970SG3 |
| 8680829700293 | Sterilization Gusset Reel 100 mm 50 mm x 100 m, 70g 2 Indicators | 868082970SG3 |
| 8680829700316 | Sterilization Gusset Reel 400 mm x 80 mm 100 m, 70g 3 Indicators | 868082970SG3 |
| 8680829700323 | Tyvek Reel 120 mm x 70 m | 868082970SG3 |
| 8680829700330 | Sterilization Pouch 170 mm x 340 mm | 868082970SG3 |
| 8680829700347 | Sterilization Flat Reel 120 mm x 200 m, 60g 2 Indicators | 868082970SG3 |
| 8680829704000 | Tyvek Reel 250 mm x 70 m | 868082970SG3 |
| 8680829707117 | Sterilization Flat Reel 450 mm x 200 m, 70g 3 Indicators | 868082970SG3 |
| 8680829707599 | Sterilization Flat Reel 500 mm x 200 m, 70g 3 Indicators | 868082970SG3 |
| 8681176422517 | Sterilization Flat Reel 450 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176422524 | Sterilization Flat Reel 450 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176422531 | Sterilization Flat Reel 500 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176422548 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 60g 2 Indicators | 86811764SW2 |
| 8681176422555 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176422562 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176422579 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176422586 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176422593 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176422609 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176422616 | Sterilization Gusset Reel 200 mm x 50 mm x 100 m 60 g 2 Indicators | 86811764SW2 |
| 8681176422623 | Düz Sterilizasyon Rulosu 500 mm x 200 m 60gr 2 Indicators | 86811764SW2 |
| 8681176422630 | Sterilization Reel 500 mm x 200 m 70gr. 2 Indicators | 86811764SW2 |
| 8681176422647 | Sterilization Reel 500 mm x 200 m 70gr. 2 Indicators | 86811764SW2 |
| 8681176422654 | Sterilization Pouches | 86811764SW2 |
| 8681176422661 | Sterilization Reel 450 mm x 200 m | 86811764SW2 |
| 8681176422678 | Sterilization Reel 500 mm x 200 m | 86811764SW2 |
| 8681176422685 | SMS Paper | 86811764SW2 |

| | | |
|---------------|--|-------------|
| 8681176422692 | Crepe Paper Sheet | 86811764SW2 |
| 8681176413027 | Sterilization Pouch 150 mm x 210 mm | 86811764SW2 |
| 8681176413034 | Sterilization Flat Reel 400 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176413041 | Sterilization Pouch 150 mm x 350 mm | 86811764SW2 |
| 8681176413058 | Sterilization Gusset Reel 350 mm x 75 mm x 100 m, 60g 2 Indicators | 86811764SW2 |
| 8681176413065 | Sterilization Pouch 90 mm x 230 mm | 86811764SW2 |
| 8681176413072 | Sterilization Gusset Reel 75 mm x 45 mm x 100 m, 60g 2 Indicators | 86811764SW2 |
| 8681176413867 | Sterilization Flat Reel 100 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176413874 | Sterilization Gusset Reel 100 mm x 50 mm x 100 m, 60g 2 Indicators | 86811764SW2 |
| 8681176413881 | Sterilization Pouches | 86811764SW2 |
| 8681176413935 | Sterilization Pouches 100 mm x 300 mm | 86811764SW2 |
| 8681176413942 | Sterilization Pouch 90 mm x 150 mm | 86811764SW2 |
| 8681176413966 | Sterilization Flat Reel 500 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176413973 | Sterilization Flat Reel 75 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176413980 | Sterilization Flat Reel 450 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176413997 | Sterilization Reel Without Print 250 mm x 200 m | 86811764SW2 |
| 8681176414000 | Sterilization Gusset Reel 120 mm x 55 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176414017 | Sterilization Pouches 150 mm x 200 mm | 86811764SW2 |
| 8681176414024 | Sterilization Pouch 150 mm x 200 mm | 86811764SW2 |
| 8681176414031 | Sterilization Pouch 200 mm x 260 mm | 86811764SW2 |
| 8681176414048 | Sterilization Flat Reel 150 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176414055 | Sterilization Pouches 150 mm x 250 mm | 86811764SW2 |
| 8681176414062 | Sterilization Pouches 150 mm x 400 mm | 86811764SW2 |
| 8681176414079 | Sterilization Reel 500 mm x 200 m | 86811764SW2 |
| 8681176414086 | Sterilization Gusset Reel 450 mm x 80 mm 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176414093 | Sterilization Flat Reel 250 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176414116 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 60g 2 Indicators | 86811764SW2 |
| 8681176414123 | Sterilization Flat Reel 400 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414130 | Tyvek Reel 300 mm x 70 m | 86811764SW2 |
| 8681176414147 | Sterilization Flat Reel 450 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176414154 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 70g 2 Indicators | 86811764SW2 |

| | | |
|---------------|---|-------------|
| 8681176414161 | Sterilization Gusset Reel 150 mm x 55 mm x 100 m, 60g 2 Indicators | 86811764SW2 |
| 8681176414178 | Sterilization Flat Reel 120 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414192 | Sterilization Flat Reel 500 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414215 | Sterilization Reel Without Print 420 mm x 200 m | 86811764SW2 |
| 8681176414222 | Sterilization Gusset Reel 150 mm x 55 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176414253 | Tyvek Reel 400 mm x 70 m | 86811764SW2 |
| 8681176414260 | Sterilization Gusset Reel 100 mm x 50 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176414277 | Sterilization Gusset Reel 400 mm x 80 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176414284 | Sterilization Flat Reel 300 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414291 | Sterilization Pouches 300 mm x 300 mm | 86811764SW2 |
| 8681176414307 | Sterilization Flat Reel 450 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414314 | Sterilization Gusset Reel 250 mm x 65 mm x 100 m, 60g 2 Indicators | 86811764SW2 |
| 8681176414321 | Sterilization Reel Without Print 400 mm x 200 m | 86811764SW2 |
| 8681176414338 | Sterilization Reel Without Print 300 mm x 200 m | 86811764SW2 |
| 8681176414345 | Sterilization Flat Reel 200 mm x 200 m , 70g 3 Indicators | 86811764SW2 |
| 8681176414352 | Sterilization Gusset Reel 120 mm x 55 mm x 100 m, 60g 2 Indicators | 86811764SW2 |
| 8681176414420 | Sterilization Flat Reel 50 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176414437 | Sterilization Pouch 200 mm x 460 mm | 86811764SW2 |
| 8681176414444 | Sterilization Gusset Reel 300 mm x 70 mm x 100 m , 70g 2 Indicators | 86811764SW2 |
| 8681176414451 | Sterilization Reel Without Print 450 mm x 200 m | 86811764SW2 |
| 8681176414468 | Sterilization Flat Reel 50 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414475 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414482 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176414499 | Tyvek Reel 250 mm x 70 m | 86811764SW2 |
| 8681176414512 | Sterilization Pouch 250 mm x 300 mm | 86811764SW2 |
| 8681176414529 | Sterilization Gusset Reel 75 mm x 45 mm x 100 m , 60g 3 Indicators | 86811764SW2 |
| 8681176414536 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414543 | Sterilization Gusset Reel 250 mm x 65 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414567 | Sterilization Flat Reel 350 mm x 200 m, 60g 2 Indicators | 86811764SW2 |

| | | |
|---------------|--|-------------|
| 8681176414581 | Sterilization Gusset Reel 350 mm x 75 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176414598 | Sterilization Pouch 150 mm x 250 mm | 86811764SW2 |
| 8681176414611 | Sterilization Reel Without Print 120 mm x 200 m | 86811764SW2 |
| 8681176414628 | Sterilization Gusset Reel 200 mm x 60 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414635 | Tyvek Reel 400 mm x 70 m | 86811764SW2 |
| 8681176414642 | Sterilization Pouch 100 mm x 200 mm | 86811764SW2 |
| 8681176414659 | Sterilization Reel 500 mm x 200 m 70gr. 2 Indicators | 86811764SW2 |
| 8681176414666 | Sterilization Gusset Reel 100 mm x 50 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176414680 | Sterilization Gusset Reel 250 mm x 65 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176414703 | Sterilization Flat Reel 400 mm x 200 m, 70g 3 Indicators | 86811764SW2 |
| 8681176414710 | Tyvek Reel 500 mm x 70 m | 86811764SW2 |
| 8681176415113 | Sterilization Flat Reel 200 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415120 | Sterilization Flat Reel 300 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176415137 | Sterilization Gusset Reel 400 mm x 80 mm x 100 m, 60g 2 Indicators | 86811764SW2 |
| 8681176415144 | Sterilization Pouch 150 mm x 340 mm | 86811764SW2 |
| 8681176415151 | Sterilization Flat Reel 420 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415281 | Sterilization Flat Reel 75 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176415298 | Sterilization Gusset Reel 75 mm x 45 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176415304 | Sterilization Pouch 250 mm x 430 mm | 86811764SW2 |
| 8681176415311 | Sterilization Gusset Reel 420 mm x 80 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415328 | Sterilization Pouches 200 mm x 30 mm | 86811764SW2 |
| 8681176415335 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415342 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415359 | Sterilization Pouches 150 mm x 300 mm | 86811764SW2 |
| 8681176415366 | Sterilization Pouches 200 mm x 400 mm | 86811764SW2 |
| 8681176415380 | Tyvek Reel 200 mm x 70 m | 86811764SW2 |
| 8681176415397 | Tyvek Reel 75 mm x 70 m | 86811764SW2 |
| 8681176415403 | Sterilization Reel Without Print 500 mm x 200 m | 86811764SW2 |
| 8681176415410 | Sterilization Reel Without Print 200 mm x 200 m | 86811764SW2 |
| 8681176415434 | Sterilization Gusset Reel 120 mm x 55 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415441 | Sterilization Reel Without Print 50 mm x 200 m | 86811764SW2 |

| | | |
|---------------|--|-------------|
| 8681176415458 | Sterilization Flat Reel 150 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415465 | Sterilization Pouches 200 mm x 200 mm | 86811764SW2 |
| 8681176415472 | Tyvek Reel 300 mm x 70 m | 86811764SW2 |
| 8681176415489 | Sterilization Gusset Reel 420 mm x 80 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176415496 | Sterilization Reel Without Print 350 mm x 200 m | 86811764SW2 |
| 8681176415502 | Tyvek Reel 100 mm x 70 m | 86811764SW2 |
| 8681176415519 | Sterilization Flat Reel 250 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415526 | Sterilization Reel 500 mm x 200 m, 70gr. 2 Indicators | 86811764SW2 |
| 8681176415533 | Sterilization Pouches 150 mm x 150 mm | 86811764SW2 |
| 8681176415540 | Sterilization Gusset Reel 150 mm x 55 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415557 | Sterilization Reel Without Print 100 mm x 200 m | 86811764SW2 |
| 8681176415571 | Sterilization Flat Reel 350 mm x 200 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415588 | Sterilization Flat Reel 420 mm x 200 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415595 | Sterilization Flat Reel 50 mm x 200 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415601 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176415618 | Sterilization Pouches 300 mm x 400 mm | 86811764SW2 |
| 8681176415625 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415632 | Tyvek Reel 100 mm x 70 m | 86811764SW2 |
| 8681176415649 | Düz Sterilizasyon Rulosu 500 mm x 200 m, 60gr 2 Indicators | 86811764SW2 |
| 8681176415656 | Tyvek Reel 50 mm x 70 m | 86811764SW2 |
| 8681176415663 | Sterilization Flat Reel 75 mm x 200 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415670 | Sterilization Flat Reel 500 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176415687 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415694 | Sterilization Gusset Reel 300 mm x 70 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415717 | Sterilization Pouch 70 mm x 230 mm | 86811764SW2 |
| 8681176415724 | Sterilization Flat Reel 350 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415731 | Sterilization Gusset Reel 350 mm x 75 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415748 | Sterilization Gusset Reel 400 mm x 80 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415755 | Sterilization Flat Reel 450 mm x 200 m, 60g 2 Indicators | 86811764SW2 |

| | | |
|---------------|--|-------------|
| 8681176415779 | Sterilization Flat Reel 250 mm x 200 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415786 | Sterilization Pouch 250 mm x 540 mm | 86811764SW2 |
| 8681176415793 | Sterilization Flat Reel 150 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176415809 | Sterilization Flat Reel 120 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176415816 | Tyvek Reel 350 mm x 70 m | 86811764SW2 |
| 8681176415823 | Sterilization Flat Reel 100 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176415830 | Sterilization Flat Reel 120 mm x 200 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415847 | Sterilization Flat Reel 300 mm x 200 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415854 | Sterilization Flat Reel 200 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176415861 | Sterilization Flat Reel 200 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176415878 | Sterilization Gusset Reel 200 mm x 50 mm x 100 m 60 g 2 Indicators | 86811764SW2 |
| 8681176415885 | Sterilization Flat Reel 500 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176415892 | Sterilization Gusset Reel 75 mm x 45 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415908 | Sterilization Flat Reel 150 mm x 200 m, 70g 3 Indicators | 86811764SW2 |
| 8681176415915 | Sterilization Reel Without Print 150 mm x 200 m | 86811764SW2 |
| 8681176415922 | Sterilization Gusset Reel 250 mm x 65 mm x 100 m, 70g 2 Indicators | 86811764SW2 |
| 8681176415939 | Sterilization Flat Reel 300 mm x 200 m, 70g 2 Indicators | 86811764SW2 |
| 8681176415946 | Tyvek Reel 75 mm x 70 m | 86811764SW2 |
| 8681176415960 | Sterilization Gusset Reel 120 mm x 55 mm x 100 m, 70g 3 Indicators | 86811764SW2 |
| 8681176420568 | Sterilization Flat Reel 75 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420575 | Tyvek Reel 350 mm x 70 m | 86811764SW2 |
| 8681176420582 | Tyvek Reel 300 mm x 70 m | 86811764SW2 |
| 8681176420599 | Tyvek Reel 400 mm x 70 m | 86811764SW2 |
| 8681176420605 | Sterilization Flat Reel 50 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420612 | Tyvek Reel 250 mm x 70 m | 86811764SW2 |
| 8681176420629 | Tyvek Reel 200 mm x 70 m | 86811764SW2 |
| 8681176420636 | Tyvek Reel 75 mm x 70 m | 86811764SW2 |
| 8681176420643 | Sterilization Flat Reel 150 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420650 | Tyvek Reel 100 mm x 70 m | 86811764SW2 |
| 8681176420667 | Tyvek Reel 50 mm x 70 m | 86811764SW2 |

| | | |
|---------------|--|-------------|
| 8681176420674 | Sterilization Flat Reel 100 mm x 200 m , 60g 3 Ind. | 86811764SW2 |
| 8681176420681 | Tyvek Reel 150 mm x 70 m | 86811764SW2 |
| 8681176420698 | Sterilization Flat Reel 350 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420704 | Sterilization Flat Reel 300 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420711 | Sterilization Flat Reel 250 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420728 | Sterilization Flat Reel 200 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420735 | Sterilization Flat Reel 400 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420773 | Sterilization Flat Reel 380 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420872 | Sterilization Flat Reel 430 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420889 | Sterilization Flat Reel 250 mm x 100 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420902 | Sterilization Flat Reel 210 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420919 | Sterilization Flat Reel 450 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420926 | Sterilization Flat Reel 500 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420933 | Sterilization Flat Reel 430 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420940 | Sterilization Flat Reel 600 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176420957 | Sterilization Flat Reel 50 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176420964 | Sterilization Flat Reel 75 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176420971 | Sterilization Flat Reel 100 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176420988 | Sterilization Flat Reel 120 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176420995 | Sterilization Flat Reel 150 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176421008 | Sterilization Flat Reel 200 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176421015 | Sterilization Flat Reel 250 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176421022 | Sterilization Flat Reel 300 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176421039 | Sterilization Flat Reel 350 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176421046 | Sterilization Flat Reel 400 mm x 100 m 60g 3 Indicators | 86811764SW2 |

| | | |
|---------------|--|--------------|
| 8681176421053 | Sterilization Flat Reel 450 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176421060 | Sterilization Flat Reel 500 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176421077 | Sterilization Flat Reel 600 mm x 100 m 60g 3 Indicators | 86811764SW2 |
| 8681176421084 | Tyvek Reel 75 mm x 100 m | 86811764SW2 |
| 8681176421091 | Tyvek Reel 100 mm x 100 m | 86811764SW2 |
| 8681176421107 | Tyvek Reel 120 mm x 100 m | 86811764SW2 |
| 8681176421114 | Tyvek Reel 150 mm x 100 m | 86811764SW2 |
| 8681176421121 | Tyvek Reel 200 mm x 100 m | 86811764SW2 |
| 8681176421138 | Tyvek Reel 250 mm x 100 m | 86811764SW2 |
| 8681176421145 | Tyvek Reel 300 mm x 100 m | 86811764SW2 |
| 8681176421152 | Tyvek Reel 350 mm x 100 m | 86811764SW2 |
| 8681176421169 | Tyvek Reel 400 mm x 100 m | 86811764SW2 |
| 8681176421176 | Tyvek Reel 450 mm x 100 m | 86811764SW2 |
| 8681176421183 | Tyvek Reel 500 mm x 100 m | 86811764SW2 |
| 8681176421336 | Sterilization Flat Pouch | 86811764SW2 |
| 8681176421343 | Sterilization Self Sealing Pouch | 86811764SW2 |
| 8681176421350 | Sterilization Gusset Pouch | 86811764SW2 |
| 8681176421367 | Tyvek Pouch | 86811764SW2 |
| 8681176421374 | Reinforced Self Sealing Pouch | 86811764SW2 |
| 8681176421381 | Reinforced Flat Pouch | 86811764SW2 |
| 8681176421442 | Crepe Paper Sheet 350 x 350 mm | 86811764SW2 |
| 8681176421459 | Crepe Paper Sheet 400 x 400 mm | 86811764SW2 |
| 8681176421466 | SMS Paper 300 x 300 mm | 86811764SW2 |
| 8681176421473 | SMS Paper 350 x 350 mm | 86811764SW2 |
| 8681176421480 | SMS Paper 400 x 400 mm | 86811764SW2 |
| 8681176421497 | SMS Paper 750 x 750 mm | 86811764SW2 |
| 8681176421572 | Tyvek Reel 125 mm x 70 m | 86811764SW2 |
| 8681176421763 | Sterilization Flat Reel 600 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176421770 | Sterilization Flat Reel 800 mm x 200 m, 60g 2 Indicators | 86811764SW2 |
| 8681176421787 | Sterilization Flat Reel 600 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8681176421794 | Sterilization Flat Reel 800 mm x 200 m, 60g 3 Indicators | 86811764SW2 |
| 8697442750017 | Sterilization Flat Reel 50 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750024 | Sterilization Flat Reel 75 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750031 | Sterilization Flat Reel 100 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750048 | Sterilization Flat Reel 120 mm x 200 m, 60g 2 Indicators | 869744275SJD |

| | | |
|---------------|--|--------------|
| 8697442750055 | Sterilization Flat Reel 150 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750062 | Sterilization Flat Reel 200 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750079 | Sterilization Flat Reel 250 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750086 | Sterilization Flat Reel 300 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750093 | Sterilization Flat Reel 350 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750109 | Sterilization Flat Reel 400 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750116 | Sterilization Flat Reel 420 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442750123 | Sterilization Gusset Reel 75 mm x 45 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442750130 | Sterilization Gusset Reel 100 mm x 50 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442750147 | Sterilization Gusset Reel 120 mm x 55 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442750154 | Sterilization Gusset Reel 150 mm x 55 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442750178 | Sterilization Gusset Reel 250 mm x 65 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442750185 | Sterilization Gusset Reel 300 mm x 70 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442750192 | Sterilization Gusset Reel 350 mm x 75 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442750208 | Sterilization Gusset Reel 400 mm x 80 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442750215 | Sterilization Gusset Reel 420 mm x 80 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442750222 | Sterilization Flat Reel 50 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750239 | Sterilization Flat Reel 75 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750246 | Sterilization Flat Reel 100 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750253 | Sterilization Flat Reel 120 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750260 | Sterilization Flat Reel 150 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750277 | Sterilization Flat Reel 200 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750284 | Sterilization Flat Reel 250 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750291 | Sterilization Flat Reel 300 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750307 | Sterilization Flat Reel 350 mm x 200 m, 60g 3 Indicators | 869744275SJD |

| | | |
|---------------|--|--------------|
| 8697442750314 | Sterilization Flat Reel 400 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750321 | Sterilization Flat Reel 420 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442750338 | Sterilization Gusset Reel 75 mm x 45 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750345 | Sterilization Gusset Reel 100 mm x 50 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750352 | Sterilization Gusset Reel 120 mm x 55 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750369 | Sterilization Gusset Reel 150 mm x 55 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750376 | Sterilization Gusset Reel 200 mm x 60 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750383 | Sterilization Gusset Reel 250 mm x 65 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750390 | Sterilization Gusset Reel 300 mm x 70 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750406 | Sterilization Gusset Reel 350 mm x 75 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750413 | Sterilization Gusset Reel 400 mm x 80 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750420 | Sterilization Gusset Reel 420 mm x 80 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442750437 | Sterilization Flat Reel 50 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750444 | Sterilization Flat Reel 75 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750451 | Sterilization Flat Reel 100 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750468 | Sterilization Flat Reel 120 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750475 | Sterilization Flat Reel 150 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750482 | Sterilization Flat Reel 200 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750499 | Sterilization Flat Reel 250 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750505 | Sterilization Flat Reel 300 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750512 | Sterilization Flat Reel 350 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750529 | Sterilization Flat Reel 400 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750536 | Sterilization Flat Reel 420 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442750543 | Sterilization Gusset Reel 75 mm x 45 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442750550 | Sterilization Gusset Reel 100 mm x 50 mm x 100 m, 70g 2 Indicators | 869744275SJD |

| | | |
|---------------|--|--------------|
| 8697442750567 | Sterilization Gusset Reel 120 mm x 55 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442750574 | Sterilization Gusset Reel 150 mm x 55 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442750598 | Sterilization Gusset Reel 250 mm x 65 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442750604 | Sterilization Gusset Reel 300 mm x 70 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442750611 | Sterilization Gusset Reel 350 mm x 75 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442750628 | Sterilization Gusset Reel 400 mm x 80 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442750635 | Sterilization Gusset Reel 420 mm x 80 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442750642 | Sterilization Flat Reel 50 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750659 | Sterilization Flat Reel 75 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750666 | Sterilization Flat Reel 100 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750673 | Sterilization Flat Reel 120 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750680 | Sterilization Flat Reel 150 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750697 | Sterilization Flat Reel 200 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750703 | Sterilization Flat Reel 250 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750710 | Sterilization Flat Reel 300 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750727 | Sterilization Flat Reel 350 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750734 | Sterilization Flat Reel 400 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750741 | Sterilization Flat Reel 420 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442750758 | Sterilization Gusset Reel 75 mm x 45 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442750765 | Sterilization Gusset Reel 100 mm x 50 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442750772 | Sterilization Gusset Reel 120 mm x 55 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442750789 | Sterilization Gusset Reel 150 mm x 55 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442750796 | 200 mm x 60 mm x 100 m, 70g 3 Ind. Sterilization Gusset Reel | 869744275SJD |
| 8697442750802 | Sterilization Gusset Reel 250 mm x 65 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442750819 | Sterilization Gusset Reel 300 mm x 70 mm x 100 m, 70g 3 Indicators | 869744275SJD |

| | | |
|---------------|--|--------------|
| 8697442750826 | Sterilization Gusset Reel 350 mm x 75 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442750833 | Sterilization Gusset Reel 400 mm x 80 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442750840 | Sterilization Gusset Reel 420 mm x 80 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442750949 | Tyvek Sterilization Reel 500 mm x 70 m | 869744275SJD |
| 8697442750970 | Crepe Paper Sheet 500 x 500 mm | 869744275SJD |
| 8697442750987 | Crepe Paper Sheet 600 x 600 mm | 869744275SJD |
| 8697442750994 | Crepe Paper Sheet 750 x 750 mm | 869744275SJD |
| 8697442751007 | Crepe Paper Sheet 900 x 900 mm | 869744275SJD |
| 8697442751014 | Crepe Paper Sheet 1000 x 1000 mm | 869744275SJD |
| 8697442751021 | Crepe Paper Sheet 1200 x 1200 mm | 869744275SJD |
| 8697442751038 | Sterilization Reel Without Print 50 mm x 200 m | 869744275SJD |
| 8697442751045 | Sterilization Reel Without Print 75 mm x 200 m | 869744275SJD |
| 8697442751052 | Sterilization Reel Without Print 100 mm x 200 m | 869744275SJD |
| 8697442751069 | Sterilization Reel Without Print 120 mm x 200 m | 869744275SJD |
| 8697442751076 | Sterilization Reel Without Print 150 mm x 200 m | 869744275SJD |
| 8697442751083 | Sterilization Reel Without Print 200 mm x 200 m | 869744275SJD |
| 8697442751090 | Sterilization Reel Without Print 250 mm x 200 m | 869744275SJD |
| 8697442751106 | Sterilization Reel Without Print 300 mm x 200 m | 869744275SJD |
| 8697442751113 | Sterilization Reel Without Print 350 mm x 200 m | 869744275SJD |
| 8697442751120 | Sterilization Reel Without Print 400 mm x 200 m | 869744275SJD |
| 8697442751137 | Sterilization Reel Without Print 420 mm x 200 m | 869744275SJD |
| 8697442754978 | Tyvek Sterilization Reel 75 mm x 70 m | 869744275SJD |
| 8697442754985 | Tyvek Sterilization Reel 100 mm x 70 m | 869744275SJD |
| 8697442755005 | Tyvek Sterilization Reel 150 mm x 70 m | 869744275SJD |
| 8697442755012 | Tyvek Sterilization Reel 200 mm x 70 m | 869744275SJD |
| 8697442755029 | Tyvek Sterilization Reel 250 mm x 70 m | 869744275SJD |
| 8697442755036 | Tyvek Sterilization Reel 300 mm x 70 m | 869744275SJD |
| 8697442755043 | Tyvek Sterilization Reel 350 mm x 70 m | 869744275SJD |
| 8697442755050 | Tyvek Sterilization Reel 400 mm x 70 m | 869744275SJD |
| 8697442755128 | Crepe Paper Sheet 150 x 200 mm | 869744275SJD |
| 8697442755685 | Tyvek Reel 50 mm x 70 m | 869744275SJD |
| 8697442755692 | Tyvek Reel 75 mm x 70 m | 869744275SJD |
| 8697442755708 | Tyvek Reel 100 mm x 70 m | 869744275SJD |
| 8697442755715 | Tyvek Reel 120 mm x 70 m | 869744275SJD |
| 8697442755722 | Tyvek Reel 150 mm x 70 m | 869744275SJD |
| 8697442755739 | Tyvek Reel 200 mm x 70 m | 869744275SJD |
| 8697442755746 | Tyvek Reel 250 mm x 70 m | 869744275SJD |
| 8697442755760 | Tyvek Reel 350 mm x 70 m | 869744275SJD |
| 8697442755777 | Tyvek Reel 400 mm x 70 m | 869744275SJD |
| 8697442757092 | Sterilization Flat Reel 450 mm x 200 m, 60g 2 Indicators | 869744275SJD |
| 8697442757108 | Sterilization Flat Reel 500 mm x 200 m, 60g 2 Indicators | 869744275SJD |

| | | |
|---------------|--|--------------|
| 8697442757115 | Sterilization Gusset Reel 450 mm x 85 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442757122 | Sterilization Gusset Reel 500 mm x 85 mm x 100 m, 60g 2 Indicators | 869744275SJD |
| 8697442757139 | Sterilization Flat Reel 450 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442757146 | Sterilization Flat Reel 500 mm x 200 m, 60g 3 Indicators | 869744275SJD |
| 8697442757153 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442757160 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 60g 3 Indicators | 869744275SJD |
| 8697442757177 | Sterilization Flat Reel 450 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442757184 | Sterilization Flat Reel 500 mm x 200 m, 70g 2 Indicators | 869744275SJD |
| 8697442757191 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442757207 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 70g 2 Indicators | 869744275SJD |
| 8697442757214 | Sterilization Flat Reel 450 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442757221 | Sterilization Flat Reel 500 mm x 200 m, 70g 3 Indicators | 869744275SJD |
| 8697442757238 | Sterilization Gusset Reel 450 mm x 80 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442757245 | Sterilization Gusset Reel 500 mm x 80 mm x 100 m, 70g 3 Indicators | 869744275SJD |
| 8697442757283 | Sterilization Reel Without Print 450 mm x 200 m | 869744275SJD |
| 8697442757290 | Sterilization Reel Without Print 500 mm x 200 m | 869744275SJD |

Additional Information:

Development, production and distribution are supported by a Quality Management System according to ISO 9001:2015 AND ISO 13485:2016 requirements. The Quality Management System has been certified by the notified body Royalcert (Certificate No: 00108/ DÖR09B and No: 00108/ DÖR13A). As proof of compliance of the products with the requirements of Medical Device Regulation 2017/745 on Medical devices, the products are considered Class I nonsterile products.

Ankara Turkey

04.08.2021

Canan ÖKTEM

General Manager

4a medical Dört A Tıp Matzemeleri
 producing health for the world Sanayi İth.İhr. Tic. Ltd. Şti.
 Balneacı Mahallesi Köyü Serpmeçeri No: 795/A
 Tel: 0312 363 50 52-53 Faks: 363 52 10 AKYUT-ANKARA
 Çukura Verki Dairesi: 313 006 8060 İşletme Sicil No: 795
 NPTC: 4a-medical@4aocli.com • orkara@4amedical.com
 Mersis No: 0313 006 0600 0018