

We Pharmaplast S.A.E., hereby declare that Pharmapore® dressing consist of a non-woven absorbent pad placed centrally on a polyester backing material, wound contact layer is a non-adherent polyethylene film and the adhesive is water-resistant, latex free, zinc oxide free, rubber free and X-ray transparent. It can also stay up to 48 hours.

Name: Nourhan Hossameldin

Position: Regulatory Affairs Specialist

Signature: *Nourhan Hossam Eldin*

Date: 16.09.2023

PHARMA PLAST S.A.E
ALEXANDRIA - EGYPT