

# 1.4 Configuration

## 1.4.1 Standard Components

SONIALVISION G4 comprises the following parts.

Unit	Name	Component
X-ray diagnostic table	Remote controlled X-ray diagnostic table ZS-200	X-ray diagnostic table body <sup>*1</sup> Control cabinet (for X-ray diagnostic table) Remote control panel Anti-scatter grid <sup>*1</sup> Tabletop mat <sup>*1,*2</sup> Foot rest <sup>*1,*2</sup> , Hand grips <sup>*1,*2</sup> , Upper hand grip <sup>*1,*2</sup> , Shoulder rests <sup>*1,*2</sup>
Digital radiography unit	Digital Radiography System DR-300	Control cabinet (for digital radiography unit) Operation cabinet FPD <sup>*1,*3</sup> Acquisition monitor, Reference monitor Touch panel <sup>*1,*4</sup> Keyboard, Mouse
High-voltage generator	Diagnostic X-ray High-Voltage Generator D150BC-40S (D150VC-40S<Not sold in the US>)	Control cabinet (for high-voltage generator) Hand switch Phototimer receiver <sup>*1,*3</sup>
X-ray tube unit	Medical X-Ray Tube Unit 0.7/1.2JG326D-265 (0.4/0.7JG326D-265< Not sold in the US>)	X-ray tube <sup>*1,*3</sup>
Collimator	Collimator R-300	Collimator <sup>*1,*3</sup>

\*1: Equipment suitable for use in patient environment

\*2: The patient's body comes into contact with these parts (applied parts).

\*3: Incorporated in the X-ray diagnostic table

\*4: Incorporated in the remote control panel