

## Mobile 2 RPW



# Mobile 2 RPW

# Chattanooga's Radial Pressure Wave technology is mobile.

The Mobile 2 RPW is compact and light enough to be easily carried while also delivering a host of features and benefits.

When there is need for a powerful yet portable Radial Pressure Wave device in the clinic or on the sports field, the Chattanooga Mobile 2 RPW is the ideal choice. A lightweight and compact design together with a quiet built-in compressor make it perfect for portable use. A user-friendly experience is at the heart of the Mobile 2 RPW's design. The simple, ergonomically laid out controls enable all of the key parameters to be easily changed, while a clear, illuminated LED display removes ambiguity when monitoring pressure, frequency and number of pulses. Pulses can be initiated via a trigger mechanism on the handpiece.

## User Profile

This Radial Pressure Wave clinical unit is ideally suited to the following applications:

- OrthopaedicsRheumatology
- Physiotherapy Sport Clubs

• Acupuncture

AestheticVeterinary

Podiatrists

Chiropractors

## Indications

Physical Medicine

- Myofascial pain
- Muscle tendon tissue activation
- Tendon attachment & site illness
- Acupuncture Shockwave Therapy

The Mobile RPW D-Actor Applicator allows quick and efficient uses of applicator transmitters for desired therapeutic effects during treatment.

#### D-ACTOR® Applicator

#### Applications:

- Myofascial trigger points
- Other disorders of tendon insertions
- Activation of muscle and connective tissue
- Acupuncture shock wave therapy
- Cellulite treatment (stages I III)
- Muscle and connective tissue tightening
- Elimination of lymphatic congestion
- Body shaping

- Regulation of muscle tone
- Medical massage



## Facts and figures


- Weight of the unit: Approx. 9 kg
- Dimensions (L x W x H): 32 x 32 x 21 cm
- Built in compressor
- Single and Continuous shocks Frequency D-actor: 1 - 21 Hz Frequency V-actor: 1 - 31 Hz
- Application pressure: 1 5 Bar Highly customisable with precise adjustments of 0.1 increments.
- Oscillating D-Actor
- V-Actor Hand piece option
- USB connectivity



#### User Interface Features

Hand piece docking port Chattanooga **₩**1 bar Ultrasound gel holder 1. Pressure screen 1 - 5 Bar 3 侣 Hz 6 2. Pulse counter screen 3. Frequency screen 4. Intensity adjustment 5. Pulse counter reset 6. Frequency adjustment

#### **D-ACTOR®** Transmitters





## Mobile 2 RPW

### Ordering Information

Part Number	Description
2910	Mobile 2 RPW

Standard Accessories (included with Mobile 2 RPW as default)

2905	Chattanooga Mobile RPW 2 Unit Only
32000	R-SW Hand Piece Applicator Set
28178	Ro40 'Beam' transmitter
29724	D20-S D-Actor <sup>®</sup> 20mm Transmitter
13-5182	Conductor <sup>™</sup> Transmission Gel 250 ml
13-00060	CD User Manual

#### Optional Accessories

	pelonal recessories
28701	Revision Kit
28178	Ro40 'Beam' transmitter
31800	9 Transmitter Set In Case, Accessory
29801	A6 Point, Transmitter
29802	T10 Finger, Transmitter
29729	C15, 15 mm Cerama-X Transmitter, Black
29728	DI15, 15 mm, Deep Impact Transmitter, Black
29726	F15 Focus-Lens 15 mm Transmitter
29724	D20-S D-Actor <sup>®</sup> 20mm Transmitter
28736	D20-T Transmitter
29725	D35-S, 35 mm, D-Actor Transmitter, Black
22375	Handpiece Cushion, Mobile 2 RPW
22376	Shaft Cushion, Mobile 2 RPW
28730	V-Actor Set
28740	V-Actor V25 Transmitter
28741	V-Actor V35 Transmitter
28745	Carrying Bag





DJO Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM www.DJOglobal.eu