

Declaration of Connectivity

Document No: MMDR-200047

Manufacturer : Shimadzu Corporation
1 Nishinokyo Kuwabara-cho, Nakagyo-ku, Kyoto 604-8511, Japan
Category : General radiographic system
Product Name : RADspeed Pro, RADspeed fit

We declare that the above-mentioned medical device, when installed, maintained and used in accordance with its installation instructions, maintenance instructions/recommendations and instructions for use is compatible with the product listed below.

<Tested Product>

Manufacturer	Vieworks
Product Name	VXvue
Software version	1.0.2.3

< General radiographic system>

System	Software version
RADspeed Pro LCD type generator	Version 2.690 or later is necessary. *Test was performed with version 2.710
RADspeed Pro Sheet Key type generator	Version 2.120 or later is necessary. *Test was performed by substituting LCD type generator.
RADspeed fit	Version 1.12 or later is necessary. (for U.S.) Version 1.11 or later is necessary. (Others) *Test was performed by substituting LCD type generator.

Test report No. : EFGQ-7185
Test result : The same connection as CANON CXDI-NE except for long view imaging was confirmed.
Document for installation : M517-E809J RADspeed Pro BK/BR series Service Manual
M503-E368H RADspeed fit Installation Manual

Signed: 2020/10/08

Shimadzu Corporation



Printed Name: Yoshiaki Miura

Title: General Manager, Global Marketing Dept.