








CODE 27557

USER'S MANUAL – DONOR CHAIR – blue



SYMBOLS:

Symbol	Symbol definition
	Carefully read the instructions for use
	Keep away from sunlight
	Store in a dry place
	Product in compliance with European Directive
	Lot number
	Manufacturer address
	Product code

I. GENERAL FEATURES

The chair is made of oval, square and round profile section metal.

The plan for use made of high-density fire-retardant class 1 IM foam is covered with coated (vinyl) fire-retardant class 1 IM fabric. The synchronized control of the two cushions is adjustable by gas lift (600 N). The backrest angle is +60° and about -5° in the Trendelenburg position; the footboard angle is -80° and about -5°.

The armrests are equipped with an articulation that guarantees a perfect angle during use; they can be adjusted in height and, when removed, they can be rotated of 180°. There are four adjustable feet ensuring optimal stability.

II. CAPACITY

130 Kg.

III. CLEANING

Use common neutral detergents on the market.

WARRANTY CONDITIONS APPLIED BY GIMA

We congratulate you for purchasing our product.

This product meets the highest quality standards in terms of material and manufacture.

The warranty is valid for 12 months from the date of supply by GIMA.

During the period of validity of the warranty, the company will take care of free repair and/or replacement of any defective parts for manufacture reasons that are well detected, with the exception of the cost of labour or any transfers, transport and packaging. The warranty doesn't include all components subject to wear.

The replacement or repair during the warranty period will not extend the duration of the warranty.

The warranty is not valid in the event of: repairs performed by unauthorized personnel or use of non-original spare parts and defects caused by negligence or misuse.

The warranty is void if the above is not respected and if the serial number (if any) has been removed, deleted or altered.

The defective products must be returned only to the retailer where the product was purchased.

Products sent to GIMA will be rejected.



GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy