

Chapter 25. Accessories

Warning

- Use the accessories specified in the Operator's Manual; using other accessories may damage the monitor, or cannot reach the performance claimed in this Manual.
- The operating and storage environment of the monitor should meet the requirements of the accessories. Please refer to the manual of the accessories for these requirements.
- Disposable accessories can only be used once, because repeated use can cause performance degradation.
- If the packaging or accessories have any sign of damage, do not use such accessories.

25.1 Standard Accessory List

Standard accessories are as follows:

No.	Name	Specifications & characteristics	Standard quantity
1	5-lead ECG cable	Adult 5-lead all-in-one clip	1
2	Electrode slice	50mm×55mm, adult, disposable	1
3	Oxygen probe	Adult finger clip	1
4	Blood pressure cuff	Adult repetitive single-tube cuff	1
	Blood pressure gas nozzle	Metal	1
	Blood pressure extension tube	3M gray	1
5	2.25K-type temperature probe	Adult body surface	1
6	2600mAh lithium battery	Detachable three-cell lithium battery (14.8V 2600mAh)	1
7	Ground wire	V0020A	1



No.	Name	Specifications & characteristics	Standard quantity
8	Power cord	3×0.75mm ² 10A 250V L=2M	1

25.2 ECG

12-pin integrated 5-lead left/right clamp ECG cable

Supported lead wires	Symbol standard	Type	Applicable objects	Length (main cable)	Length (5-lead ECG cable)
3/5-lead wire	IEC	Defibrillation type	Adult	2.7m	0.9m

25.3 SpO₂


Adult finger clip oxygen probe

Applicable objects	Applicable parts	Type	Length
Adult	Fingers, toes	Reusable	3m

25.4 NIBP

1. Blood pressure cuff

Applicable objects	Applicable parts	Type	Remark
Adult	Arm or thigh	Reusable	Single tube



2. Blood pressure gas nozzle

Operating temperature	Applicable parts	Mechanical life (connecting and disconnecting)	Remark
-30°C~+60°C	Blood pressure cuff	> 1000 times	Pull-off force of gas connector and extension tube is greater than 5kg

3. Blood pressure extension tube

Applicable objects	Type	Length	Remark
Adult, child	Reusable	3m	Used together with cuff

25.5 Temp

2.25K-type temperature probe

Applicable objects	Applicable parts	Type	Length
Adult	Body surface	Reusable	3m

