



EC Declaration of Conformity to: Medical Devices Directive 93/42/EEC

Legal Manufacturer:	BD Switzerland Sàrl Route de Crassier 17, Business Park Terre-Bonne, Batiment A4, 1262 Eysins, Switzerland
EU Authorised Representative:	Becton Dickinson Ireland Ltd. Donore Road, Drogheda Co. Louth, A92 YW26 Ireland
Manufacturing Site (s):	Plexus RO S.R.L Eugeniu Carada Street, no 2-4, Oradea, 410610, Bihor, Romania
Device Description/Family:	Alaris™ GH Syringe Pump <i>(See attached Product Schedule)</i>
EC Product Classification:	Class IIb, Annex IX, Rule 11
GMDN:	<i>13217 – Syringe Pump</i> A mains electricity (AC-powered) device designed to precisely drive the plunger of a syringe down its barrel to infuse a solution when it must be administered with a high degree of volume accuracy and rate consistency. Because of the lower flow settings and flow resolution (e.g., 0.1 ml/hr), it is especially appropriate for neonatal, infant, and critical care applications in which small volumes of concentrated drugs are to be delivered over an extended period. It can also be used to administer epidural analgesia. It will typically have internal batteries that allow the device to operate for a short period of time when no line power is available (e.g., during transport or a power outage).

We herewith declare that the product(s) listed above and detailed in the attached Product Schedule meet the provisions of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

Unofficial Copy if Printed – See Intranet for Latest Released Version	
For internal use only. This document contains confidential, proprietary information of CareFusion or one of its subsidiaries. It may not be copied or reproduced without prior written permission from CareFusion.	Page 1 of 4



Applied Standards and Directives:	<ul style="list-style-type: none"> - <i>Medical Device Directive 93/42/EEC</i> - <i>EMC Directive 2014/30/EU</i> - <i>RoHS Directive 2011/65/EU with amendment: Commission Delegated Directive 2015/863</i> - <i>Machinery Directive 2006/42/EC</i> - <i>Waste electrical and electronic equipment 2002/96/EC</i> - <i>Product Liability 85/374/EEC</i> - <i>REACH 1907/2006</i> - <i>Packaging and Packaging Waste Directive 94/62/EC</i> - <i>Battery 2006/66/EC</i> - <i>Electronic Instructions for Use of Medical Devices 207/2012</i> - <i>EN ISO 13485:2016</i> - <i>EN ISO 14971:2012</i> - <i>IEC 60601-1:2005 +A1:2012</i> - <i>IEC 60601-1-2:2014</i> - <i>IEC 60601-1-6:2010 +A1:2013</i> - <i>IEC 62366:2007 +A1:2014</i> - <i>IEC 60601-1-8:2006 +A1:2012</i> - <i>IEC 60601-2-24:2012</i> - <i>IEC 62304:2006</i> - <i>EN ISO 15223-1:2016</i> - <i>EN 1041:2008+A1:2013</i> - <i>IEC 60529:1989+A1:1999+A2:2013</i> - <i>EN 1789:2007 +A2:2014</i> - <i>ISTA-1A-2014 Edition</i>
Notified Body:	<p>BSI Group The Netherlands B.V. Say Building, John M. Keynesplein 9, 1066 EP Amsterdam. Notified Body Number: 2797</p> <p>Former Notified Body: BSI Group, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes, MK5 8PP. Notified body number 0086</p>
CE Certificate Number:	<i>Annex II (EC Certificate No. 502238)</i>
Date of issuance of original CE certificate:	16 November 2005



STED File: 009

Issue Level: 25

Signed:

Ondina Bennaim

DocuSigned by:
Ondina Bennaim
 Signer Name: Ondina Bennaim
Signing Reason: I approve this document
Signing Time: 08-Aug-2023 | 2:39:46 AM PDT
0EF48A97D3A94213AA44C0E205189A5C

Date: 08-Aug-2023

Senior Director Regulatory Affairs
International Infusion

Unofficial Copy if Printed – See Intranet for Latest Released Version

For internal use only. This document contains confidential, proprietary information of CareFusion or one of its subsidiaries. It may not be copied or reproduced without prior written permission from CareFusion.

Page
3 of 4



Product Schedule Alaris™ GH Syringe Pump

GMDN Number: 13217

Part Number	Description	EC Product Class
8002TIG03	Alaris™ GH Plus Syringe Pump	IIb
8002TIG03-G	Alaris™ GH Plus Guardrails™ Syringe Pump	IIb
8002TIG01-CHINA	Alaris™ GH Plus Syringe Pump For China	IIb

Unofficial Copy if Printed – See Intranet for Latest Released Version

For internal use only. This document contains confidential, proprietary information of CareFusion or one of its subsidiaries. It may not be copied or reproduced without prior written permission from CareFusion.

Page
4 of 4