




exeol surf optimal



Alcohol-free disinfectant detergent for surfaces and medical devices

► Ready-to-use solution:

- **Medical device application:** for cleaning and disinfection by wiping of surfaces of invasive and non-invasive reusable medical devices.
- **Biocidal product application:** for cleaning and disinfection by wiping of surfaces, furniture and heavy equipments.
- Original synergistic combination of two microbial active ingredients (SCALE effect).
- Active from 2 min. (bactericidal, yeasticidal and virucidal against enveloped viruses).
- Cleanability effect enhanced by a surfactant selected to confer a low surface tension ($ST < 27 \text{ mN/m}$).
- Broad compatibility with materials.
- Suitable for surfaces in contact with foodstuffs.



BIOCIDAL PT2/PT4

MD CLASS IIa



Microbiological properties: Bactericidal, active against MRB *Escherichia coli*, *Klebsiella pneumoniae* (Carbapenem-resistant) and *Staphylococcus aureus* (Methicillin-resistant), yeasticidal and active against *Candida auris*, active against *Aspergillus fumigatus*, tuberculocidal, virucidal against enveloped viruses, active against Norovirus and Adenovirus.


> Box 12 x 750mL spray: EXS0038

SURFACES CLEANING AND DISINFECTION



*except for surfaces in direct contact with skin, mucous membranes, foodstuffs or if a subsequent disinfection step is required.

exeol wipes optimal



Disinfectant detergent wipes for surfaces and medical devices, alcohol-free

► Pre-impregnated wipes:

- **Medical device application:** for cleaning and disinfection by wet wiping of surfaces of invasive and non-invasive reusable medical devices.
- **Biocidal product application:** for cleaning and disinfection by wet wiping of surfaces, furnitures and heavy equipments.



BIOCIDAL PT2/PT4

MD CLASS IIa



- Original synergistic combination of two microbial active ingredients (SCALE effect).
- Active from 2 min. (bactericidal, yeasticidal and virucidal against enveloped viruses).
- Cleanability effect enhanced by a surfactant selected to confer a low surface tension ($ST < 27 \text{ mN/m}$).
- Broad compatibility with materials.
- Suitable for surfaces in contact with foodstuffs.

Microbiological properties of the liquid extracted from the wipes: Bactericidal, yeasticidal, active on *Candida auris*, virucidal against enveloped viruses, active on Norovirus and Adenovirus.

> Box 12 x 100 wipes (180x200 mm): EXS0041