

AKASPRAY TÜCHER

Disinfection Wipes for Surfaces – Alcohol Based



USAGE AREAS:

- Used for surface disinfection in every place where health services are provided such as operating rooms, intensive care, dental units, doctor's office, veterinary clinic, laboratory and public areas such as home, office, factory, school etc.

PRODUCT PROPERTIES:

- Ready-to-use product.
- Does not contain phenol and aldehyde.
- Effective in 30 seconds
- Dries quickly and is pleasantly fragrant.
- High material compatibility
- Biocidal product license dated 08.06.2018 and numbered 2018/60, obtained from the Turkish Ministry of Health

COMPOSITION:

- Ethanol, 2-propanol, didecyldimethylammonium propionate

MICROBIOLOGICAL ACTIVITY:

- Bactericidal, fungicidal, virucidal, tuberculocidal in 30 seconds

USAGE:

- Sufficient amount of wipes are taken according to the size of the equipment to be wiped. It is applied by wiping method as a surface disinfectant. It is waited for 30 seconds to take effect.

PACKAGE:

QUANTITY	PACKAGING TYPE	BARCODE
100 pcs	12 x CYLINDER BOX	8680097103574



SURFACE
DISINFECTANTS

detrox[®]
Integrated Disinfection Solutions