




Product Catalog

TIANDA MEDICAL(HUAIAN)CO., LTD.
No.106 East Songjiang Road,Huaiyin Economic&Technological
Development Zone 223002 Huaian City,Jiangsu China
Tel :+86 517 84379688
Fax :+86 517 84379788
www.tianda-medical.com



Company Introduction

TianDa Medical(HuaiAn) Co Ltd.established in 2006,which is specialized in diabetic care and surgical medical products,we dedicate to developing,manufacturing and sales of medical instruments, we regard quality management as the core key,meanwhile we have passed CE,ISO13485 and FDA successfully.

Over ten years of unremitting efforts,our products have exported all over the world and widely recognized by customers,we will continue to forge ahead with our philosophical philosophy'work together,improve together'to keep improving the stability of quality and our own competitiveness,we will try our best to be the most trustworthy medical supplier in the future.

Enterprise tenet : Changes for Better ,Creating a better future

Business Philosophy : Work together,Improve together

Employing philosophy : We put talented and virtuous person in the important position; we put virtuous but talentless person in the normal position;we never use the person just have talented.

Certificate



Blood Lancets

- Material:LDPE,Stainless Needle
- Size:21g,23g,26g,28g,30g
- Package:100 pcs/box 20000 pcs/carton
- Feature:Triple bevel needle tip,automatic production
- Sterilization:Gamma Irradiation



REF NO	0121	0123	0126	0128	0130
SIZE	21G	23G	26G	28G	30G
COLOR					

Steel Lancets

- Material:Stainless Steel
- Size: S/30mm,L/40mm
- Package:200 pcs/box
20000 pcs/carton
- Feature:Tungsten steel mold,high
quality needle



REF NO	0230	0240
SIZE	30mm	40mm

Alcohol Pads

- Material:Isopropyl,nonwovens
- Size:30x65mm,30x30mm,65X56mm
45X90mm
- Package:100 pcs/box 10000 pcs/carton
- Feature:70% isopropyl,eternal use



REF NO	0301	0302	0303	0304
SIZE	30x65	30x30	65X60	45X90 (mm)



Safety Lancets PA&PAII

- Material:ABS,LDPE,Spring,Stainless Needle
- Size:21Gx2.2mm,21Gx2.4mm,23Gx1.8mm,26Gx1.8mm,28Gx1.8mm,30Gx1.8mm
- Package:100 pcs/box 2000 pcs/carton
- Feature:Safety,press to use,triple bevel needle tip,automatic assembly
- Sterilization:Gamma Irradiation



REF NO	042122	042124	042318	042618	042818	043018	PA
	052122	052124	052318	052618	052818	053018	PAII
SIZE	21Gx2.2	21Gx2.4	23Gx1.8	26Gx1.8	28Gx1.8	30Gx1.8	(mm)
COLOR							

Lancing Device

- Material:, Fine copper, Spring
- Size: 138mm
- Package:1 pcs/box 500 pcs/carton
- Feature:Continuous button,rapid rebound, 1 million times service life



REF NO	0601
SIZE	138mm



Surgical Blades



- Material:Carbon,stainless
- Size:10#.11#,12#,15#,20#,21#,22#,23#,24#
- Package:100 pcs/box 5000 pcs/carton
- Feature:automatic grinding keep sharpness, Tungsten steel mold
- ensure stable dimension
- Sterilization:Gamma Irradiation

REF NO	<u>0710</u>	<u>0711</u>	<u>0712</u>	<u>0715</u>	<u>0720</u>	<u>0721</u>	<u>0722</u>	<u>0723</u>	<u>0724</u>	carbon
	0810	0811	0812	0815	0820	0821	0822	0823	0824	stainless
SIZE	10#	11#	12#	15#	20#	21#	22#	23#	24#	

Disposable Scalpels



- Material:Carbon,stainless
- Size:10#,11#,12#,15#,20#,21#,22#,23#,24#
- Package:10 pcs/box 500 pcs/carton
- Feature:Easy for use,blades are made by our own factory
- Sterilization:Gamma Irradiation

REF NO	0910	0911	0912	0915	0920	0921	0922	0923	0924	carbon
	1010	1011	1012	1015	1020	1021	1022	1023	1024	stainless
SIZE	10#	11#	12#	15#	20#	21#	22#	23#	24#	

Mini Scalpels I

- Material:ABS,carbon,ctainless
- Size: 10#,11#,12#,15#
- Package:50 pcs/box 2000 pcs/carton
or bulk packing
- Feature:Popular in Kit,electron industry



REF NO	<u>1410</u> 1510	<u>1411</u> 1511	<u>1412</u> 1512	<u>1415</u> 1515	carbon stainless
SIZE	10#	11#	12#	15#	

Mini Scalpels II

- Material:ABS,carbon,stainless
- Size: 10#,11#,12#,15#
- Package:50 pcs/box 2000 pcs/carton
or bulk packing
- Feature:Petroleum jelly as surface coating,
castrate animals.



REF NO	<u>1610</u> 1710	<u>1611</u> 1711	<u>1612</u> 1712	<u>1615</u> 1715	carbon stainless
SIZE	10#	11#	12#	15#	

Stitch Cutter Blades

- Material:Carbon,stainless
- Size: 110mm,65mm
- Package:100 pcs/box 5000 pcs/carton
- Feature:cut the wide range of suture on the wound



REF NO	<u>1801</u> 1901	<u>1802</u> 1902	carbon Stainless
SIZE	110mm	65mm	

Skin Graft Blades

- Material:Carbon,stainless
- Size: 158mm
- Package:100 pcs/box,1000 pcs/carton
- Feature:use for skin removal or transplant during surgery



REF NO	<u>2001</u> 2101	carbon stainless
SIZE	158mm	

Umbilical Cord Clamps

- Material:PE
- Size:50mm,55mm
- Package:100 pcs/box,5000 pcs/carton
- Feature:No toxic,no allergic, one injection mold
- Sterilization:Ethylene oxygen



REF NO	1101	1102
SIZE	50mm	55mm
COLOR		

Identification Bracelet


- Material:PVC
- Size: Adult,child,Mother-Child
- Package:100 pcs/box 5000 pcs/carton
- Feature:Soft and comfortable wear,no toxic,no allergic,water resistance




REF NO	1201 1203	1202 1204
SIZE	Adult Card insert Adult Hand writing	Child Card insert Child Hand writing
COLOR		



Tourniquet




- Material:TBE
- Size: 450x25x0.6 mm
- Package:25 pcs/box 100 pcs/carton
- Feature:Soft and latex free,degradable,single use only

REF NO	1301				
SIZE	450x25x0.6				
COLOR					

Urinary Drainage Bags

- Material:PVC,ABS,PE
- Size:2000ml,2600ml with Air-Vent
- Package:250 pcs/carton
- Feature:Double hook,bed sheet clamp,
Needleless sampling port with
cap,Rigid dog House
- Sterilization:Ethylene oxygen




REF NO	2001	2002	2003	2004
Volumes	2000ml	2000ml	2000ml	2600ml
Size	T VALUE	Pull Value	Cross Value	Air Vent Type

Lancing Device



- **Material:**
ABS, Spring
- **Size:**
100 mm
- **Package:**
1pcs/box,200 pcs/carton
- **Feature:**
9 adjustable depth penetration
100,000 times life cycle
easy for use and carry
Rapid launch and rebound for use
Rapid launch and rebound for use



REF NO	SIZE (mm)	COLOR
0601	100mm	BLUE/WHITE/ROSE

How To Use



1 Remove the lancing device cover.



2 Place a lancet firmly into the carrier.



3 Twist off the protective cap of lancet.



4 Reload the cover on the lancing device.



5 Adjust the appropriate penetration depth. Slide the end barrel until you hear a click, then release.



6 Place the cap on the chosen site, then press the trigger button to fire.



7 Remove the cover again and insert the needle into the protective cap.



8 Slide the ejector to dispose of the lancet and put it on the cap for next use.