



## MacConkey Broth (Double strength) w/ Neutral Red

M539

MacConkey Broth (Double strength) w/ Neutral Red is recommended for the primary isolation of coliforms from large samples such as water and wastewater.

### Composition\*\*

Ingredients	Gms / Litre
Peptic digest of animal tissue	40.000
Lactose	20.000
Bile salts	10.000
Sodium chloride	10.000
Neutral red	0.150
Final pH ( at 25°C)	7.4±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 80.15 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Distribute into test tubes with inverted Durham tubes and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool the tubes before inoculation.

### Principle And Interpretation

MacConkey Broth is widely used as a differential medium for detection and enumeration of coliforms from a wide variety of clinical, food and water samples. Identification is based on colour change of the medium due to the presence of the indicator neutral red (1, 2). Peptic digest of animal tissue provides necessary nitrogen source. Lactose serves as the fermentable carbohydrate source. Sodium chloride maintains the osmotic balance of the cells. The selective action of these media is attributed to the presence of bile salts, which are inhibitory to most species of gram-positive bacteria. Gram-negative bacteria usually grow well on these media and are differentiated by their ability to ferment lactose. The colour change of the medium shown by lactose-fermenters is due to production of acid from lactose and a subsequent colour change of the indicator dye when the pH of the media falls below 6.8. Lactose non-fermenting strains, such as *Shigella* and *Salmonella* do not alter the appearance of the media.

MacConkey Broth (Double Strength) w/ Neutral Red M539 is recommended for the primary isolation of coliforms from large samples such as water and wastewater. The medium turns pink in case of lactose fermentors and yellow in case of lactose- non-fermenters, due to neutral red. MacConkey Broth Double Strength w/Neutral Red M539 has the same composition in double strength to that of MacConkey Broth (M007), which contains neutral red as an indicator and is considered as a standard medium for the primary isolation as well as presumptive identification of coliform-aerogenes group of organisms in food and water.

### Quality Control

#### Appearance

Light yellow to pink homogeneous free flowing powder

#### Colour and Clarity of prepared medium

Red coloured clear solution without any precipitate

#### Reaction

Reaction of 8.02% w/v aqueous solution at 25°C. pH : 7.4±0.2

#### pH

7.20-7.60

#### Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours .

#### Cultural Response

---

Organism	Inoculum (CFU)	Growth	Acid	Gas
<b>Cultural Response</b>				
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	positive reaction	positive reaction
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	positive reaction	positive reaction
<i>Salmonella Choleraesuis</i> ATCC 12011	50-100	fair - good	negative reaction	negative reaction
<i>Staphylococcus aureus</i> ATCC 25923	$\geq 10^3$	inhibited		

### Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

### Reference

1. MacConkey, 1900, The Lancet, ii:20.
2. MacConkey, 1905, J. Hyg., 5:333.

Revision : 2 / 2015

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.