Disposable Skin Stapler and Remover



Catalogue No.	Cartridge Specifications (mm)	Staple Number	Staple Diameter (mm)	Staple Remover
PFW	7.0 × 4.5	15/25/35/45/55	0.6	PQW

Disposable Clip Applier



Unit: mm

Model	Size	Clip Quanity	Length		Diameter		Clip height		Clip width	
TQ-10 10	10	10 20PCS	General length	Max. deviation	Diameter	Max. deviation	Size	Max. deviation	Size	Max. deviation
			471	±10	10	±1	7.5	±0.2	5.5	±0.12

Super interlock clips

• More secure fixation of clip to tissue provides greater ligation security.

Glare-resistant shaft and black jaw with clear clip window

• Matte gray shaft minimizes endoscope glare, while clear window allows visibility to number of clips left in the cartridge.

Safety- interlock design

• Prevents jaw closing to damage the tissue when the instruments is empty.