



| Order Information | | | | |
|--------------------|--|-----------------|----------------------------------|----------|
| ORDER CODE | DESCRIPTION | CARTRIDGE COLOR | STAPLE SIZE (INNER TO OUTER ROW) | QUANTITY |
| EGIA45CTAV* | Endo GIA™ 45 mm Curved Tip Articulating Vascular Reload | Gray | 2 mm, 2 mm, 2 mm | Pack 6 |
| EGIA30CTAVM | Endo GIA™ 30 mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology | Tan | 2 mm, 2.5 mm, 3 mm | Pack 6 |
| EGIA45CTAVM | Endo GIA™ 45 mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology | Tan | 2 mm, 2.5 mm, 3 mm | Pack 6 |
| EGIA60CTAVM | Endo GIA™ 60 mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology | Tan | 2 mm, 2.5 mm, 3 mm | Pack 6 |
| EGIA45CTAMT | Endo GIA™ 45 mm Curved Tip Articulating Medium/Thick Reload with Tri-Staple™ Technology | Purple | 3 mm, 3.5 mm, 4 mm | Pack 6 |
| EGIA60CTAMT | Endo GIA™ 60 mm Curved Tip Articulating Medium/Thick Reload with Tri-Staple™ Technology | Purple | 3 mm, 3.5 mm, 4 mm | Pack 6 |

OVERVIEW

Visibility and Maneuverability

The Endo GIA™ curved tip reload with Tri-Staple™ technology is designed with a curved tip at the distal end of the anvil to provide visibility and maneuverability around target tissues and vessels.

[Order Information >](#)

About Medtronic

[Medtronic Corporate ANZ](#)

[Privacy Statement](#)

[Terms of Use](#)

[Terms & Conditions of Supply](#)

[External Research Program](#)

Minimally Invasive Therapies

[Products](#)

[Terms & Conditions of Supply](#)

[Sitemap](#)

[Australia/New Zealand Customer Service](#)



[Contact Us](#)

Country Office

Level 3 - Building 5
Central Corporate Centre
666 Great South Road
Penrose, Auckland
New Zealand