

# "Coolview®" – CLED53 New Generation LED Examination Light

**Coolview CLED53** LED Examination Lights from Brandon Medical are the perfect light source for a wide range of procedures. Providing a high intensity, bright natural looking light, Coolview CLED53 is perfect for medical consultations and examinations in primary care and medical clinics.

#### Powerful, variable Light Intensity

*Coolview CLED53* emits a powerful 50,000 Lux of light at 0.5Metres, providing an intense, clear

light source for excellent visualization. The easily accessible and intuitive controls offer 5 levels of brightness



#### **Infection Control**

The fully enclosed light head is sealed to protect from dust and light liquid ingress.

#### Mountings

The **CLED53 Examination Light** operates on a spring balanced rigid arm offering optimum reach, positioning and stability. The Light can be specified for mounting on the Wall (with or without Backplate), Desk, "Medirail" or on a Mobile Wheeled Base

For fuller details on features and specifications please refer to the Technical Data overleaf or contact Brandon Medical

Brandon Medical Co Ltd, Elmfield Road, Morley, Leeds, LS27 OEL T: +44 (0)845 1243 666 F: +44 (0)845 1243 667 T: +44 (0)113 2777 393 F: +44 (0)113 2728 844 E: enquiries@brandon-medical.com www.brandon-medical.com



CLED53FX Rail Mounting (illustrated)

#### **Perfect Colour**

The powerful light source is matched with an selectable colour temperature feature, which can be adjusted between 3,300, 3,800 and 4,400<sup>o</sup>Kelvin. This provides the Clinician with accurate representation and visualization of the site under examination.

Coolview CLED53 offers a high CRI Ra of  $\geq$ 93 which describes the exceptionally high quality of white Light and excellent overall colour rendition, whilst the high CRI R9 value of  $\geq$ 90 can help to enhance visualization of areas where red is more predominant.







# "Coolview®" – CLED53 New Generation LED Examination Light

## **Technical Data**


|                     | CLED53  |                        | Protection Class II                                |  |
|---------------------|---|------------------------|--|--|
|                     | Degree of Protection as per IEC 60529                                   |                        | IP20 (lamphead)                                    |  |
|                     | Classification according to EU-Regulation 2017/745<br>(MDR), Article 51 |                        | Class I  |  |
| UO.                 | Electrical Safety and EMC A   | pprovals               | EN/IEC 60601-1, EN/IEC 60601-1-2, BS EN 60601-2-41 |  |
| mat                 | Type of luminaire (according  | to BS EN60601-2-41)    | Luminaire for Diagnosis                            |  |
| General Information | Light Intensity Control   |                        | Lamphead Keypad (5 step dimming)                   |  |
|                     | Lamphead On/Off switch  |                        | Spring Arm Switch                                  |  |
| Jera                | Maximum Total System Wei  | ght (kg) [approx.]     | 1.7 kg   |  |
| B                   | Total Arm Reach (mm) [appi  | ox.]                   | 850mm  |  |
|                     | Rotation About Central Axis of mount (max)                              |                        | 360°   |  |
|                     | Arm Angle bottom joint (max)  |                        | 69°  |  |
|                     | Lamphead Roll (Rotation about arm) (max)                                |                        | 300°   |  |
|                     | Mains Power Supply Lead   |                        | 1.8m lead  |  |
|                     | Rated Input Voltage   |                        | 100-240 V  |  |
| Electrical<br>Data  | Frequency Range   |                        | 50/60Hz  |  |
| Це                  | Nominal Power Consumptio  | n                      | 12-18W (15-22VA)                                   |  |
|                     | Number of Light Sources (LE   | ED Bulb)               | 7  |  |
|                     | Type of Light Source  |                        | LED  |  |
| <b>0</b> *          | Central Illuminance Ec (Light   | Intensity) @ 0.5m      | 50,000 Lux *                                       |  |
| Photometric Data*   | Light Field Diameter  | (d10 Distance of 0.5m) | 180 mm   |  |
| ţi.                 |   | (d50 Distance of 0.5m) | 90 mm  |  |
| me                  | Colour Temperature  |                        | 4400K / 3800K / 3300K                              |  |
| oto                 | Colour rendering index Ra   |                        | 93   |  |
| £.                  | Colour rendering index R9   |                        | 90   |  |
|                     | Total Irradiance Ee at maximum intensity                                |                        | <180 mW/m <sup>2</sup>                             |  |
|                     | LED Life  |                        | 50,000 Hours                                       |  |
| 100/ /              |   |                        |  |  |

\*-10% / +20%

As part of Brandon Medical policy of continuous improvement, Brandon Medical reserves the right to change specifications as and when product improvements become available. Brandon Medical Co Ltd | Elmfield Road | Morley | Leeds | LS27 0EL Tel: 0113 277 7393 Fax: 0113 272 8844 Email: enquiries@brandon.medical.com

#### Product Variants -

| CLED53FXBP | CLED53 FX Back Plate (Wall) Mounted Examination Light |
|------------|---|
| CLED53FXDP | CLED53 FX Desk (Clamp) Mounted Examination Light      |
| CLED53FXMP | CLED53 FX Mobile Examination Light                    |
| CLED53FXRP | CLED53 FX MediRail Mounted Examination Light          |
| CLED53FXWP | CLED53 FX Wall Mounted Examination Light              |





Brandon Medical Co Ltd, Elmfield Road, Morley, Leeds, LS27 OEL T: +44 (0)845 1243 666 F: +44 (0)845 1243 667 T: +44 (0)113 2777 393 F: +44 (0)113 2728 844 E: enquiries@brandon-medical.com www.brandon-medical.com

Brilliant by design

To: Whom It May Concern

Date: 06/06/2023.

#### Authorization

We, Company-manufacturer BRANDON MEDICAL COMPANY LTD., of ELMFIELD ROAD, MORLEY, LEEDS, LS27 0EL

U.K., hereby appoint following company:

F.C.P.C. "DataControl" S.R.L.,

17/6, N.Testimiteanu street,

MD-2025, Chișinau, Republic of Moldova,

to be our official Representative at the responsible authorities of the Republic of Moldova for registration, renewal,

variation of registration etc. .

Signature for and on behalf of Manufacturer, title of the authorized person, seal



Arthur Loftus Regional Director, EMEA

Brandon Medical Company Limited | Elmfield Road | Morley | Leeds | UK | LS27 0EL T:+44 113 277 7393 / 0845 1243 666 F:+44 113 272 8844 / 0845 1243 667 E: <u>enquiries@brandon-medical.com</u> W: <u>www.brandon-medical.com</u> Brandon Medical Company Limited. Registered in England & Wales No. 02827189 VAT No. GB734 591421. A subsidiary of Brandon Group Limited.

Brandon Medical is always looking for opportunities to improve our products and customer service and we greatly value any feedback you can provide: <u>feedback@brandon-medical.com</u>





# CE

#### **Brandon Medical Company Limited (BMCL)**

Elmfield Road, Morley, Leeds, West Yorkshire, LS27 0EL, UK

SRN (Non-EU Manufacturer): GB-MF-000005501

SRN (European Authorised Representative): IE-AR-000003999

Declares that under our sole responsibility, the product: C

Coolview 50 Coolview 53

(Refer to end of declaration for variant UDI-DI & Basic UDI-DI)

With intended purpose of:

#### **Examination Lamp**

Is in conformity to the general safety and performance requirements of the following EC Regulation:

- 2017/745 Medical Device Regulation Risk Classification as per Annex VIII
- Class I Self-Declaration Medical Device
- Conformity assessment route: <u>Class I</u>: EC conformity self-declaration according to MDR 2017/745 Article 19, and Annexes II and III.

The following standards were applied in full:

- EN ISO 13485:2016 Medical devices Quality management systems Requirements for regulatory purposes
- ISO 9001:2015 Quality management systems Requirements
- BS EN ISO 14971:2012 Medical devices Application of risk management to medical devices
- EN ISO 15223-1:2016 Medical devices Symbols to be used with medical device labels, labelling and information to be supplied Part 1: General requirements
- EN60601-1:2006+A1:2013 Medical electrical equipment Part 1: General requirements for basic safety and essential performance
- BS EN60601-1-2:2015 Medical electrical equipment. General requirements for basic safety and essential performance. Collateral Standard: Electromagnetic disturbances. Requirements and tests
- BS EN 60601-2-41 2009+A11:2011 Particular requirements for basic safety and essential performance of surgical luminaires and luminaires for diagnosis
- EN 60601-6:2010+A1:2015
- Waste from Electrical and Electronic Equipment (WEEE) Directive

This product has been developed, tested and verified in accordance with the international standards listed above.

| Person Responsible for<br>Regulatory Compliance (PRRC)                           | Signature | Date        |
|--|-----------|-------------|
| Annette Hemmings<br>Head of Quality Assurance<br>Brandon Medical Company Limited | A.E. A.   | 17 MAY 2022 |
| Countersignature   | Signature | Date        |
| Tom Partington<br>Product Manager<br>Brandon Medical Company Limited             | A         | 17-05-22    |

European Authorised Representative: European Healthcare & Device Solutions (Ireland) Ltd. Stratton House, Bishopstown Road, Bishopstown, Cork, Ireland



| Model variant | Model description   | UDI-DI (GTIN)  | Basic UDI-DI    |
|---------------|---|----------------|-----------------|
| CLED50FXBP    | Coolview CLED50 FX Articulated Arm Wall<br>Mounted with Backplate Examination Lamp with<br>UK Plug Lead | 05060413447439 | 506041344CLED2L |
| CLED50FXC     | Coolview CLED50FX Articulated Arm Ceiling<br>Mounted Examination Lamp (Limited Rotation)                | 05060413448511 | 506041344CLED2L |
| CLED50FXDP    | Coolview CLED50 FX Articulated Arm Desk<br>Mounted Examination Lamp with UK Plug Lead                   | 05060413447446 | 506041344CLED2L |
| CLED50FXMP    | Coolview CLED50 FX Articulated Arm Mobile<br>Mounted Examination Lamp with UK Plug Lead                 | 05060413447453 | 506041344CLED2L |
| CLED50FXRP    | Coolview CLED50 FX Articulated Arm Rail<br>Mounted Examination Lamp with UK Plug Lead                   | 05060413447460 | 506041344CLED2L |
| CLED50FXWP    | Coolview CLED50 FX Articulated Arm Wall<br>Mounted Examination Lamp with UK Plug Lead                   | 05060413447477 | 506041344CLED2L |
| CLED50SXBP    | Coolview CLED 50 SX Flexible Arm Wall Mounted<br>with Backplate Examination Lamp with UK Plug<br>Lead   | 05060413447484 | 506041344CLED2L |
| CLED50SXDP    | Coolview CLED 50 SX Flexible Arm Desk<br>Mounted with Lamp with UK Plug Lead                            | 05060413447491 | 506041344CLED2L |
| CLED50SXMP    | Coolview CLED 50 SX Flexible Arm Mobile<br>Mounted with Lamp with UK Plug Lead                          | 05060413447507 | 506041344CLED2L |
| CLED50SXRP    | Coolview CLED 50 SX Flexible Arm Rail Mounted<br>with Lamp with UK Plug Lead                            | 05060413447514 | 506041344CLED2L |
| CLED50SXWP    | Coolview CLED 50 SX Flexible Arm Wall Mounted<br>with Lamp with UK Plug Lead                            | 05060413447521 | 506041344CLED2L |
| CLED53FXBP    | Coolview CLED53 FX Articulated Arm Wall<br>Mounted with Backplate Examination Lamp with<br>UK Plug Lead | 05060413447545 | 506041344CLED2L |
| CLED53FXDP    | Coolview CLED53 FX Articulated Arm Desk<br>Mounted Examination Lamp with UK Plug Lead                   | 05060413447552 | 506041344CLED2L |
| CLED53FXMP    | Coolview CLED53 FX Articulated Arm Mobile<br>Mounted Examination Lamp with UK Plug Lead                 | 05060413447569 | 506041344CLED2L |
| CLED53FXRP    | Coolview CLED53 FX Articulated Arm Rail<br>Mounted Examination Lamp with UK Plug Lead                   | 05060413447576 | 506041344CLED2L |
| CLED53FXWP    | Coolview CLED53 FX Articulated Arm Wall<br>Mounted Examination Lamp with UK Plug Lead                   | 05060413447583 | 506041344CLED2L |
| CLED53FXC     | Coolview CLED53FX Articulated Arm Ceiling<br>Mounted Examination Lamp (Limited Rotation)                | 05060413448528 | 506041344CLED2L |





BS EN ISO 13485:2016

Certificate Number 326 Issue 3

# Certificate of Registration

*Certificate of registration of an organisation or business of assessed capability.* Business Quality Assurance International Ltd hereby grants to

# Brandon Medical Company Ltd., Elmfield Road, Morley, Leeds, West Yorkshire, LS27 OEL, UK.

(hereafter called the Business) the right and licence to be listed in the BQAI Register of Businesses of Assessed Capability in respect of the goods or services particularly described below. Goods or services shall be provided by the business, only at the address(es) given above, all in accordance with BS EN ISO 13485:2016.

The Certificate of Registration is granted subject to BQAI regulations governing the registration referred to above.

Goods or services in respect of which the business is registered: Design, manufacture, installation and servicing of: operating theatre lights, medical luminaires, and accessories. Distribution, installation and servicing of operating tables, accessories and medical pendants.





Tuchart J. Copy.

BQAI Governing Board Chairman, Prof. Michael J. Hynes

Signed:

Signed:

BQAI Chief Executive Officer, Mr. Declan O'Connor

Original date of registration: **31st August 2018** Current certification three-year cycle start date: **31st August 2021** Date of governing board endorsement: **13th September 2021** Current certificate expiry date: **31st August 2024** Date of extension of scope change approval: **13th September 2021** 

> The approval is subject to the supplier maintaining its Quality Management System to the required standard which will be monitored by BQAI. This Certificate remains the property of Business Quality Assurance International Ltd. and can be recalled at any time.



# CE

#### **Brandon Medical Company Limited (BMCL)**

Elmfield Road, Morley, Leeds, West Yorkshire, LS27 0EL, UK SRN (Non-EU Manufacturer): GB-MF-000005501

SRN (European Authorised Representative): IE-AR-000003999

Declares that under our sole responsibility, the product:

Quasar eLite Fixed Surgical Lighting Qe60 Basic UDI-DI: 506041344QE6023 Qe30 Basic UDI-DI: 506041344QE30ZP UDI-DI: Refer to appended list

With intended purpose of:

#### Surgical Luminaire

Is in conformity to the general safety and performance requirements of the following EC Regulation:

- 2017/745 Medical Device Regulation Risk Classification as per Annex VIII
- Class I Self-Declaration Medical Device
- Conformity assessment route: <u>Class I:</u> EC conformity self-declaration according to MDR 2017/745 Article 19, and Annexes II and III.

The following standards were applied in full:

- EN ISO 13485:2016 Medical devices Quality management systems Requirements for regulatory purposes
- ISO 9001:2015 Quality management systems Requirements
- BS EN ISO 14971:2012 Medical devices Application of risk management to medical devices
- EN ISO 15223-1:2016 Medical devices Symbols to be used with medical device labels, labelling and information to be supplied Part 1: General requirements
- EN60601-1:2006+A1:2013 Medical electrical equipment Part 1: General requirements for basic safety and essential performance
- BS EN60601-1-2:2015 Medical electrical equipment. General requirements for basic safety and essential performance. Collateral Standard: Electromagnetic disturbances. Requirements and tests
- BS EN 60601-2-41 2009+A11:2011 Particular requirements for basic safety and essential performance of surgical luminaires and luminaires for diagnosis
- Radio equipment directive 2014/53/EU
- BS EN 62366:2008 Medical devices. Application of usability engineering to medical devices
- BS EN1041:2008+A1:2013 Information supplied by the manufacturer of medical devices
- Waste from Electrical and Electronic Equipment (WEEE) Directive

This product has been developed, tested and verified in accordance with the international standards listed above.

| Person Responsible for<br>Regulatory Compliance (PRRC)                     | Signature  | Date                           |
|--|--|--------------------------------|
| Nigel Davill<br>Chief Technical Officer<br>Brandon Medical Company Limited | All in the second secon | 23 <sup>rd</sup> February 2023 |

European Authorised Representative: European Healthcare & Device Solutions (Ireland) Ltd. Stratton House, Bishopstown Road, Bishopstown, Cork, Ireland





| Product Code | Description   | UDI-DI (GTIN) |
|--------------|---|---------------|
| Qe3010lonc   | Quasar eLite Qe30 lamp head (160Klux), Qe10 satellite (100Klux), low ceiling, light spring arms, tandem bracket for ceiling mount.  | 5060413449631 |
| Qe3010lotc   | Quasar eLite Qe30 lamp head (160Klux), Qe10 satellite (100Klux), low ceiling, light spring arms, Titan central axis for ceiling mount.  | 5060413449648 |
| Qe303030satc | Quasar eLite Qe30 lamp head (160Klux), Qe30 satellite (160Klux), Qe30 2nd satellite (160Klux), single cardanic, AC spring arms, Titan central axis for ceiling mount.   | 5060413449051 |
| Qe3030GDsvtc | Quasar eLite Qe30 lamp (160Klux), Qe30 satellite (160Klux), 1.5-21Kg monitor arm, cabled with stops, single cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060928390541 |
| Qe3030lonc   | Quasar eLite Qe30 lamp head (160Klux), Qe30 satellite (160Klux), low ceiling, light spring arms, tandem bracket for ceiling mount.  | 5060413449655 |
| Qe3030lotc   | Quasar eLite Qe30 lamp head (160Klux), Qe30 satellite (160Klux), low ceiling, light spring arms, Titan central axis for ceiling mount.  | 5060413449662 |
| Qe3030LTJ    | Quasar eLite Qe30 lamp head (160Klux), Qe30 satellite (160Klux), single cardanic, tandem bracket for ceiling mount.   | 5060413446852 |
| Qe3030sanc   | Quasar eLite Qe30 lamp head (160Klux), Qe30 satellite (160Klux), single cardanic, AC spring arms, tandem bracket for ceiling mount.   | 5060413449679 |
| Qe3030sanc   | Quasar eLite Qe30 lamp head (160Klux), Qe30 satellite (160Klux), single cardanic, AC spring arms, tandem bracket for ceiling mount.   | 5060413449679 |
| Qe3030satc   | Quasar eLite Qe30 lamp head (160Klux), Qe30 satellite (160Klux), single cardanic, AC spring arms, Titan central axis for ceiling mount.   | 5060413449242 |
| Qe3030svtc   | Quasar eLite Qe30 lamp head (160Klux), Qe30 satellite (160Klux), single cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449686 |
| Qe3030TJ     | Quasar eLite Qe30 lamp head (160Klux), Qe30 satellite (160Klux), single cardanic, light spring arms, tandem bracket for ceiling mount.  | 5060413449556 |
| Qe30looc     | Quasar eLite Qe30 lamp head (160Klux), low ceiling, light spring arms, light horizontal for ceiling mount.  | 5060413449693 |
| Qe30lvtc     | Main Qe30 Lamp head (160kLux) - Low-Ceiling-V Spring-Arm - Titan Central Axis Ceiling Mounted   | 5060413449709 |
| Qe30svtc     | Main Qe30 Lamp head (160kLux) -Single Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449716 |
| QESUSVIC     |   | 5000415449710 |
| Qe3C10svtc   | Quasar eLite Qe30 lamp (160Klux) with camera and stops, Qe10 satellite (100Klux), single cardanic, V spring arms,<br>Titan central axis for ceiling mount.  | 5060413449723 |
| Qe3C30GDsvtc | Quasar eLite Main Qe30 Lamphead with HD Fixed Focus Camera and stop (160kLux), Satellite Qe30 Lamphead (160kLux), Satellite Accessory Arm 21kg Stopped-Single Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted | 5060928390695 |
| Qe3C30svtc   | Quasar eLite Qe30 lamp (160Klux) with camera and stops, Qe30 satellite (160Klux), single cardanic, V spring arms,<br>Titan central axis for ceiling mount.  | 5060413449730 |
| Qe3C60svtc   | Quasar eLite Qe30 lamp (160Klux) with camera and stops, Qe60 satellite (160Klux), single cardanic, V spring arms,<br>Titan central axis for ceiling mount.  | 5060413449747 |
| Qe3Csvtc     | Quasar eLite Qe30 lamp (160Klux) with camera and stops, single cardanic, V spring arms, Titan central axis for ceiling mount.   | 5060413449754 |
| Qe3H10svtc   | Main Qe30 Lamp head with HD Fixed Focus Camera 360° (160kLux) Satellite Qe10 Lamp head Lamp (100kLux) -<br>Single Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449761 |
| Qe3H30sanc   | Quasar eLite Main Qe30 Lamp head with HD Fixed Focus Camera 360° (160kLux) Satellite Qe30 Lamp head<br>(160kLux)-Single Cardanic-AC Spring-Arm-Tandem (45mm Stem) Ceiling Mounted                                       | 5060413449778 |
| Qe3H30svtc   | Main Qe30 Lamp head with HD Fixed Focus Camera 360° (160kLux) Satellite Qe30 Lamp head (160kLux) -Single<br>Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449785 |
| Qe3Hlvtc     | Main Qe30 Lamp head with HD Fixed Focus Camera 360° (160kLux) -Low-Ceiling-V Spring-Arm -Titan Central Axis<br>Ceiling Mounted  | 5060413449792 |
| Qe3Hsvtc     | Main Qe30 Lamp head with HD Fixed Focus Camera 360° (160kLux) -Single Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449808 |
| Qe603010dvtc | Main Qe60 Lamp head (160kLux) Satellite Qe30 Lamp head (160kLux) Satellite Qe10 Lamp head Lamp (100kLux)-<br>Double Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449815 |
| Qe603010svtc | Main Qe30 Lamp head with HD Fixed Focus Camera 360° (160kLux) Satellite Qe30 Lamp head (160kLux) Satellite Qe10 Lamp head Lamp (100kLux) -Single Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted              | 5060413449822 |
| Qe603030dvtc | Quasar eLite Qe60 lamp head (160Klux), Qe30 satellite (160Klux), Qe30 2nd satellite (160Klux), double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449839 |
| Qe603030svtc | Quasar eLite Qe60 lamp head (160Klux), Qe30 satellite (160Klux), Qe30 2nd satellite (160Klux), single cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449846 |
| Qe6030dvtc   | Quasar eLite Qe60 lamp head (160Klux), Qe30 satellite (160Klux), double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449006 |
| Qe6030GDdvtc | Quasar eLite Qe60 lamp (160Klux), Qe30 satellite (160Klux), 1.5-21Kg monitor arm, cabled with stops, double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060928390602 |
| Qe6030GDsvtc | Quasar eLite Qe60 lamp (160Klux), Qe30 satellite (160Klux), 1.5-21Kg monitor arm, cabled with stops, single cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060928390558 |
| Qe6030lvtc   | Quasar eLite Main Qe60 Lamp head (160kLux) Satellite Qe30 Lamp head (160kLux) -Low-Ceiling-V Spring-Arm -Titan<br>Central Axis Ceiling Mounted  | 5060413449037 |
| Qe6030svtc   | Quasar eLite Qe60 lamp head (160Klux), Qe60 satellite (160Klux), single cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449143 |
| Qe606030dvtc | Quasar eLite Qe60 lamp head (160Klux), Qe60 satellite (160Klux), Qe30 2nd satellite (160Klux), double cardanic, V   | 5060413449853 |



| Qe606030svtc         | Quasar eLite Qe60 lamp head (160Klux), Qe60 satellite (160Klux), Qe30 2nd satellite (160Klux), single cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449860 |
|----------------------|---|---------------|
| Qe606060dvtc         | Quasar eLite Qe60 lamp head (160Klux), Qe60 satellite (160Klux), Qe60 2nd satellite (160Klux), double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449877 |
| Qe606060svtc         | Quasar eLite Qe60 lamp head (160Klux), Qe60 satellite (160Klux), Qe60 2nd satellite (160Klux), single cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449884 |
| Qe6060dvtc           | Quasar eLite Qe60 lamp head (160Klux), Qe60 satellite (160Klux), double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449013 |
| Qe6060lvtc           | Quasar eLite Main Qe60 Lamp head (160kLux) Satellite Qe60 Lamp head (160kLux) -Low-Ceiling-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449167 |
| Qe6060svtc           | Quasar eLite Qe60 lamp head (160Klux), Qe60 satellite (160Klux), single cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449150 |
| Qe60dvtc             | Quasar eLite Qe60 lamp head (160Klux), double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449112 |
| Qe60GAdvtc           | Quasar eLite Qe60 lamp (160Klux), 1.5-21Kg monitor arm, 3P with 360° power connectors, double cardanic, V   | 5060928390725 |
| Qe60GDdvtc           | spring arms, Titan central axis for ceiling mount.<br>Quasar eLite Qe60 lamp (160Klux), 1.5-21Kg monitor arm, cabled with stops, double cardanic, V spring arms, Titan  | 5060928390480 |
| Occolute             | central axis for ceiling mount.<br>Quasar eLite Main Qe60 Lamp head (160kLux) -Low-Ceiling-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449136 |
| Qe60lvtc<br>Qe60svtc | Quasar eLite Qe60 lamp head (160Klux), single cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449130 |
| Qe6C30GDdvtc         | Quasar eLite Main Qe60 Lamp head with HD Camera and Stop (160kLux) Satellite Qe30 Lamp head (160kLux)<br>Satellite Accessory Arm 21kg Stopped-Double Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted  | 5060413449891 |
| Qe6C30svtc           | Quasar eLite Qe60 lamp (160Klux) with camera and stops, Qe30 satellite (160Klux), single cardanic, V spring arms,   | 5060413449075 |
| Qe6C60dvtc           | Titan central axis for ceiling mount.<br>Quasar eLite Main Qe60 Lamp head with HD Camera and Stop (160kLux) Satellite Qe60 Lamp head (160kLux)-   | 5060413449907 |
| Qe6C60svtc           | Double Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted<br>Quasar eLite Qe60 lamp (160Klux) with camera and stops, Qe60 satellite (160Klux), single cardanic, V spring arms,   | 5060413449914 |
| Qe6Clvtc             | Titan central axis for ceiling mount.<br>Quasar eLite Main Qe60 Lamp head with HD Camera and Stop (160kLux) -Low-Ceiling-V Spring-Arm -Titan Central  | 5060413449921 |
| Qe6Csvtc             | Axis Ceiling Mounted<br>Quasar eLite Qe60 lamp (160Klux) with camera and stops, single cardanic, V spring arms, Titan central axis for  | 5060413449938 |
| Qe6H3010LCCT         | ceiling mount.<br>Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) Satellite Qe30 Lamp head (160kLux) Satellite<br>Qe10 Lamp head (100kLux) -Low-Ceiling-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413448832 |
| Qe6H3030dvtc         | Quarter Lamp head (160kLux) -Cow-Cennig-V Spring-Arm - Intal Central Axis Cennig Modified<br>Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) Satellite Qe30 Lamp head (160kLux) Satellite<br>Qe30 Lamp head (160kLux)-Double Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted | 5060413449563 |
| Qe6H3030svtc         | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) Satellite Qe30 Lamp head (160kLux) Satellite   | 5060928390312 |
| Qe6H30dvtc           | Qe30 Lamp head (160kLux), single Cardanic, V Spring-Arm, Titan Central Axis Ceiling Mounted<br>Quasar eLite Qe60 lamp (160Klux) with camera, Qe30 satellite (160Klux), double cardanic, V spring arms, Titan central<br>axis for ceiling mount.   | 5060413449228 |
| Qe6H30GAdvtc         | Quasar eLite Qe60 lamp (160Klux) with camera, Qe30 satellite (160Klux), 1.5-21Kg monitor arm, 3P with 360° power connectors, double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449945 |
| Qe6H30GDdvtc         | Quasar eLite Qe60 lamp (160Klux) with camera, Qe30 satellite (160Klux), 1.5-21Kg monitor arm, cabled with stops, double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449464 |
| Qe6H30HDdvtc         | Quasar eLite Qe60 lamp (160Klux) with camera, Qe30 satellite (160Klux), 26Kg monitor arm, cabled with stops, double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449952 |
| Qe6H30lvtc           | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) Satellite Qe30 Lamp head (160kLux) -Low-Ceiling-<br>V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449235 |
| Qe6H30svtc           | Quasar eLite Qe60 lamp (160Klux) with camera, Qe30 satellite (160Klux), single cardanic, V spring arms, Titan central axis for ceiling mount.   | 5060413449020 |
| Qe6H6030dvtc         | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) Satellite Qe60 Lamp head (160kLux) Satellite Qe30 Lamp head (160kLux)-Double Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449969 |
| Qe6H6030svtc         | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) Satellite Qe60 Lamp head (160kLux) Satellite Qe30 Lamp head (160kLux)-Single Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060928390848 |
| Qe6H6060dvtc         | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) Satellite Qe60 Lamp head (160kLux) Satellite Qe60 Lamp head (160kLux)-Double Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449976 |
| Qe6H60dvtc           | Quasar eLite Qe60 lamp (160Klux) with camera, Qe60 satellite (160Klux), double cardanic, V spring arms, Titan central axis for ceiling mount.   | 5060413449068 |
| Qe6H60GAdvtc         | Quasar eLite Qe60 lamp (160Klux) with camera, Qe60 satellite (160Klux), 1.5-21Kg monitor arm, 3P with 360° power connectors, double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449983 |
| Qe6H60GDdvtc         | Quasar eLite Qe60 lamp (160Klux) with camera, Qe60 satellite (160Klux), 1.5-21Kg monitor arm, cabled with stops, double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060413449600 |
| Qe6H60HDdvtc         | Quasar eLite Qe60 lamp (160Klux) with camera, Qe60 satellite (160Klux), 26Kg monitor arm, cabled with stops, double cardanic, V spring arms, Titan central axis for ceiling mount.  | 5060928390008 |
| Qe6H60lvtc           | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) Satellite Qe60 Lamp head (160kLux) -Low-Ceiling-<br>V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449211 |
| Qe6H60svtc           | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) Satellite Qe60 Lamp head (160kLux) -Single<br>Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted  | 5060413449204 |



| Qe6Hdvtc   | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) -Double Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted                                       | 5060413449174 |
|------------|--|---------------|
| Qe6HGDdvtc | Quasar eLite Qe60 lamp (160Klux) with camera, 1.5-21Kg monitor arm, cabled with stops, double cardanic, V spring arms, Titan central axis for ceiling mount.           | 5060928390015 |
| Qe6Hlvtc   | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) -Low-Ceiling-V Spring-Arm -Titan Central Axis Ceiling Mounted   | 5060413449181 |
| Qe6Hsvtc   | Quasar eLite Main Qe60 Lamp head with HD Camera 360° (160kLux) -Single Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted                                       | 5060413449198 |
| Qe6P30dvtc | Quasar eLite Main Qe60 Lamp head with Camera Preparation (160kLux) Satellite Qe30 Lamp head (160kLux)-Double Cardanic-V Spring-Arm -Titan Central Axis Ceiling Mounted | 5060413449426 |
| Qe6P60dvtc | Quasar eLite Qe60 lamp (160Klux) with camera provision, Qe60 satellite (160Klux), double cardanic, V spring arms, Titan central axis for ceiling mount.                | 5060928390022 |
| QeGA0vtc   | Quasar eLite 1.5-21Kg monitor arm, 3P with 360° power connectors, no bracket, V spring arms, Titan central axis for ceiling mount.                                     | 5060413449440 |
| QeGA60dvtc | Quasar eLite Accessory Arm 21kg 3 Pole 360° Satellite Qe60 Lamphead (160kLux)-Double Cardanic-V Spring-Arm -<br>Titan Central Axis Ceiling Mounted                     | 5060928390770 |
| QeGC60dvtc | Quasar eLite Accessory Arm 21kg 5P + 360° video connections Satellite Qe60 Lamp head (160kLux)-Double Cardanic-<br>V Spring-Arm -Titan Central Axis Ceiling Mounted    | 5060413449570 |
| QeGDtvtc   | Quasar eLite 1.5-21Kg monitor arm, cabled with stops, VESA twin monitor bracket, V spring arms, Titan central axis for ceiling mount.                                  | 5060413449372 |
| QeGDvvtc   | Quasar eLite Accessory Arm 21kg Stopped -with Single Symposia-view VESA Monitor Bracket (fixed) -V Spring-Arm -<br>Titan Central Axis Ceiling Mounted                  | 5060413449259 |
| QeHDtvtc   | Quasar eLite 26Kg monitor arm, cabled with stops, VESA twin monitor bracket, V spring arms, Titan central axis for ceiling mount.                                      | 5060413449471 |







ISO 9001:2015 c

Certificate Number 325 Issue 2

# Certificate of Registration

Certificate of registration of an organisation or business of assessed capability.

Business Quality Assurance International Ltd hereby grants to

# Brandon Medical Company Ltd. Elmfield Road, Morley, Leeds, West Yorkshire, LS27 OEL, UK.

(hereafter called the Business) the right and licence to be listed in the BQAI Register of Businesses of Assessed Capability in respect of the goods or services particularly described below. Goods or services shall be provided by the business, only at the address(es) given above, all in accordance with ISO 9001:2015.

The Certificate of Registration is granted subject to BQAI regulations governing the registration referred to above.

Goods or services in respect of which the business is registered: Design, manufacture, installation and servicing of: operating theatre lights, medical luminaries and associated control panels, medical video systems, emergency power units and earth bonding bars. Distribution of medical pendants, operating tables, UCV canopies, IPS and UPS systems.

Signed:

BQAI Governing Board Chairman, Prof. Michael J. Hynes

Signed:

BQAI Chief Executive Officer, Mr. Declan O'Connor

Original date of registration: **31st August 2018** Current certification three-year cycle start date: **31st August 2021** Date of governing board endorsement: **31st August 2021** Current certificate expiry date: **31st August 2024** Date of extension of scope / name change approval: **N/A** 

> The approval is subject to the supplier maintaining its Quality Management System to the required standard which will be monitored by BQAI. This Certificate remains the property of Business Quality Assurance International Ltd. and can be recalled at any time.

Anexa nr. 1 La Procedurile administrative pentru notificarea dispozitivelor medicale care dețin marcajul CE

> Către Agenția Medicamentului și Dispozitivelor Medicale

#### NOTIFICARE

pentru înregistrarea dispozitivelor medicale în Registrul de stat al dispozitivelor medicale

nr. 47 din 09.06.2023

Solicitantul **FCPC** "**DataControl" S.R.L.**, cu sediul **mun. Chișinău, str. N. Testemițanu, 17/6**, tel./fax: 022 27 37 12, e-mail: <u>contact@datacontrol.md</u>, solicit înregistrarea în Registrul de stat al dispozitivelor medicale a următoarelor categorii și tipuri de dispozitive medicale pentru introducerea și punerea la dispoziție pe piață a:

#### **Product:** Examination Lamp

**Model: COOLVIEW CLED53 FX** 

CLED53FXBP CLED53FXDP CLED53FXMP CLED53FXRP CLED53FXWP CLED53FXC

Se anexează următoarele acte:

- 1. Declarație de Conformitate CE
- 2. Actul prin care producătorul își desemnează reprezentantul
- 3. Declarație pe propria răspundere.

Data 09.06.2023

A CONTROL OF A CON

Semnătura

# Tabelul de recepționare a notificării

- 1<sup>-1</sup>

(se completează de către Agenție în momentul depunerii notificării de către solicitant)

| Comentarii cu privire la                  |   |
|---|---|
| acceptul/refuzul recepționării            |   |
| acceptul/reluzurreceptionani              | A BRACK STOR  |
| notificării, inclusiv motivul refuzului   |   |
| Data/nr. de ordine atribuit notificării   | we have the second some desident a contact  |
| de către Agenție (în cazul acceptării     | Served States And Toolald He  |
|   |   |
| recepționării)                            | and the state of the |
| Numele, prenumele, funcția                |   |
| Numere, prenamere, rang, a                |   |
| persoanei responsabile de                 | All a lot and a lot a |
| recepționarea dosarului                   |   |
| Semnătura persoanei responsabile          | and the second second second second second second second  |
| Semilarula persoaner responsablic         | tionaristic and the state of the second of some second  |
| the second second in an a weather and the | min is committed and an an an and and   |

and the second

Către Agenția Medicamentului și Dispozitivelor Medicale

# DECLARAȚIE PE PROPRIE RĂSPUNDERE

Solicitantul F.C.P.C. "DataControl" S.R.L., cu sediul în mun. Chișinău, str. N. Testemițanu 17/6, tel./fax: 022 27 37 12, e-mail: contact@datacontrol.md,

declar pe proprie răspundere, cunoscînd prevederile art. **352<sup>1</sup>**, Codul Penal al Republicii Moldova cu privire la falsul în declarații, că documentele și datele furnizate pentru notificarea dispozitivului medical:

# **Product:** Examination Lamp

# Model: COOLVIEW CLED53 FX

CLED53FXBP CLED53FXDP CLED53FXMP CLED53FXRP CLED53FXWP CLED53FXC

Sunt autentice și corespund realității

Alexandru Grabazei, director



Data: 09.06.2023

Medical power and control, lighting, architectural and audio video



To: Whom It May Concern

Date: 06/06/2023.

Authorization

We, Company-manufacturer BRANDON MEDICAL COMPANY LTD., of ELMFIELD ROAD, MORLEY, LEEDS, LS27 OEL

U.K., hereby appoint following company:

F.C.P.C. "DataControl" S.R.L.,

17/6, N.Testimiteanu street,

MD-2025, Chișinau, Republic of Moldova,

to be our official Representative at the responsible authorities of the Republic of Moldova for registration, renewal,

variation of registration etc. .

Signature for and on behalf of Manufacturer, title of the authorized person, seal



Arthur Loftus Regional Director, EMEA



Brandon Medical Company Limited | Elmfield Road | Morley | Leeds | UK | LS27 0EL T:+44 113 277 7393 / 0845 1243 666 F:+44 113 272 8844 / 0845 1243 667

E: enquiries@brandon-medical.com W: www.brandon-medical.com Brandon Medical Company Limited. Registered in England & Wales No. 02827189 VAT No. GB734 591421. A subsidiary of Brandon Group Limited.

A subsidiary of brancom local periods. Brandon Medical is always looking for opportunities to improve our products and customer service and we greatly value any feedback you can provide: <u>feedback@brandon-medical.com</u>





CE

Brandon Medical Company Limited (BMCL) Elmfield Road, Morley, Leeds, West Yorkshire, LS27 0EL, UK SRN (Non-EU Manufacturer): GB-MF-000005501

SRN (European Authorised Representative): IE-AR-000003999

Declares that under our sole responsibility, the product:

**Coolview 53** (Refer to end of declaration for variant UDI-DI & Basic UDI-DI)

With intended purpose of:

#### **Examination Lamp**

**Coolview 50** 

Is in conformity to the general safety and performance requirements of the following EC Regulation:

- 2017/745 Medical Device Regulation Risk Classification as per Annex VIII
- **Class I Self-Declaration Medical Device**
- Conformity assessment route: Class I: EC conformity self-declaration according to MDR 2017/745 Article 19, and Annexes II and III.


The following standards were applied in full:

- EN ISO 13485:2016 Medical devices Quality management systems Requirements for regulatory . purposes
- ISO 9001:2015 Quality management systems Requirements
- BS EN ISO 14971:2012 Medical devices Application of risk management to medical devices
- EN ISO 15223-1:2016 Medical devices Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
- EN60601-1:2006+A1:2013 Medical electrical equipment Part 1: General requirements for basic . safety and essential performance
- BS EN60601-1-2:2015 Medical electrical equipment. General requirements for basic safety and essential performance. Collateral Standard: Electromagnetic disturbances. Requirements and tests
- BS EN 60601-2-41 2009+A11:2011 Particular requirements for basic safety and essential performance of surgical luminaires and luminaires for diagnosis DataContro
- EN 60601-6:2010+A1:2015
- Waste from Electrical and Electronic Equipment (WEEE) Directive

This product has been developed, tested and verified in accordance with the international standards lister above.

| Person Responsible for<br>Regulatory Compliance (PRRC)                           | Signature | Date        |
|--|-----------|-------------|
| Annette Hemmings<br>Head of Quality Assurance<br>Brandon Medical Company Limited | A. E. Q.  | 17 MAY 2022 |
| Countersignature   | Signature | Date        |
| Tom Partington<br>Product Manager<br>Brandon Medical Company Limited             | N         | 17-05-22    |

European Authorised Representative: European Healthcare & Device Solutions (Ireland) Ltd. Stratton House, Bishopstown Road, Bishopstown, Cork, Ireland





Brandon

ILLIANT BY DESIGH

EUDoC CL5053 Revision 3

Page 2 of 2