

17. MAI. 2021

CERTIFICAT

Prin prezentul se confirmă faptul deținerii de către **"LISMEDFARM" SRL**, IDNO **1003600113573**, a contului curent în MDL deschis la BC "EXIMBANK" S.A. Sucursala Nr. 19, codul băncii EXMMMD22, IBAN MD73EX0000000222401475MD.

Certificatul este eliberat pentru a fi prezentat după destinație.

Iurie Bețivu
Director Sucursala nr. 19



Ex. N.Jechiu
Tel.022 54 33 63

REPUBLICA



MOLDOVA

CERTIFICAT DE ÎNREGISTRARE

SOCIETATEA CU RĂSPUNDERE LIMITATĂ "LISMEDFARM"
ESTE ÎNREGISTRATĂ LA CAMERA ÎNREGISTRĂRII DE STAT

Numărul de indentificare de stat - codul fiscal

1003600113573

Data înregistrării


07.12.2001

Data eliberării

17.01.2005

Bobeica Ion, registrator de stat

*Funcția, numele, prenumele persoanei
care a eliberat certificatul*


semnătura



MD 0018019

L.Ș.

I.P. "AGENȚIA SERVICII PUBLICE"

Departamentul înregistrare și licențiere a unităților de drept

**Extras
din Registrul de stat al persoanelor juridice
nr. 115057 din 19.07.2023**



Denumirea completă: **Societatea cu Răspundere Limitată "LISMEDFARM"**

Denumirea prescurtată: **"LISMEDFARM" S.R.L.**

Forma juridică de organizare: **Societate cu răspundere limitată**

Numărul de identificare de stat și codul fiscal: **1003600113573**

Data înregistrării de stat: **07.12.2001**

Sediu: **MD-2002, Șoseaua Muncești 167/B, mun. Chișinău, Republica Moldova**

Genurile de activitate:

1. Comerțul cu ridicata al produselor farmaceutice;
2. Comerțul cu amănuntul al produselor farmaceutice și de parfumerie;
3. Comerțul cu ridicata al parfumurilor și produselor cosmetice;
4. Importul, fabricarea, comercializarea, asistența tehnică și (sau) reparația dispozitivelor medicale și (sau) a opticii;
5. Comerțul cu ridicata al produselor alimentare, băuturilor și produselor din tutun;
6. Comerțul cu amănuntul neefectuat prin magazine;
7. Intermedieri pentru vânzarea unui asortiment larg de mărfuri;
8. Comerțul cu ridicata al mărfurilor nealimentare de larg consum;
9. Activități de cercetare a pieței și de sondaj al opiniei publice;
10. Publicitate;
11. Comerț cu ridicata cu cafea, ceai, cacao și condimente;
12. Comerț cu amănuntul al articolelor medicale și ortopedice, în magazine specializate;
13. Comerț cu amănuntul al produselor cosmetice și de parfumerie, în magazine specializate;
14. Comerț cu amănuntul al produselor alimentare, băuturilor și produselor din tutun efectuat prin standuri, chioșcuri și piețe;

Capitalul social: **918180 Lei**

Administrator: **CHITIC ECATERINA IDNP 0970808153221**

Asociați:

1. **CHITIC ECATERINA (IDNP 0970808153221), partea socială 918180 Lei, ce constituie 100%**

Beneficiari efectivi: **CHITIC ECATERINA (IDNP 0970808153221)**

Prezentul extras este eliberat în temeiul art. 34 al Legii nr.220/2007 privind înregistrarea de stat a persoanelor juridice și a întreprinzătorilor individuali și confirmă datele din Registrul de stat la data de 19.07.2023

Specialist coordonator

Ludmila Ciur

tel. 022-207-838

**LISMEDFARM S.R.L.**

Șos. Muncești, 167/B, MD – 2002, mun. Chișinău, Republica Moldova
c/f: 1003600113573, TVA: 0304618, web: <https://lismedfarm.md>,
tel.: 022-80-47-98, 022-55-64-38, 022-56-94-91, e-mail: oficiu@lismedfarm.md

Număr B-012/2021 din 22.04.2021

Către: TOATE PERSOANELE INTERESATE

Lista Fondatorilor

Chitic Ecaterina Nicolai, anul nașterii: 30 decembrie 1969, cetățean al Republicii Moldova (MD),
cod personal: 0970808153221.

Cotă: 100%

Cu respect,

Directorul companiei “Lismedfarm” SRL, Ecaterina Chitic

Ex.: Vlad Chitic

Director executiv

+373 799 81 005

vlad.chitic@lismedfarm.md



GUVERNUL
REPUBLICII
MOLDOVA



SERVICIUL FISCAL DE STAT



CERTIFICAT

privind lipsa sau existența restanțelor față de bugetul public național

Nr.
Nº

Din
От

1

DATE DESPRE CONTRIBUABIL/ИНФОРМАЦИЯ О НАЛОГОПЛАТЕЛЬЩИКЕ

Codul fiscal/Numărul de identificare

Фискальный код/Идентификационный номер

Denumirea

Наименование

2

ATESTAREA LIPSEI SAU EXISTENȚEI RESTANȚELOR CONFORM DATELOR SISTEMULUI INFORMAȚIONAL AUTOMATIZAT/ПОДТВЕРЖДЕНИЕ ОТСУТСТВИЯ ИЛИ НАЛИЧИЯ НЕДОИМКИ СОГЛАСНО ДАННЫХ ИНФОРМАЦИОННОЙ АВТОМАТИЗИРОВАННОЙ СИСТЕМЫ

La data emiterii prezentului certificat restanța față de bugetul public național constituie/На дату выдачи данной справки, недоимка перед национальным публичным бюджетом составляет

lei/лей

3

VALABIL PÂNĂ LA/ДЕЙСТВИТЕЛЕН ДО



Prezentul certificat este eliberat în temeiul Art. 131, alin. (5³) din Codul fiscal nr. 1163/1997, în baza datelor furnizate de Serviciul Fiscal de Stat în Portalul Guvernamental al Antreprenorului/Сертификат выдан в соответствии со ст. 131, п. (5³) Налогового кодекса №1163/1997, на основании данных, предоставленных Государственной налоговой службой на Портале Правительства Предпринимателя

Generat și semnat de Portalul Guvernamental al Antreprenorului (<https://mcabinet.gov.md>) la

Prezentul certificat este semnat electronic în conformitate cu Legea nr. 124 din 19.05.2022 /
Сертификат подписан электронной подписью в соответствии с Законом № 124 от 19.05.2022

Certificatul este descărcat de pe Portalul Guvernamental al Antreprenorului (<https://mcabinet.gov.md>) și este semnat electronic de către posesorul acestui portal și are aceeași valoare juridică ca și documentele eliberate pe suport de hârtie de către organele cu atribuții de administrare fiscală.
Verificarea autenticității semnăturii electronice poate fi realizată la adresa: <https://msign.gov.md>.

Сертификат выгружен с Правительственного Портала Предпринимателя (<https://mcabinet.gov.md>) и подписан электронной подписью владельца портала и имеет такую же юридическую силу, как и документы, выдаваемые на бумаге органами налоговой администрации.
Проверку подлинности электронной подписи можно осуществить по адресу: <https://msign.gov.md>.

Recipisa

Respondent

Codul fiscal: 1003600113573, denumire: LISMEDFARM S.R.L.

A prezentat raportul: RSF1_21

Pentru perioada fiscala: A/2022

Data prezentarii: 15.05.2023

Marca temporală a raportului înregistrat în Sistemul de Raportare Electronică și expediat pentru procesare în Sistemul Informațional al BNS : 15.05.2023 14:06:58

Recipisa 2

Respondent

Codul fiscal: 1003600113573, denumire: LISMEDFARM S.R.L.

A prezentat raportul: RSF1_21

Pentru perioada fiscala: A/2022

Data prezentarii: 15.05.2023

Marca temporală a raportului înregistrat în Sistemul Informațional al BNS : 15.05.2023 15:10:31

Biroul Național de Statistică (BNS) a recepționat varianta electronică a raportului, expediat de DVs. Urmează verificarea și validarea raportului de către specialistul BNS pe domeniu.

SITUAȚIILE FINANCIARE

pentru perioada 01.01.2022 - 31.12.2022

Entitatea: LISMEDFARM S.R.L.

Cod CUIFO: 38183026

Cod IDNO: 1003600113573

Sediul:

MD: 2002

Raionul(municipiul): 103, DDF BOTANICA

Cod CUATM: 0110, SEC.BOTANICA

Strada: sos. Muncești, 167/B, mun. Chisinau

Activitatea principală: G4646, Comerț cu ridicata al produselor farmaceutice

Forma de proprietate: 15, Proprietatea privată

Forma organizatorico-juridică: 530, Societăți cu răspundere limitată

Date de contact:

Telefon: 022804799

WEB:

E-mail: contabil@lismedfarm.md

Numele și coordonatele al contabilului-șef: DI (dna) Alii Veronica Tel.

Numărul mediu al salariaților în perioada de gestiune: 22 persoane.

Persoanele responsabile de semnarea situațiilor financiare* Chitic Ecaterina

Unitatea de măsură: leu

BILANȚUL

Anexa 1

la 31.12.2022

Nr. cpt.	Indicatori	Cod rd.	Sold la	
			Începutul perioadei de gestiune	Sfârșitul perioadei de gestiune
1	2	3	4	5
	A C T I V			
A.	ACTIVE IMOBILIZATE			
	I. Imobilizări necorporale			
	1. Imobilizări necorporale în curs de execuție	010	0	0
	2. Imobilizări necorporale în exploatare, total	020	11931	11906
	din care:			
	2.1. concesiuni, licențe și mărci	021	0	0
	2.2. drepturi de autor și titluri de protecție	022	0	0
	2.3. programe informatice	023	11011	11906
	2.4. alte imobilizări necorporale	024	920	0
	3. Fond comercial	030	0	0
	4. Avansuri acordate pentru imobilizări necorporale	040	0	0
	Total imobilizări necorporale (rd.010 + rd.020 + rd.030 + rd.040)	050	11931	11906
	II. Imobilizări corporale			
	1. Imobilizări corporale în curs de execuție	060	0	0

	2. Terenuri	070	233698	233698
	3. Mijloace fixe, total	080	8771294	8841484
	din care:			
	3.1. clădiri	081	8006266	7716121
	3.2. construcții speciale	082	427538	270381
	3.3. mașini, utilaje și instalații tehnice	083	224939	198435
	3.4. mijloace de transport	084	58777	614026
	3.5. inventar și mobilier	085	36441	10503
	3.6. alte mijloace fixe	086	17333	32018
	4. Resurse minerale	090	0	0
	5. Active biologice imobilizate	100	0	0
	6. Investiții imobiliare	110	0	0
	7. Avansuri acordate pentru imobilizări corporale	120	0	0
	Total imobilizări corporale (rd.060 + rd.070 + rd.080 + rd.090 + rd.100 + rd.110 + rd.120)	130	9004992	9075182
	III. Investiții financiare pe termen lung			
	1. Investiții financiare pe termen lung în părți neafiliate	140	0	0
	2. Investiții financiare pe termen lung în părți afiliate, total	150	0	0
	din care:			
	2.1. acțiuni și cote de participație deținute în părțile afiliate	151	0	0
	2.2 împrumuturi acordate părților afiliate	152	0	0
	2.3 împrumuturi acordate aferente intereselor de participare	153	0	0
	2.4 alte investiții financiare	154	0	0
	Total investiții financiare pe termen lung (rd.140 + rd.150)	160	0	0
	IV. Creanțe pe termen lung și alte active imobilizate			
	1. Creanțe comerciale pe termen lung	170	0	0
	2. Creanțe ale părților afiliate pe termen lung	180	0	0
	inclusiv: creanțe aferente intereselor de participare	181	0	0
	3. Alte creanțe pe termen lung	190	0	0
	4. Cheltuieli anticipate pe termen lung	200	0	0
	5. Alte active imobilizate	210	0	0
	Total creanțe pe termen lung și alte active imobilizate (rd.170 + rd.180 + rd.190 + rd.200 + rd.210)	220	0	0
	TOTAL ACTIVE IMOBILIZATE (rd.050 + rd.130 + rd.160 + rd.220)	230	9016923	9087088
B.	ACTIVE CIRCULANTE			
	I. Stocuri			
	1. Materiale și obiecte de mică valoare și scurtă durată	240	15367	25806
	2. Active biologice circulante	250	0	0
	3. Producția în curs de execuție	260	0	0
	4. Produse și mărfuri	270	17299426	19326489
	5. Avansuri acordate pentru stocuri	280	1149545	5040292
	Total stocuri (rd.240 + rd.250 + rd.260 + rd.270 + rd.280)	290	18464338	24392587

II. Creanțe curente și alte active circulante			
1. Creanțe comerciale curente	300	8526798	15129505
2. Creanțe ale părților afiliate curente	310	0	0
inclusiv: creanțe aferente intereselor de participare	311	0	0
3. Creanțe ale bugetului	320	2946435	3461127
4. Creanțele ale personalului	330	0	31000
5. Alte creanțe curente	340	130403	33345
6. Cheltuieli anticipate curente	350	7028	14153
7. Alte active circulante	360	0	0
Total creanțe curente și alte active circulante (rd.300 + rd.310 + rd.320 + rd.330 + rd.340 + rd.350 + rd.360)	370	11610664	18669130
III. Investiții financiare curente			
1. Investiții financiare curente în părți neafiliate	380	0	0
2. Investiții financiare curente în părți afiliate, total	390	0	0
din care:			
2.1. acțiuni și cote de participație deținute în părțile afiliate	391	0	0
2.2. împrumuturi acordate părților afiliate	392	0	0
2.3. împrumuturi acordate aferente intereselor de participare	393	0	0
2.4. alte investiții financiare în părți afiliate	394	0	0
Total investiții financiare curente (rd.380 + rd.390)	400	0	0
IV. Numerar și documente bănești	410	8479693	1975580
TOTAL ACTIVE CIRCULANTE (rd.290 + rd.370 + rd.400 + rd.410)	420	38554695	45037297
TOTAL ACTIVE (rd.230 + rd.420)	430	47571618	54124385
P A S I V			
C.	CAPITAL PROPRIU		
	I. Capital social și neînregistrat		
1. Capital social	440	918181	918181
2. Capital nevărsat	450	(0)	(0)
3. Capital neînregistrat	460	0	0
4. Capital retras	470	(0)	(0)
5. Patrimoniul primit de la stat cu drept de proprietate	480	0	0
Total capital social și neînregistrat (rd.440 + rd.450 + rd.460 + rd.470 + rd.480)	490	918181	918181
II. Prime de capital	500	0	0
III. Rezerve			
1. Capital de rezervă	510	0	0
2. Rezerve statutare	520	4050	4050
3. Alte rezerve	530	0	0
Total rezerve (rd.510 + rd.520 + rd.530)	540	4050	4050
IV. Profit (pierdere)			
1. Corecții ale rezultatelor anilor precedenți	550	X	0

	2. Profit nerepartizat (pierdere neacoperită) al anilor precedenți	560	26453997	24612816
	3. Profit net (pierdere netă) al perioadei de gestiune	570	X	11187534
	4. Profit utilizat al perioadei de gestiune	580	X	(0)
	Total profit (pierdere) (rd.550 + rd.560 + rd.570 + rd.580)	590	26453997	35800350
	V. Rezerve din reevaluare	600	0	0
	VI. Alte elemente de capital propriu	610	0	0
	TOTAL CAPITAL PROPRIU (rd.490 + rd.500 + rd.540 + rd.590 + rd.600 + rd.610)	620	27376228	36722581
D.	DATORII PE TERMEN LUNG			
	1. Credite bancare pe termen lung	630	0	0
	2. Împrumuturi pe termen lung	640	7948370	8026112
	din care:			
	2.1. împrumuturi din emisiunea de obligațiuni	641	0	0
	inclusiv: împrumuturi din emisiunea de obligațiuni convertibile	642	0	0
	2.2. alte împrumuturi pe termen lung	643	7948370	8026112
	3. Datorii comerciale pe termen lung	650	0	0
	4. Datorii față de părțile afiliate pe termen lung	660	0	0
	inclusiv: datorii aferente intereselor de participare	661	0	0
	5. Avansuri primite pe termen lung	670	0	0
	6. Venituri anticipate pe termen lung	680	0	0
	7. Alte datorii pe termen lung	690	0	0
	TOTAL DATORII PE TERMEN LUNG (rd.630 + rd.640 + rd.650 + rd.660 + rd.670 + rd.680 + rd.690)	700	7948370	8026112
E.	DATORII CURENTE			
	1. Credite bancare pe termen scurt	710	0	0
	2. Împrumuturi pe termen scurt, total	720	0	0
	din care:			
	2.1. împrumuturi din emisiunea de obligațiuni	721	0	0
	inclusiv: împrumuturi din emisiunea de obligațiuni convertibile	722	0	0
	2.2. alte împrumuturi pe termen scurt	723	0	0
	3. Datorii comerciale curente	730	12152222	8938350
	4. Datorii față de părțile afiliate curente	740	0	0
	inclusiv: datorii aferente intereselor de participare	741	0	0
	5. Avansuri primite curente	750	25912	72026
	6. Datorii față de personal	760	0	0
	7. Datorii privind asigurările sociale și medicale	770	0	72075
	8. Datorii față de buget	780	67035	43974
	9. Datorii față de proprietari	790	0	0
	10. Venituri anticipate curente	800	0	0
	11. Alte datorii curente	810	1851	249267
	TOTAL DATORII CURENTE (rd.710 + rd.720 + rd.730 + rd.740 + rd.750 + rd.760 + rd.770 + rd.780 + rd.790 + rd.800 + rd.810)	820	12247020	9375692

F.	PROVIZIOANE			
	1. Provizioane pentru beneficiile angajaților	830	0	0
	2. Provizioane pentru garanții acordate cumpărătorilor/clientilor	840	0	0
	3. Provizioane pentru impozite	850	0	0
	4. Alte provizioane	860	0	0
	TOTAL PROVIZIOANE (rd.830 + rd.840 + rd.850 + rd.860)	870	0	0
	TOTAL PASIVE (rd.620 + rd.700 + rd.820 + rd.870)	880	47571618	54124385

SITUAȚIA DE PROFIT ȘI PIERDERE

de la 01.01.2022 până la 31.12.2022

Anexa 2

Indicatori	Cod rd.	Perioada de gestiune	
		precedenta	curenta
1	2	3	4
Venituri din vânzări, total	010	84401052	95300233
din care:			
venituri din vânzarea produselor și mărfurilor	011	84362297	95299733
venituri din prestarea serviciilor și executarea lucrărilor	012	38755	500
venituri din contracte de construcție	013	0	0
venituri din contracte de leasing	014	0	0
venituri din contracte de microfinanțare	015	0	0
alte venituri din vânzări	016	0	0
Costul vânzărilor, total	020	76106363	84519235
din care:			
valoarea contabilă a produselor și mărfurilor vândute	021	76106363	84519235
costul serviciilor prestate și lucrărilor executate terților	022	0	0
costuri aferente contractelor de construcție	023	0	0
costuri aferente contractelor de leasing	024	0	0
costuri aferente contractelor de microfinanțare	025	0	0
alte costuri aferente vânzărilor	026	0	0
Profit brut (pierdere brută) (rd.010 - rd.020)	030	8294689	10780998
Alte venituri din activitatea operațională	040	11159619	8316357
Cheltuieli de distribuie	050	4520275	4965557
Cheltuieli administrative	060	1101278	1259190
Alte cheltuieli din activitatea operațională	070	274308	192839
Rezultatul din activitatea operațională: profit (pierdere) (rd.030 + rd.040 - rd.050 - rd.060 - rd.070)	080	13558447	12679769
Venituri financiare, total	090	1878801	2138571
din care:			
venituri din interese de participare	091	0	0
inclusiv: veniturile obținute de la părțile afiliate	092	0	0
venituri din dobânzi	093	9577	1830
inclusiv: veniturile obținute de la părțile afiliate	094	0	0
venituri din alte investiții financiare pe termen lung	095	0	12028
inclusiv: veniturile obținute de la părțile afiliate	096	0	0

venituri aferente ajustărilor de valoare privind investițiile financiare pe termen lung și curente	097	0	0
venituri din ieșirea investițiilor financiare	098	0	0
venituri aferente diferențelor de curs valutar și de sumă	099	1869224	2124713
Cheltuieli financiare, total	100	2259094	2070177
din care:			
cheltuieli privind dobânzile	101	0	0
inclusiv: cheltuielile aferente părților afiliate	102	0	0
cheltuieli aferente ajustărilor de valoare privind investițiile financiare pe termen lung și curente	103	0	0
cheltuieli aferente ieșirii investițiilor financiare	104	0	0
cheltuieli aferente diferențelor de curs valutar și de sumă	105	2152103	1905887
Rezultatul: profit (pierdere) financiar(ă) (rd.090 - rd.100)	110	-380293	68394
Venituri cu active imobilizate și excepționale	120	7270	5310
Cheltuieli cu active imobilizate și excepționale	130	0	0
Rezultatul din operațiuni cu active imobilizate și excepționale: profit (pierdere) (rd.120 - rd.130)	140	7270	5310
Rezultatul din alte activități: profit (pierdere) (rd.110 + rd.140)	150	-373023	73704
Profit (pierdere) până la impozitare (rd.080 + rd.150)	160	13185424	12753473
Cheltuieli privind impozitul pe venit	170	1602850	1565939
Profit net (pierdere netă) al perioadei de gestiune (rd.160 - rd.170)	180	11582574	11187534

SITUAȚIA MODIFICĂRILOR CAPITALULUI PROPRIU

de la 01.01.2022 până la 31.12.2022

Anexa 3

Nr. d/o	Indicatori	Cod rd	Sold la începutul perioadei de gestiune	Majorări	Diminuări	Sold la sfârșitul perioadei de gestiune
1	2	3	4	5	6	7
I.	Capital social și neînregistrat					
	1. Capital social	010	918181	0	0	918181
	2. Capital nevărsat	020	(0)	(0)	(0)	(0)
	3. Capital neînregistrat	030	0	0	0	0
	4. Capital retras	040	(0)	(0)	(0)	(0)
	5. Patrimoniul primit de la stat cu drept de proprietate	050	0	0	0	0
	Total capital social și neînregistrat (rd.010 + rd.020 + rd.030 + rd.040 + rd.050)	060	918181	0	0	918181
II.	Prime de capital	070	0	0	0	0
III.	Rezerve					
	1. Capital de rezervă	080	0	0	0	0
	2. Rezerve statutare	090	4050	0	0	4050
	3. Alte rezerve	100	0	0	0	0
	Total rezerve (rd.080 + rd.090 + rd.100)	110	4050	0	0	4050
IV.	Profit (pierdere)					

	1. Corecții ale rezultatelor anilor precedenți	120	X	0	0	0
	2. Profit nerepartizat (pierdere neacoperită) al anilor precedenți	130	26453997	0	1841181	24612816
	3. Profit net (pierdere netă) al perioadei de gestiune	140	X	11187534	0	11187534
	4. Profit utilizat al perioadei de gestiune	150	X	(0)	(0)	(0)
	Total profit (pierdere) (rd.120 + rd.130 + rd.140 + rd.150)	160	26453997	11187534	1841181	35800350
V.	Rezerve din reevaluare	170	0	0	0	0
VI.	Alte elemente de capital propriu	180	0	0	0	0
	Total capital propriu (rd.060 + rd.070 + rd.110 + rd.160 + rd.170 + rd.180)	190	27376228	11187534	1841181	36722581

SITUAȚIA FLUXURILOR DE NUMERAR

de la 01.01.2022 până la 31.12.2022

Anexa 4

Indicatori	Cod rd	Perioada de gestiune	
		precedentă	curentă
1	2	3	4
Fluxuri de numerar din activitatea operațională			
Încasări din vânzări	010	65297079	66148952
Plăți pentru stocuri și servicii procurate	020	51775872	69150464
Plăți către angajați și organe de asigurare socială și medicală	030	2065727	2343801
Dobânzi plătite	040	0	0
Plata impozitului pe venit	050	2439641	1172166
Alte încasări	060	2187716	4428799
Alte plăți	070	3628918	3008563
Fluxul net de numerar din activitatea operațională (rd.010 - rd.020 - rd.030 - rd.040 - rd.050 + rd.060 - rd.070)	080	7574637	-5097243
Fluxuri de numerar din activitatea de investiții			
Încasări din vânzarea activelor imobilizate	090	0	0
Plăți aferente intrărilor de active imobilizate	100	339536	0
Dobânzi încasate	110	9577	1830
Dividende încasate	120	0	0
inclusiv: dividende încasate din străinătate	121	0	0
Alte încasări (plăți)	130	0	0
Fluxul net de numerar din activitatea de investiții (rd.090 - rd.100 + rd.110 + rd.120 ± rd.130)	140	-329959	1830
Fluxuri de numerar din activitatea financiară			
Încasări sub formă de credite și împrumuturi	150	0	399000
Plăți aferente rambursării creditelor și împrumuturilor	160	0	430000
Dividende plătite	170	10238950	1730709
inclusiv: dividende plătite nerezidenților	171	0	0
Încasări din operațiuni de capital	180	0	0
Alte încasări (plăți)	190	294399	1796
Fluxul net de numerar din activitatea financiară	200	-9944551	-1759913

(rd.150 - rd.160 - rd.170 + rd.180 ± rd.190)			
Fluxul net de numerar total (± rd.080 ± rd.140 ± rd.200)	210	-2699873	-6855326
Diferențe de curs valutar favorabile (nefavorabile)	220	-92161	351213
Sold de numerar la începutul perioadei de gestiune	230	11271727	8479693
Sold de numerar la sfârșitul perioadei de gestiune (± rd.210 ± rd.220 + rd.230)	240	8479693	1975580

Documente atașate - Notă explicativă (fișierul pdf)

	LISMEDFARM S.R.L. Șos. Muncești, 167/B, MD – 2002, mun. Chișinău, Republica Moldova tel.: 022-80-47-98, 022-55-64-38, 022-56-94-91, e-mail: oficiu@lismedfarm.md, web: https://lismedfarm.md, c/f: 1003600113573, TVA: 0304618, director – Ecaterina Chitic
---	--

CENTRUL PENTRU ACHIZITII PUBLICE CENTRALIZATE IN SANATATE

Declarație

Declarăm următoarele pentru LP nr. ocds-b3wdp1-MD-1685346316289 din 04.08.2022:

1. Termenul de valabilitate restant al Bunurilor la momentul livrării va constitui: pentru articolele parafarmaceutice livrate de la 01.01.2024 - 30.06.2024, va constitui:
 - a. - nu mai puțin de 60% din termenul de valabilitate inițial pentru bunurile cu o valabilitate de 2 ani și mai mulți;
 - b. - și de minim 80% din cel inițial pentru Bunurile cu o valabilitate de până la 2 ani.
2. Termenul de valabilitate restant al Bunurilor la momentul livrării, pentru articolele parafarmaceutice livrate de la 01.07.2024-31.12.2024 va constitui:
 - a. - nu mai puțin de 40% din termenul de valabilitate inițial pentru bunurile cu o valabilitate de 2 ani și mai mulți;
 - b. - și de minim 60% din cel inițial pentru Bunurile cu o valabilitate de până la 2 ani.

Cu respect,

Director executiv Lismedfarm S.R.L.

Vlad Chitic



Tianchang City Anrui Medical Equipments Co., Ltd. is one of the biggest Latex Surgical Gloves & Foley Catheter manufacturers in China, founded in 2013 and since then it is well known for its first-class products and excellent service in domestic and abroad markets.

Tianchang City Anrui Medical Equipments Co., Ltd.

Mob: +86-13635507675

WhatsApp: +86-13635507675

Tel: 0086-550-7834188

Fax: 0086-550-7834699

Mailbox: Daniel@anruimedicequipment.com

Website: www.anruimedical.com

Address: Renhe Industrial District, Tianchang City, Anhui Province, China

Add (Office in Shanghai): 910A, Building C, No.1126, Shenbin S.Road, Minhang District, Shanghai, China



安瑞医疗
ANRUI MEDICAL



TIANCHANG CITY

ANRUI MEDICAL EQUIPMENTS CO., LTD.

www.anruimedical.com



CONTENTS

COMPANY PROFILE	02
LATEX SURGICAL GLOVES.....	03
LATEX FOLEY CATHETERS	05
SILICONE FOLEY CATHETERS.....	07
GYNACOLOGICAL GLOVES	09
STERILE MALE CIRCUMCISION KITS	10

CERTIFICATES



COMPANY PROFILE

Tianchang City Anrui Medical Equipments Co., Ltd, is one of the biggest Latex Surgical Gloves & Foley Catheter manufacturers in China, founded in 2013, Located in Tianchang City, Anhui Province. We mainly produce Latex Surgical Gloves, Gynecological Gloves, Latex Foley Catheter, Silicone Foley Catheter, Operation Kits and so on. With world-class products, excellent services, complete producing and quality management team, we serve our customers very well worldwide. Sincerely look forward to building long-term and good cooperation with you in near Future !



BLOWING TEST



WATER TIGHTNESS TEST



LATEX SURGICAL GLOVES

Quality Certification: CE; ISO13485; CFDA; FSC

Size: 6.0 6.5 7.0 7.5 8.0 8.5, 9.0

Pre-Powdered / Powder Free

Material: Natural Latex

Disinfecting Type: GAMMA Sterilization

Shelf Life: 5 years


HS: 40151200

MOQ:50000 PCS

Packing: 1 Pair/Pouch;50 Pairs/Box; 10 Boxes/ Carton

Dimension: 67 X 25 X 40 CM

REF: SS0602





LATEX FOLEY CATHETERS

Quality Certification: CE, ISO 13485,

Material: Latex (100% Silicone coated)

Disinfecting Type: EO Sterilization

Shelf Life: 5 years

HS: 9018390000

MOQ: 2000 PCS

Packing: 1 Pc/Pouch;10 Pcs/Box; 50Boxes/ Carton

Dimension: 54 X 52 X 40 CM



Products Name	Size (Fr/Ch)	Balloon Capacity (ml)
2-Way	8 Fr	3 ml
	9 Fr	3.5 ml
	10 Fr	3.5 ml
	12 Fr	5-10ml
	14 Fr	5-15ml
	16 Fr	5-15ml
	18 Fr	5-15ml
	20 Fr	30 ml
	22 Fr	30 ml
	24 Fr	30 ml
	26 Fr	30 ml
	28 Fr	30 ml
	30 Fr	30 ml
	32 Fr	30 ml
3-Way	14 Fr	5-15ml
	16 Fr	5-15ml
	18 Fr	5-15ml
	20 Fr	30 ml
	22 Fr	30 ml
	24 Fr	30 ml
	26 Fr	30 ml
	28 Fr	30 ml
	30 Fr	30 ml
	32 Fr	30 ml





SILICONE FOLEY CATHETERS

Quality Certification: CE, ISO 13485,

Material: Silicone

Disinfecting Type: EO Sterilization

Shelf Life: 5 years

HS: 9018390000

MOQ: 2000 PCS

Packing: 1 Pc/Pouch;10 Pcs/Box; 20Boxes/ Carton

Dimension: 52 X 32 X 27 CM



Products Name	Size (Fr/Ch)	Balloon Capacity (ml)
2-Way	6 Fr	3ml
	8 Fr	3-5 ml
	10 Fr	3-5 ml
	12 Fr	5-10ml
	14 Fr	5-15ml
	16 Fr	5-15ml
	18 Fr	5-15ml
	20 Fr	30 ml
	22 Fr	30 ml
	24 Fr	30 ml
3-Way	14 Fr	5-15ml
	16 Fr	5-15ml
	18 Fr	5-15ml
	20 Fr	30 ml
	22 Fr	30 ml
	24 Fr	30 ml
	26 Fr	30 ml



GYNAECOLOGICAL GLOVES

Quality Certification: CE; ISO13485; CFDA; FSC

Size: S M L XL

Pre-Powdered / Powder Free

Material: Natural Latex

Disinfecting Type: GAMMA Sterilization

Shelf Life: 5 years

HS: 40151200

MOQ: 20000 PCS

Packing: 1 Pair/Pouch; 25 Pairs/Box; 4 Boxes/ Carton

Dimension: 38 X 36 X 38 CM



STERILE MALE CIRCUMCISION KITS

- Plastic Trays 1 Piece
- Stainless Steel Forceps 20cm
- Stainless Steel Forceps 14cm
- Stainless Steel Forceps 14cm Non-Toothed
- Dissecting Scissors Curved
- Stainless Steel Stitch Cutting Scissors
- Stainless Steel Needle Holder
- Retractable Scalpel Blade
- Chronic Gut 3/0 Suture
- Paraffin Gauze Swab
- Gauze Swab 10cm x 10cm
- Hypo Allergenic Tape 12.5cm x 3m
- Alcohol Pad
- Needle 21 g
- Needle 23 g
- Syringe 10 ml
- Disposable Apron
- Fenestrated Drape (Green)
- Examination Gloves
- Surgical Gloves



Technical Construction File

(File No.: AR/CE01)

Version: A/6

Surgical Gloves

According to

Medical Directive 93/42/EEC

Compiled by: Colin Lv Date: 2023.03.01

Reviewed by: Zhang Bi Wen Date: 2023.03.01

Approved by: Zhang Bi Wen Date: 2023.03.01

Tianchang City Anrui Medical Equipments Co., Ltd.,
Renhe Industrial District, 239331 Tianchang City, Anhui Province,
People's Republic of China

Index

No.	Doc. No.	Doc. Name	Ver.	Effective date
Part A				
1	AR/CE01-01	Declaration of Conformity	A/6	2024.05.26
2	AR/CE01-02	Brief Introduction	A/6	2024.05.26
3	AR/CE01-03	List of Harmonized Standard	A/6	2024.05.26
4	AR/CE01-04	Essential Requirements Checklist	A/6	2024.05.26
5	AR/CE01-05	Risk Analysis Plan	A/6	2024.05.26
6	AR/CE01-06	Product Descriptions	A/6	2024.05.26
7	AR/CE01-07	Label and Instruction for Use	A/6	2024.05.26
8	AR/CE01-08	Pre-clinical Assessment	A/6	2024.05.26
9	AR/CE01-09	Clinical Evaluation	A/6	2024.05.26
10	AR/CE01-10	Production and Quality Control	A/6	2024.05.26
11	AR/CE01-11	Usability Evaluation	A/1	2024.05.26
12	AR/CE01-12	Packaging Verification Summary	A/6	2024.05.26
13	AR/CE01-13	Sterilization Verification Summary	A/6	2024.05.26
14	AR/CE01-14	List of CE Marking Products	A/6	2024.05.26
15	AR/CE01-15	Medical Device Post Market Surveillance	A/6	2024.05.26
Part B				
14	Appendix 1	EC Authority Representative Agreement		
15	Appendix 2	Instruction for use, Labeling		
16	Appendix 3	Product type test reports		
17	Appendix 4	Product Bio-compatibility test		
18	Appendix 5	Packaging Validation Report		
19	Appendix 6	Sterilization Validation Report		
20	Appendix 7	Risk Analysis Report		
21	Appendix 8	Accelerated Aging Test Report		张必文

List of Harmonized Standard Surgical Gloves			Prepared by	Zhang Bi Yong
			Checked by	Zhang Bi Wen
Doc. No.	AR/CE01-03		Zhang Bi Wen	Zhang Bi Wen
Effective date	2024.05.26	Ver. A/6	Page No.	Page 1 of 2

List of harmonized standard

List of European Union Harmonized Standards which this product applies.

Standard/Directive	Name of document
93/42/EEC	EU Council Directive concerning medical devices
EN 455-1:2020	Medical gloves for single use - Part 1: Requirements and testing for freedom from holes
EN 62366:2015	Medical devices - Part 1: Application of usability engineering to medical devices
EN 455-3:2015	Medical gloves for single use - Part 3: Requirements and testing for biological evaluation
EN 455-4:2009	Medical gloves for single use - Part 4: Requirements and testing for shelf life determination
EN ISO 20417:2021	Medical devices - Information to be supplied by the manufacturer
EN ISO 10993-1:2018	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
EN ISO 11137-1:2015	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
EN ISO 11137-2:2015	Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose
EN ISO 11607-1:2020	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO 11607-2:2020	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
EN ISO 11737-1:2018 +A1:2021	Sterilization of health care products - Microbiological methods - Part 1: Determination of a population of microorganisms on products
EN ISO 11737-2:2020	Sterilization of health care products - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
EN ISO 14971:2019 +A11:2021	Medical devices - Application of risk management to medical devices
EN ISO 15223-1:2021	Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements

List of Harmonized Standard Surgical Gloves			Prepared by	Zhang Bi Yong
			Checked by	Zhang Bi Wen
Doc. No.	AR/CE01-03		Zhang Bi Wen	Zhang Bi Wen
Effective date	2024.05.26	Ver. A/6	Page No.	Page 2 of 2

EN 455-2:2015	Medical gloves for single use - Part 2: Requirements and testing for physical properties
ISO 10282:2014	Single-use sterile rubber surgical gloves - Specification
ISO 10993-10:2010	Biological Evaluation of Medical Device-Part 10: stimulation and allergic reaction
EN ISO 374-5:2016	Protective gloves against dangerous chemicals and micro-organisms — Part 5: Terminology and performance requirements for micro-organisms risks
EN ISO 374-1:2016	Protective gloves against dangerous chemicals and micro-organisms — Part 1: Terminology and performance requirements for chemical risks
EN ISO 374-4:2019	Protective gloves against dangerous chemicals and micro-organisms — Part 4: Determination of resistance to degradation by chemicals
EN ISO 374-2:2019	Protective gloves against dangerous chemicals and micro-organisms — Part 2: Determination of resistance to penetration



张必文.

Product Descriptions Surgical Gloves			Prepared by	Colin Lv
			Checked by	Zhang Bi Wen
Doc. No.	AR/CE01-06		Approved by	Zhang Bi Wen
Effective date	2024.05.26	Ver. A/6	Page No.	Page 1 of 3

Product Descriptions

6.1 Description of the product

Product name: Surgical Gloves


Intended use: Surgical Gloves is intended for use in surgery operation by doctors to help prevent contamination, cross-infection. It is for single-use and sterile unless the package is opened or damaged.

Materials: Natural latex

Validity: 5 years.

6.2 Structure and parameters

6.2.1 Device structure



6.2.2 Device specifications

Unit: mm

No.	Specification	Code	Width	Min Length	Min Thickness
1	6#	6	77±5	260	Smooth: 0.10 Textured: 0.13
2	6.5#	6.5	83±5	260	
3	7#	7	89±5	270	
4	7.5#	7.5	95±5	270	
5	8#	8	102±6	270	

Product Descriptions Surgical Gloves			Prepared by	Colin Lv
			Checked by	Zhang Bi Wen
Doc. No.	AR/CE01-06		Approved by	Zhang Bi Wen
Effective date	2024.05.26	Ver. A/6	Page No.	Page 2 of 3

6	8.5#	8.5	108±6	280	
7	9#	9	114±6	280	
According to powder, there are two types: pre-powdered, powder-free.					
According to finish, there are two types: textured surface over part or all of the glove, smooth surface..					
According to design, there are two types: gloves with straight fingers, gloves with fingers curved in the palmar direction.					

6.2.3 Device technical requirement

Property	Requirement
Water-tightness	No leakage shall be detected
Min. Force at break before accelerated aging /N	12.5
Min elongation at break before aging /%	700
Max force required to produce 300% elongation before accelerated aging /N	2.0
Min. Force at break after accelerated aging /N	9.5
Min elongation at break after aging /%	550

Powder Free gloves: Have a powder residual limit of 2.0 mg/piece.

6.2.4 Main materials and its suppliers

No.	Components/Materials name	Supplier name	Supplier address
1	Natural latex	Fuzhou Dongfangyuyuan Industry Co., Ltd.	Room 228, Science&Tech. Development Center, No.83 Zhujun Road, Fuzhou, China
		Zhejiang Xinwo Latex Products Co., Ltd.	Room 702, Unit 5, No.1066 Beichun Road, Yiwu City, Zhejiang, China
2	Internal packaging materials	HANGZHOU Z-SHINE MEDICAL PACKING CO.,LTD	Zhongsha Village, Dangshan Town, Xiaoshan District, Hangzhou, China
		Peak Medical Packaging Co., Ltd.	Yongda Road, Tianliao Industry Zone, Dongguan, China
3	Powder-free coating	Shanghai Longsha Trade	Room 203, No.46, Lane 700

Product Descriptions Surgical Gloves			Prepared by	Colin Lv
			Checked by	Zhang Bi Wen
Doc. No.	AR/CE01-06		Approved by	Zhang Bi Wen
Effective date	2024.05.26	Ver. A/6	Page No.	Page 3 of 3

		Development Co., Ltd.	Yongde Road, Minhang, Shanghai
4	Starch powder coating	Shanghai Dingzhao Trading Co., Ltd.	Room A130, Building 5, No.277 Huqingping Road, Minhang, Shanghai
5	E-beam sterilization	Yangzhou Yangfu Technology Co., Ltd.	Environmental Protection Technology Industrial Park, 225125 Yangzhou, PEOPLE'S REPUBLIC OF CHINA



张必文.



JUNDAMED

HUZHOU JUNDA IMPORT AND EXPORT CO., LTD.
ADD: ROOM 604,BUILDING1,NO. 567, LINGZHI EAST ROAD, ANJI COUNTY,
HUZHOU CITY, ZHEJIANG PROVINCE, CHINA
TEL: 0086 188 5726 1771

PRODUCT SPECIFICATION MD-1685346316289/34

PART I. PRODUCT IDENTIFICATION

Brand	JUNDAMED®
Code No.	GLVPPL
EAN13	4841948001423 (MD market)
Product Name	Vinyl examination gloves
Specification	Examination gloves, smooth or microtexture, non-sterile, vinyl, powdered, Size L. - Material: Polyvinylchlorid (PVC) - Unit of measure: piece - Packaging: 100 pieces - With microorganism protection properties; resistant to penetration of viruses and chemicals - Breaking resistance - Color: Clear - Ambidextrous
Characteristics	Latex protein-free, accelerator-free and DOP-free; suitable for people with type I latex allergy or type IV hypersensitivity. Smooth glove surface for a comfortable feel on patient contact.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	None
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD

Classification	Class I(Non-sterile)
Standard	EN 62366:2008, EN 1041:2008, EN ISO 15223-1:2016, EN455, EN420, EN374-2
Quality Management Systems & No.	None
	None
Production Quality Assurance System & No.	Not applicable
	Not applicable
CE Conformity Assessment Route	MDD 93/42/EEC Annex IX (Self-declaration)
Certification Body	Not applicable
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION

Size	Individual Packaging		Secondary Packaging		Outer Packaging		
	Qty	Material	Qty	Material	Qty	Dimensions	Material
L	100	Box	/	/	1000	33x25x26cm	carton

湖州骏达进出口有限公司
HUZHOU JUNDA IMPORT AND EXPORT CO., LTD



JUNDAMED

HUZHOU JUNDA IMPORT AND EXPORT CO., LTD.
ADD: ROOM 604,BUILDING1,NO. 567, LINGZHI EAST ROAD, ANJI COUNTY,
HUZHOU CITY, ZHEJIANG PROVINCE, CHINA
TEL: 0086 188 5726 1771

PRODUCT SPECIFICATION MD-1685346316289/35

PART I. PRODUCT IDENTIFICATION	
Brand	JUNDAMED®
Code No.	GLVPPM
EAN13	4841948001430 (MD market)
Product Name	Vinyl examination gloves
Specification	Examination gloves, smooth or microtexture, non-sterile, vinyl, powdered, Size M. - Material: Polyvinylchlorid (PVC) - Unit of measure: piece - Packaging: 100 pieces - With microorganism protection properties; resistant to penetration of viruses and chemicals - Breaking resistance - Color: Clear - Ambidextrous
Characteristics	Latex protein-free, accelerator-free and DOP-free; suitable for people with type I latex allergy or type IV hypersensitivity. Smooth glove surface for a comfortable feel on patient contact.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	None
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class I(Non-sterile)
Standard	EN 62366:2008, EN 1041:2008, EN ISO 15223-1:2016, EN455, EN420, EN374-2
Quality Management Systems & No.	None
	None
Production Quality Assurance System & No.	Not applicable
	Not applicable
CE Conformity Assessment Route	MDD 93/42/EEC Annex IX (Self-declaration)
Certification Body	Not applicable
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION						
Size	Individual Packaging		Secondary Packaging		Outer Packaging	
	Qty	Material	Qty	Material	Qty	Dimensions Material
M	100	Box	/	/	1000	33x25x26cm carton

湖州骏达进出口有限公司
HUZHOU JUNDA IMPORT AND EXPORT CO., LTD



JUNDAMED

HUZHOU JUNDA IMPORT AND EXPORT CO., LTD.
ADD: ROOM 604,BUILDING1,NO. 567, LINGZHI EAST ROAD, ANJI COUNTY,
HUZHOU CITY, ZHEJIANG PROVINCE, CHINA
TEL: 0086 188 5726 1771

PRODUCT SPECIFICATION MD-1685346316289/36

PART I. PRODUCT IDENTIFICATION	
Brand	JUNDAMED®
Code No.	GLVPPS
EAN13	4841948001447 (MD market)
Product Name	Vinyl examination gloves
Specification	Examination gloves, smooth or microtexture, non-sterile, vinyl, powdered, Size S. - Material: Polyvinylchlorid (PVC) - Unit of measure: piece - Packaging: 100 pieces - With microorganism protection properties; resistant to penetration of viruses and chemicals - Breaking resistance - Color: Clear - Ambidextrous
Characteristics	Latex protein-free, accelerator-free and DOP-free; suitable for people with type I latex allergy or type IV hypersensitivity. Smooth glove surface for a comfortable feel on patient contact.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	None
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class I(Non-sterile)
Standard	EN 62366:2008, EN 1041:2008, EN ISO 15223-1:2016, EN455, EN420, EN374-2
Quality Management Systems & No.	None
	None
Production Quality Assurance System & No.	Not applicable
	Not applicable
CE Conformity Assessment Route	MDD 93/42/EEC Annex IX (Self-declaration)
Certification Body	Not applicable
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION							
Size	Individual Packaging		Secondary Packaging		Outer Packaging		
	Qty	Material	Qty	Material	Qty	Dimensions	Material
S	100	Box	/	/	1000	33x25x26cm	carton

湖州骏达进出口有限公司
HUZHOU JUNDA IMPORT AND EXPORT CO., LTD



ANJI SPENQ INDUSTRIAL CO., LTD.
ADD: TANGPU ECONOMIC DEVELOPMENT ZONE, ANJI, ZHEJIANG, CHINA.
TEL: +86-572-5882803 FAX: +86-572-5882115 WWW.SPENQ.COM



PRODUCT SPECIFICATION MD-1685346316289/38

PART I. PRODUCT IDENTIFICATION	
Brand	SPENQ MEDICAL®
Code No.	SN01
EAN13	4841948000778 (MD market)
Product Name	Disposable Syringe
Specification	Disposable sterile syringes for insulin - non-removable steel / stainless steel needle, 30 G x ½ - Volume 1ml (U-100), Volume Measurement and International Units - consisting of 3 parts: transparent cylinder, piston, rubber gasket - seamlessly made - Hypoallergenic, non-pyrogenic needle, latex-free - individual packaging
Characteristics	The Barrel is transparent enough to allow easy measurement of the volume contained in the syringe and detection of air bubble. The plunger fits the inside of barrel very well, promising permitting freedom of movement.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	EO
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class II(sterile)
Standard	ISO 7864:2016, ISO 7886-1:2017, ISO 8537:2016 (indicated on secondary packing)
Quality Management Systems & No.	ISO13485 Q6 080946 0012 Rev. 02
Production Quality Assurance System & No.	CE(IIa) G2 080946 0007 Rev. 03
CE Conformity Assessment Route	MDD 93/42/EEC Annex V.3
Certification Body	TUV SUD Product Service GmbH
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION						
Size	Individual Packaging		Secondary Packaging		Outer Packaging	
	Qty	Material	Qty	Material	Qty	Dimensions Material
1ml	1	PE Bag	250	BAG	5000	64x43x51cm carton

安吉鑫洋进出口有限公司
ANJI SPENQ INDUSTRIAL CO., LTD

李强



ANJI SPENQ INDUSTRIAL CO., LTD.
ADD: TANGPU ECONOMIC DEVELOPMENT ZONE, ANJI, ZHEJIANG, CHINA.
TEL: +86-572-5882803 FAX: +86-572-5882115 WWW.SPENQ.COM



PRODUCT SPECIFICATION MD-1685346316289/39

PART I. PRODUCT IDENTIFICATION	
Brand	SPENQ MEDICAL®
Code No.	SN01
EAN13	4841948000808 (MD market)
Product Name	Disposable Syringe
Specification	Disposable sterile syringes - 3 pieces (piston, body, gasket), - Capacity of 10 ml, - 21Gx1½ needle - needle connector type Luer-Slip; - transparent - gradation in ml in ml, contrast gradation scale - stopper to prevent the piston from coming out of the syringe - 2-stage piston sealing gasket that prevents refilling of the solution during administration - Resistant to pressure - uniformly slipping the syringe plunger - non-toxic, non-pyrogenic, latex-free - packaged individually
Characteristics	The Barrel is transparent enough to allow easy measurement of the volume contained in the syringe and detection of air bubble. The plunger fits the inside of barrel very well, promising permitting freedom of movement.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	EO
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class II(sterile)
Standard	ISO 7864:2016, ISO 7886-1:2017 (indicated on secondary packing)
Quality Management Systems & No.	ISO13485 Q6 080946 0012 Rev. 02
Production Quality Assurance System & No.	CE(Ila) G2 080946 0007 Rev. 03
CE Conformity Assessment Route	MDD 93/42/EEC Annex V.3
Certification Body	TUV SUD Product Service GmbH
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION						
Size	Individual Packaging		Secondary Packaging		Outer Packaging	
	Qty	Material	Qty	Material	Qty	Dimensions
10ml	1	PE Bag	100	BAG	2000	64x43x51cm

安吉鑫洋进出口有限公司
ANJI SPENQ INDUSTRIAL CO., LTD

李福

LeLun, Life Guard!
ГГГГГГГГ ГГГГ ГГГГГГГГ



工厂地址：江苏省常州市横山桥镇

Factory Add: Hengshanqiao Town, Changzhou, Jiangsu, China

办公地址：中国常州关河东路66号九洲环宇商务广场1808室

Office Add: Room 1808, 66 East Guanhe Road, ECO Business Plaza Changzhou, China

Tel: +86-519-88168398 / 88602288 Fax: +86-519-88168398 / 88603638

E-mail: info@lelun.com Website: www.lelun.com



LELUN
MEDICAL APPLIANCES
乐伦医疗

CHANGZHOU MEDICAL APPLIANCES
GENERAL FACTORY CO.,LTD.

About US


Changzhou Medical Appliances General Factory Co., Ltd was built in 1988, it is a classic factory which has been specializing in producing disposable medical appliances in China for more than 30 years. The geographical location is close to Shanghai and convenient to import and export.

The area of the factory is 50000m² and purifying workshop is 15000m². The fixed assets are about USD 5,000,000. Our main products are Disposable syringe, Safety syringe, Auto-disable syringe, Disposable infusion set, Hernia mesh, Medical stapler, Disposable blood transfusion sets, Urine bag, IV cannula, Oxygen mask, Examination glove, Surgical glove, Vaginal Speculum etc. Our capacity can be up to more than 100,000,000 sets per year.

Through 30 years' scrupulous work, our factory has gone above and beyond in the industry, also has become one of the most stellar manufacturers with various export certificates in China. We won the title of "High quality product" by Jiangsu Province in 1991. In 2000, our factory passed through the ISO13485 and got the certificate of CE. In 2005, we got the license for export. In 2011, GMP was approved by Chinese SFDA.

Our products are not only circulating in Chinese market, but also cover more than 60 countries like Germany, USA, Canada, UK, Italy, Austria, Romania, France, Bulgaria, Ghana, Morocco, Turkey, Colombia, Mexico, Bangladesh, Peru, Chile, Sri Lanka, Azerbaijan, South Africa, Libya, Sudan, Yemen, Pakistan, Togo and so on and all have got good feedback.

We'll be more diligent in our business and keep consistent with good quality and service. We are looking forward to having more cooperation with you in the future.





Production Workshop and Equipment ADVANCED EQUIPMENTS IS CONSOLIDATE THE BASIS OF QUALITY



CONTENTS

page 01~04	Disposable Syringe Series
page 05~06	Oral / Enteral Syringe Series
page 07~10	Infusion Set Series
page 11~12	Needle Series
page 13~14	Urine / Stool Container Series
page 15	Mask Series
page 16	Medical Glove Series
page 16	Hernia Mesh Series
page 17	Components

► DISPOSABLE SYRINGE SERIES



Disposable Three-Part Syringe Luer Slip Without Needle

- Item Code: CMSSxx
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml



Disposable Three-Part Syringe Luer Slip With Needle

- Item Code: CMSSxxGx
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml



Disposable Three-Part Syringe Luer Lock Without Needle

- Item Code: CMSLxx
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml, 100ml, 200ml

Disposable Three-Part Syringe Luer Lock With Needle

- Item Code: CMSLxxGx
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml, 100ml, 200ml



Two-Part Syringe Without Needle

- Item Code: CMSTxx
- Size: 2ml, 5ml, 10ml, 20ml

Two-Part Syringe With Needle

- Item Code: CMSTxxGx
- Size: 2ml, 5ml, 10ml, 20ml

Insulin Syringe

- Item Code: CMSI
- Size: 0.3ml, 0.5ml, 1ml(U-100/U-40)
- Fixed needle size: 27g, 28g, 29g, 30g, 31g, 32g
- Fixed needle length: 6mm, 8mm, 13mm




Irrigation Syringe

- Item Code: CMSC
- Size: 20ml, 30ml, 50ml, 60ml, 100ml, 120ml, 150ml, 200ml, 250ml, 300ml, 350ml, 500ml
- 50ml/60ml with bulb
- 100ml with ring



Pump Syringe

- Item Code: CMSP
- Size: 50/60ml



Safety Syringe with retractable needle (Manually)

- Item Code: CMSA
- Size: 1ml, 2ml, 3ml, 5ml, 10ml, 20ml
- Simple hypodermic needle available
- Needle be retractable into the barrel and locked after use




Auto Disable Syringe

- Item Code: CMSAD
- Size: 0.05ml, 0.1ml, 0.3ml, 0.5ml, 1ml, 2/3ml, 5ml, 10ml
- Simple hypodermic needle available
- Plunger be locked in the barrel after use



Veterinary Syringe

- Item Code: CMSV
- Size: 5ml, 8ml, 10ml, 13ml, 20ml, 30ml, 60ml
- Main Usage: 1. Cow mastitis syringe
2. Goat & Sheep Wormer syringe
3. Oral paste syringe for horse
4. Feeding syringe for calves to supplement
5. Insecticide gel syringe



> ORAL / ENTERAL SYRINGE SERIES

Oral / Enteral Syringe with Tip Cap

- Item Code: CMOSC
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- Non-luer tip, can not fit with hypodermic needles
- Latex or latex free, medical grade PP
- Transparent barrel, colored plunger
- With tip cap for racap after use



Oral / Enteral Syringe with Adaptor

- Item Code: CMOSA
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- Non-luer tip, can not fit with hypodermic needles
- Latex or latex free, medical grade PP
- Transparent barrel, colored plunger
- With transparent adaptor for bottle, fitting the bottle neck ID 14mm~24mm



Oral / Enteral Syringe Female Reverse Luer Lock

- Item Code: CMOSLL
- Size: 1ml, 3ml, 5ml, 10ml, 20ml, 30ml, 60ml
- Feature: Enfit syringe, High transparency, accurate and clear scale marking



Oral Syringe (Luer Slip Tip) with Dispensing Tube

- Item Code: CMOSDT(LS)
- Size: 1ml, 2/3ml, 5ml, 10ml, 20ml
- Luer tip
- With protective cap and dispensing tube



Oral Syringe (Two Parts) with Adaptor

- Item Code: CMOSA-1
- Size: 2/3ml, 5ml, 10ml, 20ml
- Luer tip
- With adaptor fitting the bottle neck ID 18mm~22mm



► INFUSION SET SERIES



Disposable Infusion Set with Latex Connector

- Item Code: CMIF-1
- Type: Luer slip or luer lock
- Material: Medical grade PVC
- Specification: 1. Drip chamber with air vent
2. Medical grade PVC tube, 120~180mm
3. Rubber drug-adding port
4. Polybag or blister packing
5. With hypodermic needle or with scalp vein set
6. Drip chamber: 48mm or 54mm, 20 drips / 1ml for adult; 60 drips / 1ml for pediatric
7. Gravity use only




Disposable Infusion Set with Y-connector

- Item Code: CMIF-2
- Type: Luer slip or luer lock
- Material: Medical grade PVC
- Specification: 1. Drip chamber with air vent
2. Medical grade PVC tube, 120~180mm
3. Y drug-adding port
4. Polybag or blister packing
5. With hypodermic needle or with scalp vein set
6. Drip chamber: 48mm or 54mm, 20 drips / 1ml for adult; 60 drips / 1ml for pediatric
7. Gravity use only





Burette Infusion Set

- Item Code: CMIFB
- Size: 100ml, 150ml
- 60 drops = 1 ± 0.1 ml
- Material: medical grade PVC or PETG



Blood Transfusion Set

- Item Code: CMBT
- Product details: Spike with or without air vent
ABS or PP roller clamp
Latex bulb or Y connector
Luer slip or luer lock connector
Approximately 20 drips/ml.
With 200 micron filter and 18 G hypodermic needle



Infusion Set for Pump

- Item Code: CMIP
- Material: medical grade PVC
- Specification: 1. Small / big chamber with air vent
2. Medical grade PVC tube, 120~180mm
3. Polybag or blister packing
4. Drip chamber: 48mm or 54mm, 20 drips / 1ml
5. It can be used with the pump


NEEDLE SERIES



Hypodermic Needle

- Item code: CMHN
- Size: 16G / 18G / 19G / 20G / 21G / 22G / 23G / 24G / 25G / 26G / 27G / 28G / 29G / 30G

Size		Color		Length					
OD(mm)	Gauge	Color & Code		1/2" (13mm)	5/8" (15mm)	3/4" (20mm)	1" (25mm)	1 1/4" (32mm)	1 1/2" (38mm)
0.3	30G		Light yellow	√					
0.33	29G		Red	√					
0.36	28G		Blue-green	√					
0.4	27G		Grey	√					
0.45	26G		Brown	√	√				
0.5	25G		Orange		√	√	√		
0.55	24G		Purple			√	√		
0.6	23G		Blue				√	√	√
0.7	22G		Black				√	√	√
0.8	21G		Green				√	√	√
0.9	20G		Yellow				√		
1.1	19G		Beige				√		
1.2	18G		Pink				√		√
1.6	16G		White						√



Scalp Vein Set

- Item code: CMIV
- Size: 16G / 17G / 18G / 19G / 20G / 21G / 22G / 23G / 24G / 25G / 26G / 27G
- Luer slip or luer lock
- Butterfly wings can make better handling and attachment to skin



I.V. Cannula

- Item Code: CMCI-2
- Size: 14G / 16G / 18G / 20G / 22G / 24G / 26G



Spinal Needle

- Item Code: CMSN
- Quincke Bevel & Pencil Point
- Size: 18G / 19G / 20G / 21G / 22G / 23G / 24G / 25G / 26G / 27G
- Specification: 90mm



Special Needles

- Item Code: CMSNS



Blood Collection Needle with Wings

- Item Code: CMBC
- Size: 20G / 21G / 22G / 23G

► URINE / STOOL CONTAINER SERIES



Urine Cup / Container

- Item code: CMUC
- Specification: 30ml, 40ml, 60ml, 100ml, 120ml
- Sterile or Non-sterile
- With or without label



Stool Cup / Container with Spoon

- Item code: CMSC-1
- Specification: 30ml, 40ml, 60ml, 100ml, 120ml
- Sterile or Non-sterile
- With or without label



Specimen Container

- Item code: CMSC-2
- Specification: 40ml, 60ml, 90ml
- Sterile or Non-sterile
- With or without label



Vacuum Urine Collector - Urine Collection Cup

- Item code: CMUB-3
- Size: 100ml, 120ml
- Secure collection of urine for ECU with the vacuum tube



Urine Drainage Bag

- Item code: CMUB
- Size: 500ml, 1000ml, 1500ml, 2000ml
- Type: without outlet (CMUB-1)
with pull-push valve (CMUB-2)
with screw valve (CMUB-3)
with T/cross valve (CMUB-4)
Leg bag 500ml (CMUB-5)



CMUB-1



CMUB-2



CMUB-3



CMUB-4



CMUB-5



MASK SERIES



Nebulizer Mask

Item code: CMNM

Features:

- Made of medical PVC
- Without nose clip, face-fit shape
- Medium concentration
- Flexible connector, 360° rotation, more convenient for use
- Available with anti-crush 2m tubing
- 6cc or 10cc nebulizer cup



Oxygen Mask

Item code: CMOM

Features:

- Human body imitation design, makes a good fit with the chin, and a good air-tightness
- Transparent and clear medical grade PVC for patient comfort and visual assessment
- Adjustable nose clip which can be optional assures comfortable fit
- Available with oxygen mask without nose clip



Non-rebreather Oxygen Mask (with Reservoir Bag)

Item code: CMRM

Features:


- Mask with elastic strap
- Adjustable nose clip
- With 2m tubing
- Size: XS, S, M, L, L3, XL
- Bag: 1000ml or 600ml

Nasal Oxygen Cannula

Item code: CMNOC

Features:

- Size: Adult, Child, Infant, Neonate
- Free movement with the tube, can work with 15kpa. of airflow
- Available with standard or universal connector, and different lengths of oxygen tubing
- Double hole nares



Oxygen Connection Tube

Item code: CMOCT

Features:

- Free movement with the tube, can work with 15kpa. of airflow
- Match with 5mm star lumen tube, anti-crush
- Being made of all new material, enables it to be soft, smooth and elastic, without twist or knot
- Available with standard or universal connector, and different lengths



MEDICAL GLOVE SERIES



Latex Surgical Glove

- Item code: CMGS
- Material: Nature latex
- Specification: 6, 6.5, 7, 7.5, 8, 8.5
- Powdered and powder free
- Package: One pair/pouch, 50 pairs per box



Latex Examination Glove


- Item code: CMGE
- Material: Nature latex
- Specification: S, M, L, XL
- Grade: AQL1.5, AQL2.5, AQL4.0
- Powdered and powder free
- Package: 100 pcs per box



Nitrile Examination Glove

- Item code: CMNGE
- Material: Nitrile
- Specification: S, M, L, XL

HERNIA MESH SERIES



Item code: CMHM

Features:

- Using medical non-toxic polypropylene mesh
- The surgical trauma is small, only local anesthesia is needed, and the operation is convenient
- After the mesh plug is placed in the defect area of the groin, it can effectively disperse the abdominal pressure
- It can be cut appropriately according to the size of the hernia defect
- It can quickly promote the response of fibroblasts after surgery, and grow into the gap of the mesh to form a solid fibrous tissue wall
- Can be fixed in the human body for a long time

COMPONENTS

Disposable Safety Hypodermic Needle

- Size:
16G / 17G / 18G / 19G / 20G / 21G / 22G / 23G / 24G /
25G / 26G / 27G
- Features:
 1. Safety sheath: Safety mechanism ensures the safety
 2. Standard hub: Fit with standard luer connectors, ensure qualified medicine residues
 3. Four locking points: Totally locked, ensure safety



Drip Chamber

- Length: 48mm / 54mm
- Size: 20 drips / 1ml; 60 drips / 1ml



Combi Stopper

- Fit for both male and female luer
- Individual blister package



Syringe Filter

- Material: Hydrophilic PVDF
- Size: 0.22μm / 0.45μm
- Diameter: 13mm / 25mm
- Connector: Luer tip




Latex Connector

- Used with the infusion set
- Material: Rubber



Y Connector

- Material: PVC



Flow Regulator with Horn



- Material: PE



Adaptor

- Material: PE

CUSTOMERS FROM ALL OVER THE WORLD



Changzhou Medical Appliances General Factory Co., Ltd

Address:Hengshanqiao Town,Changzhou,Jiangsu,China


Product Description of Syringe for Single Use

Disposable syringe

General Description

The product consists of Barrel, Plunger, Piston (3 part) and injection needle.
The material of the barrel and plunger is PP; Needle material is SUS304, and the material of the piston is rubber or IR. Barrel is transparent, contrast marked by ml. Stopper prevents the piston from coming out of the syringe. Two-stage piston sealing gasket that prevents refilling of the solution during administration. Resistant to pressure. The plunger slips the syringe uniformly.
The product is sterilized by ethylene oxide. The product is sterile, non-toxic, pyrogen-free, latex-free, haemolytic reaction free.
It is for single use only. Connector type: Luer Slip. Individually packaged.
Needle connector is ECCENTRIC for 20ml.


The outward appearance of the product
Sterile hypodermic needle for single use



1. Zero graduation line
2. Graduation lines (ml)
3. Nominal capacity graduation line (ml)
4. Total graduated capacity line (ml)
5. Fiducial line
6. Finger grips
7. Nozzle cap
8. Nozzle lumen
9. Nozzle
10. Barrel
11. Piston
12. Seal
13. Plunger
14. Push-button

Changzhou Medical Appliances General Factory Co., Ltd

Address:Hengshanqiao Town,Changzhou,Jiangsu,China



1. Needle hub, made of PP
2. Needle, made of stainless steel tube
3. Protective cap of female conical fitting

Intended use

It is suitable for intradermal, subcutaneous, intramuscular, intravenous injection or extraction of liquids.

Specification

Syringe size: 1ml,2ml,3ml,5ml,10ml,20ml,30ml,50ml,60ml

Needle size:18G-26G

Instructions for use

1. Open the single packaging bag, tighten the needle seat of the injection needle, and remove the protective cover.
2. Slowly pull up the plunger and inhale the required liquid medicine or blood.
3. Subcutaneous, intramuscular or intravenous infusion of medicinal fluids to the human body.
4. The use of the product must comply with the requirements of the relevant operating specifications and relevant laws and regulations of the medical department, and it is limited to trained doctors or nurses.

Storage conditions

This product should be stored in a non-corrosive gas and well-ventilated room.

Product codes:

10ml, 3-part, needle 22Gx1 1/2" - REF: CMSS10mlG22x1 1/2"; EAN13 MD: 4841948001379

20ml, 3-part, needle 20Gx1 1/2" - REF: CMSS20mlG20x1 1/2"; EAN13 MD: 4841948001386

3ml, 3-part, needle 23Gx1" - REF: CMSS3mlG23x1"; EAN13 MD: 4841948001393

5ml, 3-part, needle 22Gx1 1/2" - REF: CMSS5mlG22x1 1/2"; EAN13 MD: 4841948000990

常州医疗器材总厂股份有限公司
CHANGZHOU MEDICAL APPLIANCES
GENERAL FACTORY CO., LTD.

陈 兵

PRODUCT SPECIFICATION MD-1685346316289/42

PART I. PRODUCT IDENTIFICATION	
Brand	SPENQ MEDICAL®
Code No.	SN01
EAN13	4841948000815 (MD market)
Product Name	Disposable Syringe
Specification	Disposable sterile syringes - 3 pieces (piston, body, gasket), - Capacity of 20 ml, - 21Gx1½ needle - needle connector type Luer-Slip; - transparent - gradation in ml in ml, contrast gradation scale - stopper to prevent the piston from coming out of the syringe - 2-stage piston sealing gasket that prevents refilling of the solution during administration - Resistant to pressure - uniformly slipping the syringe plunger - non-toxic, non-pyrogenic, latex-free - site of needle stick is ECCENTRIC - packaged individually
Characteristics	The Barrel is transparent enough to allow easy measurement of the volume contained in the syringe and detection of air bubble. The plunger fits the inside of barrel very well, promising permitting freedom of movement.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	EO
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class II(sterile)
Standard	ISO 7864:2016, ISO 7886-1:2017 (indicated on secondary packing)
Quality Management Systems & No.	ISO13485 Q6 080946 0012 Rev. 02
Production Quality Assurance System & No.	CE(Ila) G2 080946 0007 Rev. 03
CE Conformity Assessment Route	MDD 93/42/EEC Annex V.3
Certification Body	TUV SUD Product Service GmbH
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION							
Size	Individual Packaging		Secondary Packaging		Outer Packaging		
	Qty	Material	Qty	Material	Qty	Dimensions	Material
20ml	1	PE Bag	50	BAG	1200	64x43x51cm	carton

安吉鑫洋进出口有限公司
 ANJI SPENQ INDUSTRIAL CO., LTD

李强



ANJI SPENQ INDUSTRIAL CO., LTD.
ADD: TANGPU ECONOMIC DEVELOPMENT ZONE, ANJI, ZHEJIANG, CHINA.
TEL: +86-572-5882803 FAX: +86-572-5882115 WWW.SPENQ.COM



PRODUCT SPECIFICATION MD-1685346316289/44

PART I. PRODUCT IDENTIFICATION	
Brand	SPENQ MEDICAL®
Code No.	SN01
EAN13	4841948001409 (MD market)
Product Name	Disposable Syringe
Specification	Disposable sterile syringes - 3 pieces (piston, body, gasket), - Capacity of 2 ml, - 23Gx1¼ needle - needle connector type Luer-Slip; - transparent - gradation in ml in ml, contrast gradation scale - stopper to prevent the piston from coming out of the syringe - 2-stage piston sealing gasket that prevents refilling of the solution during administration - Resistant to pressure - uniformly slipping the syringe plunger - non-toxic, non-pyrogenic, latex-free - packaged individually
Characteristics	The Barrel is transparent enough to allow easy measurement of the volume contained in the syringe and detection of air bubble. The plunger fits the inside of barrel very well, promising permitting freedom of movement.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	EO
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class II(sterile)
Standard	ISO 7864:2016, ISO 7886-1:2017 (indicated on secondary packing)
Quality Management Systems & No.	ISO13485 Q6 080946 0012 Rev. 02
Production Quality Assurance System & No.	CE(IIa) G2 080946 0007 Rev. 03
CE Conformity Assessment Route	MDD 93/42/EEC Annex V.3
Certification Body	TUV SUD Product Service GmbH
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION							
Size	Individual Packaging		Secondary Packaging		Outer Packaging		
	Qty	Material	Qty	Material	Qty	Dimensions	Material
2ml	1	PE Bag	200	BAG	4200	64x43x51cm	carton

安吉鑫洋进出口有限公司
ANJI SPENQ INDUSTRIAL CO., LTD

李强



ANJI SPENQ INDUSTRIAL CO., LTD.
 ADD: TANGPU ECONOMIC DEVELOPMENT ZONE, ANJI, ZHEJIANG, CHINA.
 TEL: +86-572-5882803 FAX: +86-572-5882115 WWW.SPENQ.COM



PRODUCT SPECIFICATION MD-1685346316289/45

PART I. PRODUCT IDENTIFICATION	
Brand	SPENQ MEDICAL®
Code No.	SN01
EAN13	4841948000792 (MD market)
Product Name	Disposable Syringe
Specification	Disposable sterile syringes - 3 pieces (piston, body, gasket), - Capacity of 5 ml, - 22Gx1¼ needle - needle connector type Luer-Slip; - transparent - gradation in ml in ml, contrast gradation scale - stopper to prevent the piston from coming out of the syringe - 2-stage piston sealing gasket that prevents refilling of the solution during administration - Resistant to pressure - uniformly slipping the syringe plunger - non-toxic, non-pyrogenic, latex-free - packaged individually
Characteristics	The Barrel is transparent enough to allow easy measurement of the volume contained in the syringe and detection of air bubble. The plunger fits the inside of barrel very well, promising permitting freedom of movement.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	EO
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class II(sterile)
Standard	ISO 7864:2016, ISO 7886-1:2017 (indicated on secondary packing)
Quality Management Systems & No.	ISO13485 Q6 080946 0012 Rev. 02
Production Quality Assurance System & No.	CE(Ila) G2 080946 0007 Rev. 03
CE Conformity Assessment Route	MDD 93/42/EEC Annex V.3
Certification Body	TUV SUD Product Service GmbH
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION							
Size	Individual Packaging		Secondary Packaging		Outer Packaging		
	Qty	Material	Qty	Material	Qty	Dimensions	Material
5ml	1	PE Bag	150	BAG	3000	64x43x51cm	carton

安吉鑫洋进出口有限公司
 ANJI SPENQ INDUSTRIAL CO., LTD

李强



ANJI SPENQ INDUSTRIAL CO., LTD.
ADD: TANGPU ECONOMIC DEVELOPMENT ZONE, ANJI, ZHEJIANG, CHINA.
TEL: +86-572-5882803 FAX: +86-572-5882115 WWW.SPENQ.COM



PRODUCT SPECIFICATION MD-1685346316289/47

PART I. PRODUCT IDENTIFICATION	
Brand	SPENQ MEDICAL®
Code No.	SN04
EAN13	4841948000761 (MD market)
Product Name	Disposable Infusion Set
Specification	Disposable infusion set - With the needle inserted into the polymer bottle. Sterile. - flexible PVC tube for medical use, semirigid, resuming quickly after bending, on which is mounted the camera for counting drops with a filter that prevents the passage of impurities; - air chamber with valve and dropper that creates 20 drops of (1.0 + -0.1) g distilled water at temperature (20 + -2) ° C; - the clamping lever for the flow control; - a rubber sleeve between the tube and the metal needle for administration of medicines; - the type of connection between tube and needle – Luer-Slip - Tube length 150cm (+/-5cm) - The metallic needle of the infusion set is made of stainless steel, size 21G x 1 1/2, pyrogen-free, non-toxic. - Model: VC-A 21G x1 1/2"
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	EO
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class II(sterile)
Standard	ISO 8536-4:2019
Quality Management Systems & No.	ISO13485
Production Quality Assurance System & No.	Q6 080946 0012 Rev. 02
CE Conformity Assessment Route	CE(Ila)
Certification Body	G2 080946 0007 Rev. 03
U.S. FDA Registering No.	MDD 93/42/EEC Annex V.3
Listing Number	TUV SUD Product Service GmbH
	/
	/

PART III. PACKAGING INFORMATION							
Size	Individual Packaging		Secondary Packaging		Outer Packaging		
	Qty	Material	Qty	Material	Qty	Dimensions	Material
150cm/ 21G*11/2"	1	PE Bag	/	/	500	61x38x32.5cm	carton

安吉鑫洋进出口有限公司
ANJI SPENQ INDUSTRIAL CO., LTD

李强



ANJI SPENQ INDUSTRIAL CO., LTD.
ADD: TANGPU ECONOMIC DEVELOPMENT ZONE, ANJI, ZHEJIANG, CHINA.
TEL: +86-572-5882803 FAX: +86-572-5882115 WWW.SPENQ.COM



PRODUCT SPECIFICATION MD-1685346316289/48

PART I. PRODUCT IDENTIFICATION	
Brand	SPENQ MEDICAL®
Code No.	SN04
EAN13	4841948001416 (MD market)
Product Name	Disposable Infusion Set (no sleeve)
Specification	Disposable Infusion Set (no sleeve) - With the needle inserted into the polymer bottle. Sterile. - flexible PVC tube for medical use, semirigid, resuming quickly after bending, on which is mounted the camera for counting drops with a filter that prevents the passage of impurities; - air chamber with valve and dropper that creates 20 drops of (1.0 + -0.1) g distilled water at temperature (20 + -2) ° C; - the clamping lever for the flow control; - the type of connection of the Luer-Slip tube and needle, does not contain latex - Tube length 150cm(+/-5cm) - The metallic needle of the infusion set is made of stainless steel, size 21G x 1 1/2, pyrogen-free, non-toxic.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	EO
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class II(sterile)
Standard	ISO 8536-4:2019
Quality Management Systems & No.	ISO13485
Production Quality Assurance System & No.	Q6 080946 0012 Rev. 02
CE Conformity Assessment Route	CE(Ila)
Certification Body	G2 080946 0007 Rev. 03
U.S. FDA Registering No.	MDD 93/42/EEC Annex V.3
Listing Number	TUV SUD Product Service GmbH
	/
	/

PART III. PACKAGING INFORMATION						
Size	Individual Packaging		Secondary Packaging		Outer Packaging	
	Qty	Material	Qty	Material	Qty	Dimensions
150cm/ 21G*11/2"	1	PE Bag	/	/	500	61x88x32.5cm

安吉鑫洋进出口有限公司
ANJI SPENQ INDUSTRIAL CO., LTD


李强



ANJI SPENQ INDUSTRIAL CO., LTD.
 ADD: TANGPU ECONOMIC DEVELOPMENT ZONE, ANJI, ZHEJIANG, CHINA.
 TEL: +86-572-5882803 FAX: +86-572-5882115 WWW.SPENQ.COM



PRODUCT SPECIFICATION MD-1685346316289/49

PART I. PRODUCT IDENTIFICATION	
Brand	SPENQ MEDICAL®
Code No.	SN05
EAN13	4841948001232 (MD market)
Product Name	Disposable transfusion Set
Specification	Transfusion set - With the needle inserted into the polymer bottle. Sterile. - flexible PVC tube for medical use, semirigid, resuming quickly after bending, on which is mounted the camera for counting drops with a filter that prevents the passage of impurities; - air chamber with valve and dropper that creates 20 drops of (1.0 + -0.1) g distilled water at temperature (20 + -2) ° C; - the clamping lever for the flow control; - a rubber sleeve between the tube and the metal needle for additional medication; - the type of connection of the Luer-Slip tube and needle, - Tube length 160cm (-/+ 5cm) - The metal needle of the infusor is made of stainless steel, size 18G x 1 1/2, non-toxic, pyrogen-free.
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	EO
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class II(sterile)
Standard	ISO 1135-4:2015
Quality Management Systems & No.	ISO13485
	Q6 080946 0012 Rev. 02
Production Quality Assurance System & No.	CE(IIa)
	G2 080946 0007 Rev. 03
CE Conformity Assessment Route	MDD 93/42/EEC Annex V.3
Certification Body	TUV SUD Product Service GmbH
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION						
Size	Individual Packaging		Secondary Packaging		Outer Packaging	
	Qty	Material	Qty	Material	Qty	Dimensions
TYPE A	1	PE Bag	/	/	500	33x36x42cm
						carton

安吉鑫洋进出口有限公司
 ANJI SPENQ INDUSTRIAL CO., LTD

李强



JUNDAMED

HUZHOU JUNDA IMPORT AND EXPORT CO., LTD.
ADD: ROOM 604,BUILDING1,NO. 567, LINGZHI EAST ROAD, ANJI COUNTY,
HUZHOU CITY, ZHEJIANG PROVINCE, CHINA
TEL: 0086 188 5726 1771

PRODUCT SPECIFICATION MD-1685346316289/50

PART I. PRODUCT IDENTIFICATION	
Brand	JUNDAMED®
Code No.	GRW90L0010
EAN13	4841948001454 (MD market)
Product Name	Medical gauze
Specification	Medical non-sterile gauze 90 cm, density of 32 g / m2 - Cotton 100% - is rolled and conditioned in individual packaging without stitching and weaving defects - the fabric of the crotch does not allow loose threads to come off at the edges Physical and mechanical characteristics: - Width 90cm ±1.5cm - Length 10m - Density of 32 g / m2 - physicochemical characteristics: - hydrophilicity less than 10 seconds - whitening agent used: nontoxic, non-allergic, non-flammable - number of threads per cm2 = 24
Characteristics	Bleached, purified textile, plain weave. Non-detectable by X-ray
Applications	Medical use
Type of Product	Medical disposable
Sterilization Method	None
Storage Instruction	store in dark, dry and clean conditions, do not give heavy press
Picture	

PART II. QUALITY SYSTEM & STANDARD	
Classification	Class I(Non-sterile)
Standard	YY 0331:2006
Quality Management Systems & No.	None
Production Quality Assurance System & No.	Not applicable
CE Conformity Assessment Route	MDD 93/42/EEC Annex IX (Self-declaration)
Certification Body	Not applicable
U.S. FDA Registering No.	/
Listing Number	/

PART III. PACKAGING INFORMATION						
Size	Individual Packaging		Secondary Packaging		Outer Packaging	
	Qty	Material	Qty	Material	Qty	Dimensions Material
90cmx10m	1	Paper	-	-	200	52x42x46cm carton

湖州骏达进出口有限公司
HUZHOU JUNDA IMPORT AND EXPORT CO., LTD

**LISMEDFARM S.R.L.**

Șos. Muncești, 167/B, MD – 2002, mun. Chișinău, Republica Moldova
tel.: 022-80-47-98, 022-55-64-38, 022-56-94-91, e-mail: oficiu@lismedfarm.md,
web: <https://lismedfarm.md>, c/f: 1003600113573, TVA: 0304618, director – Ecaterina Chitic

Act de primire-predare a bunurilor
din data de 14.08.2023

Locul: mun. Chișinău.

Temeiul: mostrele bunurilor oferite pentru LP ocds-b3wdp1-MD-1685346316289 din 04.08.2023 privind "Achiziționarea centralizată a articolelor parafarmaceutice conform necesităților instituțiilor medico-sanitare publice (IMSP) pentru anul 2024".

La momentul dat vă predăm:

Nr. lot	Denumirea bunurilor	U/M	Cantitate
14	Manuși chirurgicale, sterile, din latex, cu pudra, N 6 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China	pereche	2
15	Manuși chirurgicale, sterile, din latex, cu pudra, N 7 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China	pereche	2
16	Manuși chirurgicale, sterile, din latex, cu pudra, N 7.5 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China	pereche	2
17	Manuși chirurgicale, sterile, din latex, cu pudra, N 8 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China	pereche	2
18	Manuși chirurgicale, sterile, din latex, cu pudra, N 8,5 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China	pereche	2
19	Manuși chirurgicale, sterile, din latex, cu pudra, N 9 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China <u>Vă rugăm să acceptați și să evaluați mostre de alte mărimi din loturi 14-19, pentru că produsele sunt identice cu excepția mărimii.</u>	pereche	0
20	Manuși chirurgicale, sterile, din latex, Fara pudra, N 6 REF: SS0602	pereche	2



LISMEDFARM S.R.L.

Șos. Muncești, 167/B, MD – 2002, mun. Chișinău, Republica Moldova
tel.: 022-80-47-98, 022-55-64-38, 022-56-94-91, e-mail: oficiu@lismedfarm.md,
web: <https://lismedfarm.md>, c/f: 1003600113573, TVA: 0304618, director – Ecaterina Chitic

	Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China		
21	Manuși chirurgicale, sterile, din latex, Fara pudra, N 7 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China	pereche	2
22	Manuși chirurgicale, sterile, din latex, Fara pudra, N 7,5 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China	pereche	2
23	Manuși chirurgicale, sterile, din latex, Fara pudra, N 8.0 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China	pereche	2
24	Manuși chirurgicale, sterile, din latex, Fara pudra, N 8.5 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China	pereche	2
25	Manuși chirurgicale, sterile, din latex, Fara pudra, N 9 REF: SS0602 Producător: Tianchang City Anrui Medical Equipments Co., Ltd., China <u>Vă rugăm să acceptați și să evaluați mostre de alte mărimi din loturi 20-24, pentru că produsele sunt identice cu excepția mărimii.</u>	pereche	0
34	Mănuși pentru examinare, netede, vinil cu pudră, Nesterile, L REF: GLVPPL Producător: Huzhou Junda Import & Export Co. Ltd., China	bucată	2
35	Mănuși pentru examinare, netede, vinil cu pudră, Nesterile, M REF: GLVPPM Producător: Huzhou Junda Import & Export Co. Ltd., China	bucată	2
36	Mănuși pentru examinare, netede, vinil cu pudră, Nesterile, S REF: GLVPPS Producător: Huzhou Junda Import & Export Co. Ltd., China	bucată	2
38	Seringa sterilă cu ac nedetasabil, pentru insulina, ac 30 G x ½ REF: SN01 Producător: Anji Spenq Industrial Co., Ltd., China	bucată	2
39	Seringa, cu ac, 10ml, 3 compon, ac 21Gx1½, sterila REF: SN01	bucată	2

**LISMEDFARM S.R.L.**

Șos. Muncești, 167/B, MD – 2002, mun. Chișinău, Republica Moldova
tel.: 022-80-47-98, 022-55-64-38, 022-56-94-91, e-mail: oficiu@lismedfarm.md,
web: <https://lismedfarm.md>, c/f: 1003600113573, TVA: 0304618, director – Ecaterina Chitic

	Producător: Anji Spenq Industrial Co., Ltd., China		
40	Seringa, cu ac, 10ml, 3 compon, ac 22Gx1 1/2", sterila REF: CMSS10mlG22x1 1/2" Producător: Changzhou Medical Appliances General Factory Co., Ltd, China	bucată	2
41	Seringa, cu ac, 20ml, 3 compon, ac 20Gx1 1/2", sterila, REF: CMSS20mlG20x1 1/2" Producător: Changzhou Medical Appliances General Factory Co., Ltd, China	bucată	2
42	Seringa, cu ac, 20ml, 3 compon, ac 21Gx1 1/2", sterila, REF: SN01 Producător: Anji Spenq Industrial Co., Ltd., China	bucată	2
43	Seringa, cu ac, 3ml, 3 compon, ac 23Gx1, sterila REF: CMSS3mlG23x1" Producător: Changzhou Medical Appliances General Factory Co., Ltd, China	bucată	2
44	Seringa, cu ac, 2ml, 3 compon, ac 23Gx1 1/4", sterila, REF: SN01 Producător: Anji Spenq Industrial Co., Ltd., China	bucată	2
45	Seringa, cu ac, 5ml, 3 compon, ac 22Gx1 1/4", sterila, REF: SN01 Producător: Anji Spenq Industrial Co., Ltd., China	bucată	2
46	Seringa, cu ac, 5ml, 3 compon, ac 22Gx1 1/2", sterila REF: CMSS5mlG22x1 1/2" Producător: Changzhou Medical Appliances General Factory Co., Ltd, China	bucată	2
47	Sisteme de perfuzie a soluțiilor cu ac metalic, L-tub-150cm REF: SN04 Producător: Anji Spenq Industrial Co., Ltd., China	bucată	2
48	Sisteme de perfuzie a soluțiilor cu ac metalic, L- tub-150cm (fără manșon) REF: SN04 Producător: Anji Spenq Industrial Co., Ltd., China	bucată	2
49	Sisteme de transfuzie a sîngelui, cu ac polimer	bucată	2

**LISMEDFARM S.R.L.**

Șos. Muncești, 167/B, MD - 2002, mun. Chișinău, Republica Moldova
tel.: 022-80-47-98, 022-55-64-38, 022-56-94-91, e-mail: oficiu@lismedfarm.md,
web: <https://lismedfarm.md>, c/f: 1003600113573, TVA: 0304618, director - Ecaterina Chitic

	REF: SN05 Producător: Anji Spenq Industrial Co., Ltd., China		
50	Tifon medical nesteril, 90 cm, densitatea min. 32 g/m2 REF: GRW90L0010 Producător: Huzhou Junda Import & Export Co. Ltd., China <u>Mostra nu are marcarea CE, pentru că producerea produsului cu marcarea CE s-a început la 04.08.2023 (conform declarației de conformitate CE atașate la notificare AMDM pentru acest lot care va fi prezentată ca răspuns la solicitare CAPCS nr. Rg02-2232 din 11.08.2023). Livrarea va fi cu marcarea CE conform machetei de ambalaj primar predat.</u>	metri	20 (2 rulouri de 10 metri)

Împreună cu mostre se predau machete de ambalaj primar, care trebuie să fie referite ca ambalaje care vor fi folosite pentru livrări în cazul atribuirii contractului sub propria răspundere Lismedfarm SRL. Ambalajele mostrelor predate sunt doar orientative. Pentru referință trebuie să fie folosite machete de ambalaje primare predate.

Am predat mostrele (din partea Lismedfarm S.R.L.):

Administratoarea Ecaterina Chitic

(Nume, Prenume Funcția Semnătura, Stampila)

Am primit mostrele (din partea Centrului pentru Achiziții Publice Centralizate în Sănătate):

14.08.2023

Seferina

(Nume, Prenume Funcția Semnătura)

↓ ÎNLĂȚURAȚI INVELIȘUL ↓
PENTRU A DESCHIDE

LOT 14-19

MĂNUȘI CHIRURGICALE DIN LATEX

Se conformă standardului EN455
Sterile dacă ambalajul nu este deschis
Nu folosiți dacă sunteți alergici la latex



EN ISO 374-1:2016



TYPE C

EN ISO 374-5:2016



VIRUS



O PERECHE

Cu Pudră
Mărime: 7.5

REF

SS0602

LOT

XXXXXXXX

STERILE R

XXXX-XX

EXP

XXXX-XX




Tianchang City Anrui Medical Equipments Co. Ltd.
Renhe Industrial District, 239331 Tangchang City,
Anhui Province, People's Republic of China

EC REP

SUNGO EUROPE B.V.
24, Olympisch Stadion, 1076DE Amsterdam,
The Netherlands

Importator, Reprezentant autorizat în RM:

Lismedfarm S.R.L.
167/B, sos. Muncesti, MD-2002 Chisinau,
mun. Chisinau, Republica Moldova
Tel: +373 22 569491, +373 22 569490



↓ ÎNLĂȚURAȚI INVELIȘUL ↓
PENTRU A DESCHIDE

LOT 19-25

MĂNUȘI CHIRURGICALE DIN LATEX

Se conformă standardului EN455
Sterile dacă ambalajul nu este deschis
Nu folosiți dacă sunteți alergici la latex



EN ISO 374-1:2016



TYPE C

EN ISO 374-5:2016



VIRUS



O PERECHE

Fără Pudră
Mărime: 7.5

REF

SS0602

LOT

XXXXXXXX

STERILE R

XXXX-XX

EXP

XXXX-XX




Tianchang City Anrui Medical Equipments Co. Ltd.
Renhe Industrial District, 239331 Tangchang City,
Anhui Province, People's Republic of China

EC REP

SUNGO EUROPE B.V.
24, Olympisch Stadion, 1076DE Amsterdam,
The Netherlands

Importator, Reprezentant autorizat în RM:

Lismedfarm S.R.L.
167/B, sos. Muncesti, MD-2002 Chisinau,
mun. Chisinau, Republica Moldova
Tel: +373 22 569491, +373 22 569490



4 841948 000891

LOT 34

Size: 10cmx7.5cm

Color: ■ 100% Pantone 301C

MĂNUȘI P/U EXAM. DIN VINIL

Cantitate: 100 bucăți

REF: GLVPPL

Specificație:

- Mărime „L“, netede, pudrate
- Rezistență la rupere
- Protecție împotriva microorganismelor
- Rezistente la penetrarea virusilor și a substanțelor chimice
- Standarde: EN 455; EN 420; EN 374-2



Importator, Reprezentant autorizat în RM:

Lismedfarm SRL

167/B, șos. Muncești, MD-2002 Chișinău

Contacte: +37322569490

Web: <https://lismedfarm.md>

Producător:

Huzhou Junda Import & Export Co. Ltd.

Huzhou City, Zhejiang Province, China



LOT: MM-YYYY

EXP: MM-YYYY

LOT 35

Size: 10cmx7.5cm

Color: ■ 100% Pantone 301C

MĂNUȘI P/U EXAM. DIN VINIL

Cantitate: 100 bucăți

REF: GLVPPM

Specificație:

- Mărime „M“, netede, pudrate
- Rezistență la rupere
- Protecție împotriva microorganismelor
- Rezistente la penetrarea virusilor și a substanțelor chimice
- Standarde: EN 455; EN 420; EN 374-2

Importator, Reprezentant autorizat în RM:

Lismedfarm SRL

167/B, șos. Muncești, MD-2002 Chișinău

Contacte: +37322569490

Web: <https://lismedfarm.md>

Producător:

Huzhou Junda Import & Export Co. Ltd.

Huzhou City, Zhejiang Province, China



LOT: MM-YYYY
EXP: MM-YYYY

LOT 36

Size: 10cmx7.5cm

Color: ■ 100% Pantone 301C

MĂNUȘI P/U EXAM. DIN VINIL

Cantitate: 100 bucăți

REF: GLVPPS

Specificație:

- Mărime „S“, netede, pudrate
- Rezistență la rupere
- Protecție împotriva microorganismelor
- Rezistente la penetrarea virusilor și a substanțelor chimice
- Standarde: EN 455; EN 420; EN 374-2

Importator, Reprezentant autorizat în RM:

Lismedfarm SRL

167/B, șos. Muncești, MD-2002 Chișinău

Contacte: +37322569490

Web: <https://lismedfarm.md>

Producător:

Huzhou Junda Import & Export Co. Ltd.

Huzhou City, Zhejiang Province, China



LOT: MM-YYYY
EXP: MM-YYYY



Seringă sterilă cu ac nedetașabil pentru insulină 1ml (U-100)
Стерильный шприц с неснимаемой иглой для инсулина 30G x 1/2"

Producător: Anji SPENQ Industrial Co., Ltd., China
 Importator, Reprez. autoriz. în RM: Lismedfarm SRL
 Adresa: 167/B, sos. Muncești, MD-2002 Chișinău
 tel.: +373-22-94-90, +373-22-56-94-91
 Nr. înregistrării: DM000154790



STERILE EO

REF: SN01
 LOT: XXXXXX
 EXP: XX-XXXX





Seringă sterilă cu ac (3 componente)
Стерильный шприц с иглой (3 компонента)

10ml
21G x 1 1/2"

Producător: Anji SPENQ Industrial Co., Ltd., China
 Importator, Reprez. autoriz. în RM: Lismedfarm SRL
 Adresa: 167/B, sos. Muncești, MD-2002 Chișinău
 tel.: +373-22-94-90, +373-22-56-94-91
 Nr. înregistrării: Dm000154772



STERILE EO

REF: SN01
 LOT: XXXXXX
 EXP: XX-XXXX



LOT 40



Producător: Changzhou Medical Appliances
General Factory, China
Importator, Reprez. autoriz în RM: Lismedfarm SRL
Adresa: 167/B, sos. Muncești, MD-2002 Chișinău
tel.: +373-22-94-90, +373-22-56-94-91

Seringă sterilă cu ac (3 componente)
Стерильный шприц с иглой (3 компонента)



STERILE EO

REF: CMSS10ml
G22x1 1/2"
LOT: XXXXXX
EXP: XX-XXXX

10ml
22G x 1 1/2"



LOT 41



Producător: Changzhou Medical Appliances
General Factory, China
Importator, Reprez. autoriz. în RM: Lismedfarm SRL
Adresa: 167/B, sos. Muncești, MD-2002 Chișinău
tel.: +373-22-94-90, +373-22-56-94-91

Seringă sterilă cu ac (3 componente)
Стерильный шприц с иглой (3 компонента)



STERILE EO

REF: CMSS20ml
G20x1 1/2"
LOT: XXXXXX
EXP: XX-XXXX

20ml
20G x 1 1/2"



LOT 42



Seringă sterilă cu ac (3 componente)
Стерильный шприц с иглой (3 компонента)

20ml
21G x 1 1/2"

Producător: Anji SPENQ Industrial Co., Ltd., China
Importator, Reprez. autoriz. în RM: Lismedfarm SRL
Adresa: 167/B, sos. Muncești, MD-2002 Chișinău
tel.: +373-22-94-90, +373-22-56-94-91
Nr. înregistrării: DM000154774



STERILE EO

REF: SN01
LOT: XXXXXX
EXP: XX-XXXX



LOT 43



Producător: Changzhou Medical Appliances
General Factory, China
Importator, Reprez. autoriz. în RM: Lismedfarm SRL
Adresa: 167/B, sos. Muncești, MD-2002 Chișinău
tel.: +373-22-94-90, +373-22-56-94-91

Seringă sterilă cu ac (3 componente)
Стерильный шприц с иглой (3 компонента)



STERILE EO

REF: CMSS3ml
G23x1"
LOT: XXXXXX
EXP: XX-XXXX

3ml
23G x 1



LOT 44



Seringă sterilă cu ac (3 componente)
Стерильный шприц с иглой (3 компонента)

2ml
23G x 1 1/4"

Producător: Anji SPENQ Industrial Co., Ltd., China
Importator, Reprez. autoriz. în RM: Lismedfarm SRL
Adresa: 167/B, sos. Muncești, MD-2002 Chișinău
tel.: +373-22-94-90, +373-22-56-94-91
Nr. înregistrării: DM000154768



STERILE EO

REF: SN01
LOT: XXXXXX
EXP: XX-XXXX



LOT 45



Producător: Anji SPENQ Industrial Co., Ltd., China
Importator, Reprez. autoriz. în RM: Lismedfarm SRL
Adresa: 167/B, sos. Muncești, MD-2002 Chișinău
tel.: +373-22-94-90, +373-22-56-94-91
Nr. înregistrării: DM000154778

Seringă sterilă cu ac (3 componente)
Стерильный шприц с иглой (3 компонента)

CE 0123

STERILE EO

REF: SN01
MFD: XXXXXX
EXP: XX-XXXX

5ml
22G x 1 1/4"



LOT 46



Producător: Changzhou Medical Appliances
General Factory, China
Importator, Reprez. autoriz. în RM: Lismedfarm SRL
Adresa: 167/B, sos. Muncești, MD-2002 Chișinău
tel.: +373-22-94-90, +373-22-56-94-91


Seringă sterilă cu ac (3 componente)
Стерильный шприц с иглой (3 компонента)



STERILE EO

REF: CMSS5ml
G22x1 1/2"
LOT: XXXXXX
EXP: XX-XXXX

5ml
22G x 1 1/2"



LOT 47 SISTEM DE PERFUZIE A SOLUȚIILOR CU AC METALIC


Model: VC-A 21G x1 ½", 20 picături,
conector Luer Slip, tub 150cm(±5cm), **cu manson**



DESCHIDE AICI



De folosit imediat după desfacerea ambalajului. Nu folosiți dacă ambalajul este deschis sau deteriorat. Sistemul se folosește pentru perfuzii clinice de medicamente sau de alimentare parenterală. Închideți clema de control. Scoateți protectorul de pe vârful de plastic și introduceți vârful în orificiul de evacuare al recipientului pentru soluție. Suspendați recipientul. Umpleți camera de picurare prin strângere și eliberare. Nu umpleți complet. Scoateți protectorul de la capătul luer atașați ferm acul. Ținând acul în poziție verticală, deschideți clema și lăsați soluția să deplaseze complet aerul în tub și ac. Închideți clema. Efectuați puncția venoasă și reglați fluxul prin deschiderea clemei. A se păstra într-un loc răcoros uscat, curat.



Producator: Anji SPENQ Industrial Co., Ltd., China
EC REP: SUNGO Europe B.V., The Netherlands

Importator, Reprezentant autorizat in RM:
Lismedfarm SRL, sos. Muncesti 167/B, 2002 Chisinau
Contacte: +373-22-56-94-91; <https://lismedfarm.md>

REF: SN04
ISO 8536-4:2019
LOT: XXXXX
EXP: XX-XXXX

STERILE EO



LOT 48 SISTEM DE PERFUZIE A SOLUȚIILOR CU AC METALIC


Model: VC-A 21G x1 ½", 20 picături,
conector Luer Slip, tub 150cm(±5cm), **fără** manson



De folosit imediat după desfacerea ambalajului. Nu folosiți dacă ambalajul este deschis sau deteriorat. Sistemul se folosește pentru perfuzii clinice de medicamente sau de alimentare parenterală. Închideți clema de control. Scoateți protectorul de pe vârful de plastic și introduceți vârful în orificiul de evacuare al recipientului pentru soluție. Suspendați recipientul. Umpleți camera de picurare prin strângere și eliberare. Nu umpleți complet. Scoateți protectorul de la capătul luer atașați ferm acul. Ținând acul în poziție verticală, deschideți clema și lăsați soluția să deplaseze complet aerul în tub și ac. Închideți clema. Efectuați puncția venoasă și reglați fluxul prin deschiderea clemei. A se păstra într-un loc răcoros uscat, curat.

Producator: Anji SPENQ Industrial Co., Ltd., China
EC REP: SUNGO Europe B.V., The Netherlands

Importator, Reprezentant autorizat in RM:
Lismedfarm SRL, sos. Muncesti 167/B, 2002 Chisinau
Contacte: +373-22-56-94-91; <https://lismedfarm.md>



DESCHIDE AICI

STERILE EO



REF: SN04
ISO 8536-4:2019
LOT: XXXXX
EXP: XX-XXXX

LOT 49 SISTEM DE TRANSFUZIE A SÂNGELUI CU AC METALIC

Model: ac 18G x 1 1/2", 20 picături,
conector Luer Slip, tub 160cm (±5cm)




DESCHIDE AICI



De folosit imediat după desfacerea ambalajului. Nu folosiți dacă ambalajul este deschis sau deteriorat. Sistemul se folosește pentru transfuzii clinice de sânge. Scoateți capacul de protecție, închideți regulatorul de debit și introduceți vârful în punga de sânge. Fixați racordul conic cu acul pentru a vă asigura că nu există scurgeri de sânge. Eliberați regulatorul de debit, sângele curge în tub și evacuează complet aerul, închideți regulatorul de debit. Operați puncția venoasă normală, deschideți regulatorul de flux și faceți transfuzii de sânge după cum este necesar. A se păstra într-un loc răcoros uscat, curat.

Producator: Anji SPENQ Industrial Co., Ltd., China
EC REP: SUNGO Europe B.V., The Netherlands

Importator, Reprezentant autorizat în RM:
Lismedfarm SRL, sos. Muncești 167/B, 2002 Chișinău
Contacte: +373-22-56-94-91; <https://lismedfarm.md>



STERILE EO



REF: SN05
ISO 1135-4:2015
LOT: XXXXX
EXP: XX-XXXX

LOT 50

Size: 10cmx7.5cm

Color: ■ 100% Pantone 301C

TIFON MEDICAL

Cantitate: 1 rolă

REF: GRW90L0010

Specificație:

- 100% bumbac
- Mărime 90cm±1.5cm x 10m
- Densitate 32g/m², 24 fire pe cm²
- Hidrofilie pînă la 10 secunde
- Standard: YY 0331-2006 equiv. cu EN 14079:2003

Importator, Reprezentant în RM:

Lismedfarm SRL

șos. Muncești 167/B, Chișinău 2002

Contacte: +37322569490

Web: <https://lismedfarm.md>

Producător:

Huzhou Junda Import & Export Co. Ltd.

Huzhou City, Zhejiang Province, China



LOT: MM-YYYY

EXP: MM-YYYY



中国认可
国际互认
检测
TESTING
CNAS L0604

scan to see the report



QDHL2206008796MD-01

Test Report

Report No.: QDHL2206008796MD-01

Applicant: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Manufacturer: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Sample Description: STERILE LATEX SURGICAL GLOVES

Sample Type: 7.5 PRE POWDERED

Test Type: SUBMITTED BY CLIENT

SGS does not verify authenticity of any Brand/Trademark of products. Buyers must check if the product is genuine with the Brand/Trademark owner directly.

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Statement

1. The report is considered invalidated if no testing seal of SGS.
2. The test report cannot be reproduced in any way, except in full content, without prior approval in writing by the laboratory.
3. The copy report is considered invalidated if no red testing seal of SGS.
4. The report is considered invalidated if no approval signature.
5. The report is considered invalidated if altered.
6. Should you have any queries or objection to the test report, please contact us in writing within 15 days after receiving the report.
7. The results shown in this test report refer only to the sample(s) tested.
8. The sample information and sample(s) was/were submitted and certified by the client, SGS quoted the information with no responsibility as to the accuracy, adequacy and/or completeness.

Company Name:

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Company Address:

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China.

Tel : +86(0)532-68999187

Fax : +86(0)532-80991952

Zip : 266101



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgs.com www.sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008796MD-01

Page 1 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Report Home Page

Sample Description	STERILE LATEX SURGICAL GLOVES		Sample No.	QDHL2206008796MD-01
	Submitted by Applicant (✓)	Random sampling (—)		
Brand	/		Type/Specifications	7.5 PRE POWDERED
Applicant	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Test Type	SUBMITTED BY CLIENT
Applicant Address	RENHE INDUSTRIAL DISTRICT, TIANCHANG CITY, ANHUI PROVINCE, CHINA		Lot No.	220310
Manufacturer	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Inspection Report No.	—
Inspected Company	—		Manufacture Date	2022.03
Sampled By	—		Sample Quantity	100PAIRS
Sampled Location	—		Lot Quantity	100000PCS
Sampled Date	—		Testing Location	QINGDAO
Sample Received Date	JUN.23,2022		Testing Period	JUN.23,2022 TO JUN.29,2022
Test items	Water tightness test for detection of holes			
Testing Accordance	EN 455-1:2020 Medical Gloves for Single Use – Part 1: Requirements and Testing for Freedom from Holes Clause 5			



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgs.com.cn sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008796MD-01

Page 2 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test conclusion	This report only provides the test results and individual assessment, details please refer to the following page(s). Issue date: JUL.05,2022
Remark	1) “——”=Not Applicable, “/”=Blank. 2) The declaration of conformity is only based on the actual value of laboratory activity, measurement uncertainty of the results not take into account. 3) Unless otherwise stated the results shown in this test report refer only to the sample(s) tested. This document cannot be used for publicity, without prior written approval of the SGS. 4) This report is to supersede test Report No.: QDHL2206008796MD, Date: JUN.29,2022. The original report shall become invalid as of the date of issuance of this report. 5) The marked▲part is the modification information.

Approver/Date:


JUL.05,2022



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgs.com.cn sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008796MD-01

Page 3 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Results

No.	Test Items	Standard /Clause	Requirement	Test Result	Assessment	Remark
1	Water tightness test for detection of holes	EN 455-1: 2020 Clause 5	Sample quantity: 200pcs AQL: 0.65 Ac: 3 Re: 4	Found: 0	Pass	/



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgs.com.cn sgs.china@sgs.com

Member of the SGS Group (SGS SA)

Test Report Picture Page

Picture and Description ▲



Product information is provided by applicant without verification of authentication of the brand.



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008796MD-01

Page 5 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Sample Description
Color: Rubber color; Size: 7.5 Pre Powdered; Material: 100% natural rubber latex; Shelf Life Time: 5 Years; Package Type: Sterilization Unit Packaging▲; Number of Contents: 1 pair▲ ; Surgical Gloves; Textured surface; Powdered Gloves
Type/Specifications/Others
/

End of Report



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgsgroup.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

scan to see the report



QDHL2206008797MD

Test Report

Report No.: QDHL2206008797MD

Applicant: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Manufacturer: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Sample Description: STERILE LATEX SURGICAL GLOVES

Sample Type: 7.5 PRE POWDERED

Test Type: SUBMITTED BY CLIENT

SGS does not verify authenticity of any Brand/Trademark of products. Buyers must check if the product is genuine with the Brand/Trademark owner directly.

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Statement

1. The report is considered invalidated if no testing seal of SGS.
2. The test report cannot be reproduced in any way, except in full content, without prior approval in writing by the laboratory.
3. The copy report is considered invalidated if no red testing seal of SGS.
4. The report is considered invalidated if no approval signature.
5. The report is considered invalidated if altered.
6. Should you have any queries or objection to the test report, please contact us in writing within 15 days after receiving the report.
7. The results shown in this test report refer only to the sample(s) tested.
8. The sample information and sample(s) was/were submitted and certified by the client, SGS quoted the information with no responsibility as to the accuracy, adequacy and/or completeness.

Company Name:

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Company Address:

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China.

Tel : +86(0)532-68999187

Fax : +86(0)532-80991952

Zip : 266101



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgsgroup.com.cn sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008797MD

Page 1 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Report Home Page

Sample Description	STERILE LATEX SURGICAL GLOVES		Sample No.	QDHL2206008797MD
	Submitted by Applicant (✓)	Random sampling (—)		
Brand	/		Type/Specifications	7.5 PRE POWDERED
Applicant	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Test Type	SUBMITTED BY CLIENT
Applicant Address	RENHE INDUSTRIAL DISTRICT, TIANCHANG CITY, ANHUI PROVINCE, CHINA		Lot No.	220310
Manufacturer	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Inspection Report No.	—
Inspected Company	—		Manufacture Date	2022.03
Sampled By	—		Sample Quantity	30PAIRS
Sampled Location	—		Lot Quantity	100000PCS
Sampled Date	—		Testing Location	QINGDAO
Sample Received Date	JUN.23,2022		Testing Period	JUN.23,2022 TO JUL.06,2022
Test items	Dimensions (Length, Width), Strength (Force at break, Force at break after challenge testing)			
Testing Accordance	EN 455-2:2015 Medical Gloves for Single Use—Part 2: Requirements and Testing for Physical Properties clause 4.2, 4.3, 5.2, 5.3			



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008797MD

Page 2 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test conclusion	This report only provides the test results and individual assessment, details please refer to the following page(s). Issue date: JUL.06,2022
Remark	1) “——”=Not Applicable, “/”=Blank. 2) The declaration of conformity is only based on the actual value of laboratory activity, measurement uncertainty of the results not take into account. 3) Unless otherwise stated the results shown in this test report refer only to the sample(s) tested. This document cannot be used for publicity, without prior written approval of the SGS.

Approver/Date:

JUL.06,2022



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgsgroup.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008797MD

Page 3 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Results

No.	Test Items	Standard/ Clause	Requirement	Test Result	Assessment	Remark	
1	Dimensions						
1.1	Length	EN 455-2: 2015 Clause 4.2	Median value: Size 7.5: ≥270mm	Sample quantity: 13pcs	See appendix 1 for details	Pass	/
1.2	Width	EN 455-2: 2015 Clause 4.3	Median value: Size 7.5: 95±5mm			Pass	/
2	Strength						
2.1	Force at break	EN 455-2: 2015 Clause 5.2	Median value: a): ≥9.0N	Sample quantity: 13pcs	See appendix 2 for details	Pass	/
2.2	Force at break after challenge testing	EN 455-2: 2015 Clause 5.3	Median value: a): ≥9.0N			Pass	/



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgsgroup.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008797MD

Page 4 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Appendix 1: Dimensions

Size	7.5	
Specimen No.	Length (mm)	Width (mm)
1	283	100
2	281	100
3	281	100
4	284	100
5	287	100
6	288	100
7	285	100
8	282	99
9	285	99
10	282	99
11	284	100
12	283	99
13	288	99
Standard requirement	≥270	95±5
Median value	284	100

Appendix 2: Strength

Size: 7.5			
Before aging		After aging	
Specimen No.	Force at break (N)	Specimen No.	Force at break (N)
1	18.1	1	17.7
2	18.7	2	17.9
3	18.2	3	18.4
4	20.7	4	19.2
5	18.9	5	17.2
6	18.1	6	18.6
7	19.1	7	18.0
8	18.3	8	17.6
9	18.9	9	19.1
10	18.1	10	15.6
11	21.2	11	18.8
12	18.5	12	18.3
13	18.5	13	17.7
Standard requirement	≥9.0	Standard requirement	≥9.0
Median value	18.5	Median value	18.0



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)

Test Report Picture Page

Picture and Description



Product information is provided by applicant without verification of authentication of the brand.



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008797MD

Page 6 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Sample Description
Color: Rubber color; Size: 7.5 Pre Powdered; Material: 100% natural rubber latex; Shelf Life Time: 5 Years; Package Type: Sterilization Unit Packaging / Multi-unit packaging; Number of Contents: 1 pair / 50 pairs; Surgical Gloves; Textured surface; Powdered Gloves
Type/Specifications/Others
/

End of Report



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgsgroup.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

scan to see the report



QDHL2206008798MD-01

Test Report

Report No.: QDHL2206008798MD-01

Applicant: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Manufacturer: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Sample Description: STERILE LATEX SURGICAL GLOVES

Sample Type: 7.5 PRE POWDERED

Test Type: SUBMITTED BY CLIENT

SGS does not verify authenticity of any Brand/Trademark of products. Buyers must check if the product is genuine with the Brand/Trademark owner directly.

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgsgroup.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Statement

1. The report is considered invalidated if no testing seal of SGS.
2. The test report cannot be reproduced in any way, except in full content, without prior approval in writing by the laboratory.
3. The copy report is considered invalidated if no red testing seal of SGS.
4. The report is considered invalidated if no approval signature.
5. The report is considered invalidated if altered.
6. Should you have any queries or objection to the test report, please contact us in writing within 15 days after receiving the report.
7. The results shown in this test report refer only to the sample(s) tested.
8. The sample information and sample(s) was/were submitted and certified by the client, SGS quoted the information with no responsibility as to the accuracy, adequacy and/or completeness.

Company Name:

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Company Address:

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China.

Tel : +86(0)532-68999187

Fax : +86(0)532-80991952

Zip : 266101



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgsgroup.com.cn sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008798MD-01

Page 1 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Report Home Page

Sample Description	STERILE LATEX SURGICAL GLOVES		Sample No.	QDHL2206008798MD-01
	Submitted by Applicant (✓)	Random sampling (——)		
Brand	/		Type/Specifications	7.5 PRE POWDERED
Applicant	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Test Type	SUBMITTED BY CLIENT
Applicant Address	RENHE INDUSTRIAL DISTRICT, TIANCHANG CITY, ANHUI PROVINCE, CHINA		Lot No.	220310
Manufacturer	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Inspection Report No.	——
Inspected Company	——		Manufacture Date	2022.03
Sampled By	——		Sample Quantity	20PAIRS
Sampled Location	——		Lot Quantity	100000PCS
Sampled Date	——		Testing Location	QINGDAO
Sample Received Date	JUN.23,2022		Testing Period	JUN.23,2022 TO JUN.29,2022
Test items	Powder			
Testing Accordance	EN 455-3:2015 Medical Gloves for Single Use—Part 3: Requirements and Testing for Biological Evaluation clause 5.2			
Test conclusion	This report only provides the test results and individual assessment, details please refer to the following page(s).			
				Issue date: JUL.05,2022



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008798MD-01

Page 2 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Remark	1) “——”=Not Applicable, “/”=Blank. 2) The declaration of conformity is only based on the actual value of laboratory activity, measurement uncertainty of the results not take into account. 3) Unless otherwise stated the results shown in this test report refer only to the sample(s) tested. This document cannot be used for publicity, without prior written approval of the SGS. 4) This report is to supersede test report No.: QDHL2206008798MD, Date: JUN.29,2022. The original report shall become invalid as of the date of issuance of this report. 5) The marked ▲part is the modification information.
--------	---

Approver/Date:

Jessica Cao
JUL.05,2022



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgsgroup.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008798MD-01

Page 3 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Results

No.	Test Items	Standard /Clause	Requirement	Test Result	Assessment	Remark
1	Powder	EN 455-3: 2015 Clause 5.2 EN ISO 21171: 2006 Method A	Sample quantity: 2pcs	0.9mg per glove	Record only	/



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)

Test Report Picture Page

Picture and Description ▲



Product information is provided by applicant without verification of authentication of the brand.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008798MD-01

Page 5 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Sample Description
Color: Rubber color; Size: 7.5 Pre Powdered; Material: 100% natural rubber latex; Shelf Life Time: 5 Years; Package Type: Sterilization Unit Packaging▲; Number of Contents: 1 pair▲; Textured surface; Surgical Gloves; Powdered Gloves
Type/Specifications/Others
/

End of Report



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

scan to see the report



QDHL2206008612MD-01

Test Report

Report No.: QDHL2206008612MD-01

Applicant: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Manufacturer: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Sample Description: STERILE LATEX SURGICAL GLOVES

Sample Type: 7.0 POWDER FREE

Test Type: SUBMITTED BY CLIENT

SGS does not verify authenticity of any Brand/Trademark of products. Buyers must check if the product is genuine with the Brand/Trademark owner directly.

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Statement

1. The report is considered invalidated if no testing seal of SGS.
2. The test report cannot be reproduced in any way, except in full content, without prior approval in writing by the laboratory.
3. The copy report is considered invalidated if no red testing seal of SGS.
4. The report is considered invalidated if no approval signature.
5. The report is considered invalidated if altered.
6. Should you have any queries or objection to the test report, please contact us in writing within 15 days after receiving the report.
7. The results shown in this test report refer only to the sample(s) tested.
8. The sample information and sample(s) was/were submitted and certified by the client, SGS quoted the information with no responsibility as to the accuracy, adequacy and/or completeness.

Company Name:

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Company Address:

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China.

Tel : +86(0)532-68999187

Fax : +86(0)532-80991952

Zip : 266101



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgsgroup.com.cn sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008612MD-01

Page 1 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Report Home Page

Sample Description	STERILE LATEX SURGICAL GLOVES		Sample No.	QDHL2206008612MD-01
	Submitted by Applicant (✓)	Random sampling (—)		
Brand	/		Type/Specifications	7.0 POWDER FREE
Applicant	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Test Type	SUBMITTED BY CLIENT
Applicant Address	RENHE INDUSTRIAL DISTRICT, TIANCHANG CITY, ANHUI PROVINCE, CHINA		Lot No.	220310
Manufacturer	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Inspection Report No.	—
Inspected Company	—		Manufacture Date	2022.03
Sampled By	—		Sample Quantity	100PAIRS
Sampled Location	—		Lot Quantity	100000PCS
Sampled Date	—		Testing Location	QINGDAO
Sample Received Date	JUN.23,2022		Testing Period	JUN.23,2022 TO JUN.29,2022
Test items	Water tightness test for detection of holes			
Testing Accordance	EN 455-1:2020 Medical Gloves for Single Use – Part 1: Requirements and Testing for Freedom from Holes Clause 5			



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008612MD-01

Page 2 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test conclusion	This report only provides the test results and individual assessment, details please refer to the following page(s). Issue date: JUL.05,2022
Remark	1) “——”=Not Applicable, “/”=Blank. 2) The declaration of conformity is only based on the actual value of laboratory activity, measurement uncertainty of the results not take into account. 3) Unless otherwise stated the results shown in this test report refer only to the sample(s) tested. This document cannot be used for publicity, without prior written approval of the SGS. 4) This report is to supersede test report No.: QDHL2206008612MD, Date: JUN.29,2022. The original report shall become invalid as of the date of issuance of this report. 5) The marked▲part is the modification information.

Approver/Date:

JUL.05,2022



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008612MD-01

Page 3 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Results

No.	Test Items	Standard /Clause	Requirement	Test Result	Assessment	Remark
1	Water tightness test for detection of holes	EN 455-1: 2020 Clause 5	Sample quantity: 200pcs AQL: 0.65 Ac: 3 Re: 4	Found: 0	Pass	/



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888


www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)

Test Report Picture Page

Picture and Description▲



Product information is provided by applicant without verification of authentication of the brand.



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008612MD-01

Page 5 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Sample Description
Color: Rubber color; Size: 7.0 Powder Free; Material: 100% natural rubber latex; Shelf Life Time: 5 Years; Package Type: Sterilization Unit Packaging▲; Number of Contents: 1 pair▲; Surgical Gloves; Textured surface; Powder-Free Gloves
Type/Specifications/Others
/

End of Report



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgsgroup.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

scan to see the report



QDHL2206008794MD

Test Report

Report No.: QDHL2206008794MD

Applicant: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Manufacturer: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Sample Description: STERILE LATEX SURGICAL GLOVES

Sample Type: 7.0 POWDER FREE

Test Type: SUBMITTED BY CLIENT

SGS does not verify authenticity of any Brand/Trademark of products. Buyers must check if the product is genuine with the Brand/Trademark owner directly.

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgsgroup.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Statement

1. The report is considered invalidated if no testing seal of SGS.
2. The test report cannot be reproduced in any way, except in full content, without prior approval in writing by the laboratory.
3. The copy report is considered invalidated if no red testing seal of SGS.
4. The report is considered invalidated if no approval signature.
5. The report is considered invalidated if altered.
6. Should you have any queries or objection to the test report, please contact us in writing within 15 days after receiving the report.
7. The results shown in this test report refer only to the sample(s) tested.
8. The sample information and sample(s) was/were submitted and certified by the client, SGS quoted the information with no responsibility as to the accuracy, adequacy and/or completeness.

Company Name:

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Company Address:

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China.

Tel : +86(0)532-68999187

Fax : +86(0)532-80991952

Zip : 266101



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgsgroup.com.cn sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008794MD

Page 1 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Report Home Page

Sample Description	STERILE LATEX SURGICAL GLOVES		Sample No.	QDHL2206008794MD
	Submitted by Applicant (✓)	Random sampling (—)		
Brand	/		Type/Specifications	7.0 POWDER FREE
Applicant	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Test Type	SUBMITTED BY CLIENT
Applicant Address	RENHE INDUSTRIAL DISTRICT, TIANCHANG CITY, ANHUI PROVINCE, CHINA		Lot No.	220310
Manufacturer	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Inspection Report No.	—
Inspected Company	—		Manufacture Date	2022.03
Sampled By	—		Sample Quantity	30PAIRS
Sampled Location	—		Lot Quantity	100000PCS
Sampled Date	—		Testing Location	QINGDAO
Sample Received Date	JUN.23,2022		Testing Period	JUN.23,2022 TO JUL.06,2022
Test items	Dimensions (Length, Width), Strength (Force at break, Force at break after challenge testing)			
Testing Accordance	EN 455-2:2015 Medical Gloves for Single Use—Part 2: Requirements and Testing for Physical Properties clause 4.2, 4.3, 5.2, 5.3			



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008794MD

Page 2 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test conclusion	This report only provides the test results and individual assessment, details please refer to the following page(s). Issue date: JUL.06,2022
Remark	1) “——”=Not Applicable, “/”=Blank. 2) The declaration of conformity is only based on the actual value of laboratory activity, measurement uncertainty of the results not take into account. 3) Unless otherwise stated the results shown in this test report refer only to the sample(s) tested. This document cannot be used for publicity, without prior written approval of the SGS.

Approver/Date:

JUL.06,2022



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgs.com.cn sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008794MD

Page 3 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Results

No.	Test Items	Standard/ Clause	Requirement		Test Result	Assessment	Remark
1	Dimensions						
1.1	Length	EN 455-2: 2015 Clause 4.2	Median value: Size 7.0: ≥270mm	Sample quantity: 13pcs	See appendix 1 for details	Pass	/
1.2	Width	EN 455-2: 2015 Clause 4.3	Median value: Size 7.0: 89±5mm			Pass	/
2	Strength						
2.1	Force at break	EN 455-2: 2015 Clause 5.2	Median value: a): ≥9.0N	Sample quantity: 13pcs	See appendix 2 for details	Pass	/
2.2	Force at break after challenge testing	EN 455-2: 2015 Clause 5.3	Median value: a): ≥9.0N			Pass	/



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008794MD

Page 4 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Appendix 1: Dimensions

Size	7.0	
Specimen No.	Length (mm)	Width (mm)
1	291	92
2	289	93
3	286	93
4	284	92
5	282	93
6	290	93
7	291	92
8	285	92
9	290	93
10	285	92
11	284	93
12	285	93
13	291	93
Standard requirement	≥270	89±5
Median value	286	93

Appendix 2: Strength

Size: 7.0			
Before aging		After aging	
Specimen No.	Force at break (N)	Specimen No.	Force at break (N)
1	19.3	1	21.6
2	23.0	2	19.4
3	20.5	3	18.6
4	21.9	4	18.2
5	20.4	5	19.1
6	21.3	6	18.6
7	21.6	7	18.8
8	18.2	8	20.2
9	21.4	9	19.6
10	21.7	10	20.1
11	21.0	11	20.3
12	21.6	12	19.4
13	22.2	13	21.3
Standard requirement	≥9.0	Standard requirement	≥9.0
Median value	21.4	Median value	19.4



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)

Test Report Picture Page

Picture and Description



Product information is provided by applicant without verification of authentication of the brand.



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008794MD

Page 6 of 6

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Sample Description
Color: Rubber color; Size: 7.0 Powder Free; Material: 100% natural rubber latex; Shelf Life Time: 5 Years; Package Type: Sterilization Unit Packaging / Multi-unit packaging; Number of Contents: 1 pair / 50 pairs; Surgical Gloves; Textured surface; Powder-Free Gloves
Type/Specifications/Others
/

End of Report



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

scan to see the report



QDHL2206008795MD-01

Test Report

Report No.: QDHL2206008795MD-01

Applicant: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Manufacturer: TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD

Sample Description: STERILE LATEX SURGICAL GLOVES

Sample Type: 7.0 POWDER FREE

Test Type: SUBMITTED BY CLIENT

SGS does not verify authenticity of any Brand/Trademark of products. Buyers must check if the product is genuine with the Brand/Trademark owner directly.

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Statement

1. The report is considered invalidated if no testing seal of SGS.
2. The test report cannot be reproduced in any way, except in full content, without prior approval in writing by the laboratory.
3. The copy report is considered invalidated if no red testing seal of SGS.
4. The report is considered invalidated if no approval signature.
5. The report is considered invalidated if altered.
6. Should you have any queries or objection to the test report, please contact us in writing within 15 days after receiving the report.
7. The results shown in this test report refer only to the sample(s) tested.
8. The sample information and sample(s) was/were submitted and certified by the client, SGS quoted the information with no responsibility as to the accuracy, adequacy and/or completeness.

Company Name:

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Company Address:

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China.

Tel : +86(0)532-68999187

Fax : +86(0)532-80991952

Zip : 266101



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101 t (86-532) 68999888 www.sgsgroup.com.cn sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008795MD-01

Page 1 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Report Home Page

Sample Description	STERILE LATEX SURGICAL GLOVES		Sample No.	QDHL2206008795MD-01
	Submitted by Applicant (✓)	Random sampling (—)		
Brand	/		Type/Specifications	7.0 POWDER FREE
Applicant	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Test Type	SUBMITTED BY CLIENT
Applicant Address	RENHE INDUSTRIAL DISTRICT, TIANCHANG CITY, ANHUI PROVINCE, CHINA		Lot No.	220310
Manufacturer	TIANCHANG CITY ANRUI MEDICAL EQUIPMENTS CO.,LTD		Inspection Report No.	—
Inspected Company	—		Manufacture Date	2022.03
Sampled By	—		Sample Quantity	20PAIRS
Sampled Location	—		Lot Quantity	100000PCS
Sampled Date	—		Testing Location	QINGDAO
Sample Received Date	JUN.23,2022		Testing Period	JUN.23,2022 TO JUN.29,2022
Test items	Powder (Powder-free gloves)			
Testing Accordance	EN 455-3:2015 Medical Gloves for Single Use—Part 3: Requirements and Testing for Biological Evaluation clause 4.4			
Test conclusion	This report only provides the test results and individual assessment, details please refer to the following page(s).			

Issue date: JUL.05,2022



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008795MD-01

Page 2 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Remark	1) “—”=Not Applicable, “/”=Blank. 2) The declaration of conformity is only based on the actual value of laboratory activity, measurement uncertainty of the results not take into account. 3) Unless otherwise stated the results shown in this test report refer only to the sample(s) tested. This document cannot be used for publicity, without prior written approval of the SGS. 4) This report is to supersede test report No.: QDHL2206008795MD, DATE: JUN.29,2022. The original report shall become invalid as of the date of issuance of this report. 5) The marked▲part is the modification information.
--------	---

Approver/Date:

Jessica Bao

JUL.05,2022



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008795MD-01

Page 3 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Test Results

No.	Test Items	Standard/ Clause	Requirement	Test Result	Assess ment	Remark
1	Powder (Powder-free gloves)	EN 455- 3:2015 Clause 5.2 EN ISO 21171:2006 Method C	Sample quantity: 6 pcs Average: ≤2mg per glove	0.2mg per glove	Pass	/



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.
Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgs.com

sgs.china@sgs.com

Member of the SGS Group (SGS SA)

Test Report Picture Page

Picture and Description ▲



Product information is provided by applicant without verification of authentication of the brand.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com



中国认可
国际互认
检测
TESTING
CNAS L0604

Report No.: QDHL2206008795MD-01

Page 5 of 5

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Sample Description
Color: Rubber color; Size: 7.0 Powder Free; Material: 100% natural rubber latex; Shelf Life Time: 5 Years; Package Type: Sterilization Unit Packaging [▲] ; Number of Contents: 1 pair [▲] ; Textured surface; Surgical Gloves; Powder-Free Gloves
Type/Specifications/Others
/

End of Report



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 83071443, or email: CN.Doccheck@sgs.com

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

SGS Center, No.143, Zhuzhou Road, Laoshan District, Qingdao, Shandong, China 266101

t (86-532) 68999888

www.sgsgroup.com.cn

sgs.china@sgs.com

Member of the SGS Group (SGS SA)



HUZHOU JUNDA IMPORT AND EXPORT CO., LTD.

ADD: ROOM 604, BUILDING 1, NO. 567, LINGZHI EAST ROAD, ANJI COUNTY, HUZHOU CITY, ZHEJIANG PROVINCE, CHINA

Tel: 0086 188 5726 1771

湖州骏达进出口有限公司
HUZHOU JUNDA IMPORT & EXPORT CO., LTD
纱布卷质量检验单
CERTIFICATE OF ANALYSIS OF ABSORBENT GAUZE ROLL

规格 (SPECIFICATION)		卷数 (QUANTITY)	箱数 (CARTONS)
Absorbent Gauze Roll 90CM X 10M REF: GRW90L0010		100Roll	1
序号 NO.	检验项目 Testing Items	标准要求 Standard	检验结果 Results
1	纤维性状 Fibre Identification	A.纤维长度不少于 10mm 长。 A.Each Fibre is not less than 10mm.	OK
		B.纤维放在碘氯化锌的溶液中会变紫。 B.Treat with iodinated zinc chloride solution. The fibres become violet.	OK
		C.0.1G 纤维放入 10ml 氯化锌的溶液中加热至 40℃, 并放置 2.5 个小时不会褪色。 C.To 0.1g add 10ml of zinc chloride solution, heat to 40℃ and stand 2.5 hours. The fibres don't dissolve.	OK
2	酸碱度 Acidity & Alkalinity	100ml 的供试液中加入酚酞指示剂不显粉红色, 加溴甲酚紫指示剂不显黄色。 PH value = 7.0	7.0
3	外来物质 Foreign Fibres	显微镜观察, 除了可能存在的极少数孤立的外来物质外, 都是棉质纤维。 Exam from microscope, it is seen to consist exclusively of typical cotton fibres, except that occasionally a few isolated foreign fibres may be present.	OK
4	荧光物 Fluorescence	365nm 紫外灯下观察只允许显微棕紫色荧光和少数黄色颗粒, 除少数分离纤维外, 无强蓝色荧光。 No intense blue fluorescence.	OK
5	密度 Thread Count	纬线: 每 100MM 含 69-76 线, 经线: 每 10CM 含 40-45 线。 Warp: 69-76 threads per 100mm, Weft: 40-45 threads for 10cm.	74threads/ 42 threads
6	平方克重 Mass per square meter	纱布每平方米重量不少于 32 克。 Not less than 32g/m2 for gauze.	32.5G/M2



HUZHOU JUNDA IMPORT AND EXPORT CO., LTD.

ADD:ROOM 604,BUILDING1,NO. 567, LINGZHI EAST ROAD, ANJI COUNTY ,HUZHOU CITY ,ZHEJIANG PROVINCE, CHINA

Tel: 0086 188 5726 1771

7	最小断裂强度 Minimum Breaking Load	纬线: 每 50MM 不少于 70 牛顿, 经线: 每 50MM 不少于 40 牛顿。 Warp: Not less than 70N for per 50mm, Weft: Not less than 40N for per 100mm	90N/56N
8	基本尺寸 Basic Dimension	长度偏为±50CM, 宽度偏差±2CM。 More or less 50cm for length, more or less 2cm for width.	OK
9	吸水时间 Sinking Time	5g 试样沉入液面以下所需时间在 10s 内。 Less than 10s for 5g cotton sinking into water completely.	5s
10	醚中可溶物 Ether-soluble substances	100ml 的供试液中遗留残渣小于 0.5% 。 Not more than 0.5% for 100ml testing water.	0.45%
11	表面活性物 Surface-active Substances	供试液中表面活性物质泡沫不超过 2mm。 Less than 2mm for the height of froth.	1.8mm
12	水中可溶物 Water-soluble substances	100ml 的供试液中遗留残渣小于 0.5% 。 No more than 0.5% for 100ml testing water.	0.43%
13	淀粉和糊精 Starch and dextrin	供试液中添加淀粉和糊精不会显示任何蓝色、 紫色、浅红色或褐色。 The solution shall not show any blue,violet,reddish or brownish color.	OK
14	可提取变色物质 Extractable colouring matter	取 10G 放置过滤器, 用 50ml 酒精慢慢提取得 到的萃取物质与参考颜色 GY6 接近。 In a narrow percolator, slowly extract 10G with alcohol R until 50ml of extract is obtained, and it's not more intensely colored than reference solution GY6.	OK
15	干燥失重 Loss on drying	2g 试样干燥至恒重, 质量损失不大于 8.0% Not more than 8% of lose for 2g cotton.	7.1%
16	易氧化物 Easily oxidized substances	40ml 的供试液中加高锰酸钾红色在 5min 内 不完全消失。 Add the KMnO4 into 40ml testing water,the red color will be not disappeared within 5min.	OK
17	吸水量 Water-holding capacity	每克试样的吸水量不少于 23g。 Not less than 23g water per gram cotton.	24.5



HUZHOU JUNDA IMPORT AND EXPORT CO., LTD.

ADD:ROOM 604,BUILDING1,NO. 567, LINGZHI EAST ROAD, ANJI COUNTY ,HUZHOU CITY ,ZHEJIANG
PROVINCE, CHINA
Tel: 0086 188 5726 1771

18	炽灼残渣 Sulfated Ash	2g 试样遗留残渣不大于 0.4% 。 Not more than 0.4% for the 2g cotton.	0.35%
19	白度 Whiteness	不低于 80 度。 ≥80 PATE	85
质检科意见 Comments From Inspection Section		经质检科确认，产品质量符合出口欧标。 The above results complied with exporting China standard YY 0331-2006 (equivalent to EN 14079:2003).	
备注 (Remarks)		Batch No.: 20230818D Mfg. Date: 2023-08 Exp. Date: 2028-08	

检验 Inspector: Chenlan 核校 Corrector: Wang Jingjing
批准 Ratifier: Wang Shen 日期 Date: 2023-08

湖州骏达进出口有限公司
HUZHOU JUNDA IMPORT AND EXPORT CO., LTD.