

1.1 Applications

■ Intended Use

The ceiling type x-ray tube support CH-200/CHU-200 is used to hold the X-ray tube and collimator to position the X-ray tube unit and determine the radiography area after locating the radiography position of patient when combined with an X-ray high voltage generator, an X-ray tube unit, an X-ray collimator, and if necessary, an X-ray radiography stand, an X-ray radiography table, and digital radiography system.

■ Indications for Use

This unit is used for radiographic examinations of whole body except mammography when combined with an X-ray high voltage generator, an X-ray tube unit, an X-ray collimator, and if necessary, an X-ray radiography stand, an X-ray radiography table, and digital radiography system.

This unit must only be operated by qualified personnel, such as radiography technicians or those with equivalent qualifications.

This system is used for total patient population including pediatric examination.



CAUTION

- **Do NOT use the system for any other purpose than indicated above.**
- **Do NOT connect the system to any other equipment, either by electrical or mechanical means, or modify the system.**



NOTE

Tube Support CH-200 and CHU-200 are hereafter referred to as "Tube Support" or "equipment".