


Europclip™ Titanium Haemostatic Clips


Key Features

- Chevron shape design for **„tip-first” closing action** allowing for secure and precise positioning, even on delicate structures while reducing the risk of tissue slippage
- **Unique inner cross-serration** to maximize contact surface ensuring secure but less traumatic vessel sealing
- **Heart-shaped profile** for reliable vessel occlusion also reducing the risk of scissoring




... more Features

- Medical grade titanium for maximized biocompatibility
- Soft-load cartridge design for convenient loading action
- 100% made in Germany
- Color coded cartridges for convenient clip size and applier identification
- Sterile packed for immediate use
- Multiple packaging configurations
- Adhesive at the back of the cartridge for easy fixation providing precise handling



Pistol Handle Appliers for Endo Surgery



	Diameter			Color coding	Item no.
	Ø 5 mm	Ø 10 mm	Ø 12.5 mm		
Endo Appliers for small Europclips™, pistol handle, rotatable, 330 mm	•			Yellow	13-1613-105AS
Endo Appliers for medium Europclips™, pistol handle, rotatable, 330 mm		•		Blue	13-1610AS
Endo Appliers for medium-large Europclips™, pistol handle, rotatable, 330 mm		•		Green	13-1615AS
Endo Appliers for large Europclips™, pistol handle, rotatable, 330 mm			•	Orange	13-1613-106AS

Townsend Handle Appliers for Endo Surgery



	Diameter			Color coding	Item no.
	Ø 5 mm	Ø 10 mm	Ø 12.5 mm		
Endo Appliers for small Europclips™, townsend handle, rotatable, 330 mm	•			Yellow	13-1613-105A
Endo Appliers for medium Europclips™, townsend handle, rotatable, 330 mm		•		Blue	13-1610A
Endo Appliers for medium-large Europclips™, townsend h., rotatable, 330 mm		•		Green	13-1615A
Endo Appliers for large Europclips™, townsend handle, rotatable, 330 mm			•	Orange	13-1613-106A

	Color coding of item	Item no.
Open surgery appliers for medium Europclips™, straight jaws, 19 cm / 8 "	●	AI-321526DD
Open surgery appliers for medium Europclips™, straight jaws, 28 cm / 11 "	●	AI-321522DD
Open surgery appliers for medium Europclips™, right-angled jaws, 19 cm / 8 "	●	AI-321520DD
Open surgery appliers for medium Europclips™, right-angled jaws, 28 cm / 11 "	●	AI-321528DD
Endo Appliers for medium Europclips™, pistol handle, rotatable, Ø 10 mm, 330 mm	●	13-1610AS
Endo Appliers for medium Europclips™, townsend handle, rotatable, Ø 10 mm, 330 mm	●	13-1610A
Endo Appliers for medium Europclips™, townsend handle, rot., demountable, Ø 10 mm, 330 mm	●	13-1610AZ
Medium-large Europclips™, 120 clips / box	●	AI74CML6
Medium-large Europclips™, 180 clips / box	●	AI74CML6-30C
Medium-large Europclips™, 576 clips / box	●	AI74CML6/24
Open surgery appliers for medium-large Europclips™, straight jaws, 19 cm / 8 "	●	AI-321540DD
Open surgery appliers for medium-large Europclips™, straight jaws, 28 cm / 11 "	●	AI-321542DD
Open surgery appliers for medium-large Europclips™, right-angled jaws, 19 cm / 8 "	●	AI-321544DD
Open surgery appliers for medium-large Europclips™, right-angled jaws, 28 cm / 11 "	●	AI-321546DD
Endo Appliers for medium-large Europclips™, pistol handle, rotatable, Ø 10 mm, 330 mm	●	13-1615AS
Endo Appliers for medium-large Europclips™, townsend handle, rotatable, Ø 10 mm, 330 mm	●	13-1615A
Endo Appliers for medium-large Europclips™, townsend handle, rot., demount., Ø 10 mm, 330 mm	●	13-1615AZ
Large Europclips™, 120 clips / box	●	AI74CL6
Large Europclips™, 180 clips / box	●	AI74CL6-30C
Large Europclips™, 576 clips / box	●	AI74CL6/24
Open surgery appliers for large Europclips™, straight jaws, 19 cm / 8 "	●	AI-321550DD
Open surgery appliers for large Europclips™, straight jaws, 28 cm / 11 "	●	AI-321552DD
Open surgery appliers for large Europclips™, right-angled jaws, 19 cm / 8 "	●	AI-321554DD
Open surgery appliers for large Europclips™, right-angled jaws, 28 cm / 11 "	●	AI-321556DD
Endo Appliers for large Europclips™, pistol handle, rotatable, Ø 12.5 mm, 330 mm	●	13-1613-106AS
Endo Appliers for large Europclips™, townsend handle, rotatable, Ø 12.5 mm, 330 mm	●	13-1613-106A
Endo Appliers for large Europclips™, townsend handle, rot., demountable, Ø 12.5 mm, 330 mm	●	13-1613-106AZ