

PAGINA 0

**CREATA PENTRU A
PUTEA SEMNA
ELECTRONIC
DOCUMENTUL**

Re-entry Catheter

The Enteer™ Re-entry System, consisting of the catheter and guidewire, gives you intuitive control to reliably target the true lumen from the subintimal channel above or below the knee.

The system requires no capital equipment. It's designed to be nothing less than a precise extension of your own expert hand.

Flat-shaped self-orienting balloon

The Enteer™ Catheter's unique balloon design self-orientes towards the true lumen within the subintimal space when inflated.

180° and offset exit ports

Offset exit ports are located on each side of the device, allowing the Enteer™ guidewire to re-enter the correct port into the true lumen.

OTW 0.014" and 0.018" guidewire compatible

Allows for flexibility during your case and minimizes guidewire exchanges.

Order information

Product catalogue number	Balloon size (W x H x L mm)	Working length (cm)	Guidewire compatibility (inch)	Crossing profile (max inch)	Sheath compatibility
ENB-375-20-135	3.75 x 1.5 x 20	135	≤ 0.018	0.066	5 F
ENB-275-20-150	2.75 x 1.0 x 20	150	≤ 0.018	0.066	5 F

Product catalogue number	Product	Description
ENW-FX-014-300	Enteer™ Guidewire	0.014" x 300 xm flexible
ENW-SD-014-300	Enteer™ Guidewire	0.014" x 300 cm standard
ENW-SF-014-300	Enteer™ Guidewire	0.014" x 300 cm stiff

