








# Original KLEINSASSER Instruments

for Endolaryngeal Microsurgery

KLEINSASSER Instruments, working length 23 cm, for use with Handle 8597



- |   |        |  |
|---|--------|--|
|   | 8595 A |  |
|    | 8595 A | <b>Knife</b> , oval, straight                                |
|    | 8595 B | <b>Same</b> , angled 45°                                     |
|    | 8595 C | <b>Same</b> , sickle-shaped, curved, pointed                 |
|    | 8595 D | <b>Same</b> , straight cut, pointed                          |
|    | 8595 E | <b>Same</b> , golf club-shaped, curved, round                |
|    | 8595 F | <b>Same</b> , round, cuts vertically                         |
|  | 8596 A | <b>Hook</b> , blunt, with probe end                          |
|  | 8596 B | <b>Hook</b> , sharp  |
|  | 8596 C | <b>Hook</b> , blunt, angled 90°, with probe end, length 5 mm |
|  | 8596 E | <b>Needle</b> , curved to right, working length 23 cm        |
|  | 8596 F | <b>Same</b> , curved to left                                 |
|  | 8596 H | <b>Suction Elevator</b> , working length 23 cm               |
|  | 8596 T | <b>Knot Tier</b>   |

**Lot 41**

8597 KLEINSASSER **Handle**, for use with 8595 A – 8596 T, 8655 A – K, 8693 A/B



- |  |         |   |
|--|---------|---|
|  | 8596 W  |   |
|  | 8596 W  | <b>Cotton Applicator</b> , straight, working length 25 cm                                 |
|  | 8596 WJ | <b>Cotton Applicator</b> , for endolaryngeal microsurgery, straight, working length 21 cm |
|  | 8596 WK | <b>Same</b> , working length 18 cm  |

1-942