

TECHNICAL DATA SHEET

AUTOCLAVABLE BIOHAZARD BAG

PRODUCT CODE	:	OT
RELATED STANDARDS	:	Autoclavable Biohazard Bag meets the requirements of the Council Directive 94/62/EC on Packaging.
INTENDED USE	:	Autoclavable Biohazard Bags are used for the sterilization of medical wastes before disposal.
OPERATION CONDITIONS	:	Autoclavable Biohazard Bags are designed for use in steam sterilization methods. The recommended operation temperature is max. 141 °C. Other sterilization conditions should be determined by the end user regarding to material to be sterilized. During sterilization, bag should be left open for an effective sterilization.

PRODUCT SPECIFICATIONS

PROPERTIES	METHOD	UNIT	VALUE	
THICKNESS (one side)	ASTM D 374	micron	50 ± 5	
TENSILE STRENGTH	ASTM D 882	N/15 mm	MD	22,4
			CD	14,4
ELONGATION AT BREAK	ASTM D 882	%	MD	640
			CD	510

PACKAGING : Autoclavable Biohazard Bags are packaged as the listed below.

<u>Product Code</u>	<u>Dimensions</u>	<u>Unit</u>
OT5370	53 cm x 70 cm	300 pcs/box
OT6389	63 cm x 89 cm	200 pcs/box
OT87110	87 cm x 110 cm	120 pcs/box



PMS Steri-test

PMS MEDİKAL Amb. San. Tic Ltd. Şti
pms@pmsmedikal.com