

E.R.C.P. Catheters

Single Use

Features:

- "Pre-Curved" → Easier Positioning to the Papilla
- Metal Tips/ Metal Rings → Maximum X-Ray Visualization
- Guide Wire Compatibility → Contrast Injection while Guide Wire in Place
- Color Markings → Macroscopically Controlled Cannulation Depth
- Removable Stylet (1.2m) → Kink Free Cannulation
- Single Use → Safety and Convenience

Benefits:

Order Information:

"Pre-Curved"

Catalog No.	Description	Outer Diameter	Guide Wire
GER-01-18-035	Standard Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-02-18-035	5 mm Tapered Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-03-18-021	10 mm Tapered Tip "Pre-Curved"	1.8 mm	.021" (0.89 mm)
GER-04-18-035	Metal Cone Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-05-18-035	Round Metal Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-06-18-035	"Bullet Type" Metal Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)
GER-07-18-035	Metal Ball Tip "Pre-Curved"	1.8 mm	.035" (0.89 mm)

Manufactured by:

 Medi-Globe

Medi-Globe GmbH
 Medi-Globe-Strasse 1-5
 83101 Achenmühle, Germany
 Phone: ++49 8032 973 379
 Fax: ++49 8032 973 339
 e-mail: sales@medi-globe.de
 www.medi-globe.de

Distributed by:



Ver. 2.0 - 10/07