


Product-Information-Sheet

Polypectomy, EMR and ESD

Cold and Hot
Snares



- For polyp ligation prior to polypectomy
- Prevents postpolypectomy bleeding
- Preloaded loop, ready to use with handle
- Nylon-loop ligation keeps tissue damage to a minimum
- Opening diameter of 30 mm enables ligation of large polyps



Article number	Article Name	Quantity value	Min. working channel ø value	Working length value	Loop ø	Area Compatibility
N5382130	HX-400U-30	5	 2.8	2300 mm	30 mm	