



iSyringe

Conventional

Conventional syringes are available with and without an attached needle. Syringe-needle combinations feature a clear barrel with a bold scale.


Auto disable

Auto disable syringes are distinguished by the fact that they cannot be reused: following injection, the syringe is automatically blocked, thus becoming unusable.




3 ml	Ordinary Syringe	23G x 1 1/4
5 ml	Ordinary Syringe	21G x 1 1/2
10 ml	Ordinary Syringe	18G x 1 1/2
20 ml	Ordinary Syringe	18G x 1 1/2
50 ml	Ordinary Syringe	18G x 1 1/2

Other needle sizes are available




Luer slip tip

Our luer slip tip provides a frictionfit connection that requires you to push and twist the syringe tip into the needle hub.



Luer Lock tip

Our Luer-Lok tip is generally used for injections requiring a secure connection of the syringe to another device.



Eccentric luer slip tip

The eccentric luer slip tip is generally used for venipunctures and aspiration of fluids, allowing close proximity to the skin.

