

Anexa tehnica, lot 6

Parametri solicitati	Parametri oferiti
LOT 6 Set de diagnostic care contine un amestec de enzime, tampon de reactie si un amplificator de detectie pentru reactiile One Step RT-PCR cantitativ	
set compatibil cu tehnologia TaqMan contine 3 componente: 1) 2X RT-PCR Buffer 2) 25X RT-PCR Enzyme Mix 3) Nuclease-free Water utilizeaza enzima AmpliTaq Gold ADN Polimeraza este inclus in componenta sa ROX refference dye set compatibil cu instrumente Real Time: Applied Biosystems 7500 Real-Time PCR System, Applied Biosystems 7900HT Fast Real-Time PCR System, Stratagene Mx3000P System, Cepheid SmartCycler System La livrare produsul sa fie insotit de certificat de calitate si conformitate.	set compatibil cu tehnologia TaqMan contine 3 componente: 1) 2X RT-PCR Buffer 2) 25X RT-PCR Enzyme Mix 3) Nuclease-free Water utilizeaza enzima AmpliTaq Gold ADN Polimeraza este inclus in componenta sa ROX refference dye set compatibil cu instrumente Real Time: Applied Biosystems 7500 Real-Time PCR System, Applied Biosystems 7900HT Fast Real-Time PCR System, Stratagene Mx3000P System, Cepheid SmartCycler System La livrare produsul va fi insotit de certificat de calitate(analiza) si conformitate. AM1005, AgPath-ID™ One-Step RT-PCR Reagents, 100 teste
LOT 6.1 Amplificator de detectie (Detection Enhancer, for real-time PCR)	
1 tub x 190 µL stocare -5 / -30°C	1 tub x 190 µL, stocare -5 / -30°C A44810, Detection Enhancer, for real-time PCR, 100 teste
LOT 6.2 Primerii pentru H5, H7, N1	
Pentru H5: Sense primer H5LH1 5'-ACA TAT GAC TAC CCA CAR TAT TCA G-3'; Antisense primer H5RH1 5'-AGA CCA GCT AYC ATG ATT GC-3'; Probe H5PRO 5'-FAM- TCW ACA GTG GCG AGT TCC CTA GCA-TAMRA-3'; Pentru H7: Sense primer LH6H7 5'-GGC CAG TAT TAG AAA CAA CAC CTA TGA-3'; Antisense primer RH4H7 5'-GCC CCG AAG CTA AAC CAA AGT AT-3'; Probe H7pro11 5'-FAM-CCG CTG CTT AGT TTG ACT GGG TCA ATC T-BHQ1-3'; Pentru N1: N1 Forw 5'-GTT TGA GTC TGT TGC TTG GTC -3' N1 Rev 5'- TGA TAG TGT CTG TTA TTA TGC C -3' Sonda N1 5' VIC – TTG TAT TTC AAT ACA GCC AC- MGB -3'	Pentru H5: Sense primer H5LH1 5'-ACA TAT GAC TAC CCA CAR TAT TCA G-3'; Antisense primer H5RH1 5'-AGA CCA GCT AYC ATG ATT GC-3'; Probe H5PRO 5'-FAM- TCW ACA GTG GCG AGT TCC CTA GCA-TAMRA-3'; Pentru H7: Sense primer LH6H7 5'-GGC CAG TAT TAG AAA CAA CAC CTA TGA-3'; Antisense primer RH4H7 5'-GCC CCG AAG CTA AAC CAA AGT AT-3'; Probe H7pro11 5'-FAM-CCG CTG CTT AGT TTG ACT GGG TCA ATC T-BHQ1-3'; Pentru N1: N1 Forw 5'-GTT TGA GTC TGT TGC TTG GTC -3' N1 Rev 5'- TGA TAG TGT CTG TTA TTA TGC C -3' Sonda N1 5' VIC – TTG TAT TTC AAT ACA GCC AC- MGB -3' 4331348, Custom TaqMan™ Gene Expression Assay, FAM 4331348, Custom TaqMan™ Gene Expression Assay, FAM 4448508, Custom TaqMan™ Gene Expression Assay, VIC