


Amplatz Renal Dilation Set

Hydrophilic coated

Used for progressive dilation of the nephrostomy tract prior to percutaneous stone removal. Radiopaque dilators.



Art. No.	Size (Fr/CH)
ARDS 8/30	8 to 30

Amplatz Renal Dilators

Hydrophilic coated

Used for progressive dilation of the nephrostomy tract prior to percutaneous stone removal. Radiopaque dilators.

Single

Art. No.	Size (Fr/CH)	Length (cm)
RD 18	18	30
RD 20	20	30
RD 22	22	30
RD 24	24	30
RD 26	26	30
RD 28	28	30
RD 30	30	30

Amplatz Sheath

Hydrophilic coated

Used during renal dilatation to provide atraumatic working tracks after removal of the dilator.



Art. No.	Size (Fr/CH)	Length (cm)
AM 10	10	30
AM 24	24	17
AM 26	26	17
AM 28	28	17
AM 30	30	17