

For use with 24/26 Fr. Resectoscope Sheaths

The cutting loops are delivered with a wire diameter of 0.35 mm.

Loops with 30 as the last digit of the order number indicate a wire diameter of 0.30 mm.



27040 GP1

Distal Tip	24/26 Fr.	Instrument Description	for use in the	
			Bladder	Prostate

Cutting Loops, color code: yellow

	27040 GP1	Cutting Loop , bipolar	●	●
	27040 GD1	Cutting Loop , bipolar, small	●	●
	27040 BL1	Coagulation Electrode , bipolar, pointed	●	●

Special Bladder Cutting Loops, color code: yellow/orange

	27040 GP130	Cutting Loop , bipolar, diameter 0.30 mm	●	–
	27040GP140	Cutting Loop, bipolar, diameter 0.40 mm, 24/26 Fr., color code: yellow		
	27040 JB1	Cutting Loop , bipolar, longitudinal	●	–
	27040 JB130	Cutting Loop , bipolar, longitudinal, diameter 0.30 mm	●	–

Vaporization Electrode, color code: yellow

	27040 NB	Vaporization Electrode HALF MOON® , bipolar, ball-shaped	●	●
--	----------	---	---	---



280

280

Protection Tube, for sterilization and storage of electrodes, loops, curettes and knives

High Frequency Cords see page 73

High Frequency Surgical Units see chapter 15, UNITS