




XN System(Primary/Revision)

Femoral Condyle (Unit:mm)

Femoral Condyle	#2	#3	#4	#5	#6	#7	
ML	60.5	63	65	67.5	70.5	75	
AP	54	56	59	62	65	68	

Tibial Insert (Unit:mm)

Tibial Insert	#2/#3/#4	#5/#6	#7	
ML	68	75	84	
AP	45	51	56	
Thickness	9、11、13、15、17、19			

Tibial Tray (Unit:mm)

Tibial Tray	#2	#3	#4	#5	#6	#7	
ML	65	68	70	72	75	79	
AP	43	45	47	49	51	53	

Straight Stem (Unit:mm)

Diameter	Φ9	Φ10	Φ11	Φ12	Φ13	Φ14	Φ15	Φ16	
Length	60、100、130								

Stem with Offset (Unit:mm)

Diameter	Φ9	Φ10	Φ11	Φ12	Φ13	Φ14	Φ15	Φ16	
Length	100、130								

Adjustable Wedge 1.(Unit:mm)

Adjustable Wedge	#2	#3	#4	#5	#6	#7	
Distal Wedge	5	5	5	5	5	5	
	10	10	10	10	10	10	
Posterior Condyle Wedge	5	5	5	5	5	5	
	10	10	10	10	10	10	

Adjustable Wedge 2.(Unit:mm)

Adjustable Wedge 2	#2	#3	#4	#5	#6	#7	
Tibial Half Wedge	5 10	5 10	5 10	5 10	5 10	5 10	
Tibial Half Wedge with Angle	16°	16°	16°	16°	16°	16°	
Tibial Total Wedge with Angel	7°	7°	7°	7°	7°	7°	

The Indications of XN Revision System:**Case 1, RA severe knee abduction with deformity.**

Preoperative:

Postoperative:

**Case 2, The prosthesis is loosening after 9 years of postoperation of TKA.**

Preoperative:

Postoperative:

**Case 3, Osteolysis after postoperation of TKA,the prosthesis is loosening.**

Preoperative:

Postoperative:

