


# Product Catalog





#### TIANDA MEDICAL(HUAIAN)CO., LTD.

No.106 East Songjiang Road, Huaiyin Economic&Technological Development Zone 223002 Huaian City, Jiangsu China

Tel:+86 517 84379688
Fax:+86 517 84379788







#### Company Introduction

TianDa Medical(HuaiAn) Co Ltd.established in 2006, which is specialized in diabetic care and surgical medical products, we dedicate to developing, manufacturing and sales of medical instruments, we regard quality management as the core key, meanwhile we have passed CE, ISO13485 and FDA successfully.

Over ten years of unremitting efforts, our products have exported all over the world and widely recognized by customers, we will continue to forge ahead with our philosophical philosophy'work together, improve together'to keep improving the stability of quality and our own competitiveness, we will try our best to be the most trustworthy medical supplier in the future.

Enterprise tenet: Changes for Better, Creating a better future

Business Philosophy: Work together, Improve together

**Employing philosophy:** We put talented and virtuous person in the important position; we put virtuous but talentless person in the normal position; we never use the person just have talented.

#### Certificate










#### **Blood Lancets**

- Material:LDPE,Stainless Needle
- Size:21g,23g,26g,28g,30g
- Package:100 pcs/box 20000 pcs/carton
- Feature: Triple bevel needle tip, automatic production
- Sterilization:Gamma Irradiation



REF NO	0121	0123	0126	0128	0130
SIZE	21G	23G	26G	28G	30G
COLOR					

#### Steel Lancets

Material:Stainless Steel

• Size: S/30mm,L/40mm

Package:200 pcs/box
 20000 pcs/carton

Feature:Tungsten steel mold,high
 quality needle



REF NO	0230	0240	
SIZE	30mm	40mm	

#### Alcohol Pads

- Material:Isopropyl,nonwovens
- Size:30x65mm,30x30mm,65X56mm
   45X90mm
- Package:100 pcs/box 10000 pcs/carton
- Feature:70% isopropyl, eternal use



REF NO	0301	0302	0303	0304
SIZE	30x65	30x30	65X60	45X90 (mm)

### Safety Lancets PA&PAII

- Material: ABS, LDPE, Spring, Stainless Needle
- Size:21Gx2.2mm,21Gx2.4mm,23Gx1.8mm,26Gx1.8mm,28Gx1.8mm, 30Gx1.8mm
- Package:100 pcs/box 2000 pcs/carton
- Feature:Safety, press to use, triple bevel needle tip, automatic assembly
- Sterilization:Gamma Irradiation



REF NO	042122 052122	042124 052124	042318 052318	042618 052618	042818 052818	043018 053018	
SIZE	21Gx2.2	21Gx2.4	23Gx1.8	26Gx1.8	28Gx1.8	30Gx1.8	(mm)
COLOR							

# Lancing Device

Material:, Fine copper, Spring

• Size: 138mm

Package:1 pcs/box 500 pcs/carton

 Feature:Continuous button,rapid rebound, 1 million times service life



REF NO	0601	
SIZE	138mm	

# Surgical Blades



- Material: Carbon, stainless
- Size:10#.11#,12#,15#,20#,21#,22#,23#,24#
- Package:100 pcs/box 5000 pcs/carton
- Feature: automatic grinding keep sharpness, Tungsten steel mold
- ensure stable dimension
- Sterilization:Gamma Irradiation

REF NO	0710	0711	0712	0715	0720	0721	0722	0723	0724 carbon
	0810	0811	0812	0815	0820	0821	0822	0823	0824 stainless
SIZE	10#	11#	12#	15#	20#	21#	22#	23#	24#

# Disposable Scalpels



- Material: Carbon, stainless
- Size:10#,11#,12#,15#,20#,21#,22#,23#,24#
- Package:10 pcs/box 500 pcs/carton
- Feature: Easy for use, blades are made by our own factory
- Sterilization:Gamma Irradiation

REF NO	0910	0911	0912	0915	0920	0921	0922	0923	0924 carbon
	1010	1011	1012	1015	1020	1021	1022	1023	1024 stainless
SIZE	10#	11#	12#	15#	20#	21#	22#	23#	24#

### Mini Scalpels I

Material:ABS,carbon,ctainless

• Size: 10#,11#,12#,15#

Package:50 pcs/box 2000 pcs/carton
 or bulk packing

Feature: Popular in Kit, electron industry



REF NO	$\frac{1410}{1510}$				carbon stainless	
SIZE	10#	11#	12#	15#		

### Mini Scalpels II

Material:ABS, carbon, stainless

• Size: 10#,11#,12#,15#

Package:50 pcs/box 2000 pcs/carton
 or bulk packing

 Feature: Petroleum jelly as surface coating, castrate animals.



REF NO	$\frac{1610}{1710}$	$\frac{1611}{1711}$			carbon stainless	
SIZE	10#	11#	12#	15#		

#### Stitch Cutter Blades


Material:Carbon,stainless

• Size: 110mm,65mm

Package:100 pcs/box 5000 pcs/carton

Feature:cut the wide range of

suture on the wound




### Skin Graft Blades

Material: Carbon, stainless

• Size: 158mm

Package:100 pcs/box,1000 pcs/carton

 Feature:use for skin removal or transplant during surgery



# Umbilical Cord Clamps

- Material:PE
- Size:50mm,55mm
- Package:100 pcs/box,5000 pcs/carton
- Feature:No toxic,no allergic, one injection mold
- Sterilization: Ethylene oxygen



REF NO	1101	1102	
SIZE	50mm	55mm	
COLOR			


### Identification Bracelet

- Material:PVC
- Size: Adult, child, Mother-Child
- Package:100 pcs/box 5000 pcs/carton
- Feature:Soft and comfortable wear,no toxic,no allergic,water resistance



REF NO	1201 1203	1202 1204	
SIZE	Adult Card insert Adult Hand writing	Child Card insert Child Hand writing	
COLOR			

# Tourniquet




- Material:TBE
- Size: 450x25x0.6 mm
- Package:25 pcs/box 100 pcs/carton
- Feature:Soft and latex free, degradable, single use only

REF NO	1301
SIZE	450x25x0.6
COLOR	

# Urinary Drainage Bags

- Material:PVC,ABS,PE
- Size:2000ml,2600ml with Air-Vent
- Package:250 pcs/carton
- Feature:Double hook,bed sheet clamp,
   Needleless sampling port with cap,Rigid dog House
- Sterilization: Ethylene oxygen





REF NO	2001	2002	2003	2004
Volumes	2000ml	2000ml	2000ml	2600ml
Size	TVALUE	Pull Value	Cross Value	Air Vent Type

TIANDA MEDICAL(HUAIAN) CO, LTD

0

#### Lancing Device

• Material: ABS, Spring

- Size: 100 mm
- Package: 1pcs/box,200 pcs/carton
- Feature:

9 adjustable depth penetration 100,000 times life cycle easy for use and carry Rapid launch and rebound for use Rapid launch and rebound for use

SIZE (mm)	COLOR	
100mm	BLUE/WHITE/ROSE	
	A side	



#### How To Use








Place a lancet firmly into thec arrier.







Reload the cover on the lancing device.



Adjust the appropriate penetration depth Slide the end barrel until hear a click then release.



Place the cap on the chosen site then press the trigger button to fire.



Remove the cover again and insert the needle into the protective cap.



Slide the ejector to dispose of the lancet and put on the cap for next use.