

MEDICAL QUALITY MANAGEMENT SYSTEM CERTIFICATE

Universal GmbH

This certificate is granted to the organization,

Zeyni Medikal Tekstil Ins. San. Ve Tic. Ltd. Sti.

Yunus Emre Mah. 12. Yildirim Sok. No 17 Yildirim-Bursa/Turkey

by review of SA2-5076 numbered report for the scope

Manufacture And Sales Of Disposable Surgical Sterile Gowns, Drapes And Custom Packs, Non-Sterile Medical Textile Products

to certify that a management system in accordance with standard's clauses is established and being implemented

DIN EN ISO 13485:2016

Certificate No: MDMS 0520 006861

Original Certification Date: 22.04.2020

Issue / Revised Date: 13.04.2022

Expiry Date: 21.04.2023

Certification Period: 3 years (3rd year)





Universal GmbH



QUALITY MANAGEMENT SYSTEM CERTIFICATE

Universal GmbH

This certificate is granted to the organization,

Zeyni Medikal Tekstil Ins. San. Ve Tic. Ltd. Sti.

Yunus Emre Mah. 12. Yildirim Sok. No 17 Yildirim-Bursa/Turkey

by review of SA2-5076 numbered report for the scope

Manufacture And Sales Of Disposable Surgical Sterile Gowns, Drapes And Custom Packs, Non-Sterile Medical Textile Products

to certify that a management system in accordance with standard's clauses is established and being implemented

DIN EN ISO 9001:2015

Certificate No: QMS 0520 006861

Original Certification Date: 22.04.2020

Issue / Revised Date: 13.04.2022

Expiry Date: 21.04.2023

Certification Period: 3 years (3rd year)











Notified Body 1023
INSTITUTE FOR TESTING AND CERTIFICATION, Inc.,
třída Tomáše Bati 299, Louky, 763 02 Zlín, Czech Republic

EC Certificate - Production Quality Assurance No. 09 0816 QS/NB

The quality system of manufacturer

Zeyni Medikal Tekstil Ins. San. Ve Tic. Ltd. Sti.
Yunus Emre Mah. 12. Yıldırım Sok. No 17 Yıldırım –BURSA, Turkey

has been certified as meeting the requirements of

Directive 93/42/EEC

on medical devices, Annex V

for the following product category(ies):

Sterile Disposable Surgical Gowns and Drapes and Drape Sets

The Notified Body No. 1023 declares that the aforementioned manufacturer has implemented a quality assurance system for manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex V. The Notified Body has audited this system with limitation to those aspects of manufacture concerned with securing and maintaining sterile conditions. This part of quality assurance system conforms to the requirements of this Directive and is subjected to periodical surveillance.

Valid from: 2020-06-23 Valid until: 2022-11-24 First Issued: 2009-11-26

Revision: f

Date: 2020-06-23

TUTTEROAD

Mgr. Jiří Heš

Representative of the Notified Body No. 1023

Zeyni Medikal Tekstill Ins. San. Ve Tic. Ltd. Sti Yunus Emre Mah. 12. Yildirim Sok. N° 17 Yildirim -BURSA, Turkey

Product(s):

Name:

Drape sets

Trade name(s):

Model(s):

See table below

Class:

GMDN:

47783

Drape set types	Code number
General surgery drape set	101-001-S
Laparoscopy drape set	201-001-S
Minor surgery drape set	102-001-S
Dressing surgery drape set	103-001-S
Spinal surgery drape set	301-001-S
Craniotomy surgery drape set	401-001-S
Shunt surgery drape set	401-002-S
Valve Remplacement surgery drape set	501-001-S
Cardiovascular Thoracic drape set	202-001-S
Anjyo drape set	601-001-S
By-pass drape set	602-001-S
Ophthalmic drape set 1	701-001-S
Ophthalmic drape set 2	701-002-S
Caesarean drape set	801-001-S
O.P.U surgery drape set	802-001-S
Delivery drape set	802-002-S
Dental drape set	803-001-S
Biopsie drape set	804-001-S
Cystoscopy drape set	901-001-S

GENERAL DIRECTOR:

Baki Conard

07.06.22

STAMP & SIGNATURE

ZEYNI MEDI TEKS.İNŞ. SAN VE TİZ Yunusemre Mh. 12, Yıldırım Sk. No. Yıldırım V.D.: 998 141 1400 Tizlaret Mersis No: 0998141 40002

Zeyni Medikal Tekstill Ins. San. Ve Tic. Ltd. Sti Yunus Emre Mah. 12. Yildirim Sok. N° 17 Yildirim -BURSA, Turkey

Drape set types	Code number
Percutaneous drape set	901-002-S
T.U.R drape set	901-003-S
Ear-nose-throat drape set	902-001-S
Vertebra surgery drape set	903-001-S
Arthroscopy drape set	904-001-S
Extremity drape set	904-002-S
Hip surgery drape set	904-003-S
Rinoplasty set	905-001-S
Vertebra set	906-001-S
I.V.F Drape set	907-001-S
E.T set	907-002-S
Abdominal drape set	908-001-S
Hand surgery set	904-004-S
Chemotherapy drape set	909-001-S

Product(s):

Name:

GOWNS

Trade name(s):

Model(s):

See table below

Class:

Is

GMDN:

35778

Cown Types	Code number	
Gown Types	100-001-S	
Standard Surgery Gowns	200-001-S	
Drotocted Surgery Gowns		

GENERAL DIRECTOR:

Batr Conavor

09.06,22

STAMP & SIGNATURE

EKS. INS. /SAN VE/FIC.L nusemre Mh. Q. Wildryn S. Vo: 1/1 Yıldırım V. (J. 1990 V. 1990 Lizaret No.

ZEYNI

MEDICAL

For medical devices, regulated in accordance with EC Directive 93/42/EEC

EC Declaration of Conformity

11042020-UB

Manufacturer : ZEYNİ MEDİKAL TEKSTIL INSAAT SAN. VE TIC. LTD. IST.

Address : Yunus Emre Mah. 12. Yildirim Sok. No: 17 Yildirim-Bursa / TURKEY

Trade name : Zeyni Medical

Product Name: STERILE AND DISPOSABLE SURGICAL SUPPORTS, COVERS AND COVER SETS

STERILE AND DISPOSABLE SURGICAL LINERS, DRAWER AND LID SETS

(The list of products is given in appendix-1.)

Classification: Class 1-Sterile (with sterile function)

Conformity assessment according to Annex VII and the annexes of the applied directives (Medical Device Directive 93/42/EEC; Class 1-as sterile device)

Applied harmonized standards: ISO 11135:2014, ISO 14644-3 2019, 13795 1 2019 ISO 14971 2012, 150 15223-1:2016, ISO 10993-7:2008, ISO 14937 2009

Notified body: ITC (Institute for testing and Certification, a s] Trida T Bati 299,76421 Zlin-Loky/Czech Republic

Notified body number: CE 1023

EC certificate: 090816QS/NB

Validity of the CE certificate: 24.11.2022

GENERAL DIRECTOR:

Baki Conaucr

07.06,22

STAMP & SIGNATURE

S.INS. SAN XX Samre Mh. 12. Yıldının Sa. No. 171 Vidinim V.D.: 998 144 400 / Vy 200 Mersis No. 1998 147 400 / V **ZEYNI**

MEDICAL

Zeyni Medikal Tekstil İnş San. and Tic. Ltd. Şti hereby declares that the above mentioned

products are in compliance with all applicable provisions of Directive 93/42/EEC. The product is

safe under foreseeable and reasonably foreseeable conditions of storage and use.

The company has put in place measures to ensure that all products of the type mentioned above

are safe and meet the essential requirements of Directive 93/42/EEC

The company has put in place a systematic procedure to review experience gained with the

devices in the post-production phase and to apply appropriate means for any necessary

corrective action and keeps this systematic procedure up to date. The company undertakes to

report to the competent authority any malfunction or alteration in the properties, performance

or inadequacy of the product's instructions for use that could lead to the death of the patient or

serious damage to the patient's health, as well as the or medical reasons that led to the

systematic recall of the product by the manufacturer.

If the device is modified without the consent of the undersigned, this statement becomes invalid

with respect to the modified product.

Regulated in accordance with EC Directive 93/42/EEC concerning medical devices

EC DECLARATION OF CONFORMITY

11042020-UB

GENERAL DIRECTOR:

Bati Canavar

07,06,22

STAMP & SIGNATURE



"HYGIENE FOR HEALTH,
DISPOSABLE STERILE PRODUCTS FOR HYGIENE"







Our company was established in Bursa. Our product range includes sterile and non-sterile products especially surgical gowns, surgical drape and other products to be applicable in clinical application. Our customers include hospitals, clinics, pharmacies, retailers and reputable industrial companies.

Our products are sterilized as controlled in our owner site with EO unit. As a manufacturer and exporter of medical devices Zeyni Medical operates in accordance with the following regulations and standards:

- The European Medical Device Directive 93/42/EEC (2007/47/EC)
- The Quality Management Standards ISO 9001 and
- The Medical Management Standards ISO 13485

Zeyni Medical Ltd holds a Quality System Certificates and CE certificates. Besides Turkey, Zeyni Medical products are also available in Europe, South America, Middle East and Africa area.

We are committed to be at the forefront of technological innovation in the healthcare industry. As our valued customers, you can rest assured that you will be access to the very best & latest technology in the field of medical disposables.

Hakkımızda

Bursa'da kurulmuş olan firmamız, steril ve non-steril ürünler özellikle cerrahi önlükler, cerrahi örtüler ve klinik uygulamalarında kullanılan diğer ürünleri üretmektedir. Müşterilerimiz hastaneler, klinikler, eczaneler, perakendeciler ve endüstri firmalarıdır.

EO sterilizasyonun kendi tesislerinde kurduğu sistemle kontrollü olarak gerçekleşmektedir. ZEYNİ MEDİKAL, Medikal ürünler üreticisi ve ihracatçısı olarak aşağıdaki yasa ve standartlar doğrultusunda faaliyet göstermektedir.

- Avrupa tarafından yayınlanmış Medikal Cihaz Direktifi 93/42/EEC
- Kalite Yönetim Sistem Standardı ISO 9001 ve
- Medikal Kalite Yönetim Sistem Standardı ISO 13485.

Zeyni Medikal, yukarıda adı geçen kalite yönetim sistem sertifikaları ve CE sertifikasına sahip bir firmadır.







BASIC PRODUCT INFORMATION / TEMEL ÜRÜN BİLGİSİ

RESISTANCE

Two and three layer materials, used for gowns and drapes production, are laminated with plastic film that works as a perfect barrier against body fluids and microorganism penetration.

LINTING

The surface is resistant to abrasion, material does not release any fibres.

ABSORBENCY

The drape surface is made from absorbent viscose with high capacity of absorbency 250 ml water/sqm.

ADHESIVE ZONES

Sides of drapes and apertures are lined with adhesive zones, which allow simple effective and safe fixation of the drape to patient body.

ACCESSORIES

Following elements may be integrated in the drape upon request:

- Incision drape
- Complements for tubes and cables fixation
- Equipment pouch
- Fluid collection bags, etc.

FOLDING AND PACKING

Special folding technique enables that drapes and integrated accesories are applied without any contact with sterile and non-sterilee surfaces.

FUNCTIONAL PROPERTIES

- Each part is anatomically shaped
- Materials perfectly adhere to the surface
- All components are ranked into logic sets which include all that is necessary for a surgery preparation
- Always ready to use









USED MATERIALS / KULLANILAN MATERYALLER



1- TWO LAYER LAMINATE 54 gsm

İKİ KATLI NONWOVEN 54 gr/m²

Two layer laminate, upper layer is made of absorbent viscose 28 g/sqm, bottom layer is bacterial barrier polyethylene foil 25µm. Application:

- Drapes & Reinforcement Zones
- EquipmentCovers

İki kat lamisyonlu, üst katı emici medikal nonwoven 28gr/m², alt katı 25 mikron bakteri bariyerli polietilen.

2- TWO LAYER SPUNBOND LAMINATE 54 gsm İKİ KATLI NONWOVEN 54 gr/m²

Non-woven material composed of two layers, upper layer is absorbent spunbond 30 g/sqm, bottom layer polyethylene foil 25µm.

Application:

- Drapes & Reinforcement Zones
- EquipmentCovers

İki kattan oluşmuş nonwoven materyal, üst katı emici medikal spunbond 30 gr/m², alt katı 25 mikron bakteri bariyerli polietilendir.

3- THREE LAYER LAMINATE 70 gsm ÜÇ KATLI NONWOVEN 70 gr/m²

Three layer laminate, upper layer is made of absorbent viscose 28 g/sqm, middle layer is made of polyethylene foil 25µm and bottom layer is made of absorbent non-woven textile. Application:

Drapes & Reinforcement Zones


Üç kat laminasyonlu, üst katı emici medikal nonwoven 28 gr/m. orta katı 25 mikron bakteri bariyerli polietilen ve alt katı emici nonwoven materyal

4- SMMS 50 gsm

SMMS fabrics is composed of two layers of spunbonded fabrics and two layer of Meltblown fabrics. With good tensile strength, high filtration and remarkable softness. It prevent and isolate dust particle, alcohol, blood, bacteria and virus. Application:

- Gowns&Drapes
- Sterilization Wraps

İki kat medikal nonwoven ve iki kat meltblown'dan oluşur. İyi çekme mukavemeti sayesinde yüksek filtrasyon ve dikkate değer bir yumuşaklık sağlar. Toz partikül, kan, alkol, bakteri ve virüslere karşı koruma ve izolasyon sağlar.













5- SPUNBOND 50 gsm

Spunbond material which doesn't show resistance to liquid is used in the production of non-sterile gowns and textiles.

Spunbond, sıvı direnci olmayan, non-steril tekstil ürünleri ve önlük imalatı yapılan bir materyaldir.



6- POLYETHYLENE FOIL 60µm (PE 60 micron)

A turquoise very drapeable, micro embossed PE coex film, Impermeable to fluids, non toxic, very high strength. Application:

- Back Table Cover&Mayo Stand Cover
- EquipmentCovers

Turkuaz renginde, geçirimsiz, toksik olmayan, yüksek dirençli bir üründür.



It has high absorbency and lint free. It is using for drying the hands of surgeons before operation.

Yüksek emici ve hava bırakmama özelliklerine sahip olup, ameliyat öncesi doktorların el kurulaması işleminde kullanılır.



8- STERILIZATION WRAP 60 gsm STERİLİZASYON BOHÇASI 60 gr/m²

It is used for covering and keeping of drape packs during the sterilization process, and it supplies Hygienic security and tight packing.

Örtü setlerini sterizasyon prosesi süresince korumayı sağlar, ayrıca hijyen ve sıkı bir paketleme sağlar.





INDEX / İÇİNDEKİLER

STERILE PRODUCTS	STERİL ÜRÜNLER	PAGE SAYFA
STERILE SURGICAL GOWNS	STERİL CERRAHİ ÖNLÜK	6-7
STERILE SURGICAL DRAPES	STERİL CERRAHİ ÖRTÜLER	8-9
GENERAL SURGERY PACKS	GENEL CERRAHİ SETLERİ	10-11
ORTHOPEDIC PACKS	ORTOPEDI SETLERI	12-14
CARDIOVASCULAR PACKS	KARDIOVASKÜLER SETLER	15-16
OBSTETRICS&GYNEGOLOGY PACKS	KADIN DOĞUM&JİNEKOLOJİK SETLER	17-18
UROLOGY PACKS	ÜROLOJİ SETLERİ	19-20
CARDIOLOGY PACKS	KARDİYOLOJİ SETLERİ	21
E.N.T PACKS	K.B.B SETLERİ	22
NEUROSURGERY & BIOPSY PACKS	NÖROŞİRURJİ & BİYOPSİ SETLERİ	23-24
OPHTHMOLOGY & DENTISTRY PACKS	GÖZ-IMPLANT CERRAHİ SETLER	25
CHEMOTHERAPY&PLASTIC SURGERY PACKS	KEMOTERAPI-PLASTIK CERRAHI SETLERI	26
PAEDIATRIC PACKS	PEDIATRIK SETLER	27
NON-STERILE PRODUCTS	STERİL OLMAYAN ÜRÜNLER	
NON-STERILE GOWNS	STERİL OLMAYAN ÖNLÜKLER	28
DISPOSABLE SLIPPERS	TEK KULLANIMLIK TERLİK	29
COLONOSCOPY PANTS/SHORTS	KOLONOSKOPİ PANTOLONU/ŞORTU	29
STERILIZATION WRAPS	STERİLİZASYON BOHÇALARI	30
BED SHEETS	YATAK ÖRTÜLERİ	30



Surgical Gown Cerrahi Önlük















"Hygiene for health, Disposable Sterile Products for hygiene"

SURGICAL GOWN CLASSIC



- Made from non-woven textile hydrophobic material SMMS that protects against fluid penetration
- Gowns are fixed with four ties system
- Sleeves are finished with knitted elastic cuffs
- Surgical gowns are a tie-round type which enables to keep the sterility of surgeons back side
- There are three types of surgical gown classic; standard, reinforced, full reinforced
- Reinforced and full reinforced types have front side and sleeves with impermeable material that assures perfect protection of surgeon during major surgeries
- Fulfil requirements of EN 13795 standard







6
















"Hygiene for health, Disposable Sterile Products for hygiene"

SURGICAL GOWN COMFORT

- Made from non-woven textile triplex repellent material against fluid penetration
- Gowns are fixed with four ties system
- Sleeves are finished with knitted elastic cuffs
- Surgical gowns are a tie-round type, which enables to keep the sterility of surgeons back side
- Fully reinforced
- Fulfil requirements of EN 13795 standard










Plain Drapes

Material	Size (cm)
sms	50x50
	50x70
	75x75
	75x90
	150x180
	200x240
bilaminate	50x50
	50x70
	75x75
	75x90
	150x180
	200x240
trilaminate	50x50
	50x70
	75x75
	75x90
	150x180
	200x240

Surgical Drapes with Adhesive Tape

Material	Size (cm)
sms	50x50
	50x70
	75x75
	75x90
	150x180
	200x240
bilaminate	50x50
	50x70
	75x75
	75x90
	150x180
	200x240
trilaminate	50x50
	50x70
	75x75
	75x90
	150x180
	200x240







STERILE SURGICAL DRAPES / STERIL CERRAHI ÖRTÜLER

Fenestrated Drapes

Material	Size (cm)	
bilaminate	50x70	7
	50x70	10
	75x90	7
	75x90	10
	75x110	10
	100×150 50×50	6x8 6x8
	50x70	6x8
	160×200	12
trilaminate	50x70	7
	50x70	10
	75x90	7
	75x90	10
	75x110	10
	100x150	6x8
	50x50	6x8
	50x70	6x8

Fenestrated Drapes with Adhesive Tape

Material	Size (cm)	
bilaminate	50x70	7
	50x70	10
	75x90	7
	75x90	10
	75x110	10
	90x110 120x150	7 10
	45x75	7
	240x280	10
trilaminate	50x70	7
	50x70	10
	75x90	7
	75x90	10
	75x110	10
	90x110	7
	120x150	10
	45x75	7













General Surgery Sets Genel Cerrahi Setleri














"Hygiene for health, Disposable Sterile Products for hygiene"

GENERAL SURGERY PACK / GENEL CERRAHİ ÖRTÜ SETİ


	Product Name Ürün Adı		Pcs. Adet
1.	Anesthesia Drape / Anestezi Örtüsü	150x240	1 Pc.
2.	Foot Drape with Adhesive Tape / Yapışkan Bantlı Ayak Örtüsü	150x180	1 Pc.
3.	Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
4.	Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
5.	Drape with Adhesive Tape / Yapışkan Bantlı Örtü	75x100	2 Pc.
6.	Reinforced Surgical Gown / Korumalı Cerrahi Önlük		2 Pc.
7.	Medical Towel / Medikal Havlu	40x40	2 Pc.
8.	Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



LAPAROSCOPY PACK / LAPAROSKOPİ SETİ

U

	Product Name Ürün Adı	Size Ebat	Pcs. Adet
1.	Laparascopy Drape / Laparoskopi Örtüsü	200x300	1 Pc.
2.	Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
3.	Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
4.	Medical Towel / Medikal Havlu	40x40	2 Pc.
5.	Reinforced Surgical Gown / Korumalı Önlük		2 Pc.
6.	Sponge / Spanç	7.5x7.5	10 Pc.
7.	Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.





General Surgery Sets Genel Cerrahi Seti















"Hygiene for health, Disposable Sterile Products for hygiene"

MINOR SURGERY PACK / MINÖR CERRAHI SETI



Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Minor Surgery Drape / Minör Cerrahi Örtüsü	75x75	1 Pc.
2. Medical Towel / Medikal Havlu	40x40	1 Pc.
3. Back Table Cover / Alet Masa Örtüsü	40x50	1 Pc.
4. Sterilization Wrap / Sterilizasyon Bohçası	50x50	1 Pc.

11 DRESSING PACK / PANSUMAN SETI



	Product Name Ürün Adı	Size Ebat	Pcs. Adet
1.	Fenestrated Drape with Adhesive Tape/Yapışkan Bantlı Delikli Örtü	60x60	1 Pc.
2.	250 cc Solution Cup / Solüsyon Kabı		1 Pc.
3.	Betadine Sponge / Boyama Süngeri		1 Pc.
4.	Universal Bag / Universal Torba	30x35	1 Pc.
5.	Sponge / Spanç	7.5x7.5	5 Pc.
6.	Gloves / Eldiven		1 Pair
7.	Medical Towel / Medikal Havlu	40x40	1 Pc.
8.	Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



Orthopedic Packs Ortopedi Setleri













"Hygiene for health, Disposable Sterile Products for hygiene"


ARTHROSCOPY PACK / ARTROSKOPİ ÖRTÜ SETİ

Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Arthroscopy Drape with LCP / Poşlu Artroskopi Örtüsü	200x300	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
3. Leggings / Bacak Kılıfı	40x75	2 Pc.
4. Medical Towel / Medikal HavluXL	40x40	2 Pc.
5. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		2 Pc.
6. Skin Marker & Ruler / Cilt Kalemi ve Cetveli		1 Pc.
7. Elastic Bandage / Elastik Bandaj	15x150	1 Pc.
8. Universal Bag / Universal Torba	30x35	2 Pc.
9. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x30	2 Pc.
10. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



ARTHROSCOPY PACK WITHOUT POUCH / POSSUZ ARTROSKOPİ SETİ


Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Arthroscopy Drape / Artroskopi Örtüsü	200x310	1 Pc.
2. Leggings / Bacak Kılıfı	40x75	2 Pc.
3. Elastic Bandage / Elastik Bandaj	15x150	1 Pc.
4. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x50	1 Pc.
5. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.












Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Extremity Knee Drape / Ekstremite Diz Örtüsü	200x310	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü XL	150x200	1 Pc.
3. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		2 Pc.
4. Medical Towel / Medikal Havlu	40x40	2 Pc.
5. Skin Marker & Ruler / Cilt Kalemi ve Cetveli		1 Pc.
6. Universal Bag / Universal Torba	30x35	2 Pc.
7. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.

13 HIP PACK / KALÇA ÖRTÜ SETİ



Product Name Ürün Adı		Pcs. Adet
1. Drape with U Aperture / U Yarıklı Örtü	200x240	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü XL	150x200	1 Pc.
3. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		2 Pc.
4. Medical Towel / Medikal Havlu	40x40	2 Pc.
5. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x50	2 Pc.
6. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	150x240	1 Pc.
7. Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
8. Leggings / Bacak Kılıfı	40x75	8 Pc.
9. Universal Bag / Universal Torba	30x35	2 Pc.
10. Elastic Bandage / Elastik Bandaj	10x150	1 Pc.
11. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



Orthopedic Packs Ortopedi Setleri














"Hygiene for health, Disposable Sterile Products for hygiene"


VERTEBRA PACK / VERTEBRA ÖRTÜ SETİ

Product Name	Size	Pcs.
1. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	150x240	2 Pc.
2. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	50x120	2 Pc.
3. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	100x120	2 Pc.
4. Back Table Cover / Alet Masa Örtüsü	150x200	2 Pc.
5. Universal Bag / Universal Torba	30x35	1 Pc.
6. Skin Marker & Ruler / Cilt Kalemi ve Cetveli XL		1 Pc.
7. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		2 Pc.
8. Medical Towel / Medikal Havlu	40x40	2 Pc.
9. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x50	2 Pc.
0. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



HAND/FOOT DRAPE PACK / EL-AYAK CERRAHİ ÖRTÜ SETİ

Product Name Ürün Adı		Pcs. Adet
1. Hand/Foot Drape / El-Ayak Örtüsü	200x300	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
3. Medical Towel / Medikal Havlu	40x40	3 Pc.
4. Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
5. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x50	1 Pc.
6. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.





Cardiovascular Packs Kardiyovasküler Örtü Setleri



UNICERT ISO 13485:2016




UNICERT ISO 9001:2015



€1023

"Hygiene for health, Disposable Sterile Products for hygiene"

BY PASS DRAPE PACK / BY PASS ÖRTÜ SETİ



Product Name Ürün Adı		Pcs. Adet
1. Bypass Drape / Bypass Örtüsü	240x340	1 Pc.
2. Perineum Drape with U Aperture / U Yarıklı Perine Örtüsü	150x240	1 Pc.
3. Back Table Cover / Alet Masa Örtüsü	150x200	2 Pc.
4. Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
5. Leggings / Ayak Kılıfı XXL		Pair.
6. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		3 Pc.
7. Medical Towel / Medikal Havlu	40x40	3 Pc.
8. Universal Bag / Universal Torba	30x35	2 Pc.
9. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x50	2 Pc.
10. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.

15 CARDIOVASCULAR PACK / KARDIOVASKÜLER ÖRTÜ SETI



Product Name Ürün Adı		Pcs. Adet
1. C.V.C Drape / K.V.C Örtüsü	200x300	1 Pc.
2. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	75x90	3 Pc.
3. Back Table Cover / Alet Masa Örtüsü	150x200	2 Pc.
4. Mayo Stand Cover / Mayo Masa Örtüsü XXL	78x145	1 Pc.
5. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		3 Pc.
6. Medical Towel / Medikal Havlu	40x40	3 Pc.
7. Universal Bag / Universal Torba	30x35	2 Pc.
8. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



Cardiovascular Packs Kardiyovasküler Örtü Setleri














"Hygiene for health, Disposable Sterile Products for hygiene"


VALVE REPLACEMENT PACK / KAPAK DEĞİŞİM SETİ

	Product Name Ürün Adı		Pcs. Adet
1.	Valve Replacement Drape / Kapak Örtüsü	240x340	1 Pc.
2.	Back Table Cover / Alet Masa Örtüsü	150x200	2 Pc.
3.	Mayo Stand Cover / Mayo Masa Örtüsü XXL	78x145	1 Pc.
4.	Reinforced Surgical Gown / Korumalı Cerrahi Önlük		3 Pc.
5.	Medical Towel / Medikal Havlu	40x40	3 Pc.
6.	Universal Bag / Universal Torba	35x35	2 Pc.
7.	Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x50	2 Pc.
8.	Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



THORACIC PACK / TORAKS ÖRTÜ SETİ

Product Name Ürün Adı		Pcs. Adet
1. Thoracic Drape / Toraks Örtüsü	250x340	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
3. Mayo Stand Cover / Mayo Masa Örtüsü XXL	78x145	1 Pc.
4. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		3 Pc.
5. Medical Towel / Medikal Havlu	40x40	6 Pc.
6. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	100x120	1 Pc.
7. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x50	1 Pc.
8. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.





Obstetrics & Gynecology Surgery Packs Kadın Dogum & Jinekoloji Setleri



UNICERT ISO 13485:2016




UNICERT ISO 9001:2015



€1023


"Hygiene for health, Disposable Sterile Products for hygiene"

CAESAREAN SECTION PACK / SEZERYAN ÖRTÜ SETİ



Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Section Drape L.T.T / Sezeryan Örtüsü L.T.T	200x300	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
3. Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
4. Standard Surgical Gown XL / Standart Cerrahi Önlük XL		2 Pc.
5. Medical Towel / Medikal Havlu	40x40	2 Pc.
6. Baby Drape / Bebek Örtüsü	75x75	1 Pc.
7. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.

17 O.P.U DRAPE PACK / O.P.U ÖRTÜ SETİ



Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. O.P.U Drape with hollow / Delikli O.P.U Örtüsü	130x150	1 Pc.
2. Plain Drape / Düz Örtü	75x100	1 Pc.
3. Transparent Drape / Şeffaf Örtü	60x60	1 Pc.
4. Back Table Cover / Alet Masa Örtüsü	100x150	1 Pc.
5. Standard Surgical Gown XL / Standart Cerrahi Önlük XL		2 Pc.
6. Medical Towe I/ Medikal Havlu	40x40	2 Pc.
7. Leg Cover / Bacak Kılıfı	50x110	2 Pc.
8. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



Obstetrics & Gynecology Surgery Packs Kadın Dogum & Jinekoloji Setleri














"Hygiene for health, Disposable Sterile Products for hygiene"


DELIVERY DRAPE PACK / DOĞUM ÖRTÜ SETİ


Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Abdominal Drape / Abdominal Örtü	75x150	1 Pc.
2. Plain Drape / Düz Örtü	75x100	1 Pc.
3. Back Table Cover / Alet Masa Örtüsü	100x150	1 Pc.
4. Reinforced Surgical Gown XL / Korumalı Cerrahi Önlük XL		2 Pc.
5. Medical Towel / Medikal Havlu	40x40	2 Pc.
6. Leg Cover / Bacak Kılıfı	50x110	2 Pc.
7. Baby Drape / Bebek Örtüsü	75x75	1 Pc.
8. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



E.T DRAPE PACK / E.T ÖRTÜ SETİ

Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Perforated Drape With Adhesive Tape/DelikliYapışkan Bantlı Örtü	80x80	1 Pc.
2. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	50x90	1 Pc.
3. Sponge / Spanç	7,5x7,5	10 Pc.
4. Sterilization Wrap / Sterilizasyon Bohçası	75x75	1 Pc.

















Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. TUR Drape / TUR Örtüsü	180x200	1 Pc.
2. Plain Drape / Düz Örtü	100x120	1 Pc.
3. Back Table Cover / Alet Masa Örtüsü	100x150	1 Pc.
4. Reinforced Surgical Gown XL / Korumalı Cerrahi Önlük XL		2 Pc.
5. Medical Towel / Medikal Havlu	40x40	2 Pc.
6. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.

19 UROLOGY DRAPE PACK / ÜROLOJİ ÖRTÜ SETİ



Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Urology Drape / Üroloji Örtüsü	100x150	1 Pc.
2. Plain Drape / Düz Örtü	100x120	1 Pc.
3. Back Table Cover Standart / Alet Masa Örtüsü	150x200	1 Pc.
4. Surgical Gown XL / Standart Cerrahi Önlük XL		2 Pc.
5. Medical Towel / Medikal Havlu	40x40	2 Pc.
6. Leg Cover / Bacak Kılıfı	60x115	2 Pc.
7. Sterilization Wrap / Sterilizasyon Bohçası	75x75	1 Pc.



Urology Drape Packs Üroloji Örtü Setleri














"Hygiene for health, Disposable Sterile Products for hygiene"

CYSTOSCOPY DRAPE PACK / SISTOPKOPI ÖRTÜ SETİ


Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Cystoscopy Drape / Sistoskopi Örtüsü	75x150	1 Pc.
2. Plain Drape / Düz Örtü	100x120	1 Pc.
3. Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
4. Standart Surgical Gown XL / Standart Cerrahi Önlük XL		2 Pc.
5. Medical Towel / Medikal Havlu	40x40	2 Pc.
6. Leg Cover / Bacak Kılıfı	60x115	2 Pc.
7. Sterilization Wrap / Sterilizasyon Bohçası	75x75	1 Pc.




PERCUTANEOUS DRAPE PACK / PERKÜTAN ÖRTÜ SETİ

		١
4	U	J

Product Name Ürün Adı	Size Ebat	Pcs. Adet
I. Percutaneous Drape / Perkütan Örtüsü	150x300	1 Pc.
2. Camera Cover / Kamera Kılıfı		1 Pc.
3. Scopy Cover / C Skopi Kılıfı		1 Pc.
Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
5. Reinforced Surgical Gown XL / Korumalı Cerrahi Önlük XL		2 Pc.
6. Medical Towel / Medikal Havlu	40x40	2 Pc.
7. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x30	2 Pc.
3. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.







Cardiology Packs Kardiyoloji Örtü Setleri















"Hygiene for health, Disposable Sterile Products for hygiene"

ANGIOGRAPHY PACK / ANJIOGRAFI ÖRTÜ SETİ -I-



Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Angiography Drape / Anjio Örtüsü	180x300	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü XL	100x150	1 Pc.
3. Standard Surgical Gown / Standart Cerrahi Önlük		2 Pc.
4. Medical Towel / Medikal Havlu	40x40	2 Pc.
5. Sponge / Spanç	7.5x7.5	10 Pc.
6. Fluoroscopy Cover 86cm / Fluoroskopi Kılıfı 86cm		1 Pc.
7. Solution Bowl 200cc / Solüsyon Kabı 200cc		2 Pc.
8. Tray / Tepsi		1 Pc.
9. Betadine Sponge / Boyama Süngeri		1 Pc.
10. Opaque Set / Opak Set		1 Pc.
11. Scalpel Tip No:11/ Bistüri Ucu No:11		1 Pc.
12. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.

21 ANGIOGRAPHY PACK / ANJİOGRAFİ ÖRTÜ SETİ -II-



Product Name Ürün Adı		Pcs. Adet
1. Angiography Drape / Anjio Örtüsü	150x300	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü XL	100x150	1 Pc.
3. Standard Surgical Gown / Standart Cerrahi Önlük		2 Pc.
4. Medical Towel / Medikal Havlu	40x40	2 Pc.
5. Fluoroscopy Cover 86cm / Fluoroskopi Kılıfı 86cm		1 Pc.
6. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



E.N.T Surgery Packs KBB Örtü Setleri














"Hygiene for health, Disposable Sterile Products for hygiene"


EAR-NOSE-THROAT PACK / KBB ÖRTÜ SETİ

Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Splint E.N.T U Drape / K.B.B U Örtüsü	150x250	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü Standard	100x150	1 Pc.
3. Surgical Gown / Standart Cerrahi Önlük		2 Pc.
4. Medical Towel / Medikal Havlu	40x40	2 Pc.
5. Fenestrated Drape with Adhesive Tape / Bantlı Delikli Örtü	100x100	1 Pc.
6. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.





THYROID PACK / TROID ÖRTÜ SETİ 22

Product Name Ürün Adı	Size Ebat	Pcs. Adet
. Thyroid Drape / Troid Örtüsü	240x320	1 Pc.
. Back Table Cover / Alet Masa Örtüsü	100x150	1 Pc.
s. Standard Surgical Gown / Standart Cerrahi Önlük		2 Pc.
. Medical Towel / Medikal Havlu	40x40	2 Pc.
. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x30	2 Pc.
. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.











Neurosurgery Packs Beyin Cerrahi Örtü Setleri








UNICERT ISO 9001:2015



€1023


"Hygiene for health, Disposable Sterile Products for hygiene"

SPINAL (LAMINECTOMY) PACK / SPINAL (LAMINEKTOMI) ÖRTÜ SETİ



Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Spinal Vertebra Drape / Spinal Vertebra Örtüsü	160x300	1 Pc.
2. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	100x120	4 Pc.
3. Back Table Cover / Alet Masa Örtüsü	100x150	1 Pc.
4. Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
5. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		3 Pc.
6. Medical Towel / Medikal Havlu	40x40	4 Pc.
7. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x30	4 Pc.
8. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.

23 CRANIOTOMY PACK / KRANIOTOMI ÖRTÜ SETİ



Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Craniotomy Drape / Kraniotomi Örtüsü	220x350	1 Pc.
2. Plain Drape / Düz Örtü	150x250	1 Pc.
3. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	100x100	4 Pc.
4. Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
5. Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
6. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		3 Pc.
7. Medical Towel / Medikal Havlu	40x40	4 Pc.
8. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



Neurosurgery & Biopsy Packs Beyin Cerrahi & Biyopsi Örtü Setleri














"Hygiene for health, Disposable Sterile Products for hygiene"

SHUNT PACK / SHUNT ÖRTÜ SETİ

Product Name Ürün Adı		Pcs. Adet
1. Shunt Drape / Shunt Örtüsü 100x120	150x300	1 Pc.
2. Drape with Adhesive Tape / Yapışkan Bantlı Örtü		4 Pc.
3. Plain Drape / Düz Örtü	150x180	1 Pc.
4. Anesthesia Drape / Anestezi Örtüsü	150x240	1 Pc.
5. Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
6. Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
7. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		3 Pc.
8. Medical Towel / Medikal Havlu	40x40	4 Pc.
9. Sterilization Wrap / Sterilizasyon Bohcası	100x100	1 Pc.



BIOPSY PACK / BİYOPSİ ÖRTÜ SETİ


	\mathbf{A}
2	4

Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Biopsy Drape / Biyopsi Örtüsü	120x150	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü	100x150	1 Pc.
3. Standard Surgical Gown / Standart Cerrahi Önlük		2 Pc.
4. Medical Towel / Medikal Havlu	40x40	2 Pc.
5. Sponge / Spanç	7.5x7.5	10 Pc.
6. Betadine Sponge / Boyama Süngeri		1 Pc.
7. Solution Bowl 500cc / Solüsyon Kabı 500cc		1 Pc.
8. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.






OPTHALMIC DRAPE PACK / GÖZ ÖRTÜ SETİ



Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Opthalmic Drape with Pouch / Poşlu Göz Örtüsü	100x130	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü	100x150	1 Pc.
3. Standard Surgical Gown / Standart Cerrahi Önlük		2 Pc.
4. Medical Towel / Medikal Havlu	40x40	2 Pc.
5. Plain Drape / Düz Örtü	150x170	1 Pc.
6. Sterilization Wrap / Sterilizasyon Bohçası	75x75	1 Pc.

25 IMPLANT PACK / İMPLANT ÖRTÜ SETİ



Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Patient Upper Drape / Hasta Üst Örtüsü	75x90	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü	75x75	1 Pc.
3. Reflector Drape / Reflektör Örtüsü	15x15	2 Pc.
4. Transparent Hose Cover / Şeffaf Hortum Kılıfı		2 Pc.
5. Sponge / Spanç	7.5x7.5	5 Pc.
6. Medical Towel / Medikal Havlu	40x40	1 Pc.
7. Universal Bag / Universal Torba	30x35	1 Pc.
8. Sterilization Wrap / Sterilizasyon Bohçası	75x75	1 Pc.



Chemotherapy & Plastic Surgery Packs Kemoterapi & Plastik Cerrahi Setleri














"Hygiene for health, Disposable Sterile Products for hygiene"


CHEMOTHERAPY PACK / KEMOTERAPİ ÖRTÜ SETİ


Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Chemotherapy Drape / Kemoterapi Örtüsü	50x70	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü	50x100	1 Pc.
3. Chemotherapy Gown / Kemoterapi Önlüğü		1 Pc.
4. Chemotherapy Gloves / Kemoterapi Eldiveni		1 Pair/ <i>Çift</i>
5. Protective Glasses / Koruyucu Gözlük		1 Pc.
6. Mask / Ventilli		1 Pc.
7. Bonnet / Bone		1 Pc.
8. Medical Towel / Medikal Havlu	40x40	2 Pc.
9. Medical Waste Bag / Tıbbi Atık Torbası		1 Pc.
10. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



PLASTIC SURGERY PACKS / PLASTİK CERRAHİ ÖRTÜ SETLERİ

Product Name Ürün Adı	Size Ebat	Pcs. Adet
1. Adhesive Drape / Yapışkan Bantlı Örtü	150x170	1 Pc.
2. Back Table Cover / Alet Masa Örtüsü	100x150	1 Pc.
3. Anesthesia Drape / Anestezi Örtüsü	150x240	1 Pc.
4. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	75x90	2 Pc.
5. Sterilization Wrap / Sterilizasyon Bohçası	75x75	1 Pc.







Paediatric Packs Pediatrik Setler












€1023


"Hygiene for health, Disposable Sterile Products for hygiene"

PEDIATRIC UNIVERSAL PACK / PEDIATRIK GENEL CERRAHI SETI



Product Name Ürün Adı		Pcs. Adet
1. Anesthesia Drape / Anestezi Örtüsü	150x240	1 Pc.
2. Foot Drape with Adhesive Tape / Yapışkan Bantlı Ayak Örtüsü	175x175	1 Pc.
3. Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
4. Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
5. Drape with Adhesive Tape / Yapışkan Bantlı Örtü	75x75	2 Pc.
6. Reinforced Surgical Gown / Korumalı Cerrahi Önlük		2 Pc.
7. Medical Towel / Medikal Havlu	33x34	2 Pc.
8. Nonwoven Op-Tape / Nonwoven Yapışkan Bant	10x50	1 Pc.
9. Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.

27 PAEDIATRIC LAPAROSCOPY PACK / PEDIATRIK LAPAROSKOPI SETI



	Product Name Ürün Adı		Pcs. Adet
1.	Paediatric Laparascopy Drape (10x10 Ins. Film) / Pediatrik Laparoskopi Örtüsü (10x10 Ins. Film)	150x180	1 Pc.
2.	Back Table Cover / Alet Masa Örtüsü	150x200	1 Pc.
3.	Mayo Stand Cover / Mayo Masa Örtüsü	78x145	1 Pc.
4.	Medical Towel / Medikal Havlu	40x40	2 Pc.
5.	Reinforced Surgical Gown / Korumalı Önlük		2 Pc.
6.	Sponge / Spanç	7.5x7.5	10 Pc.
7.	Sterilization Wrap / Sterilizasyon Bohçası	100x100	1 Pc.



NON STERILE GOWNS / STERIL OLMAYAN ÖNLÜKLER





Patient Gown	М
(Long Sleeve)	L
(====,	XL
	XXL

Patient Gown (Short Sleeve)	М
	L
	XL
	XXL



Chemotherapy Gown	М
	L
	XL
	XXL




Urology Gown	М
	L
	XL
	XXL

28



COLONOSCOPY PANTS/SHORTS / KOLONOSKOPİ PANTOLONU / ŞORTU






29

Colonescent Dente	L
Colonoscopy Pants	XL
Colonoscopy Shorts	L
	XL



STERILIZATION WRAP / STERILIZASYON BOHÇALARI



Sterilization Wrap	50x50
Otorinization wrap	75x75
	90x90
	100x100
	120x120
	130x150

BED SHEETS / YATAK ÖRTÜSÜ

SMS	80x220
	150x220
Trilaminate Green	80x220
	150x220
Trilaminate Blue	80x220
	150x220
Patient Bed Set	Bed Sheet+Pillow Cover+ Blanket Cover

30















Yunusemre Mh. 12. Yıldırım Sk. No:17 Yıldırım / Bursa Tel.: 0224. 369 32 32 Faks: 0224. 369 32 36 E-mail: info@zeynimedical.com.tr





This Certificate has been awarded to

ZEYNİ MEDİKAL TEKSTİL İNŞAAT SANAYİ VE TİCARET LİMİTED ŞİRKETİ

YUNUSEMRE MAHALLESİ 12. YILDIRIM SOKAK NO:17
YILDIRIM / BURSA / TÜRKİYE

In recognition of the organization's Management System which complies with

EN 13795-1:2019

The scope of activities covered by this certificate is defined below

MANUFACTURING, SALES OF DISPOSABLE, STERILE AND NON-STERILE SURGICAL GOWNS, INSULATED VISITOR GOWNS AND COVERS

TEK KULLANIMLIK, STERİL VE STERİL OLMAYAN CERRAHİ ÖNLÜK, İZOLASYONLU ZİYARETÇİ ÖNLÜĞÜ VE ÖRTÜLERİN İMALATI, SATIŞI

Certificate Number: SISTUREN082020161 Date of Issue of Original Certificate: 21.08.2020 Date of Issue of latest certificate: 23.09.2021

Expiry Date: 20.08.2022

Re-certification Due on: 21.07.2023

Managing Director

Note: This certificate is valid only if produced with the continuation letter after the surveillance is carried out successfully.

The Organization's documentation and Implementation has been reviewed and found to comply with the relevant standard rules. This certificate of Compliance is based on the evaluation of the mentioned scope given above. Organization is responsible for maintaining the responsibilities of the relevant standard rules. Any significant changes in the scope of the certification or standard referred above render this certificate invalid

Corporate office(SIS): Unit No. 514, 5th Floor, Vipul Business Park, Sohna Road, Sector-48, Gurgaon-122018, Haryana, India. International office(SIS):- URB. Santa Ana Cal. German, Scherieber 276, San Isidro, Lima, Peru 15047. Email us :-support@siscertifications.com, info@siscertifications.co.in. Call:- +91 99105 01396, + 91 96430 73391. The status of this certificate can be verified on "http://www.siscertifications.com/verify/" Web:- http://www.siscertifications.co.in, www.siscertifications.com/









Page 1 of 6

Version No / Date: 1.0

ZEYNİ MEDİKAL TEKSTİL İNŞ.SAN. VE TİC. LTD. ŞTİ.

YUNUSEMRE MAH. 12 YILDIRIM SOKAK NO:17 YILDIRIM/BURSA

Report No.	111699085
Buyer	/
Test Item. :	Nonwoven fabric,3-layer, SMS
Item No.:	/
Colour Name. :	/

Condition at delivery. : Samples tested as received.

Test Scope.: Parameters selected by customer

Test Specification: Determination of hydrostatic pressure

Applicant's Provided Care Instruction/Label:

Sample Receiving date: 2021-10-18 (p.m)

Testing Period: 2021-10-21 to 2021-10-25

Test Result: No Comment

For and on behalf of

/ 2014-07-02

TÜV Rheinland Uluslararası Standartlar Sertifikasyon ve Denetim A.Ş

Tomris Hasançebi / Customer Relations Manager

Abdullah Akil / Physical Laboratory Manager

TÜV Rheinland Uluslararası Standartlar Sertifikasyon ve Denetim A.Ş.

Kozyatagi Mah. Saniye Ermutlu Sok. No:12

Colakoglu Plaza B Blok 34742 Kadikoy Istanbul,

Tel. +902166653200, Fax +902166653299, e-mail: info@tr.tuv.com

TÜVRheinland® Precisely Right.

Report No.: 111699085

Date: 10.25.2021

Page 2 of 6



TÜVRheinland®
Precisely Right.

Report No.: 111699085

Date: 10.25.2021

Material List: Page 3 of 6

Material No.	Material	Color	Location
M001	Textile	-	Textile
M002	Textile + coating	-	Textile w /coating

TÜVRheinland®

Report No.: 111699085

Date: 10.25.2021

Conclusion: Page 4 of 6

TEST PROPERTY M001 M002

Hydrostatic Pressure Test # #

TÜVRheinland®
Precisely Right.

Report No.: 111699085

Date: 10.25.2021

1.Hydrostatic Pressure Test Page 5 of 6

Test method : EN 20811:1993

Face Side: Original

Test Conditions: 65 cm2 Test Speed:65 cmH2O/dak

M001 Requirement

Pressure

Average 65 cm H2O 65 cm H2O -

- END -

Report No.: 111699085

Date: 10.25.2021



Page 6 of 6





1. PRODUCT

ANGIOGRAPHY SET

2. CODE

Model 601-001-S.

3. PRODUCT CONFIRMITIES

ISO 9001, ISO 13485, EN 13795, ISO 11135:2014, ISO 14644-3:2019, ISO 14791:2012, ISO 14937: 2009, ISO 10993-7:2008, EN 20811

4. CONTENT

Angiography Drape 240 CM X 370 CM 1 Pc.

Adhesive Drape 50 CM X 50 CM 1 Pc.

C-Arm Cover 140 CM 1 Pc.

Superabsorbent Drape 35 CM X 50 CM 1 Pc.

Instrument Table Field 150 CM X 190 CM 1 Pc.

5. PACKING

The product is individual packaged. On the packaging is the inscription with the name and address of the manufacturer, the product model, description, dimensions, date of manufacture, validity date, lot number, etc.

6. STORAGE

0°C - + 30°C, Avoid Sunlight, Avoid Humidity


7. Validity

Sterile EO. 3 years from the date of production





Content

















SCIENTIFIC AND TECHNOLOGICAL RESEARCH COUNCIL OF TURKEY MARMARA RESEARCH CENTER GENETIC ENGINEERING and BIOTECHNOLOGY INSTITUTE

P.K.21,41470 GEBZE–KOCAELI T02626772000 F02626412309 http://www.mam.gov.tr

CERTIFICATE OF ANALYSIS

(Industrial Technical Support Service)

Report no :16563500-125.05-89 /4654

Report date :23.08.2017

Applicant : ZEYNİ MEDİKAL TEKSTİL İNŞAAT SAN. VE TİC. LTD. ŞTİ.

Address : YUNUSEMRE MAH. 12. YILDIRIM SK. NO:17 YILDIRIM BURSA

Subject :CYTOTOXICITY TEST CARRIED OUT for 'SURGERY DRAPE' IN THE

SCOPEOF BIOCOMPATIBILITY TESTS

The results included in this report are related to only the sample analyzed.

Approved by

Assoc. Prof. Dr. FatimaYÜCEL GMBE Industrial Services Officer

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited.

The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).

Page 1/6









Report no :16563500-125.05-89/4654

:ZEYNİ MEDİKAL TEKSTİL İNŞAAT SAN. VE TİC. LTD.ŞTİ **Applicant** : YUNUSEMRE MAH. 12. YILDIRIM SK. NO:17 YILDIRIM BURSA Address of Applicant

Sample : Standardized sample

Sample number : 3

Delivery type of sample : By Cargo

Situation in delivery time: It was provided in closed sterilization

packets under sterile conditions

Sample registration no at

Expiry date: 05/2022

Institution: 17/70-GMBE

Acceptance date and hour: 11.08.2017

Analysis date :14.08.2017-17.08.2017

Witness sample ()Return to customer (x)Witness sample is ()Witness sample was not taken information: available

1-Samples

The standardized 3 samples, defined as 'Surgery Set', were analyzed for cytotoxicity test upon the application of ZEYNİ MEDİKAL TEKSTİL İNŞAAT SAN. VE TİC. LTD. ŞTİ. dated 14.08.2017 and numbered 3940.

Sample	Characteristics	Item
Sterile Surgical Drapes	The product is a set consisted of medical textile products in various sizes and forms used for various purposes during surgeries. 1- Drape (100x150cm)	3
	Production Date: 05/2017, STERILE EO LOT 000001, Latex Free, Single Use Only	

Table 1.Product tested.

Desc		

Signatories:

53192 (signature)

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, exceptPublic Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited.

The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).

Page 2/6

P.K.21,41470 GEBZE-KOCAELÏ



Report no:16563500-125.05-89/4654

 Sample	Specifications	Item
Surgery Set		3
	The state of the s	

Table1 (contd.). Product tested.

Descriptions:

Signatories:

53192

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited.

The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).

Page 3/6



Report no:16563500-125.05-89/4654

Cytotoxicity Test

The cytotoxicity tests were carried out in accordance with 'Biological evaluation of medical products: ISO 10993-5: 2009 Tests for in vitro cytotoxicity' standards.

Beginning date of test:

Expiry date of test:

Sample description: The sample characteristics are as explained in the Part 1. The samples were provided by 'ZEYNİ MEDİKAL TEKSTİL İNŞAAT SAN. VE TİC. LTD. ŞTİ.' firm. The tests werecarried out by taking parts from surgery set selected as representing the product.

Explanation of cell strain used and its reasons: L929 rat cell strain was used. It was selected since it is one of the cell strains recommended by ISO 10993-5 and appropriate for representation of mammals system.

Explanation of cell strain used and its reasons: L929 rat cell strain was used. It was selected since it is one of the cell strains recommended by ISO 10993-5 and appropriate for representation of mammals system.

Name and batch number of the firm provided the feed used and the serum and antibiotics added: DMEM/F12 (Sigma cat # D0547-10X-1, lot #SLBH4393) + 10% Fetal bovine serum (Gibco Life Sciences Cat # 10270-106, lot #41G2943K)+ penicilinle streptomycin (Gibco Life Sciences, Cat # 15140-122, lot # 1688254) were used as feed lot. Test method: Direct contact method.

Rational: The analyses of toxic effect of solid matters, dissolved in liquid, as a result of contact with cell.

Direct contact protocol: The 'Surgical Drape' was provided as sterile and not subjected to a sterilization process once again. The 'Surgical Drape' was cut to the extent of 1/10 of total surface area of cells and put into middle of pits over the cells.

Test protocol: L929 cells were counted and implanted in &' pits as 4 X 10⁵ cell/pit and incubated for 24 hours at 37°C and 5% CO₂. The samples prepared as explained above and the controls were immediately added on cells and incubated for 24 hours at 37°C and 5% CO₂. The samples were contacted with cells for 24 hours and then observed microscopically. Following this process, WST-1 agent at 1:50 rate was added onto the cells and the color formation was waited for 2 hours. 100 μ liquid from each pit was taken and poured into 3 pits in 96 sizes so that the absorbance measure was realized at microplate reader, at 450 nm and 650 nm reference wave length for viability test.

Cytotoxicity Measure method: WST-1 cell viability analyses (Colorimetric)

Ratio	ona	al:	The o	consiste	ent and	sensitive	measu	ırement	of cell	viabilit	y.
	-			,	,	<u> </u>	•	•	,	•	

Rational: The consist	ent and sensitive measuremen	t of cell viability.
Descriptions:		
Signatories:	52966	53077

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited.

The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).

Page



Report no: 16563500-125.05-89/4654

Negative, positive and other controls:

Controls:

Control 1- The plates on which any material was not put but only the cells were implanted

Control 2- The sterile and dry whatmann paper cut in 1 cm X 1 cm sizes.

Negative Control: Whatmann paper saturated with PBS (25_{ul}) and cut in 1cm X 1cm sizes.

Positive Control: Whatmann paper saturated with EMS (Ethyl methanesulfonate, Merck Millipore, 820774) (25ul) and cut in 1cm X 1cm sizes.

RESULT:

Cell response and other observations:

Quantitative Evaluation

It was observed that there was no cell direct under the sample. However, the same result was observed in PBS saturated whatmann paper samples. A zonal toxic effect was not observed around the sample. This observation was reflected in qualitative evaluation.

	Degree		Degree
Sample		Control 1, Empty	0
Surgery Drape	0	Control 2, dry Whatmann paper	0
		Negative Control, PBS saturated Whatmann	0
		paper	
		Positive Control, EMS saturated Whatmann paper	4

Table 2. Quantitative Evaluation.

Degree	Reactivity	Situation
0	None	Clear intracytoplasmic granuls, no any cell lysis, no impact on cell growth.
1	Less	The rate of cells, which are become rounded, representing poor adhesion and/or having no intracytoplasmic granules, showing morphological changes and partly lysis, is below 20%.
2	Slight	The rate of cells, which are become rounded and having no intracytoplasmic garnuls, is below 50%, there is no extensive cell lysis; growth inhibition is below 50%.
3	Moderate	Less than 70% of cells are rounded or showing lysis and growth inhibition is not more than 50%.
4	Serious	All or almost all of the cells are showing lysis.

Table 3. Evaluation criteria and scoring for cytotoxicity test.

Descriptions:		
Signatories:	52966	53077

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited.

The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).



Report no: 16563500-125.05-89/4654

According to the evaluation criteria given in the Table 3, L929 cell culture was exposed to 'Surgical Drape' for 24 hours in order to examine its potential for causing to cell death according to the control and the cell viability was found as %98.11 ±1.31 according to the control as a result of ISO 10993-5 cytotoxicity test carried out in the scope of biocompatibility.




Figure 1. The viability analyses, carried out by using extracts from samples called 'Surgical Drape', are shown in the Figure.

Cytotoxicity of 'Surgical Drape' sample; The absorbance values of samples were normalized by

using Negative Control 2 absorbance values in parallel with the samples for 100% viability.

The data, obtained for each sample, was provided as a result of triplicated process in tests of 3 extracts from 3 randomly selected samples. EVĂLUATION:

No cytotoxic effect according to the viability test carried out by using 'Surgical Drape' extraction method.

Descriptions:		
Signatories:	52966	53077
		y your

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited.

The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).



SCIENTIFIC AND TECHNOLOGICAL RESEARCH COUNCIL OF TURKEY MARMARA RESEARCH CENTER GENETIC ENGINEERING and BIOTECHNOLOGY INSTITUTE

P.K.21,41470 GEBZE–KOCAELİ T02626772000 F02626412309 http://www.mam.gov.tr

CERTIFICATE of ANALYSIS-R

(Industrial Technical Support Service)

Report no : 16563500-125.05-64/3882-2

Report date : 22.06.2017

Applicant : ZEYNİ MEDİKAL TEKSTİL İNŞAAT SAN. VE TİC. LTD. ŞTİ.

Address : YUNUSEMRE MAH. 12. YILDIRIM SK. NO:17 YILDIRIM BURSA

Subject : SENSITIZATION TEST CARRIED OUT FOR 'SURGERY SET' IN THE SCOPE

OF BIOCOMPATIBILITY TESTS

The results included in this report are related to only the sample analyzed.

Approved by

Assoc. Prof. Dr. Fatima YÜCEL

GMBE Industrial Services Officer

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited. The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).

Page 1/6









Report no : 16563500-125.05- 64 / 3882-2

Applicant : ZEYNİ MEDİKAL TEKSTİL İNŞAAT SAN. VE TİC. LTD. ŞTİ. Address of Applicant : YUNUSEMRE MAH. 12. YILDIRIM SK. NO:17 YILDIRIM BURSA

Sample : Standardized sample

Sample number 3

Delivery type of sample: By Cargo

Situation in delivery time: It was provided in closed

sterilization packets under sterile conditions.

Acceptance date and hour 16.05.2017

Sample registration no at

Institution: 17/35-GMBE

Expiry date

Analysis date : 07.06.2017 – 21.06.2017

: 05/2022

Witness sample information: () Return to customer (x) Witness sample is available () Witness sample was not taken

1- Samples:

The standardized 3 samples, defined as 'Surgery Set', were analyzed for sensitization tests upon the application of ZEYNİ MEDİKAL TEKSTİL İNŞAAT SAN. VE TİC. LTD. ŞTİ. dated 16/05/2017 and numbered 2669.

Table 1. The tested product.

Sample	Characteristics	Item
Sterile Surgical Drapes	The product is a set consisted of medical textile products in various sizes and forms used for various purposes during surgeries. 1- Drape (100x150cm)	3
	Production date: 05/2017, STERILE EO LOT 000001, Latex Free, Single Use Only	

Descriptions:

Signatories:

53192

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited.

The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).

Page 2/6



Report no: 16563500-125.05- 64 / 3882-2

Sample	Specifications	Item
Surgery Set	The state of the s	3

Table 1 (contd.). The tested product.

Descriptions:

Signatories:

53192

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited. The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute

archive). Page 3/6



Report no: 16563500-125.05- 64 / 3882-2

2- Skin Sensitization Test

Sensitization test was carried out in accordance with 'ISO 10993-10: 2010 tests for irritation and delayed-type hypersensitivity' standard protocol.

The extract was provided by applying 6 cm²/ml surface area/volume rate in compliance with the form and structure of the product. For this purpose, the incubation for 72 hours at 37°C was applied. The sensitization test was carried out by using adult female subjects of guinea pig (Cavia porcellus) family weighted between 300-500 gr. As stated in the document titled ISO 10993-10:2010, the tests were carried out with 0,1 ml subcutaneous use of material to be tested. The topical application was made to the region to which subcutaneous injection (intradermal induction phase) is not applied as left region of the animal in the 7th day and rigt region of the animal in the 14th day of the test. The application plan administered on experimental animals is shown in Figure 1.




Figure 1.

- 1- Head of experimental animal.
- A- The test regions treated by mixing Freund's Complete Anjuvant (FCA) and serum physiological solution at the rate of 50:50.
- B- The test regions treated by using only the test material.
- C- The test regions treated by mixing the sample applied in region A and the test material applied in region B at the rate of 50:50.
- D- The test material was made as its topical application is 0.3 ml to the intracapsular region.

One pair of 0.1 ml injection was made to the left and right regions of the animal during the applications in A, B and C regions.

In the region D, it was applied to the left region (D₁) in the 7^{th} day and the right topical region (D₂) in the 14^{th} day.

Descriptions:
Signatories:

53006

53192

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited.

The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).

Page 4/6



Report no: 16563500-125.05-64 / 3882-2

Negative Control

Negative control was comparatively carried out in 2 different regions in 2 different applications (Figure 2).




Figure 2.

- 1- Head of experimental animal.
- B- Serum Physiological 0.1 ml.
- C- Freund's Complete Anjuvant (FCA) and serum physiological solution at the rate of 50:50 was mixed and applied.
- D-0.3 ml serum physiological was applied to the topical regions.

The test materials were applied to experimental animals as shown in Figure 2 one day after they were shaved to provide an application field and to control animals as shown in Figure 3. All applications were made as 0.1 ml to subcutaneous. The regions were not closed in anyway after application. In the topical application, the test material for experimental animals and 0.3 ml serum physiological for control animals were applied to the skin, the application regions were bandaged with gauze bandage and all application regions were wrapped with elastic bandage after application. The gauze bandages were contacted with the regions for 48 hours. At the end of application duration, the bandages were removed and the reactions on the skin were noted. The second topical application was carried out after 7 days and the same experimental processes were followed.

Descriptions: Signatories: 53006 53192

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited. The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Institute archive).

Page 5/6



Report no: 16563500-125.05- 64 / 3882-2

Applied test material

10 animals were test material and 5 animals for control were used in the application. Total 15 animals were used in this test since there was only one test material.

Test Material: Surgical Drape

Reaction	
No visible change	
Clear or patched rash	
Moderate or confluent rash	
Intensive rash or pocks	

Table 2. Evaluation criteria and scoring.

Evaluation Average

Result
0.4
0,1

Table 3. Average score values.

Result

As stated in the test carried out for test and control samples, the observations were scored by taking into account the evaluation and scoring criteria given in the Table 2. The rash was seen on the skin of animals in the group to which 'Surgical Drape' extract was applied as a result of evaluation. The sensitization score was found as 0.4 (Table 3) as a result of observations. A significant weight loss and visible negative effects in general health condition were not observed in experimental animals. According to the results obtained, it was determined that the material tested did not have any sensitive (sensitive to matter) effect according to the evaluation criteria and protocol stated in ISO 10993-10:2010 document

document.	
Descriptions:	
Signatories:	
53006	53192
Junt 28	

This report and its results may not be partly or wholly replicated or published for commercial or advertisement purposes. Also, any firm, except Public Enterprises, may not use this report in any legal transactions.

The analysis with (*) mark in the report were accredited.

The unsigned analysis reports are not valid.

This report is consisted of 6 pages and prepared as 2 original copies (1 original for customer and 1 original for Page 6/6 Institute archive).