



Hospital Care & OPM

MEA Product Summary

Discofix® C - The stopcock and manifold that resist stress cracking



- Allows convenient access for injections and aspirations by rotating adapter
- Lets you feel the proper setting by tactile response of every 45° turn
- Resists the hazardous effects of pharmaceuticals
- Maintains a leakproof connection even during long-term application
- Prevents leakage, contamination and air embolism











Stopcocks Drug Resistant	Color	Rotating adapter	Proximal Stopcock	Extension (cm)	Sales Unit/ pcs.	Ref. no.	PVC-free	DEHP-, Latex-free	
Discofix® C 3-way stopcocks					100	16494C			
					100	16495C			
					100	16496C			
					100	16497C			
					10	50	16500C	—	
					10	50	16501C	—	
					25	50	16520C	—	
					50	50	16540C	—	
					75	50	16551C	—	
				100	50	16560C	—		
Discofix® C 3-way stopcocks with Safeflow					100	16494CSF			
					100	16495CSF			
				10	50	16500CSF	—		
				10	50	16501CSF	—		
				25	50	16520CSF	—		
				50	50	16540CSF	—		
				75	50	16551CSF	—		
				100	50	16560CSF	—		
Discofix® C 3-gang manifolds					50	16600C			
					50	16605C			
				180	50	16611C	—		
				150	50	16610C	—		
				150	50	16615C	—		
Discofix® C 5-gang manifolds					40	16608C			
					40	16609C			
				150	40	16612C	—		
				150	40	16613C	—		
				180	40	16614C	—		

Vasofix® Safety – Shielded IV Catheter with injection port

Protects you and your healthcare team against needlestick injuries and related infections.



- Complete portfolio for every vein condition, puncture site and patient age in two material types (FEP + PUR)
- Easy puncture at a wide range of angles because of optimized back cut bevel
- Minimum effort of catheter insertion
- Self-activating clip-technology => covers needle tip automatically after use
- Hygienic injection of medication via injection port without an extra needle or risk of needlestick injury
- PVC-free, latex-free

Color	Size		Vasofix® Safety		Stylet 
	G	Ø mm	FEP	PUR (Certon®)	
 yellow	24 x 3/4"	0.7 x 19	–	4269071S-01	–
 blue	22 x 1"	0.9 x 25	4268091S-01	4269098S-01	4215095
 pink/white	20 x 1"	1.1 x 25	–	4269217S-01	4215117
 pink	20 x 1 1/4"	1.1 x 33	4268113S-01	4269110S-01	4219104
 green/white	18 x 1 1/4"	1.3 x 33	4268334S-01	4269330S-01	4219139
 green	18 x 1 3/4"	1.3 x 45	4268130S-01	4269136S-01	4219120
 white	17 x 1 3/4"	1.5 x 45	4268156S-01	4269152S-01	4215150
 grey	16 x 2"	1.7 x 50	4268172S-01	4269179S-01	4219171
 orange	14 x 2"	2.2 x 50	4268210S-01	4269225S-01	4219201

Original Perfusor®-Tubings – For use with syringe pumps



- Made of PVC or PE
- Minimized priming volumes
- Luer-Lock fitting
- DEHP-free

Original accessories – for system safety

Article	Priming Volume	Inner Ø (mm)	Length cm	Sales Unit/ pcs.	Ref. no.	Latex free	PVC free
Original Perfusor®-Tubing	0.88	1.5	50	100	8255172	●	–
Original Perfusor®-Tubing	2.65	1.5	150	100	8722960	●	–
Original Perfusor®-Tubing	3.53	1.5	200	100	8722862	●	–
Original Perfusor®-Tubing	4.41	1.5	250	100	8255490	●	–
Original Perfusor®-Tubing	5.3	1.5	300	100	8255253	●	–
Original Perfusor®-Tubing with Safsite®	2.65	1.5	150	100	8722820	●	–
Original Perfusor®-Tubing 0.2 µm disc-filter	3.53	1.5	200	100	8723001	●	–
Original Perfusor®-Tubing PCA	2.94	1.5	150	50	8726019	●	–
Original Perfusor®-Tubing	0.5	0.9	75	100	8722870N	●	–
Original Perfusor®-Tubing	1.0	0.9	150	100	8255504N	●	–
Original Perfusor®-Tubing	1.3	0.9	200	100	874591N	●	–
Original Perfusor®-Tubing	0.39	1.0	50	100	8255059	●	●
Original Perfusor®-Tubing	0.79	1.0	100	100	8255067	●	●
Original Perfusor®-Tubing	1.18	1.0	150	100	8722935	●	●
Original Perfusor®-Tubing	1.57	1.0	200	100	8723060	●	●
Original Perfusor®-Tubing	1.96	1.0	250	100	8272565	●	●
Original Perfusor®-Tubing, orange	1.18	1.0	150	100	8723017	●	●
Original Perfusor®-Tubing, yellow	1.57	1.0	200	100	4185838	●	●
Original Perfusor®-Tubing with slide clamp	1.0	0.9	150	100	8722940	●	–
Original Perfusor®-Tubing with slide clamp	1.3	0.9	200	100	8722941	●	–

B. Braun Melsungen AG | 34209 Melsungen | Germany
Tel +49(0) 5661 71-0 | www.bbraun.com