


ADJUSTMENT OF DRAWING

BLOOD TUBING SET.

Form No. BL-PI-020-02



REMARK
GAMMA STERILE (CAUTION)
4.3x6.8 (65)
3.5x5.5 (65)
2.6x4.1 (65)

Change symbol "7" to "-".
A290R/V801R → A290R-V801R

Packing Method	
Item	GAMMA Sterilization
Individual Package	CE GAMMA Package
Outer Box	CE GAMMA
Attach sticker	Attach sticker barcode on outer box
IFU	CE GAMMA
Packing Q'ty	24 Sets.
Special Labelling M/C Name	FOR B. BRAUN DIALOG DN
PRIMING VOLUME A290R-V801R	151 ML

ESTABLISHED DATE	DRAWING	CHECK	APPROVE
10 Jun.'11	Ms. Yanisa	5 Apr. 17	24 Apr. 2017 SEE BACK SIDE FOR APPROVED BY NIPRO

RECORD OF CHANGE				
REVISED DATE	REVISION	APPR.NO.	DATE OF RECEIPT	FILE NO.
05 Apr.'14	5	DC-BL-14D-008		
09 Jun.'16	6	DC-BL-16F-009		REGISTERED
31 Mar.'17	7	DC-BL-17D-011		