

# **CRDLIGHT** **PROTECTIVE EQUIPMENT**

## **PRODUCT CATALOG**

Sales Manager : \_\_\_\_\_

Phone Number : \_\_\_\_\_

E-mail Address: \_\_\_\_\_

**Jiangmen Yanyang Trading Co., Ltd.**

# Company Introduction

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.



Jiangmen Yanyang Trading Co., Ltd. was established in 2008. The company is headquartered in Jiangmen City, Guangdong Province. We are a modern hightech enterprise integrating sales, testing and service.

Committed to the benchmark of the Chinese nursing brand industry chain, in addition to a number of new technologies and perfect quality management systems leading the international industry standards, the company is also equipped with a technical management team and quality inspection team of 50 people, and has 10 automatic production lines, not only ensures the output of products, but also guarantees the high quality standards of products

**营业执照**

(副本) (副本号:1-1)

统一社会信用代码  
91440704MA52768M32

名称 江门岩羊工贸有限公司  
类型 有限责任公司(自然人投资或控股)  
法定代表人 胡杰  
经营范围 研发、生产、加工、销售: LED、驱动电源、开关电源、电子整流器、智能化设备系列产品及配件、电器产品附件、低压电器产品、照明电器附件、电子产品、电气元器件、家用电器、五金、塑料制品、医疗器械、消毒产品(不含危险化学品); 医疗器械技术研发; 生产、销售; 机械设备、橡胶制品、日用品、口罩; 货物或技术进出口(国家禁止或涉及行政审批的货物和技术进出口除外)。(依法须经批准的项目, 经相关部门批准后方可开展经营活动。)

注册资本 人民币伍拾万元  
成立日期 2018年08月30日  
营业期限 长期  
住所 江门市江海区信义路18号2幢4层自编1号(信息申报制、一址多照)(经营场所: 江门市江海区信义路18号5幢5层自编A04(信息申报制))(一照多址)

登记机关 江门市江海区市场监督管理局  
2020年5月12日

扫描二维码登录“国家企业信用信息公示系统”了解更多登记、备案、许可、监管信息。

国家企业信用信息公示系统网址: <http://www.gsxt.gov.cn>


市场主体应当于每年1月1日至6月30日通过国家企业信用信息公示系统报送公示年度报告

国家市场监督管理总局监制



# Production Factory

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.



**Daily output: 500,000PCS**





# Testing Laboratory

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.



**Gas Exchange  
Pressure Tester**



**Resistance Tester**



**Ventilation  
Resistance & Pressure  
Difference Tester**



**Flame Retardant  
Tester**



**Synthetic Blood  
Penetration Tester**



**Surface Moisture  
Resistance Tester**



**Particulate Filtration  
Efficiency Tester**



**PFE Tester**



**Electronic Tensile Tester**

## Filtering Half Mask

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.


# YY0525

## FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



# CE 2163



**EN 149:2001 + A1:2009**  
**FFP2 NR**  
**Standard**



Breathable



Skin friendly



5 Layers



## Filtering Half Mask

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

# YY0525

## FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



### YY0525

25pcs/color box; color box size: 12.5x9.5x12.5cm;

40boxes/carton; carton size: 51.5x49x26.5cm;

20GP: 352ctns;

40GP: 768ctns;

40HQ: 960ctns.



# Filtering Half Mask

CRDLIGHT  
Jiangmen Yanyang Trading Co., Ltd.

# YY0525 FILTERING HALF MASK



## EU TYPE EXAMINATION CERTIFICATE

Certificate No: 2163-PPE-834

Respiratory protective devices, filtering half masks to protect against particles manufactured by

**JIANGMEN YANYANG TRADING CO., LTD.**

No.1, 4th Floor, Building 2, No. 18 Xinyi Road, Jianghai District, Jiangmen City, Guangdong Province, China

are tested and evaluated according to

### EN 149:2001 + A1:2009 Respiratory Protective Devices - Filtering Half Masks to Protect Against Particles - Requirements, Testing, Marking

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

#### Product Definition

Brand Name: CRDLIGHT Model: YY0525

Filtering half mask

Classification: FFP2 NR

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with:

- Issuing an appropriate EU Declaration of Conformity according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.
- Ongoing successful performance in fulfillment of the requirements set out in Personal Protective Equipment Regulation (EU) 2016/425 and harmonised standards, ensured by assessments based on Annex 7 (Module C2) or Annex 8 (Module D) of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on 25/06/2020 and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.



Suat KACMAZ  
UNIVERSAL CERTIFICATION  
Director



#### Oil: Paraffin Oil Testing

No. of Sample	Paraffin Oil Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result
-	3.0		
-	3.4		
-	3.3	FFP1 ≤ 20 %	
-	3.5		
-	3.5	FFP2 ≤ 6 %	
-	3.4		
-	4.7	FFP3 ≤ 1 %	
-	4.1		
-	4.5		

anical Strength  
ature Conditioning  
ceived, original  
ated wearing treatment

Practical Performance report, the likelihood of mask materials in contact with the skin causing irritation or other not reported.

of ple	Visual inspection	Requirements in accordance with EN 149:2001 + A1:2009	Result
	Burn for 1s	Filtering half mask shall not burn or not	Passed
	Burn for 1s	continue to burn for more than 5 s after removal from the flame	
	Burn for 2s		Filtering half masks fulfill requirements of the standard

ceived, original  
ature Conditioning  
be inhalation air:

CO <sub>2</sub> content of the inhalation air [%] by volume	An average CO <sub>2</sub> content of the inhalation air	Requirements in accordance with EN 149:2001 + A1:2009	Result
0.2			
0.3			
0.3	0.3 [%]	CO <sub>2</sub> content of the inhalation air shall not exceed an average of 1.0% by volume	Passed

ceived, original

Performance and TIL test reports no adverse effects have been reported for donning and remove of the mask also the s that the ear loops / head harness are capable of holding the mask firmly enough.

Performance report, no adverse effects were reported for the field of vision availability when the mask is worn.

odel under inspection have no valves.

lation

e figures gathered for 9 different samples 3 as received, 3 with temperature conditioning and 3 simulated wearing ies with the limits given in the standard for FFP1, FFP2 and FFP3 classes. This is valid for inhalation results for 30 iose at 160 L/min.





**Filtering Half Mask**

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

# YYC1028

## COLORED FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON




**CE 2163**

**EN 149:2001 + A1:2009**

**FFP2 NR**

**Standard**



Breathable



Skin friendly



5 Layers



## Filtering Half Mask

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

# YYC1028

## COLORED FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



### YYC1028

25pcs/color box; color box size: 12.5x9.5x12.5cm;

40boxes/carton; carton size: 51.5x49x26.5cm;

20GP: 352ctns;

40GP: 768ctns;

40HQ: 960ctns.



# Filtering Half Mask

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

# YYC1028

## COLORED FILTERING HALF MASK



UNIVERSAL



### EU TYPE EXAMINATION CERTIFICATE

Certificate No: 2163-PPE-1732

Respiratory protective devices, filtering half masks to protect against particles manufactured by

**JIANGMEN YANYANG TRADING CO., LTD**

No. 1, 4th Floor, Building 2, No.18 Xinyi Road, Jianghai District, Jiangmen City,  
Guangdong Province, China

are tested and evaluated according to

#### EN 149:2001 + A1:2009 Respiratory Protective Devices - Filtering Half Masks to Protect Against Particles - Requirements, Testing, Marking

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

#### Product Definition

Single shift use particle filtering half mask for protection against solid and liquid aerosols, is a folding type, 5 layers cotton and polypropylene fabrics, without valve, fitted with ear loops, with internal nose clip and inside sponge strip.

Brand Name: CRDLIGHT Model: YYC1028 Classification: FFP2 NR  
Model have white, black, grey, pink, and blue versions

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with:

- Issuing an appropriate EU Declaration of Conformity according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.
- Ongoing successful performance in fulfilment of the requirements set out in Personal Protective Equipment Regulation (EU) 2016/425 and harmonised standards, ensured by assessments based on Annex 7 (Module C2) or Annex 8 (Module D) of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on 04/12/2020 and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.



*Suat KACMAZ*  
UNIVERSAL CERTIFICATION  
Director



#### at: Sodium Chloride Testing

No. of sample	Sodium Chloride Testing: 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result
-	0,1	FFP1 ≤ 20 %	Filtering half masks fulfill the requirements of the standard EN 149:2001 + A1:2009 given in 7.9.2 in range of the first, and second protection classes. FFP1, FFP2, FFP3
-	0,2		
-	0,1		
-	0,2	FFP2 ≤ 6 %	
-	0,1		
-	0,2	FFP3 ≤ 1 %	
-	0,6		
-	0,5		
-	0,6		

Initial Strength 95 L/min = 1,6 dm³/min

Initial Strength  
ature Conditioning  
ceived, original  
ated wearing treatment  
it: Paraffin Oil Testing

No. of Sample	Paraffin Oil Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result
-	0,2	FFP1 ≤ 20 %	Filtering half masks fulfill the requirements of the standard EN 149:2001 + A1:2009 given in 7.9.2 in range of the first and second protection classes. FFP1, FFP2, FFP3
-	0,3		
-	0,2		
-	0,3	FFP2 ≤ 6 %	
-	0,2		
-	0,7	FFP3 ≤ 1 %	
-	0,6		
-	0,8		

Initial Strength  
ature Conditioning  
ceived, original  
ated wearing treatment

Practical Performance report, the likelihood of mask materials in contact with the skin causing irritation or other not reported.

No. of Sample	Visual inspection	Requirements in accordance with EN 149:2001 + A1:2009	Result
-	0,4s	Filtering half mask shall not burn or not continue to burn for more than 5 s after removal from the flame	Passed
-	0,4s		
-	0,5s		
-	0,4s		

ceived, original  
ature Conditioning  
be inhalation air:

CO <sub>2</sub> content of the inhalation air [%] by volume	As average CO <sub>2</sub> content of the inhalation air	Requirements in accordance with EN 149:2001 + A1:2009	Result
0,6827	0,68	CO <sub>2</sub> content of the inhalation air shall not exceed an average of 1,0% by volume	Passed
0,6813			
0,6810			

ceived, original

Performance and TIL test reports no adverse effects have been reported for donning and remove of the mask also the that the head harnesses are capable of holding the mask firmly enough.

7.14	FIELD OF VISION: In accordance with EN 149:2001 + A1:2009, no adverse effects were reported for the field of vision availability when the mask is worn.
Article 7.16	<p><b>Breathing Resistance: Inhalation</b></p> <p>The overall evaluation in the figures gathered for 9 different samples 3 as received, 3 with temperature conditioning and 3 Simulated wearing conditioned, complies with the limits given in the standard for FFP1, FFP2 and FFP3 classes. This is valid for inhalation results for 30 L/min, 95 L/min and exhalation at 160 L/min.</p> <p>Passed.</p>

UFR-383 12.12.2018 Rev.01





## Filtering Half Mask

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

# YY1029

## FFP3 FILTERING HALF MASK


Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



# FFP3

**99% filtration efficiency**

**CE 2163**



**EN 149:2001 + A1:2009**

# FFP3 NR

**Standard**



Breathable



Skin friendly



5 Layers



## Filtering Half Mask

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

# YY1029

## FFP3 FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



### YY1029

25pcs/color box; color box size: 16x11x12.5cm;  
40boxes/carton; carton size: 56.5x51.5x33.5cm;  
20GP: 240ctns;  
40GP: 480ctns;  
40HQ: 720ctns.



# Glove Introductions

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

**COLOR:** **CLEAR** / **BLACK** / **BLUE**



**SIZES:** **S** / **M** / **L** / **XL**

## Suggestion S Size

Palm of your hand is 7.1-8cm wide

## Suggestion M Size

The palm of your hand is 8.1-9cm wide


## Suggestion L Size

The palm of your hand is 9.1-10cm wide

## Suggestion XL Size

The palm of your hand is 10.1-11cm wide

Most women use size S and those with larger hands use size M  
Most men use size M and those with larger hands use size L





# Manufacturing

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.





## Vinyl Gloves

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

# YYPC01

## PVC EXAMINATION GLOVE (Clear)

# YYPC02

## PVC EXAMINATION GLOVE (Blue)

SIZE: S / M / L / XL

- NON-STERILE
- POWDER-FREE
- AMBIDEXTROUS
- SINGLE USE
- LATEX-FREE

### YYPC01/02

100pcs/color box (23x12.5x6.5cm)

10boxes/carton (35x27x25cm)

20GP: 1,024ctns;

40GP: 2,112ctns;

40HQ: 2,640ctns.





## Nitrile Gloves

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

# YYNG01

## NITRILE EXAMINATION GLOVE (Black)



# YYNG02

## NITRILE EXAMINATION GLOVE (Blue)

SIZE: S / M / L / XL



- NON-STERILE
- POWDER-FREE
- AMBIDEXTROUS
- SINGLE USE
- LATEX-FREE

### YYNG01/02

100pcs/color box (23x12.5x6.5cm)

10boxes/carton (35x27x25cm)

20GP: 1,024ctns;

40GP: 2,112ctns;

40HQ: 2,640ctns.





## Blended Gloves

**CRDLIGHT**  
Jiangmen Yanyang Trading Co., Ltd.

# YYBG01

## PVC/NITRILE BLENDED GLOVE (Black)



# YYBG02

## PVC/NITRILE BLENDED GLOVE (Black)

SIZE: S / M / L / XL



- NON-STERILE
- POWDER-FREE
- AMBIDEXTROUS
- SINGLE USE
- LATEX-FREE

### YYBG01/02

100pcs/color box (23x12.5x6.5cm)

10boxes/carton (35x27x25cm)

20GP: 1,024ctns;

40GP: 2,112ctns;

40HQ: 2,640ctns.

