

Radifocus® Guidewire M Stiff Type



from ACCESS
to CLOSURE




TERUMO
INTERVENTIONAL
SYSTEMS

Radifocus® Guidewire M Stiff Type

Guide Wire

Indicated for use for standard angiography in coronary and peripheral territories, catheterization, embolization, peripheral angioplasty and stenting (with or without an exchange wire) in tight, stenotic, tortuous and narrow vessels.



More on terumo-europe.com

TERUMO INTERVENTIONAL
SYSTEMS

Characteristics

Characteristics

- ◆ Reinforced superelastic Nitinol core:
 - ◆ Good torque control and stiffness
 - ◆ Adaptable shape memory and great flexibility for controlled navigation
 - ◆ Kink resistance for swift device exchange
- ◆ Radiopaque polyurethane jacket for good fluoroscopic visibility and controlled navigation
- ◆ Highly flexible, tapered tip for atraumatic and controlled navigation through the vasculature
- ◆ Hydrophilic "M" polymer coating for frictionless endovascular steerability and trackability
- ◆ One-piece construction with one-to-one torque transmission for predictable navigation
- ◆ Rounded, atraumatic tip for the prevention of vessel trauma and low friction guidewire insertion

Radifocus® Guidewire M Stiff Type
Guide Wire

 **TERUMO** INTERVENTIONAL
SYSTEMS

Specifications

Specifications

General specifications

Core Material	Stiff Tapered Nitinol
Hydrophilic Coating	"M" Polymer
Radiopaque Jacket	Polyurethane Layer Containing Tungsten

Item specifications

Guidewire Diameter - Outer Diameter	Length	Flexible Length	Distal Curve	Code
0.018 in 0.46 mm	180 cm	30 mm	Angled	RF*PA18183M
0.02 in 0.51 mm	150 cm	30 mm	Angled	RF-PA18153M
0.02 in 0.51 mm	260 cm	30 mm	Angled	RF-PA18263M
0.020 in 0.51 mm	260 cm	30 mm	Angled	RF*PA18263M
0.025 in 0.64 mm	150 cm	30 mm	Angled	RF-PA25153M
0.025 in 0.64 mm	150 cm	30 mm	Angled	RF*PA25153M
0.025 in 0.64 mm	150 cm	30 mm	Straight	RF-PS25153M
0.025 in 0.64 mm	180 cm	30 mm	Angled	RF-PA25183M
0.025 in 0.64 mm	260 cm	30 mm	Angled	RF*PA25263M
0.025 in 0.64 mm	260 cm	30 mm	Angled	RF-PA25263M
0.035 in 0.89 mm	80 cm	30 mm	Angled	RF*PA35083M
0.035 in 0.89 mm	80 cm	30 mm	Angled	RF-PA35083M
0.035 in 0.89 mm	80 cm	30 mm	Straight	RF-PS35083M
0.035 in 0.89 mm	150 cm	10 mm	Angled	RF-PA35151M

Continued on next page

Item specifications (Continued)

Guidewire Diameter - Outer Diameter	Length	Flexible Length	Core Material	Code
0.035 in 0.89 mm	150 cm	30 mm	Angled	RF*PA35153M
0.035 in 0.89 mm	180 cm	30 mm	Angled	RF*PA35183M
0.035 in 0.89 mm	180 cm	30 mm	J Angled (1.5 mm)	RF*HR35183M
0.035 in 0.89 mm	180 cm	30 mm	J Angled (2.0 mm)	RF*HN35183M
0.035 in 0.89 mm	180 cm	30 mm	Straight	RF-PS35183M
0.035 in 0.89 mm	180 cm	80 mm	Angled	RF*PA35188M
0.035 in 0.89 mm	260 cm	30 mm	Angled	RF*PA35263M
0.035 in 0.89 mm	260 cm	30 mm	Angled	RF-PA35263M
0.035 in 0.89 mm	260 cm	30 mm	Straight	RF-PS35263M
0.035 in 0.89 mm	300 cm	30 mm	Angled	RF*PA35303M
0.035 in 0.89 mm	300 cm	30 mm	Angled	RF-PA35303M
0.035 in 0.89 mm	300 cm	30 mm	J Angled (1.5 mm)	RF*HR35303M
0.038 in 0.97 mm	80 cm	30 mm	Straight	RF-PS38083M
0.038 in 0.97 mm	150 cm	30 mm	Angled	RF*PA38153M

Continued on next page

Radifocus® Guidewire M Stiff Type
Guide Wire

TERUMO INTERVENTIONAL
SYSTEMS

Item specifications (Continued)

Guidewire Diameter - Outer Diameter	Length	Flexible Length	Core Material	Code
0.038 in 0.97 mm	150 cm	30 mm	Angled	RF-PA38153M
0.038 in 0.97 mm	150 cm	30 mm	Straight	RF-PS38153M
0.038 in 0.97 mm	180 cm	30 mm	Angled	RF*PA38183M
0.038 in 0.97 mm	180 cm	30 mm	Angled	RF-PA38183M
0.038 in 0.97 mm	180 cm	30 mm	Straight	RF-PS38183M
0.038 in 0.97 mm	260 cm	30 mm	Angled	RF-PA38263M

Continued on next page

Radifocus® Guidewire M Stiff Type
Guide Wire

 **TERUMO** INTERVENTIONAL
SYSTEMS