

1. Filtru cu filtrare eficientă și cu returnarea optimală a umidității și căldurii




BACTERIAL-VIRAL HME FILTERS

- Transparent
- Low dead space
- Gas sampling port available
- Straight and angled choices
- Low resistance
- Latex-free
- Inspected and approved by third party experts

[Download PDF](#)

Bacterial-Viral HME

Filter with paper HME



Code	40820
Sterile Code	40820S
Qty/Box	45
Tidal Volume(ml)	150-1500
Dead Space(ml)	53
Bacterial-Viral Efficiency	>99,999%
Viral Efficiency	>99,9%
Resistance to Flow	30 lt/min 10 mm H ₂ O 60 lt/min 22 mm H ₂ O 90 lt/min 36 mm H ₂ O
Weight (gr)	30
Fittings	22mmM-15/22 mmMF
Humidification efficiency	37,3 mg/l H ₂ O (500 ml tid. vol)


2.Filtru p/u utilizare în anestezie sau TI



BACTERIAL-VIRAL FILTERS

- Transparent, easy to see the surface
- High efficiency with low resistance electrostatic filters
- Effective to prevent microbial cross contamination
- Remarkably efficient heat and moisture exchangers
- 15 mm/22 mm connectors conforming to ISO Standards
- Retainable cap
- Inspected and approved by third party experts

[Download PDF](#)

 Bacterial-Viral Filter	Code	40920
	Sterile Code	40920S
	Qty/Box	50
	Tidal Volume(ml)	150-1500
	Dead Space(ml)	33
	Bacterial-Viral Efficiency	>99,999%
	Viral Efficiency	>99,9%
	Resistance to Flow	30 lt/min 8 mm H ² O 60 lt/min 17 mm H ² O 90 lt/min 27 mm H ² O
	Weight (gr)	19
	Fittings	22mmM-15/22 mmMF